

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™


## Main index

Page 1 / 1 - Index: 0.0.0


[Overview of all documents](#)

No	Subject	Rev. date	Document
<b>1</b>	<b>Detection of beta lactamases</b>		
	Front page and preface	22-03-2010	1.0.1
	Range of Neo-Sensitabs™ and Diatabs™ for detection of resistance mechanisms	22-03-2010	1.0.2
	Extended spectrum beta lactamases (ESBL), ESBL+AmpC Screen kit and ESBL Confirm ID kit	22-03-2010	1.1.0
	ESBL Quality Control		1.1.1
	Inducible cephalosporinases or AmpC beta-lactamases		1.2.0
	Testing/Reporting		1.2.1
	Plasmid mediated AmpC beta lactamases. AmpC Confirm ID kit	22-03-2010	1.3.0
	Inhibitor resistant TEM beta-lactamases (IRT)		1.4.0
	Carbapenemases. KPC+ MBL Confirm ID kit	22-03-2010	1.5.0
	Detection of susceptible beta lactamases in one strain	22-03-2010	1.6.0
	Detection of beta-lactam resistance phenotypes in Enterobacteriaceae	22-03-2010	1.7.0
	Detection of beta-lactam resistance phenotypes in non-fermenters		1.8.0
<b>2</b>	<b>Screening of 16S + rRNA methylases</b>		
	Screening of 16S rRNA methylases (HLR to Aminoglycosides)	08-12-2009	2.1.0
<b>3</b>	<b>Screening of plasmid mediated quinolone resistance</b>		
	Screening of plasmid-mediated quinolone resistance		3.1.0
<b>4</b>	<b>Detection of hypermutable strains</b>		
	Detection of hypermutable strains		4.1.0
<b>5</b>	<b>Detection of heteroresistant strains</b>		
	Detection of heteroresistant strains		5.1.0
<b>6</b>	<b>Efflux pump inhibitors and NMP Diatabs</b>		
	Efflux pump inhibitors and NMP Diatabs	22-03-2010	6.1.0
<b>7</b>	<b>Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar</b>		
	Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar	22-03-2010	7.1.0
	Detection of VISA, GISA and hVISA using the Neo-Sensitabs prediffusion method	22-03-2010	7.2.0
	Detection of vancomycin resistant enterococci (VRE)	22-03-2010	7.3.0
	Detection of staphylococci and enterococci non-susceptible to daptomycin using the Neo-Sensitabs prediffusion method	22-03-2010	7.4.0
	Detection of colistin resistant P.aeruginosa, A. baumannii and Enterobacteriaceae	22-03-2010	7.5.0
<b>8</b>	<b>Resistance mechanisms and automated methods</b>		
	Resistance mechanisms and automated methods	04-12-2009	8.1.0
<b>9</b>	<b>Detection of resistance mechanisms</b>		
	Detection of resistance mechanisms (general)	04-12-2009	9.1.0
<b>10</b>	<b>Intrinsic (natural) resistance</b>		
	Intrinsic (natural) resistance	22-03-2010	10.1.0

# Detection of resistance mechanisms



E.coli AmpC positive



K. pneumoniae KPC positive

**using Neo-Sensitabs™  
and Diatabs™**

**2010**

## **PREFACE**

### **Detection of Resistance Mechanisms 2010**

Laboratories must be able to identify resistant microorganisms, trends in resistance, emerging resistance patterns among clinically relevant bacteria. In order to recognize resistance phenotypes, a suitable number of antimicrobial agents must be used in susceptibility testing, as well as specific agents (indicator drugs) implying mechanisms conferring resistances that may be less obvious in current tests.

With this approach, it would be possible for the microbiologist to:

- a) detect mechanisms of resistance, including low level expression
- b) modify clinical classifications that are inconsistent with the inferred resistance mechanism and
- c) inference of susceptibility results for antimicrobials that are not included in the antibiogram.

The Users Guide for Detection of Resistance Mechanisms includes a series of newly developed tests, particularly useful for the detection of beta-lactamases : ESBL, AmpC, Carbapenemases as well as many other enzymes such as 16S rRNA methylases, Plasmid mediated quinolone resistance, and has been written by J.B.Casals on behalf of Rosco Diagnostica.

Besides, the booklet contains information on the technique used to detect HVISA, VISA strains, using Vancomycin and Teicoplanin Neo-sensitabs by the prediffusion method. Technique for detecting Daptomycin non-susceptible staphylococci and enterococci using Daptomycin Neo-Sensitabs and the prediffusion method.

As well as the detection of colistin resistant *P. aeruginosa*, *A. baumannii* and Enterobacteriaceae using Colistin 10 ug Neo-Sensitabs by the prediffusion method.

Tests for detecting Hypermutable strains and Heteroresistant strains, as well as Efflux pump inhibitors are described too.

Finally, a chapter on resistance mechanisms and automated methods, could be interesting lecture for laboratories relying only on automated methods. It is obvious, that these laboratories will have good help in using the tests mentioned in this booklet to complement their results.

The Users Guide is available at our website [www.rosco.dk](http://www.rosco.dk) and updated information will continuously be included.

ROSCO DIAGNOSTICA is welcoming any feedback and questions from users directly ([info@rosco.dk](mailto:info@rosco.dk)) or through our representatives.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 1 / 2 - Document: 1.0.2

## Detection of beta lactamases

Range of Neo-Sensitabs™ and Diatabs™ for detection of resistance mechanisms

### Neo-Sensitabs™/Diatabs™ Range

Beta-lactamases	Code		Code
<b>ESBL screening</b> Cefpodoxime 10 µg Cefpodoxime+Clavulanate 10+1 µg	CPD10 CPD+C		
<b>ESBL</b> Cefotaxime 5 µg Ceftazidime 10 µg Cefotaxime+Clavulanate 5+1 µg Ceftazidime+Clavulanate 10+1 µg	CTX.5 CAZ10 CT5+C CZ10C	<b>ESBL Confirm ID kit</b> Cefotaxime 30 µg Cefotaxime+Clavulanate Ceftazidime 30 µg Ceftazidime+Clavulanate Cefepime 30 µg Cefepime+Clavulanate	CTX30 CTX+C CAZ30 CAZ+C FEP30 FEP+C
<b>AmpC</b> Cefotaxime 30 µg Cefotaxime+Boronic acid Cefotaxime+Cloxacillin Ceftazidime 30 µg Ceftazidime+Boronic acid Ceftazidime+Cloxacillin Boronic acid Diatabs Cloxacillin Diatabs	CTX30 CTXBO CTXCX CAZ30 CAZBO CAZCX BORON CLOXA	<b>AmpC Confirm ID kit</b> Cefotaxime 30 µg Cefotaxime+Cloxacillin Ceftazidime 30 µg Ceftazidime+Cloxacillin	CTX30 CTXCX CAZ30 CAZCX
		<b>ESBL+AmpC Screen kit</b> Cefotaxime 30 µg Cefotaxime+Clavulanate Cefotaxime+Cloxacillin Cefotaxime+Clavulanate+Cloxacillin	CTX30 CTX+C CTXCX CTXCC
<b>Carbapenemases</b> Meropenem 10 µg Meropenem + Boronic Meropenem + Cloxacillin Meropenem+ Dipicolinic Dipicolinic acid Diatabs Cloxacillin Diatabs Boronic and Diatabs	MRP10 MRPBO MRPCX MRPDP D.P.A CLOXA BORON	<b>KPC+MBL Confirm ID kit</b> Meropenem 10 µg Meropenem+Boronic Meropenem+Cloxacillin Meropenem+Dipicolinic	MRP10 MRPBO MRPCX MRPDP
<b>Metallo-β-lactamases (MBL)</b> Meropenem 10 µg Meropenem+Dipicolinic Imipenem 10 µg Imipenem+EDTA	MRP10 MRPDP IMI10 IMI0E		

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 2 / 2 - Document: 1.0.2

## Detection of beta lactamases

## Range of Neo-Sensitabs™ and Diatabs™ for detection of resistance mechanisms

Others	Code
<b><u>Screening 16S rRNA methylases</u></b> Amikacin 30 µg	AMI30 or AMIKA
<b><u>Screening of plasmid-mediated quinolone resistance</u></b> Nalidixic acid 30 µg Ciprofloxacin 1 µg Norfloxacin 10 µg Ceftazidime 30 µg Tobramycin 10 µg	NAL30 CIPR1 NORFX CAZ30 TOB10

	Code
<b><u>Detection of hypermutable strains</u></b> Fosfomycin Rifampicin 30 µg Ceftazidime 30 µg Ciprofloxacin 1 µg Imipenem 10 µg Meropenem 10 µg Tobramycin 10 µg	FO200 or FOSFO RIFAM CAZ30 CIPR1 IMI10 MRP10 TOB10
<b><u>Efflux pumps inhibitors</u></b> NMP Diatabs	NMP

Prediffusion method	Code
<b><u>hVISA/VISA</u></b> Vancomycin 30 µg Teicoplanin 30 µg	VAN30 TPN30
<b><u>Vancomycin resistant enterococci</u></b> Vancomycin 30 µg Teicoplanin 30 µg	VAN30 TPN30
<b><u>Daptomycin non-susceptible staphylococci/enterococci</u></b> Daptomycin (+Ca) 30 µg	DAPCa
<b><u>Colistin resistance in P.aeruginosa, A. baumannii, Enterobacteriaceae</u></b> Colistin 10 µg	CO.10

## Detection of beta lactamases

### Extended spectrum beta lactamases (ESBL), ESBL+AmpC Screen kit and ESBL Confirm ID kit

#### Extended-Spectrum Beta-Lactamases (ESBL)

##### Screening and Confirmatory Tests for Extended-Spectrum Beta-Lactamases (ESBL)

Transferable plasmid-mediated beta-lactamases that produce resistance towards third generation cephalosporins and monobactams (e.g. aztreonam) have been described in strains of *Klebsiella pneumoniae*, *K. oxytoca*, *E. coli* and other Enterobacteriaceae. These enzymes are classified as extended-spectrum beta-lactamases (ESBL) and they have been implicated in clinical resistance to monobactams and broad-spectrum cephalosporins such as ceftazidime (CAZ), cefotaxime (CTX), and ceftriaxone (CTR).


Some ESBLs confer high-level resistance to these beta-lactams and are easily detected as resistant (or intermediate) by disk (tablet) diffusion testing. But the ESBL may provide low levels of resistance (MIC 1-2 µg/ml) to monobactams and third generation cephalosporins that can be easily overlooked by routine susceptibility methods and current interpretative criteria (1). These latter isolates may not reach current CLSI breakpoints for resistance, yet can be clinically resistant to beta-lactam therapy (2).

Since some ESBLs are more active on CAZ, while others are more active on CTX, the choice of cephalosporins tested can also affect the ability of laboratories to detect resistant strains (3). Most ESBLs are inhibited by clavulanic acid, tazobactam or sulbactam and can be readily detected by the double-disk (tablet) synergy test (4).

#### Double disk (tablet) synergy test

Inoculate a Mueller-Hinton plate as for susceptibility testing and apply Ceftriaxone (CTR) Neo-Sensitabs, Cefotaxime (CTX) Neo-Sensitabs, Ceftazidime (CAZ) Neo-Sensitabs, Cefepime (FEP) Neo-Sensitabs and Aztreonam (AZT) Neo-Sensitabs **at approximately 20 mm** (30 mm from tablet center to tablet center) from a tablet containing Amoxycillin+Clavulanate Neo-Sensitabs (AMC) using a dispenser. Incubate overnight at 35 °C.

**Extension of the zone of inhibition** (synergism) towards the tablet containing AM+CL, indicates the presence of an extended spectrum beta-lactamase (ESBL).



**Plate 1.1.0-a.** *Klebsiella pneumoniae* (ATCC 700603) producing extended-spectrum beta-lactamases (ESBL). Note the synergy between cefotaxime Neo-Sensitabs (CTX), ceftazidim Neo-Sensitabs (CAZ) and Amoxycillin+Clavulanate Neo-Sensitabs (AMC).

Another possibility of screening for ESBL is the use of lower MIC break-points for ceftazidime and aztreonam. Livermore et al (5) showed that most ESBL producers were resistant to CAZ at 2 µg/ml and AZTRM at 1 µg/ml. The corresponding zones with Neo-Sensitabs using McFarland 0.5 inoculum are 24 mm (CAZ and CTR) and 27 mm (AZT). As a consequence, *Klebsiella* spp, *E. coli* and *Salmonella* spp. showing zones < 24 mm with Ceftazidime, Cefepime and/or Ceftriaxone Neo-Sensitabs and/or ≤ 27 mm with Aztreonam and/or Cefotaxime Neo-Sensitabs, should be suspected of ESBL production. The CLSI has adopted practically all the same MIC breakpoints.

Cefpodoxime 10 µg may also be used in the screening of ESBL. Zones < 17 mm should be suspected of strains with ESBL production (18). Recently, the CLSI changed their Cefpodoxime screening breakpoints for ESBL from ≥ 2 to ≥ 8 µg/ml (18).

## Detection of beta lactamases

### Extended spectrum beta lactamases (ESBL), ESBL+AmpC Screen kit and ESBL Confirm ID kit

In a study comparing several ESBL screening methods, Vercauteren et al. (10) found that the double tablet synergy test using Neo-Sensitabs detected 96.9 % of ESBL producers while the E-test ESBL Screen detected 81.2 %.

De Gheldre et al showed that synergism between ceftazidime and cefepime with clavulanate (Neo-Sensitabs) was very useful to detect ESBL in *Enterobacter aerogenes* from Belgian hospitals (13).

Rodriguez-Villalobos et al. (20) and Fluit et al. (21) showed that the double disk (Neo-Sensitabs) synergy test has a higher sensitivity for the detection of ESBL than all combination disks (Oxoid, E-test).

Florijn et al (16) conclude that the use of ceftazidime, ceftriaxone and amoxycillin+clavulanate as Neo-Sensitabs is a cheap and reliable method for detection of *E. coli*, *Klebsiella* spp. and *P. mirabilis* isolates suspected of carrying ESBL. It performs better in a routine setting than the E-test, which often yields a result that cannot be interpreted.

Pitout et al. (22) recommend the use of cefepime and piperacillin+tazobactam when testing against strains with high level expression of AmpC beta-lactamases (*E. cloacae*, *E. aerogenes*, *C. freundii*, *S. marcescens*). Synergism between cefepime and tazobactam indicates presence of an ESBL.

Synergism between Ticarcillin + Clavulanate and aztreonam/ceftazidime/cefepime permit the detection of ESBL producing strains of *Ps. aeruginosa* ( SFM 2001 ). These strains show currently no zone around Ceftazidime Neo-Sensitabs (14).

### Confirmatory Tests for ESBL (combination disks)

The CLSI (9) recommends the use of Ceftazidime in combination with Clavulanic acid: Ceftazidime+Clavulanate Neo-Sensitabs, as a phenotypic confirmatory test for the presence of ESBL. Perform the antibiogram using Mueller Hinton Agar and McFarland 0.5 inoculum. Test both Ceftazidime+Clavulanate, Cefepime+Clavulanate and Ceftazidime/ Cefepime Neo-Sensitabs.

An increase in zone diameter of  $\geq 5$  mm for the combination Ceftazidime+Clavulanate or Cefepime+Clavulanate compared to Ceftazidime/ Cefepime alone is confirmatory of the presence of an ESBL.

Steward et al (12) showed that synergism between cefepime and clavulanate (Cefepime + Clavulanate Neo-Sensitabs) is very useful to detect ESBL in *Klebsiella pneumoniae*, differentiating strains producing ESBL (synergy between cefepime and clavulanate) from strains producing Amp C or hyperproducers of beta-lactamase.

*Enterobacter*, *Serratia*, *Morganella morganii*, *Providencia*, *Citrobacter freundii* and *Pseudomonas aeruginosa* produce chromosomally encoded inducible **Amp C** beta-lactamase. High level expression of Amp C may prevent the recognition of ESBL. Cefepime is practically not affected by Amp C and consequently Cefepime Neo-Sensitabs should be included as an ESBL screening agent when testing *Enterobacter*, *Serratia* ect. Synergism between Amox-Clav and Cefepime will indicate ESBL production (11,12,13,17,19). Strains with Cefepime zones < 24 mm should be suspected of ESBL production.

Recently Schwaber et al (32) found that the Vitek 2, Advanced Expert System identified the ESBL phenotype in only 62.5 % isolates of *Enterobacter* spp. and erroneously reported cephalosporin susceptibility in 28 %.

Cefepime+Clavulanate (and cefepime) Neo-Sensitabs should be used in the confirmatory tests for ESBL, because they are effective in detecting ESBL in strains of *Klebsiella*, *E. coli* etc. that may produce Amp C or are hyperproducers of beta-lactamase (31).

## Detection of beta lactamases

### Extended spectrum beta lactamases (ESBL), ESBL+AmpC Screen kit and ESBL Confirm ID kit

ESBLs can be obscured by the chromosomal AmpC cephalosporinase in *P. aeruginosa* (30). Cloxacillin 500 µg or Boronic acid Diatabs can be used to inhibit AmpC, for example by prediffusing (1 hour) one of these compounds on the agar before inoculation and before adding the antibiotic tablets (Neo-Sensitabs), placed on the same spots.

With *Klebsiella oxytoca*, synergism between Amoxycillin+Clavulanate (AMC) and Aztreonam or Ceftriaxone but not with ceftazidime indicates the presence of hyperproduction of K-1 chromosomal beta-lactamase (but **negative** for ESBL). Strains producing ESBL show synergism between AMC and ceftazidime (Use Ceftazidime+Clavulanate).

The use of cefotaxime, ceftriaxone, cefepime, aztreonam with AMC may result in **false positive** results for ESBL in *Klebsiella oxytoca* (Vitek,41).

The emergence of ESBL in *Salmonellae* merits concern. They cause frequently neonatal meningitis in many developing countries and are often already resistant to ampicillin and chloramphenicol (7).

Karas et al (8) reports clinical failure due to ESBL, in spite of the organism being susceptible with disk diffusion and MIC test (CTX MIC 0.75 µg/ml). The double disk diffusion test indicated the presence of an ESBL, but the test was first performed when therapy with cefotaxime was stopped, due to treatment failure.

The laboratory report should indicate that ESBL-producing strains may be resistant clinically to all penicillins, cephalosporins and aztreonam (9).

For serious systemic infections, even if the isolate appears susceptible to Amoxycillin+Clavulanate, Ticarcillin+ Clavulanate or Piperacillin+Tazobactam, do not report it as susceptible, because resistant mutants may be selected during therapy.

For Q.C. use *Klebsiella pneumoniae* ATCC 700603: zone of Ceftazidime+Clavulanate and Cefepime+Clavulanate is  $\geq 5$  mm larger than Ceftazidime/Cefepime Neo-Sensitabs (see document **1.1.1**)

## Detection of ESBLs using Neo-Sensitabs™

### ENTEROBACTERIACEAE

At present there are no CLSI guidelines available for detection of ESBL in other species than *E.coli*, *K.pneumoniae*, *P.mirabilis* and *Salmonella* spp.

ESBLs are also found in other Enterobacteriaceae such as *Enterobacter* spp. etc (22)

Strains showing cefotaxime and/or ceftazidime MICs  $\geq 1$  µg/ml, showing reduced susceptibility to amoxicillin + clavulanate should be tested further for the presence of ESBLs.

## Detection of beta lactamases

Extended spectrum beta lactamases (ESBL), ESBL+AmpC Screen kit and ESBL Confirm ID kit

### Procedure 1 (Double disk synergy test)

Mueller Hinton agar plates are inoculated with the strain to be tested and Neo-Sensitabs applied onto the agar. Cefotaxime, Ceftazidime and Cefepime Neo-Sensitabs at a distance of 15-20 mm (edge to edge) from Amoxycillin + Clavulanate Neo-Sensitabs.

### Procedure 2 (Combination disks synergy test). ESBL Confirm ID kit



K. pneumoniae ATCC 700603, ESBL positive

Use Cefotaxime, Ceftazidime and Cefepime and their combinations with Clavulanate (ESBL Confirm ID kit).

### Interpretation

A key hole or ghost zone between Amoxycillin + Clavulanate and any of Cefotaxime, Ceftazidime or Cefepime Neo-Sensitabs indicates the presence of an ESBL.

When using the combination disks, a  $\geq 5$  mm larger zone for any of the combinations compared to the corresponding single antimicrobial indicates the presence of an ESBL.

Cefpodoxime and Cefpodoxime + Clavulanate may be used for screening purposes.

*Klebsiella oxytoca* hyperproducing K-1 beta-lactamase may show a false positive result (potentiation of cefotaxime and /or cefepime). Only when the strain is resistant to ceftazidime and shows synergism between ceftazidime and clavulanate should it be reported as ESBL positive.

## Detection of beta lactamases

### Extended spectrum beta lactamases (ESBL), ESBL+AmpC Screen kit and ESBL Confirm ID kit

#### ESBL + AmpC beta-lactamases

Current susceptibility tests are not accurate enough for surveillance of ESBL and AmpC producing pathogens. Laboratories testing for ESBL and AmpC need to be aware that some of these organisms may test false positive for ESBL, using CLSI methodology (40).

#### Combined disk test (ESBL+AmpC). ESBL+AmpC Screen kit

Apply one of each: Cefotaxime (CTX30), Cefotaxime+Clavulanate (CTX+C), Cefotaxime+Cloxacillin (CTXCX) and Cefotaxime+Clavulanate+Cloxacillin (CTXCC) Neo-Sensitabs on a MH plate inoculated (McFarland 0.5) with the strain to be tested.

#### Interpretation

		<b>Cefotaxime CTX30</b>	<b>Cefotaxime+Clav. CTX+C</b>	<b>Cefotaxime+Cloxa. CTXCX</b>
<b>ESBL</b>	CTX+C or CTXCC	≥ 5 mm -	- <4 mm	- ≥ 5 mm
<b>AmpC</b>	CTXCX or CTXCC	≥ 5 mm -	- ≥ 5 mm	- <4 mm
<b>ESBL+AmpC</b>	CTX+C <u>and</u> CTXCC	<4 mm -	- ≥ 5 mm	- ≥ 5 mm

In strains possessing both chromosomal (Enterobacter, Cit. Freundii etc. (44) or plasmidic AmpC beta lactamases and ESBLs, Boronic acid or Cloxacillin are used as inhibitors of the AmpC beta lactamase.

#### Combined disk test (KPC+ESBL)

For detection of ESBLs in clinical isolates of KPC carbapenemase-possessing Enterobacteriaceae, the following method is used (43):

Inoculate (McF 0.5) the MH agar plate with the strain to be tested and add:

1 Cefotaxime+Boronic acid and 1 Cefotaxime+Clavulanate+Boronic acid, 1 Ceftazidime+Boronic acid and 1 Ceftazidime+Clavulanate+Boronic acid.

#### Interpretation

Cefotaxime+Boronic+Clav. zone ≥ 5 mm than Cefotaxime+Boronic and/or Ceftazidime+Boronic+Clav. zone ≥ 5 mm than Ceftazidime+Boronic indicates the presence of an ESBL.

#### Detection of ESBLs in different strains

The presence of ESBLs may be masked by the overexpression of AmpC beta-lactamases or by the induction of AmpC beta-lactamase by clavulanate used in synergy tests. ESBLs may be confused with enzymes such as *K. oxytoca* chromosomal β-lactamase (K1). Laboratory staff must be aware of the increasing array of different resistance mechanisms and phenotypes.

## Detection of beta lactamases

## Extended spectrum beta lactamases (ESBL), ESBL+AmpC Screen kit and ESBL Confirm ID kit

### NON-FERMENTERS

Here are particularly *P. aeruginosa* and *A. baumannii* that may possess several types of beta-lactamases. Non-fermenters showing reduced susceptibility to ceftazidime and/or cefepime and/or aztreonam should be tested for the presence of ESBLs.

#### **Procedure**

Apply Ceftazidime, Cefepime and Aztreonam Neo-Sensitabs. At a distance of approx. 15 mm (edge to edge) apply Ticarcillin + Clavulanate Neo-Sensitabs. Separately apply Ceftazidime + Clavulanate and Cefepime + Clavulanate Neo-Sensitabs.

#### **Interpretation**

A key-hole zone or ghost zone between Ticarcillin + Clavulanate and any of Ceftazidime, Cefepime or Aztreonam Neo-Sensitabs indicates the presence of an ESBL.  
With the combination disks a  $\geq 5$  mm larger zone for Ceftazidime + Clavulanate and/or Cefepime + Clavulanate compared to the single antimicrobials indicates the presence of an ESBL.  
The prediffusion procedure with Boronic acid may also be used (30) when ESBLs can be obscured by the chromosomal AmpC cephalosporinase in *P. aeruginosa*.

### **1) E. coli, Klebsiella spp., Salmonella spp., Proteus mirabilis, Shigella sonnei:**

ESBL	Synergism Ceftazidime+Clavulanate and/or Cefepime+Clavulanate. Cefoxitin S.
VEB-1 (ESBL) (25)	Synergism Ceftazidime+Clavulanate and/or Cefepime+Clavulanate. Synergism Imipenem and Ceftazidime and/or Cefoxitin and Ceftazidime.
AmpC plasmid (no ESBL)	No synergism Ceftazidime+Clavulanate and Cefepime+Clavulanate. Cefoxitin R, Ceftazidime R. Synergism Cefotaxime/Ceftazidime and Cloxacillin and Cefotaxime or Cefotaxime and Boronic acid.
DHA (Induc. plasmid AmpC)	Antagonism Clavulanate (AMC) and 3rd generation. Synergism Cefotaxime/Ceftazidime and Cloxacillin or Boronic acid.
DHA+ ESBL (33)	Antagonism Clavulanate (AMC) and 3rd gen. cephalosporins (DHA). Synergism Tazobactam (Piperacillin+Tazobactam) and Ceftazidime/Cefepime.
Amp C + ESBL	Synergism Cefepime+Clavulanate: ESBL (31).Cefoxitin R, Ceftazidime R. Synergism Cefotaxime or Ceftazidime and Cloxacillin or Boronic acid: AmpC.
DHA + ESBL	Synergism Ceftazidime + Clavulanate and/or Cefepime + Clavulanate: ESBL. Antagonism Clavulanate (AMC) and 3rd gen. cephalosporins: DHA/ACT-1.
ESBL + Metallo-beta-lactamases (24,28)	Synergism Aztreonam + Clavulanate (Amoxicillin+Clavulanate): ESBL. Synergism Imipenem+EDTA: metallo-beta-lactamases. Synergism Meropenem and DPA: metallo-beta-lactamases.

## Detection of beta lactamases

## Extended spectrum beta lactamases (ESBL), ESBL+AmpC Screen kit and ESBL Confirm ID kit

ESBL + 16S rRNA methylases (26)	Synergism Ceftazidime+Clavulanate and/or Cefepime+Clavulanate. Cefoxitin S. No zone with Amikacin, Gentamicin, Tobramycin Neo-Sensitabs.
ESBL+carbapenemases (not MBL)	Synergism Aztreonam +Clav. (Amox+Clav) Ertapenem I/R, metallo-β-lactamases neg. Positive Hodge test.

### 2) High level K-1 (*Klebsiella oxytoca*):

(no ESBL)	No synergism Ceftazidime+Clavulanate. Ceftazidime S. Synergy is currently observed with third gen. cephalosporins and Clavulanate as well as with Cefepime+Clavulanate.
K-1 + ESBL (23)	Synergism Ceftazidime+Clavulanate. Ceftazidime I/R.

### 3a) *Enterobacter* spp., *Serratia* spp., *Providencia rettgeri*, *Citrobacter freundii*:

ESBL	Synergism Cefepime and Clavulanate (Amoxycillin+Clavulanate) and/or Ceftazidime and Clavulanate.
ESBL + 16S rRNA methylases	Synergism Cefepime and Clavulanate and/or Ceftazidime +Clavulanate. No zone with Amikacin, Gentamicin, Tobramycin Neo-Sensitabs.
ESBL + Metallo-beta-lactamases	Synergism Aztreonam (or Cefepime)+Clavulanate: ESBL Synergism Imipenem + EDTA and/or Meropenem + DPA: Metallo-beta-lactamases.

### 3b) *Morganella morganii*

ESBL	Synergism Cefepime and Tazobactam (Piperazillin+Tazobactam). Synergism Sulbactam (Ampicillin+Sulbactam) and Ceftazidime or Cefotaxime.
------	---

### 4) *Pseudomonas aeruginosa*:

ESBL	Ceftazidime S or R. Ticarcillin resistant. Synergism Ceftazidime+Clavulanate and/or Cefepime+Clavulanate. Synergism Aztreonam/Ceftazidime or Cefepime with Ticarcillin+Clavulanate
ESBL (PER-1)	Piperacillin S, Ceftazidime R. Synergy as above.
VEB-1 (ESBL)	Synergism between Imipenem and Ceftazidime (or Cefepime) in the presence of Cloxacillin. Synergism Ceftazidime+Clavulanate and/or Cefepime+Clavulanate.
OXA-18 (ESBL)	Ceftazidime R, Ticarcillin R, Aztreonam R, Meropenem I. Synergism Ceftazidime/Cefepime with Ticarcillin+Clavulanate.
ESBL + Metallo-beta-lactamases	Synergism Aztreonam (or Cefepime)+Clavulanate: ESBL Synergism Imipenem + EDTA and/or Meropenem + DPA: Metallo-beta-lactamases.

## Detection of beta lactamases

## Extended spectrum beta lactamases (ESBL), ESBL+AmpC Screen kit and ESBL Confirm ID kit

### 5) *Acinetobacter* spp. Use Double Disk Synergy Test

ESBL	Synergism Ceftazidime/Cefepime and Amoxycillin/Clavulanate. Synergism Ceftazidime/Cefepime and Ticarcillin+Clavulanate
PER-1 (ESBL)	Cephalosporins R, Aminoglycosides R. Synergy TIC+Clav, Synergy PIP+TAZO.
VEB-1 (ESBL)	Synergism Cefepime and Ticarcillin+Clavulanate (distance 15 mm). Best at 30 °C in the presence of Cloxacillin or Boronic acid.

Beceiro et al (37) have shown that the double disk synergy test gives the best results with *Acinetobacter* spp. due to *Acinetobacter*'s intrinsic susceptibility to clavulanic acid.

### 6) *Achromobacter xylosoxidans*

VEB-1 (ESBL)	Synergy between Ceftazidime and Clavulanate.
--------------	--

### 7) *Haemophilus influenzae*:

ESBL (27)	Compare Cefpodoxime and Cefpodoxime+Clavulanate. Zones: $\geq 5$ mm larger with the combination.
-----------	---

#### Note:

Although the lower cephalosporins MIC breakpoints of EUCAST and those forthcoming from CLSI are better able to separate ESBL pos/neg populations of *E.coli* and *Klebsiella* spp, significant numbers of ESBLpositive isolates remain below these breakpoints. Confirmation testing of ESBL is necessary (42).

#### **References:**

- 1) Sader H.S. et al: Prevalence of important pathogens and the antimicrobial activity of parenteral drugs at numerous Medical Centers in the U.S.. *Diagn. Microbiol. Infect. Dis.*, **20**, 203-208, 1994.
- 2) NCCLS: Performance Standards for Antimicrobial Disk Susceptibility Tests. 6th Ed., **M2-A6**, section 6.6., page 14-15, 1997.
- 3) Rasheed J.K. et al: Evolution of extended-spectrum beta-lactams resistance (SHV-8) in a strain of *E.coli* during multiple episodes of bacteremia. *Antimicrob. Agents Chemother.*, **41**, 647-653, 1997.
- 4) Casals J.B., Pringler N.: Detection in the routine laboratory of new plasmid mediated broadspectrum beta lactamases in Enterobacteriaceae. Abst 603 7th Mediterr. Congr. Chemother., Barcelona, 1990.
- 5) Livermore D.M. et al: Antibiotic resistance and production of extended spectrum betalactamases amongst *Klebsiella* spp. from intensive care units in Europe. *J. Antimicrob. Chemother.*, **38**, 409-424, 1996.
- 6) Cormican M.G. et al: Detection of extended-spectrum beta lactamases (ESBL) producing strains by the E-test ESBL Screen. *J. Clin. Microbiol.*, **34**, 1880-1884, 1996.
- 7) Wahaboglu H. et al: Resistance to extended-spectrum cephalosporins caused by PER-1 beta lactamase in *Salm.typhimurium* from Istanbul, Turkey. *J. Med. Microbiol.*, **43**, 294-299, 1995.
- 8) Karas J.A. et al: Treatment failure due to extended spectrum beta lactamase. *J. Antimicrob. Chemother.*, **37**, 203-204, 1996.
- 9) CLSI. Performance standards for Antimicrobial Susceptibility Testing. 15th Inf. Suppl. **M 100-S15**, 2005.
- 10) Vercauteren E. et al: Comparison of screening methods for detection of Extended-Spectrum Beta Lactamases and their prevalence among blood isolates of *E. coli* and *Klebsiella* spp. in a Belgian Teaching Hospital. *J. Clin. Microbiol.* **35**, 2191-2197, 1997.
- 11) Thomson K.S.: Controversies about extended-spectrum and AmpC beta lactamases. *Emerg. Infect. Dis.* 7, March/April 2001.
- 12) Steward C.D. et al: Characterization of clinical isolates of *Kl. pneumoniae* from 19 laboratories using the NCCLS ESBL detection methods. *J. Clin. Microbiol.* **39**, 2864-2872, 2001.

## Detection of beta lactamases

## Extended spectrum beta lactamases (ESBL), ESBL+AmpC Screen kit and ESBL Confirm ID kit

- 13) de Gheldre Y. et al: National epidemiologic survey of *Enterobacter aerogenes* Belgian hospitals from 1996 to 1998. J. Clin. Microbiol. **39**, 889-896, 2001.
- 14) Livermore D.M.: Detection of beta-lactamase-mediated resistance. J. Antimicrob. Chemother. **48**, Suppl S1, 59-64, 2001.
- 15) Kamile Rasheed J. et al: Characterization of ESBL Reference Strain, *Kl. pneumoniae* ATCC 700603, which produces the novel enzyme SHV-18. Antimicrob. Ag. Chemother. **44**, 2382-2388, 2000.
- 16) Florijn A. et al: Comparison of E-test and double disk diffusion test for detection of extended spectrum Beta-lactamases. Eur. J. Clin. Microbiol. Infect. Dis., **21**, 241-243, 2002.
- 17) Tzelepi E. et al: Detection of ESBL in clinical isolates of *Enterobacter cloacae* and *E. aerogenes*. J. Clin. Microbiol., **38**, 542-6, 2000.
- 18) Oliver A. et al: Mechanisms of decreased susceptibility to Cefpodoxime in *E. coli*. Antimicrob. Ag. Chemother., **46**, 3829-3836, 2002.
- 19) Bolmström A.: Cefepime ± Clavulanic acid in a E-test configuration for investigating non-determinable ESBL results for NCCLS criteria. 42nd ICAAC, Presentation D-527, 2002.
- 20) Rodriguez-Villalobos H. et al: Evaluation of the combined disk method for the detection of ESBL in Enterobacteriaceae. ICAAC 2003, presentation D-202.
- 21) Fluit A. et al: Comparison of phenotypic tests for the detection of ESBLs and AmpC beta-lactamases. ICAAC 2003, presentation D-199.
- 22) Pitout J.D.D. et al: Modification of the double-disk test for detection of Enterobacteriaceae producing ESBL and AmpC beta lactamases. J. Clin. Microbiol., **44**, 3933-5, 2003.
- 23) Decré D. et al: Outbreak of multiresistant *K. oxytoca* involving strains with ESBLs and strains with extended-spectrum activity of the chromosomal  $\beta$ -lactamase. J. Antimicrob. Chemother., **54**, 881-8, 2004.
- 24) Lartigue M.F. et al: First detection of a carbapenem-hydrolyzing metalloenzyme in an Enterobacteriaceae isolate in France. Antimicrob. Ag. Chemother., **48**, 4929-30, 2004.
- 25) Kim J-Y et al: Nosocomial outbreak by *P. mirabilis* producing extended-spectrum  $\beta$ -lactamase VEB-1 in a Korean University Hospital. J. Antimicrob. Chemother., **54**, 1144-7, 2004.
- 26) Yan J-J et al: Plasmid mediated 16S rRNA methylases conferring high-level aminoglycoside resistance in *E. coli* and *K. pneumoniae* isolates from two Taiwanese hospitals. Antimicrob. Chemother., **54**, 1007-12, 2004.
- 27) Tristram S.G. et al: Disc diffusion-based screening tests for ESBL's in *H. influenzae*. J. Antimicrob. Chemother., **55**, 570-3, 2005.
- 28) Lincopan N. et al: First isolation of metallo- $\beta$ -lactamase-producing multiresistant *K. pneumoniae* from a patient in Brazil. J. Clin. Microbiol., **43**, 516-9, 2005.
- 29) Pitout J.D.D. et al: Emergence of Enterobacteriaceae producing ESBLs in the community. J. Antimicrob. Chemother., **56**, 52-59, 2005.
- 30) Lee Seungok et al: Prevalence of Ambler class A and D beta-lactamase among clinical isolates of *P. aeruginosa* in Korea. J. Antimicrob. Chemother., **56**, 122-27, 2005.
- 31) Song Wonkeun et al: Failure of cefepime therapy in treatment of *K. pneumoniae* bacteremia. J. Clin. Microbiol., **43**, 4891-94, 2005.
- 32) Schwaben J.J. et al: Utility of the Vitek 2 Advanced Expert System for identification of ESBL production in *Enterobacter* spp. J. Clin. Microbiol., **44**, 241-3, 2006.
- 33) Muratani T. et al: Emergence and prevalence of beta-lactamase producing *K. pneumoniae* resistant to cepheims in Japan. Intl. J. Antimicrob. Ag., **27**, 491-9, 2006.
- 34) Wonkeun Song et al: Detection of ESBLs by using Boronic acid as an AmpC  $\beta$ -lactamase inhibitor in clinical isolates of *K. pneumoniae* and *E. coli*. J. Clin. Microbiol., **45**, 1180-84, 2007.
- 35) Wonkeun Song et al: Clonal spread of both oxyimino-cephalosporin- and cefoxitin-resistant *K. pneumoniae* isolates coproducing SHV-2a and DHA-1  $\beta$ -lactamase at a burns intensive unit. Int. J. Antimicrob. Agents, **28**, 520-24, 2006.
- 36) Naiemi N et al: Widely distributed and predominant CTX-M ESBLs in Amsterdam, The Netherlands. J. Clin. Microbiol., **44**, 3012-3014, 2006.
- 37) Beceiro A. et al: False ESBL detection in *Acinetobacter* spp, due to intrinsic susceptibility to clavulanic acid. J. Antimicrob. Chemother **61**, 301-8, 2008.
- 38) Sook Hoon Jeong et al: Boronic acid disk tests for identification of ESBL production in clinical isolates of Enterobacteriaceae producing chromosomal AmpC beta-lactamases. Int. J. Antimicrob. Agents **31**, 467-471, 2008.

## Detection of beta lactamases

### Extended spectrum beta lactamases (ESBL), ESBL+AmpC Screen kit and ESBL Confirm ID kit

- 39) Moland E.S et al: Newer beta-lactamases; clinical and laboratory implications. Part I Clin. Microbiol. Newsletter **30**, 71-85, 2008.
- 40) Hanson N. D. et al: Community based reservoirs reveals the presence of CTX-M imported AmpC and OXA-30 beta-lactamases in urine isolates of K. pneumoniae and E.coli in a U.S community. Antimicrob. Chemother, **52**, 3814-16, 2008.
- 41) Granier S. A et al: False susceptibility of K. oxytoca to some extended-spectrum cephalosporins. J. Antimicrob. Chemother **50**, 303-304, 2002.
- 42) Badal R. et al: Evaluation of EUCAST and CLSI cephalosporins breakpoints versus ESBL pos and ESBL neg. E.coli, K.pneumoniae and K. oxytoca in the SMART study. Abstract Presentation E-204, 49<sup>th</sup> ICAAC 2009.
- 43) Tsakris A. et al: Use of Boronic acid disk tests to detect ESBL in clinical isolates of KPC carbapenemase-possessing Enterobacteriaceae. J. Clin. Microbiol. **47**, 3420-3426, 2009.
- 44) Yano H. et al: Detection of ESBL by using Boronic acid disk in clinical isolates of Enterobacter spp producing chromosomal AmpC beta lactamases.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Detection of beta lactamases

Page 1 / 1 - Document: 1.1.1

ESBL Quality Control

### ESBL Quality Control

*Klebsiella pneumoniae* ATCC 700603

#### Original

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm	MIC µg/ml (15)
Aztreonam	30 µg	AZTRM	11-19	64
Cefepime	30 µg	CFEPM	21-26	1
Cefepime+Clavulanate	30+10 µg	CP+CL	≥5 mm larger	≤ 0.12
Cefotaxime	30 µg	CFTAX	20-26	8
Ceftazidime	30 µg	CEZDI	11-19	32
Ceftazidime-Clavulanate	30+10 µg	CZ+CL	≥5 mm larger	1
Ceftriaxone	30 µg	CETRX	16-24	16

#### CLSI/EUCAST potency Neo-Sensitabs™

CLSI/EUCAST potency NEO-SENSITABS™	POTENCY	CODE	Zone diameter in mm	MIC µg/ml (15)
Cefotaxime	30 µg	CTX30	17-25	8
Cefotaxime-Clavulanate	30+10 µg	CTX+C	≥5 mm larger	-
Ceftazidime	30 µg	CAZ30	10-18	32
Ceftazidime-Clavulanate	30+10 µg	CAZ+C	≥5 mm larger	1
Cefepime	30 µg	FEP30	19-25	1
Cefepime-Clavulanate	30+10 µg	FEP+C	≥5 mm larger	≤0.12
Cefpodoxime	10 µg	CPD10	9-16	-
Cefpodoxime-Clavulanate	10+1 µg	CPD+C	≥5 mm larger	-

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 1 / 1 - Document: 1.2.0

## Detection of beta lactamases

## Inducible cephalosporinases or AmpC beta-lactamases

### Inducible Cephalosporinases or AmpC Beta-lactamases

Inducible cephalosporinases or AmpC beta-lactamases are produced by *Enterobacter cloacae*, *E. aerogenes*, *Serratia marcescens*, *Citrobacter freundii*, *Hafnia alvei*, *Providencia stuartii* and *Morganella morganii*, and they are inhibited by aztreonam, but not by clavulanic acid, sulbactam or tazobactam. Resistant mutants with high beta-lactamase activity are present at a high frequency. As a result therapy with cephalosporins (except fourth generation agents) and monobactams may fail because of selection of such mutants.

The tablet approximation test is useful to demonstrate the presence of inducible cephalosporinases, during routine antibiogram testing.

Neo-Sensitabs containing an inducer, e.g. ceftazidime (or imipenem) and indicators such as piperacillin+tazobactam, cefotaxime or ceftazidime are placed approx. 20-25 mm apart center to center. A wider spacing (30 mm) may be preferable for e.g. *M. morganii* and *Providencia* spp.

Following overnight incubation at 35 °C in air, the presence of an inducible beta-lactamase is indicated by the blunting of the zone of inhibition around the indicator drug (piperacillin+tazobactam, cefotaxime/ceftazidime) adjacent to the inducer (ceftazidime/imipenem).

Dunne et al (1) have shown that the combination Imipenem and Piperacillin + Tazobactam has the highest sensitivity (97.1 %) followed by Imipenem and Ceftazidime (94.2 %).

The result should be reported as R (resistant) for penicillins (except temocillin), penicillin/inhibitor combinations, cephalosporins (except ceftazidime and ceftazidime), cephamycins and monobactams, irrespective of the size of the inhibition zone.



Plate 1.2.0-a



Plate 1.2.0-b

Demonstration of the presence of inducible beta-lactamases in *Enterobacter cloacae* (ATCC 13047). Note the flattened edges of Cefotaxime Neo-Sensitabs (CFTAX) and Ceftazidime Neo-Sensitabs (CEZDI) zones adjacent to Cefoxitin Neo-Sensitabs (CFOXT, Plate 1.2.0-a) and Imipenem Neo-Sensitabs (IMIPM, Plate 1.2.0-b), respectively.

### References:

1) Dunne W.M. et al: Use of several inducer and substrate antibiotic combinations in a disk approximation assay format to screen for AmpC induction in patient isolates of *P. aeruginosa*, *Enterobacter* spp., *Citrobacter* spp. and *Serratia* spp. J.C.M., **43**, 5945-9, 2005.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Detection of beta lactamases

### Testing / Reporting of Susceptibility to Beta-lactams against Enterobacteriaceae and Non-fermenters

Use the table below for testing/reporting of susceptibility to beta-lactams against Enterobacteriaceae and non-fermenters causing serious infections when inducible beta-lactamases are present <sup>a)</sup>:

	AMP	AMC	CXM	CTX	CTR	CAZ	FEP	CFO	AZT	IMI	MRP	TIM	PI+TZ
<i>E. aerogenes/cloacae</i> <i>C. freundii</i> / <i>S. marcescens</i>	R	R	R	R	R	R	T	R	R	T	T	R	R
<i>Prov. stuartii/rettgeri</i> <i>Morg. morganii</i>	R	R	R	R	R	R	T	R	T	T	T	R	T
<i>P. vulgaris/penneri</i>	R	T	R	R	R	T	T	T	R	T	T	T	T
<i>Klebsiella oxytoca</i> <sup>e)</sup>	R	T	R	R	R	T	T	T	R	T	T	T	T
<i>Hafnia alvei</i>	R	R	R	R	R	R	T	R	R	T	T	R	R
Enterobacteriaceae with ESBL (no inducible β-lactamases)	R	T <sup>b)</sup>	R	R	R	R	R	R	R	T	T	T <sup>b)</sup>	T <sup>b)</sup>
Enterobacteriaceae with inducible β-lactamases and ESBL	R	R	R	R	R	R	R	R	R	T	T	R	R
<i>Aeromonas</i> with A2 <sup>c)</sup> (most <i>A. sobria</i> )	R	R	T	T	T	T	T	T	T	R	R	R	R
<i>Aeromonas</i> with A 1 and A 2	R	R	R	R	R	R	T	R	T	R	R	R	R
<i>Ps. aeruginosa</i> <i>Burkholderia</i> spp.	R	R	R	R	R	T	T	R	T	T	T	T	T
<i>S. maltophilia</i>	R	R	R	R	R	T	R	R	R	R	R	T	R
<i>A. baumannii</i>	R	T <sup>d)</sup>	R	R	R	T	T	R	R	T	T	T	T

AMC, Amoxycillin+Clavulanate; AMP, Ampicillin; AZT, Aztreonam; CXM, Cefuroxime; CTR, Ceftriaxone; CAZ, Ceftazidime; FEP, Cefepime; CFO, Cefoxitin; CTX, Cefotaxime; IMI, Imipenem; MRP, Meropenem; PI+TZ, Piperacillin+Tazobactam; TIM, Ticarcillin+Clavulanate.

A1 Inducible cephalosporinase, the enzyme is usually found in *A. hydrophila* and *A. caviae*. These species are considered resistant to cephalosporins and cephamycins.

A2 Penicillinase/carbapenemase that hydrolyses imipenem and meropenem. The expression may be heterogeneous.

a) R, the microorganism is resistant and may possess a resistance mechanism not always detected by the diffusion method.

T, these antimicrobials may be used for testing.

b) Test isolates from urine only. Isolates from other sites are considered resistant.

c) Test also for Cefazolin. *A. caviae* does not possess a carbapenemase (A2) and can be tested against imipenem and meropenem.

d) Use Ampicillin+Sulbactam.

e) *K. oxytoca* producing a K-1 enzyme are susceptible to ceftazidime. They may show synergism between other third gen. cephalosporins and amoxicillin+clavulanate and may be mistaken as ESBL producers.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Detection of beta lactamases

Page 2 / 2 - Document: 1.2.1

Testing/Reporting

Adopted from CDS 2005 with modifications.

### References:

- 1) Sanders C.C. and Sanders W.E.: Beta-lactam resistance in gram negative bacteria: global trends and clinical impact. Clin. Infect. Dis., **15**, 824-39, 1992.
- 2) Livermore D.M. et al.: Detection of beta-lactamase mediated resistance. J. Antimicrob. Chemother., **48**, Suppl. S1, 59-64, 2001.

ORIGINAL ROSCO DOCUMENT

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 1 / 8 - Document: 1.3.0

## Detection of beta lactamases

## Plasmid mediated AmpC beta lactamases. AmpC Confirm ID kit

### Plasmid-mediated AmpC Beta-lactamases

Plasmid-mediated beta-lactamases represent a new threat, since they confer resistance to aminopenicillins, carboxypenicillins, ureidopenicillins, although they are generally susceptible *in vitro* to mecillinam and/or temocillin.

The enzymes provide resistance to third generation cephalosporins and ceftazidime. The enzymes are also active against aztreonam although for some strains the aztreonam MICs are in the susceptible range. Susceptibility to ceftazidime is little affected (inoculum effect) and the carbapenems are not affected. The enzymes are not affected by beta-lactamase-inhibitors, except for CMY-8 and CMY-9 that are inactivated by tazobactam.

Their expression is generally constitutive, nevertheless inducible plasmid AmpC (ACT-1, DHA-1, DHA-2, CFE-1, CMY-13) have been reported (6).

Plasmid-mediated AmpC beta-lactamases have been found most frequently in species naturally negative for AmpC, such as *K. pneumoniae*, *E. coli*, *K. oxytoca*, *Salmonella* and *P. mirabilis*. Recently they were also found in *Enterobacter* spp. (2).

The strains with plasmid-mediated AmpC show resistance to ceftazidime (MIC > 16 µg/ml) and ceftazidime (MIC > 32 µg/ml) corresponding to zones of inhibition < 16 mm (McF. 0.5).

Strains with plasmid-mediated AmpC do not show antagonism between ceftazidime and 3rd generation cephalosporins (are not inducible), while inducible plasmid-mediated AmpC (ACT-1, DHA-1, DHA-2, CMY-13) show antagonism between ceftazidime (or imipenem) and third generation cephalosporins.

Isolated that coproduce an ESBL and a plasmid mediated AmpC beta-lactamase may yield a positive confirmatory test for ESBL using ceftazidime and ceftazidime+clavulanate (synergism).

### Characteristics of AmpC beta-lactamases:

	Chromosomally mediated AmpC (partially derepressed AmpC mutants)	Plasmid-mediated AmpC (derepressed AmpC mutants)	Inducible plasmid-mediated AmpC (ACT-1, DHA-1, DHA-2, CFE-1, CMY-13)	ESAC in <i>P. aeruginosa</i> (21)
Ceftazidime+Clavulanate and/or Ceftazidime+Clavulanate	No synergism	No synergism (except MOX-1, MOX-2)	No synergism	No synergism
Ceftazidime, Imipenem or Amoxicillin+Clavulanate (7)	Ceftazidime R (zone < 16 mm) <b>Antagonism</b> with 3rd gen. cepha.	Ceftazidime R (zone < 16 mm) <b>No antagonism</b> with 3rd gen cepha.	Ceftazidime R (zone < 16 mm) <b>Antagonism</b> with 3rd gen. cepha.	<b>No antagonism</b> Imipenem / 3rd gen cepha.
Ceftazidime Ceftazidime	S → R S	R (zone < 20 mm) S	S → R S	Ceftazidime R Ceftazidime I/R
Cloxacillin	Synergism Cloxacillin-ceftazidime	Synergism Cloxacillin-ceftazidime Cloxacillin-ceftazidime	Synergism Cloxacillin-ceftazidime Cloxacillin-ceftazidime	Synergism Cloxacillin-ceftazidime Cloxacillin-Carbapenems
Boronic acid synergy	Ceftazidime-ceftazidime	Ceftazidime-ceftazidime	Ceftazidime-ceftazidime	Carbapenems

*Enterobacter* spp., *C. freundii*, *M. morgani*, *Hafnia alvei*, *Providencia* spp., *Proteus* indole positive and *Serratia marcescens*, all produce an inducible chromosomal AmpC beta lactamase, which is not inhibited by clavulanate. There may be seen an antagonism between amoxicillin and clavulanate (smaller zone with the combination that with amoxicillin alone) due to the presence of the inducible beta-lactamase.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Detection of beta lactamases

## Plasmid mediated AmpC beta lactamases. AmpC Confirm ID kit

All these strains should be reported as resistant to ampicillin/amoxycillin and to amoxycillin+clavulanate (except *P. vulgaris*).

Using an Amoxycillin+Clavulanate disc (Neo-Sensitabs) better performance is obtained due to the dual action of clavulanic acid: 1) induces expression of inducible plasmid mediated AmpC beta-lactamases (antagonism with 3rd gen. cephalosporins) and 2) permits the detection of an ESBL by enlarging inhibition zones of 3rd gen. cephalosporins (synergism) (7). In the presence of an ESBL + an inducible plasmid-mediated AmpC, both antagonism and synergism can be detected in the same plate (7).

### Differentiation of AmpC beta-lactamases in *E. coli*

Mirelis et al (9) and Aragon (15) found a simple phenotypic method for the differentiation between plasmid-mediated and chromosomal AmpC-β-lactamases in *E. coli* and *P. mirabilis* using Cloxacillin Diatabs and by visual examination of the antibiogram plates. The presence of scattered colonies located near the edge of the zone of inhibition of cefoxitin, cefotaxime, ceftazidime and aztreonam indicated the presence of plasmid-mediated AmpC beta-lactamases.

Cloxacillin or Boronic acid alone do not distinguish between chromosomal or plasmidic AmpC beta-lactamases.

### *E. coli* (AMC I/R, Cefotaxime I/R, Ceftazidime I/R Cefoxitin: most I/R)

	Plasmid AmpC	Inducible plasmid AmpC	Chromosomal AmpC hyperprod.	Chromosomal ESAC (14, 18)
<b>Cloxacillin</b>	Synergy with ceftazidime or cefotaxime	Synergy with ceftazidime or cefotaxime	Synergy with ceftazidime or cefotaxime	Synergy
<b>Boronic acid</b>	Synergy with Ceftazidime and/or Cefotaxime	Synergy with Ceftazidime and/or Cefotaxime	Synergy with Ceftazidime. and/or Cefotaxime	Synergy
<b>Cefoxitin Imipenem</b>	No antagonism with 3 <sup>rd</sup> gen. cephalosporins	Antagonism with 3 <sup>rd</sup> gen. cephalosporins	No antagonism	No antagonism
<b>Antibiogram</b>	Scattered colonies (resistant mutants) near the edge of the zone of cefoxitin, cefotaxime, ceftazidime and aztreonam	Scattered colonies	Well defined edge of zone	
<b>Cefepime MIC µg/ml</b>	Cefepime MIC ≤1 µg/ml	Cefepime MIC ≤1 µg/ml	Cefepime MIC ≤1 µg/ml	Cefepime MIC 1-8 µg/ml (zone < 26mm) Ceftazidime R Cefoxitin R

- ESAC=extended spectrum AmpC (14)
- The same procedure will be appropriate for *K. pneumoniae* and *P. mirabilis* strains.
- Chromosomal AmpC's of *E. coli* are not inducible.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Detection of beta lactamases

## Plasmid mediated AmpC beta lactamases. AmpC Confirm ID kit

Plasmid-mediated AmpC beta-lactamases MIC µg/ml

Beta-lactamases	FOX	CAZ	AZT	FEP	IMP	MRP	Microorganisms
AAC-1	4-8	≥32	1	0.25	0.125-0.8	0.03	E. coli, K. pneumoniae, P. mirabilis, Salmonella, C. freundii
ACT-1	> 256	4 - > 128	4 → > 128	≤ 0.06-8	1	·	E. coli, K. pneumoniae, E. cloacae (inducible)
BIL-1 (CMY-2)	R	> 16	4-16	1	0.5	0.06	E. coli
CMY-1	256	4-128	32	0.25-4	≤ 0.5	0.06	K. pneumoniae
CMY-2	32-256	32-128	16-64	0.5-4	≤ 0.5	0.06	E. coli, Salmonella, K. pneumoniae
CMY-3	128	64	32-256	1	0.25	0.03	P. mirabilis
CMY-4	8 - > 256	8-256	0.5-32	0.06-4	0.25	0.125	E. coli, Salmonella, P. mirabilis
CMY-5	R	256	64	-	0.5-1	-	K. oxytoca
CMY-6	256	256	64	0.5	0.25	0.06	E. coli
CMY-7	R	> 32	·	I	0.25	< 2	E. coli, Salmonella
CMY-8	> 256	32-64	·	·	0.25-0.5	-	K. pneumoniae
CMY-9	> 128	128	8	0.25	0.5	0.06	E. coli
CMY-10	> 128	8-64	4-128	0.12-0.5	0.25-0.5	≤ 0.125	E. coli, E. aerogenes, K. pneumoniae
CMY-11	> 256	256	128	·	·	·	E. coli
CMY-12	256	128	8-32	16	0.25-4	0.5	P. mirabilis
CMY-13	512	256	64	1	0.25	≤ 0.03	E. coli (inducible)
CMY-14	128	128	16-32	0.5-32	0.25-2	0.06	P. mirabilis
CMY-15	512	128	8-32	0.25-8	0.25-16	4	P. mirabilis
CMY-16	≥ 32	≥ 32	1	2	2	0.05	P. mirabilis / Synergy TAZO-FEP
CMY-19	≥ 128	> 128	16	4	0.25	≤ 0.06	K. pneumoniae (8)
CMY-20	≥ 128	> 128	>64	4	-	1	E. coli
CMY-21	> 64	64	32	0.5	0.25	0.06	E. coli
CMY-29	≥ 256	128	·	1	0.5	0.12	E.coli
CMY-30	128	32	128	5	0.5	0.12	E.coli
ESAC (14)			1-8				E.coli
ESAC (18)	> 256	> 256	≥32	32			E.coli
CFE-1	R	64	8	0.25	0.25	·	E. coli (inducible)
DHA-1	128-512	8-64	1-16	≤ 0.125-2	≤ 0.125 - 0.5	·	E. coli, K. pneumoniae, Salmonella, P. mirabilis (inducible)
DHA-2	16	8	2	0.03	0.25	·	K. pneumoniae (inducible)
FOX-1	128	8	1	1	0.25	≤ 0.03	E. coli, K. pneumoniae
FOX-2	256	32	2	0.13	0.5	0.03	E. coli
FOX-3	64	16	1	≤ 0.06	0.12	-	E. coli, K. oxytoca
FOX-4	> 512	> 128	64	2	0.5	0.12	E. coli
FOX-5	512	128	8-16	0.5	0.5	-	K. pneumoniae, E. coli
FOX-7	-	-	-	-	-	-	E. coli, K. pneumoniae, E. cloacae
LAT-1	64-256	> 128	64	1	0.25-2	0.06	K. pneumoniae
LAT-2 (CMY-2)	256	> 256	64-256	-	-	-	E. coli, K. pneumoniae, E. aerogenes
LAT-3 (CMY-6)	256	128	64	0.5	0.25	0.06	E. coli
LAT-4 (CMY-1)	64-256	8-256	8-128	0.125-1	0.25	0.125	E. coli
MIR-1	≥ 256	128	128	1	1	0.125	E. cloacae, K. pneumoniae
MOX-1	R	16 → R	16	-	0.5	-	K. pneumoniae
MOX-2	≥ 128	4-256	·	0.25-4	≤ 0.125	-	K. pneumoniae

Gupta et al. (3) describes isolation of multiresistant Salmonella, with plasmid-mediated AmpC beta-lactamase, from cattle and humans in the USA.

Three cases of invasive infections caused by Salmonella enterica serotype cholerasuis found in Taiwan (5). The strains were resistant to ciprofloxacin (mutations gyrA and ParC) and to ceftioxone (presence of plasmid-mediated CMY-2 beta-lactamase).

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 4 / 8 - Document: 1.3.0

## Detection of beta lactamases

### Plasmid mediated AmpC beta lactamases. AmpC Confirm ID kit

Recent studies (4) show that the outcome of cephalosporin treatment in serious infections due to AmpC beta-lactamase producing *K. pneumonia* isolates was poor. A standard test for detection of plasmid-mediated AmpC beta-lactamases is needed. Emergence of cefepime-hydrolyzing CMY-19 in Japan (8).

#### Detection of Plasmid-mediated-AmpC beta-lactamases



The capability to detect AmpC is important to improve the clinical management of infections and provide sound epidemiological data. Reduced susceptibility to cefoxitin in the Enterobacteriaceae may be an indicator of AmpC activity, but it should be confirmed by other tests. Laboratories should be able to recognize AmpC derepressed strains and those with plasmid AmpC. Guidelines from the CLSI are not yet available for detection of bacteria with AmpC beta-lactamases.

#### Screening

Derepressed/plasmid AmpC should be suspected when we see:

- Resistance to 3rd generation cephalosporins – NOT Cefepime.
- Resistance to Cefoxitin (inhibition zone < 16 mm).
- No cephalosporin / Clav. synergism.
- I / R to Amoxycillin + Clav.
- AmpC derepressed *Serratia* are S to ceftazidime.
- *Providencia*, *Morganella* and *Serratia* inducible & derepressed may appear S / I to cefoxitin.
- Strains producing AAC-1 beta-lactamase are susceptible to cefoxitin

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 5 / 8 - Document: 1.3.0

## Detection of beta lactamases

Plasmid mediated AmpC beta lactamases.  
AmpC Confirm ID kit

### Confirmation

#### 1. Combined disk test. AmpC Confirm ID kit

Apply Cefotaxime, Cefotaxime+Boronic, Ceftazidime, Ceftazidime+Boronic, Cefotaxime+Cloxacillin and Ceftazidime+ Cloxacillin in an inoculated plate.



E.coli ATCC FN9414 AmpC positive

### Interpretation

A Cefotaxime+Cloxacillin inhibition zone  $\geq 5$  mm than Cefotaxime alone and/or a Ceftazidime+Cloxacillin zone  $\geq 5$  mm than Ceftazidime alone indicates the presence of an AmpC.

A Cefotaxime+Boronic inhibition zone  $\geq 5$  mm than Cefotaxime alone and/or a Ceftazidime+Boronic inhibition zone  $\geq 5$  mm than Ceftazidime alone indicates the presence of an AmpC.


Inhibition zone  $\geq 5$  mm than Ceftazidime alone, indicates the presence of AmpC.

### 2. Double disk synergy test

Apply one Cefotaxime (CTX) and one Ceftazidime (CAZ) Neo-Sensitabs on an inoculated MH agar plate. In between apply one Boronic Acid Diatabs (BOR) at a distance of approx. 10 mm (edge to edge). If the strain is totally resistant to the Cephalosporins combination, the distance should be reduced to 5 mm.

Apply one Ceftazidime (CAZ) and one Cefotaxime (CTX) Neo-Sensitabs. In between at a distance of 5-10 mm edge to edge, apply one Cloxacillin Neo-Sensitabs.

Instead of Cefotaxime and Ceftazidime their combinations with clavulanate may be used for the synergy test (19).



E.coli ATCC FN9414 AmpC positive

### Interpretation

A keyhole or ghost zone (synergism) between Boronic Acid and any of Cefotaxime or Ceftazidime indicates the presence of an AmpC beta-lactamase.

A keyhole or ghost zone between Cloxacillin and Ceftazidime and/or Cefotaxime indicates the presence of an AmpC beta-lactamase.

Plasmid mediated AmpC differ from chromosomal AmpC in being uninducible (few exceptions). Strains producing inducible plasmid AmpC beta-lactamases (ACT-1, DHA-1, DHA-2, CFE-1, CMY-13) will show antagonism (distorted zone) between Cefoxitin or Imipenem and 3rd generation cephalosporins.

Strains of Klebsiella spp, Salmonella spp and P. mirabilis showing synergism with Boronic Acid and/or Cloxacillin possess presumptively plasmid mediated AmpC beta-lactamases.

The method cannot distinguish between chromosomal and plasmid mediated AmpC beta-lactamases in E. coli, but the test is useful to select strains for further analysis. Plasmid mediated are often multiresistant and may show scattered colonies near the edge of the zone of third gen. cephalosporins and aztreonam disks.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 7 / 8 - Document: 1.3.0

## Detection of beta lactamases

Plasmid mediated AmpC beta lactamases.  
AmpC Confirm ID kit

### Inducible phenotype

The inducible phenotype is identified by a tablet approximation test, using Imipenem or Cefoxitin against 3rd generation cephalosporins (distance 15 mm from edge to edge).

Distorted zones indicate the presence of an inducible AmpC beta-lactamase.

Treatment with 3rd generation cephalosporins should be avoided in severe Enterobacter, C. freundii, Serratia and Morganella infections except in UTI, because of risk for selection of cephalosporin-resistance during therapy.

### AmpC + ESBL. ESBL +AmpC Confirm ID kit

Screening criterion for ESBL presence among AmpC-producing Enterobacter, C. freundii and Serratia marcescens is Cefepime MIC > 1 ug/ml (inhibition zone < 26 mm).

High level expression of AmpC may prevent recognition of an ESBL. Use of Cefepime is more reliable to detect these strains because high AmpC production has little effect on cefepime activity.

### Combined disk test (ESBL+AmpC)

Apply one of each:

- A) Cefotaxime (CTX 30)
- B) Cefotaxime+Clavulanate (CTX+C)
- C) Cefotaxime+Cloxacillin (CTXCX)
- D) Cefotaxime+Clavulanate+Cloxacillin (CTXCX) on the inoculated MH plate.

### Interpretation

		<b>Cefotaxime CTX30</b>	<b>Cefotaxime+Clav. CTX+C</b>	<b>Cefotaxime+Cloxa. CTXCX</b>
<b>ESBL</b>	CTX+C or CTXCC	≥ 5 mm -	- <4 mm	- ≥ 5 mm
<b>AmpC</b>	CTXCX or CTXCC	≥ 5 mm -	- ≥ 5 mm	- <4 mm
<b>ESBL+AmpC</b>	CTX+C <u>and</u> CTXCC	<4 mm -	- ≥ 5 mm	- ≥ 5 mm

Neither ESBL or AmpC: All zones within 2 mm of each other CTXCC= Cefotaxime+Clavulanate+Cloxacillin.

### Detection of ESAC in P.aeruginosa (10)

ESAC in P. aeruginosa constitutes a favorable background for the selection of carbapenem-resistant strains. P. aeruginosa isolates being I/R to Imipenem and Ceftazidime should be tested for ESAC.

Place 1 Cloxacillin Diatabs between Imipenem and Cefepime Neo-Sensitabs (distance 8 mm from edge to edge). Synergy between Cloxacillin and Imipenem and/or Cefepime indicates the presence of and ESAC (extended spectrum AmpC)

### **Please note:**

KPC and other class A carbapenemases (Sme, IMI, GES, NmcA) may give a positive synergy test with Boronic acid. They do not show synergism when using Cloxacillin (differentiation from AmpC).

K. oxytoca hyperproducing K-1 enzyme may give a positive synergy test with Boronic acid, but they show susceptibility to Ceftazidime and are therefore easily detected.

## Detection of beta lactamases

## Plasmid mediated AmpC beta lactamases. AmpC Confirm ID kit

### References:

- 1)Smith Moland E. et al: Occurrence of newer beta-lactamases in *Klebsiella pneumoniae* isolates from 24 US Hospitals. Antimicrob. Agents Chemother., **46**, 3837-3842, 2002.
- 2)Philippon A. et al: Plasmid-determined AmpC-type beta-lactamases. Antimicrob. Agents Chemother., **46**, 1-11, 2002.
- 3)Gupta et al: Emergence of multidrug resistant *Salmonella enterica* serotype Newport resistant to expanded-spectrum cephalosporins in the US. J. Infect. Dis., **188**, 1707-16, 2003.
- 4)Hyunjeon Pai et al: Epidemiology and clinical features of bloodstream infections caused by AmpC type beta-lactamase-producing *K. pneumoniae*. Antimicrob. Agents Chemother., **48**, 3720-28, 2004.
- 5)Wen-Chien Ko et al: A new therapeutic challenge for old pathogens: community acquired invasive infections caused by ceftriaxone resistant *Salmonella enterica* serotype Choleraesuis. CID, **40**, 315-8, 2005
- 6)Nakano R et al: CFE-1, a novel plasmid-encoded beta-lactamase with an ampR gene originating from *C. freundii*. Antimicrob. Agents Chemother., **48**, 1151-8, 2004.
- 7)Kyungwon Lee et al: Evaluation of phenotypic screening methods for detecting plasmid mediated AmpC beta-lactamases-producing isolates of *E. coli* and *K. pneumoniae*. Diagn. Microbiol. Infect. Dis., **53**, 319-323, 2005.
- 8)Wachino Jun-ichi et al: Horizontal transfer of blaCMY-bearing plasmids among clinical *E. coli* and *K. pneumoniae* isolates and emergence of cefepime-hydrolyzing CMY-19, AAC, **50**, 534-541, 2006.
- 9)Mirelis B. et al: A simple phenotypic method for the differentiation between acquired and chromosomal AmpC beta-lactamases in *E. coli*. Enferm. Infect. Microbiol. Clin., **24**, 370-2, 2006.
- 10)Wonkeun Song et al: Use of boronic acid methods to detect the combined expression of plasmid mediated AmpC beta-lactamases and ESBLs in clinical isolates of *Klebsiella* spp., *Salmonella* spp. and *P. mirabilis*. Diagn. Microbiol. Infect. Dis., **57**, 315-18, 2007.
- 11)Prof. P. Nordmann: Evaluation de tests phenotypiques de detection de cephalosporinases integrant l'utilisation des disques de cloxacilline et d'acide boronique. Oct. 2006. Internal study.
- 12)Ruppé E. et al: First detection of the AmpC Class C1 AmpC beta-lactamase in *Citrobacter freundii* by a new, simple double-disk synergy test. J. Clin. Microbiol., **44**, 4204-4207, 2006.
- 13)Merlino J. et al: Boronic acid test for AmpC detection in *E.coli* and *K.pneumoniae* isolates. Study at Concord Hospital, Sydney (Australia) November 2007.
- 14)Mammeri H. et al: Molecular characterization of AmpC-producing *E.coli* clinical isolates recovered in a french hospital. J. Antimicrob. Chemother **61** 498-503,2008.
- 15)Aragon L. M. et al: Increase in  $\beta$ -lactam resistant *P.mirabilis* strains due to CTX-M and CMY-type as well as new VEB-inhibitor-resistant TEM-type  $\beta$ -lactamases. J. Antimicrob. Chemother **61**, 1029-32, 2008.
- 16)Ak. S. et al: Detection of plasmid mediated AmpC beta lactamases in clinical isolates by 3 different methods. Poster P-875, ECCMID, Barcelona, april 2008.
- 17)Thean Yen Tan et al: Evaluation of screening methods to detect plasmid mediated AmpC in *E.coli*, *K.pneumoniae* and *P. mirabilis*. Antimicrob. Agents Chemother **53**, 146-149, 2009.
- 18)Doi Y. et al: Reduced susceptibility to cefepime among *E.coli* clinical isolates producing novel variants of CMY-2 beta-lactamase. Antimicrob. Agents Chemother **53**, 3159-3164, 2009.
- 19)Adler H. et al: Detection of plasmid mediated AmpC beta-lactamases in Enterobacteriaceae lacking chromosomal AmpC using the Rosco double disc synergy test and the E-test. Poster ECCMID 2009.
- 20)Grigis A. et al: Detection of strains of ESBL and AmpC producing Enterobacteriaceae using synergy tests with Boronic acid and Cloxacillin (italian) 2009.
- 21) Rodriguez-Martinez JM et al: Extended-spectrum cephalosporinases in *P.aeruginosa*. Antimicrob. Agents Chemother **53**, 1766-71, 2009.
- 22) Ortega M. et al: Presumptive detection of AmpC beta-lactamases in *E.coli* and *Klebsiella* spp from UTI (spanish) presentation 85, XII, SEIMC Congress, Sevilla 2009.

### Inhibitor Resistant TEM Beta-lactamases (IRT)

Strains with this phenotype give patterns of antibiotic resistance similar to TEM 1 or 2 or SHV 1 beta-lactamases, but they are resistant to amoxicillin + clavulanate. IRT are found mainly in *E. coli* and *Klebsiella pneumoniae*.

Are R (resistant) to amoxicillin + clavulanate. Zone diameter for Amoxicillin+Clavulanate Neo-Sensitabs < 17 mm. Are S (generally susceptible) to cephalosporins: cefazolin, cefoxitin, cefotaxime.

#### References:

- 1) Reguera J.A. et al: Factors determining resistance to beta-lactams combined with beta-lactamase inhibitors. J. Antimicrob. Chemother. **27**, 569-75, 1991.
- 2) Lemozy. J. et al: First characterization of inhibitor resistant TEM (IRT) beta-lactamases in *Klebsiella pneumoniae* strains. Antimicrob. Agents Chemother., **33**, 2580-2, 1995.
- 3) Dubois V. et al: Molecular and biochemical characterization of SHV-56, a novel inhibitor resistant beta-lactamase from *K. pneumoniae*. Antimicrob. Agents Chemother **52**, 3792-94, 2008.
- 4) Miro E. et al: Prevalence of clinical isolates of *E.coli* producing IRT beta lactamases at a University hospital in Barcelona, Spain, over a 3-year period. Antimicrob. Agents Chemother **46**, 3991-94, 2002.
- 5) Venkataraman L. et al: Inhibitor resistant TEM (IRT) beta lactamases among amoxicillin/clavulanate resistant *E.coli* (ACREC). Presentation C2-4187, 48<sup>th</sup> ICAAC, 2008.
- 6) Badaracco A. et al: Variety of beta lactamases produced by amoxicillin/clavulanate resistant enterobacteria isolated in Buenos Aires, Argentina. Presentation C2-4186, 48<sup>th</sup> ICAAC, 2008.
- 7) Perez-Moreno M.O. et al: *K. pneumoniae* with reduced susceptibility to amoxicillin/clavulanate: coexistence of isolates producing IRT-11 and multiresistant qnrS2 positives producing OXA-1 associated to class 1 integrons. Presentation 155, XIII Reunion SEIMC, 2009.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Detection of beta lactamases

Page 1 / 13 - Document: 1.5.0

Carbapenemases. KPC+ MBL Confirm ID kit

### Carbapenemases

Carbapenemases are beta-lactamases that significantly hydrolyze at least imipenem and/or meropenem. Carbapenemases involved in acquired resistance are of Ambler classes A, B and D. They may be plasmid or chromosomally encoded.

Because several of these carbapenemases confer only reduced susceptibility to carbapenems in Enterobacteriaceae, they may remain underestimated, because they are not detected in the laboratory. Acquired carbapenemases are increasingly reported worldwide and consequently it is important to be able to detect them in the laboratory.

For many isolates with carbapenemases, the MICs of carbapenems are around the susceptible breakpoint making resistance difficult to detect - particularly with automated systems. Therefore special zone breakpoints are needed in first line screening

**Enterobacteriaceae** with reduced susceptibility to Imipenem 10 µg (zone < 23 mm or MIC > 1 µg/ml) or Meropenem 10 µg (inhibition zone < 25 mm or MIC > 0.5 µg/ml) or Ertapenem (zone ≤ 22 mm) on Mueller-Hinton Agar with McFarland 0.5 inoculum, should be suspected of possessing carbapenemases

**P. aeruginosa** with inhibition zones Imipenem 10 µg ( < 22 mm) or Meropenem 10 µg ( < 26 mm) should be suspected of possessing carbapenemase. Most isolates with KPC and GES enzymes are highly resistant to Ceftazidime. Ertapenem Neo-Sensitabs is the most sensitive indicator for possible carbapenemase, but in approximately 20% of cases other resistance mechanisms are involved (confirmation of carbapenemase with Modified Hodge Test is necessary). It is important to recognize small resistant colonies growing inside the Ertapenem disk zone.

### Carbapenemases classification (1)

Ambler classification	Enzymes	MICs µg / ml				Inhibited by		
		3rd gen cepha	AZT	IMP	MRP	CLAV	EDTA	Boronic acid
A	NmcA	S	4	≥ 16	2-8	± wk	no	yes
	Sme-1 to Sme-3	S	4-64	≥ 16	0.25-8	± wk	no	yes
	IMI-1 to IMI-2	S	S	≥ 64	4-32	+	no	yes
	KPC-1 to KPC-4	≥ 32	≥ 64	4→16	4→16	+ or wk	no	yes
	GES-2 to GES-5	≥ 32	16→R	0.25→16	0.5-16	+ or 0	no	yes
B Metallo-beta-lactamases	IMP 1-16	≥ 32	S→R	0.5-128	0.25→R	no	yes	no
	VIM 1-12	≥ 64	S→R	1→R	0.5→R	no	yes	no
	SPM-1	≥ 256	4	R	R	no	yes	no
	GIM-1	16-32	8-16	> 8	> 8	no	yes	no
	SIM-1	≥ 256	128	8-16	16	no	yes	no
D Oxacillinases	OXA 23-27	> 256	> 256	4-64	4-128	± wk	no	no
	OXA 40-48	S→R	S→R	2-64	0.25-64	wk	no	no
	OXA 54-55	S	S	4	0.25	wk	no	no
	OXA-60	S	R	0.5	2	no	no	no
	OXA-58	4-128	≥ 32	3-32	2→64	no	no	no

wk = weak

### References:

- 1) Nordmann P et al: Emerging carbapenemases in gram-negative aerobes. Clin Microbiol. Infect **8**, 321-331, 2002.

### Detection of acquired carbapenemases Ambler classes A and D

Class A carbapenemases are penicillinases with greater activity against imipenem than meropenem and they also give resistance to penicillins, cephalosporins and aztreonam.

Boronic acid in an inhibitor of class A carbapenemases and consequently synergy with meropenem or imipenem, is the best method to detect these enzymes (26,27,28,22).

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Detection of beta lactamases

## Carbapenemases. KPC+ MBL Confirm ID kit

Clavulanate is an inhibitor of class A carbapenemases and therefore synergy with imipenem may be useful to detect these enzymes (1,2,3,4,5).

The KPC family of enzymes confer greater resistance to third gen cephalosporins than to carbapenems (3,5).

KPC possessing *Enterobacter* spp. and *K. pneumoniae* were reported as falsely susceptible to carbapenems using automated systems (Vitek). MIC microdilution using standard inocula of  $10^4$  or  $10^5$  CFU/ml did not detect carbapenem resistance, while diffusion methods (E-test) using inocula of  $10^8$  CFU/ml detected resistance (5,7,12,18).

*K. pneumoniae* intermediate or resistant to ertapenem or meropenem should be considered resistant to all carbapenems (7). KPC possessing *E.coli* was identified in nine patients in New York. Three of the isolates possessed also ESBL: CTX M15 (19).

Pasteran et al (20) found that Boronic acid disks could be used to detect carbapenemases of type 2f (Class A) in Enterobacteriaceae. Class A producing strains showed synergy between Imipenem and Boronic acid disks (distance from edge to edge 6 mm). Strains showing zones of inhibition  $\leq 21$  mm with Imipenem 10 µg disks were screened with this test.

Carbapenemase IMI-2 is the first inducible and plasmid-encoded carbapenemase.

Please note that KPC detection may require screening multiple colonies, because carbapenemase susceptible strains may co-exist with resistant (21).

Class D carbapenemases correspond to the enzymes classified as OXA-types (oxacillinase activity). They hydrolyze imipenem and meropenem weakly and do not hydrolyze third gen cephalosporins and aztreonam (although MICs against the later drugs are often increased due to the presence of other beta-lactamases).

Clavulanate is a progressive inhibitor of most OXA carbapenemases, but not all. The synergy test (clavulanate and imipenem) may have value for the detection of these enzymes.

Clinicians should be aware of the potential for clinical failure (Class D, OXA-55 carbapenemase) when imipenem is used for treatment of serious infections caused by *S. algae* (9).

Fernandez et al report heteroresistance carbapenems in *Ac. baumannii* associated with Imipenems MIC 4-16 µg and the presence of OXA-58 (17).

Yilmaz et al (16) report oxacillinases (OXA-48) in 21 Enterobacteriaceae, mainly *K.pneumoniae*, but also in *E.coli* and *Enterobacter cloacae/aerogenes* in Turkey, and warns that oxacillinases (carbapenemases) are spreading in Enterobacteriaceae.

Castanheira et al (25) report the clonal dissemination of OXA-24 and OXA-58 producing *A. baumannii* in Houston, Texas.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Detection of beta lactamases

### Acquired carbapenemases Ambler class A and D

Ambler class.	Enzymes	MICs µg/ml					IMIPENEM + CLAV (synergy)	Mero-penem+ Boronic synergy	Organisms	Genetic Location
		PIPER	3rd gen cepha	AZTRM	IMI	MEROP				
A	NmcA	S	S (0.25-2)	4	≥ 16	2-8	± wk	yes	<i>E. cloacae</i>	Chromosomal
	Sme-1 to Sme-3	S	0.25-0.5	4-64	≥ 26	0.25-8	± wk	yes	<i>S. marcescens</i>	Chromosomal
	IMI-1	> 256	1-2	8	> 64	4	+	yes	<i>E. cloacae</i>	Chromosomal
	IMI-2	16→128	0.1-2	4-8	64	4-32	+	yes	<i>E. asburiae</i>	Plasmid
	KPC-1	> 128	≥ 32	> 64	16	16	+		<i>K. pneumoniae</i>	Plasmid
	<b>KPC-2</b>	≥ 64	≥ 8	> 16	8→16	≥ 16	+		<i>K. pneumoniae</i> /oxytoca. <i>Raoultella</i> (39) <i>Salmonella</i> <i>Enterobacter</i>	Plasmid
	<b>KPC-2</b>	.	≥ 256	> 256	256	256	(+)	yes	<i>P. aeruginosa</i>	Plasmid/Chrom.
	<b>KPC-3</b>	256	256	> 256	> 4	> 4	(+ +)		<i>K. pneumoniae</i>	Plasmid
	KPC-4				> 16	> 16	(+)		<i>Enterobacter</i>	
	KPC-5/6								<i>E. coli</i> , <i>P. mirabilis</i> , <i>Citrobacter</i> <i>Enterobacter</i>	
	<b>GES-2</b>	128	≥ 32	16	4→16	4-16	+	yes	<i>E. cloacae</i>	Plasmid, integron
D	<b>GES-3</b>	128	64→256	64	0.25	0.5	+	yes	<i>P. aeruginosa</i>	Plasmid
	GES-4	128	R	R	8	8	(+)	yes	<i>K. pneumoniae</i>	CEFOX R
	GES-5	R	R	R	8-32	8-32	+	yes	<i>K. pneumoniae</i>	Integron
	GES-11	R	> 256	> 256	4	8		yes	<i>P. aeruginosa</i> , <i>K. pneumoniae</i> <i>A. baumannii</i>	
	<b>OXA-23 to OXA-27</b>	> 256	> 256	> 256	4-64	4→128	± wk	no	<i>A. baumannii</i>	Chromosomal ± integron
	<b>OXA-40</b>	R	4→128	4→128	> 32	≥ 32	wk	no	<i>Ac. haemolyticus</i>	Plasmid
	<b>OXA-48</b>	8→R	S→R	S→R	2→64	0.25→64	wk	no	<i>K. pneumoniae</i>	Plasmid
	OXA-54	32	S	S	1	0.12	wk	no	<i>E. coli</i>	Not integron
	OXA-55	S	S	S	1-4	0.25	no wk	no	<i>Sh. putrefaciens</i>	Chromosomal
	<b>OXA-58</b>	256	4-128	≥ 32	2-32	2→64	no	no	<i>Sh. algae</i> (9)	Plasmid
	OXA-60	S	S	R	0.5	2	no	no	<i>A. baumannii</i>	Chromosomal
	OXA-62	S→R	S→R	S→R	2→64	64→128	no	no	<i>R. pickettii</i> <i>Pandorea</i> (10) <i>pnomenusa</i>	Chromosomal
D (8)	<b>OXA-23, 27, 49</b> (subgroup 1)			>16	8→32	8→>32		no	<i>Ac. baumannii</i>	Plasmid (only 23)
	<b>OXA-24, 25, 26, 40</b> (subgroup 2)	-	> 256		> 128	> 128	-	no	<i>Ac. baumannii</i>	Chromosomal
	OXA-51 + OXA-64-66, 68-71, 78-82-107 OXA-51-like (subgroup 3)				≥ 1	≥ 1		no	<i>Ac. baumannii</i>	Chromosomal plasmid
	<b>OXA-58</b> (subgroup 4)	R	R	> 16	4/16			no	<i>Ac. baumannii</i>	Plasmid (only 58)
	<b>OXA-143</b>	R	FEP4	.	32	32	-	no	<i>Ac. baumannii</i>	Plasmid

Bold = involved in outbreaks

**Procedure for KPC carbapenemases detection (Class A enzymes)**

Isolates giving negative metallo-beta-lactamase tests, may produce other carbapenemases. The most current are KPC enzymes isolated from *Enterobacteriaceae* (*K. pneumoniae*, *E. coli*, *Enterobacter spp.*, *P.mirabilis*) particularly *K. pneumoniae*, but also Sme, IMI, GES and Nunc A are found.

To detect these strains in rectal swab screening samples, direct plating on McConkey agar in the presence of Ertapenem Neo-Sensitabs and Imipenem Neo-Sensitabs may be useful.

Place one Boronic Acid Diatabs between one Ertapenem and one Imipenem Neo-Sensitabs (distance 6 mm from edge to edge).

Place one Cloxacillin Diatabs between Ertapenem and Imipenem Neo-Sensitabs - (6 mm from edge to edge). Perform Modified Hodge Test (MHT) with Ertapenem and Meropenem Neo-Sensitabs.


**Interpretation (Double disk synergy test)**

The following results will presumably indicate the presence of a KPC beta-lactamase:

- a) Negative metallo-beta-lactamases tests.
- b) Positive synergy test between Boronic Acid and the carbapenems (one or both).
- c) Negative synergy test between Cloxacillin and the Carbapenems (11)
- d) Positive synergy test between clavulanate (AMC) and carbapenems (one or both). Not always easy to see. Although isolates with ESBL + impermeability may give false positive results.
- e) Positive Modified Hodge Test.
- f) Sme, IMI, GES and Nunc A will show the same results as KPC, but the mentioned enzymes result in smaller zones around Imipenem compared to Ertapenem. With KPC enzymes zones around Imipenem and Ertapenem are similar.

**Combined disk test. KPC + MBL Confirm ID kit**

Apply Meropenem, Meropenem+DPA, Meropenem+Boronic, Meropenem+Cloxacillin on an inoculated plate.



K.pneumoniae PHA3 CL5761 KPC positive

**Interpretation (combined test)**

A Meropenem + Boronic inhibition zone  $\geq 5$  mm then Meropenem, Meropenem+DPA and Meropenem+Cloxacillin indicates a presence of a KPC enzyme (or other class A). Meropenem+Boronic and Meropenem+Cloxacillin inhibition zones  $\geq 5$  mm, than Meropenem and Meropenem+DPA indicates AmpC hyperproduction + porin loss, or efflux (30).

A Meropenem +DPA inhibition zone  $\geq 5$  mm than Meropenem, indicates the presence of a metallo- $\beta$ -lactamases (MBL).

**Please note**

Test only ertapenem-resistant strains. Ertapenem susceptible strains may provide a false positive result with Boronic Acid.

**Conclusion**

Reduced susceptibility to ertapenem, synergy between Boronic Acid and the carbapenems, and no synergy between Cloxacillin and the carbapenems is clearly indicative of KPC enzyme being present (or other class A enzymes). Isolates producing high level AmpC + impermeability can be detected by synergy between Cloxacillin and the carbapenems (11). Isolates producing ESBL + impermeability will show synergy between AMC and the carbapenems or cephalosporins.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Detection of beta lactamases

Page 6 / 13 - Document: 1.5.0

Carbapenemases. KPC+ MBL Confirm ID kit

### Procedure for Oxacillinase detection (Class D enzymes)

Strains producing oxacillinases will currently show zones of inhibition < 22 mm with Ertapenem and/or <25 mm with Meropenem Neo-Sensitabs. Most are resistant to Aztreonam. These enzymes are mainly found in *Acinetobacter baumannii* but also in Enterobacteriaceae (*K. pneumoniae*, *Enterobacter*) and *P. aeruginosa* although these are rare.

### **Interpretation**

The following results will presumably indicate the presence of oxacillinases:

- a) Negative metallo-beta-lactamase tests.
- b) Negative synergy test between Boronic acid/Cloxacillin and the carbapenems. (one or both).
- c) Negative (or weak positive) synergy test between clavulanate (AMC) and carbapenems (one or both)
- d) Positive Modified Hodge Test.

### **Summary detection of carbapenemases**

<b>Meropenem</b>	<b>MRP+DPA</b>	<b>MRP+BOR</b>	<b>MRP+Cloxa</b>
Metallo- $\beta$ -lactamases	Synergy	No synergy	No synergy
KPC	No synergy	Synergy	No synergy
AmpC impermeability	No synergy	Synergy	Synergy
Oxacillinases	No synergy	No synergy	No synergy

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 7 / 13 - Document: 1.5.0

## Detection of beta lactamases

## Carbapenemases. KPC+ MBL Confirm ID kit


### **Modified Hodge Test**

Is used to determine if resistance to carbapenems is caused by a carbapenemase. A MH agar plate (or a McConkey plate) is inoculated with the susceptible strain E. coli ATCC 25922 (Mc Farland 0.5, diluted 1/10) as for disk diffusion.

When testing Enterobacteriaceae, one Ertapenem Neo-Sensitabs and one Meropenem Neo-Sensitabs are applied onto the plate approx.30 mm apart from each other. For non-fermenters one Imipenem Neo-Sensitabs and one Meropenem Neo-Sensitabs are applied.

A suspension of the microorganism to be tested for carbapenemase is adjusted to Mc Farland 0.5 standard and a loop is used to make a heavy streak passing through the two carbapenem disks. Two more streaks are placed perpendicularly making a cross.

Thereafter incubation for 18-24 hours at 35-37 C. Alteration in the shape (indentation) of the zones of inhibition around the test organism is considered indicative of the presence of a carbapenemase (figure).



K.pneumoniae KPC positive

### **Limitations:**

- Not reliable for detection of SME from *S. marcescens*
- *P. mirabilis* swarming may give lecture problems.

### **References:**

- 1) Pottumarthy S et al: NmCA carbapenem-hydrolyzing enzyme in *E. cloacae* in North America. *Emerg Infect Dis* **9**, 999-1002, 2003.
- 2) Yigit H et al: Novel carbapenem-hydrolyzing beta-lactamase KPC-1 from a carbapenem-resistant strain of *K. pneumoniae*. *Antimicrob. Ag. Chemother.* **45**, 1151-61, 2001.
- 3) Smith-Moland E et al: Plasmid-mediated carbapenem-hydrolyzing beta-lactamase KPC-2 in *K. pneumoniae* isolates. *J. Antimicrob. Chemother.* **51**, 711-14, 2003.
- 4) Aubron C et al: Carbapenemase-producing Enterobacteriaceae, U.S. rivers. *Emerg. Infect. Dis.* **11**, 260-4, 2005.
- 5) Bratu S. et al: Detection of KPC carbapenem-hydrolyzing enzymes in *Enterobacter* spp from Brooklyn, N.Y. *Antimicrob. Ag. Chemother.* **49**, 776-8, 2005.
- 6) Lopez-Otsoa F et al: Endemic carbapenem resistance associated with OXA-40 carbapenemase among *A. baumannii* isolates from a hospital in Northern Spain. *J. Clin. Microbiol* **40**, 4741-3, 2002.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Detection of beta lactamases

Page 8 / 13 - Document: 1.5.0

Carbapenemases. KPC+ MBL Confirm ID kit

- 7) Bratu S. et al: Emergence of KPC-possessing *K. pneumoniae* in Brooklyn, N.Y.: epidemiology and recommendations for detection. *Antimicrob. Ag. Chemother.*, **49**, 3018-20, 2005.
- 8) Brown S., Amyes S.: OXA beta-lactamases in *Acinetobacter*: the story so far. *JAC*, **57**, 1-3, 2006.
- 9) Dong-Min K. et al: Treatment failure due to emergence of resistance to carbapenem during therapy for *Shewanella alga bacteremia*. *J.C.M.*, **44**, 1172-4, 2006.
- 10) Schneider I. et al: Novel carbapenem-hydrolyzing oxacillinase OXA-62 from *Pandora* sp. *Antimicrob. Ag. Chemother.*, **50**, 1330-35, 2006.
- 11) Villegas M.V. et al: First identification of *P. aeruginosa* isolates producing a KPC-type carbapenem hydrolyzing  $\beta$ -lactamase. *Antimicrob. Ag. Chemother.*, **51**, 1553-55, 2007.
- 12) Tenover F.C. et al: Carbapenem resistance in *K. pneumoniae* not detected by automated susceptibility testing. *Emerg. Inf. Dis.*, **12**, 1209-1213, 2006.
- 13) Sundin D.R.: Rapid screening for KPCs. Presentation D 1555, ICAAC september 2007, Chicago, USA.
- 14) Queenan A.M. et al: Carbapenemases: the versatile  $\beta$ -lactamases. *Clin. Microbiol. Reviews*, **20**, 440-458, 2007.
- 15) Lee K. et al: Effect of Oxgall on the Imipenem disk Hodge test in screening MBL- producing gram-negative bacilli. ECCMID, poster P-891, Barcelona 2008.
- 16) Yilmaz M. et al: Warning! Oxallinase-mediated carbapenem-resistance spreading in Enterobacteriaceae. Poster P1517, ECCMID. Barcelona, April 2008.
- 17) Fernandez F. et al: Association between heteroresistance to carbapenems in *A. baumannii*, production of beta-lactamases (AmpC and OXA). XIII Congress SEIMC presentation 12, Madrid, May 2008.
- 18) Marchaim D. et al: Isolation of Imipenem-R *Enterobacter* spp: emergence of KPC-2 carbapenemase, molecular characterisation, epidemiology and outcomes. *Antimicrob. Ag. Chemother* **52**, 1413-18, 2008.
- 19) Urban C. et al: Carbapenem-resistant *E.coli* harboring *K. pneumoniae* carbapenemase beta-lactamases associated with long term care facilities. *CID* **46** e127-e130, 2008.
- 20) Pasteran F. et al: Detection of carbapenemases of group 2f. Promising additive uses of the Boronic acid disk test. Poster 371-20836 XI, Argentinian Congress of Microbiology 2007.
- 21) Kohlerschmidt D.J. et al: *Klebsiella pneumoniae* KPC detection may require screening multiple colonies for accurate carbapenem resistance assessment. ASM, presentation C-028, June 2008.
- 22) Pasteran F. G. et al: *Klebsiella pneumoniae* KPC-2, Buenos Aires, Argentina. *Emerg. Infect. Dis.* **14**, 1178-1180, 2008.
- 23) Ya-Gang Chen et al: In vivo development of carbapenem resistance in clinical isolates of *E. aerogenes*, producing multiple beta-lactamases. *Int. J. Antimicrob. Ag.* **32**, 302-307, 2008.
- 24) Cuzon G. et al: Plasmid encoded carbapenem hydrolyzing  $\beta$ -lactamase OXA-48 in an Imipenem susceptible strain from Belgium. *Antimicrob. Ag. Chemother*, **52**, 3463-4, 2008.
- 25) Castanheira M. et al: Emergence and clonal dissemination of OXA-24 and OXA-48 producing *A.baumannii* strains in Houston, Texas: report from the SENTRY-ASP. *J. Clin. Microbiol.* **46**, 3179-80, 2008.
- 26) Doi Y. et al: Simple disk based method for detection of *K.pneumoniae* carbapenemase by use of Boronic acid compound. *J. Clin Microbiol.* **46**, 4076-83. 2008.
- 27) Giske C. et al: Personal communication 2009
- 28) Tsakris A. et al: First occurrence of KPC-2 possessing *K. pneumoniae* in a Greek hospital and recommendation for detection with Boronic acid disk tests. *J. Antimicrob. Chemother* **62**, 1257-60, 2008.
- 29) Pournaras S. et al: Clonal spread of KPC-2 carbapenemase producing *K. pneumoniae* in Greece. *JAC*, **64**. 348-52, 2009.
- 30) Smith Moland E. et al: Concerns about KPC screening and confirmatory tests. Presentation D-729, 49<sup>th</sup> ICAAC 2009.
- 31) Endimiani A. et al: Evaluation of the 2009 CLSI criteria for determining carbapenem susceptibility of *K. pneumoniae*. Presentation D-722, 49<sup>th</sup> ICAAC, 2009.
- 32) Anderson K. et al: Characterization of Enterobacteriaceae with a false positive Modified Hodge Test. Presentation D-719, 49<sup>th</sup> ICAAC, 2009.
- 33) Landman D. et al: Contribution of *OmpK 36* to carbapenem susceptibility in KPC-producing *K. pneumoniae*. *J. Med. Microbiol.* **58**, 1303-1308, 2009.
- 34) Kristo et al: Investigation of meropenem heteroresistance in KPC-2 producing *K.pneumoniae*. Presentation E-183. 49<sup>th</sup> ICAAC. 2009.

- 35) Tsakris A. et al: Evaluation of Boronic acid disk tests for differentiating KPC-possessing *K. pneumoniae* isolates in the clinical laboratory. *J. Clin. Microbiol.* **47**, 362-67, 2009.
- 36) Pasteran F. et al: Sensitive screening tests for suspected class A carbapenemase production in species of Enterobacteriaceae. *J. Clin. Microbiol.* **47**, 1631-39, 2009.
- 37) Weisenberg S.A. et al: Clinical outcomes of patients with *K. pneumoniae* KPC producing, after treatment with imipenem or meropenem. *Diagn. Micr. Infect. Dis.* **64**, 233-35, 2009.
- 38) Nicoletti A.G. et al: *K. pneumoniae* non-susceptible to carbapenems in Brazil. What are the mechanisms behind? Presentation C2-663, 49<sup>th</sup> ICAAC, 2009.
- 39) Castanheira M et al: First descriptions of KPC in *Raoultella* spp (*R. planticola* and *R. ornitholytica*): Report from the Sentry Antimicrobial Surveillance program. *J. Clin. Microbiol.* **47**, 4129-30, 2009.

#### **Detection of acquired Metallo-beta-lactamases (MBL)**

The worldwide spread of acquired metallo-beta-lactamases (MBL) in gram-negative aerobes is of great concern. MBL production in clinical isolates of key gram-negatives: *P. aeruginosa*, *E. cloacae*, *S. marcescens* and *K. pneumoniae* should be carefully monitored (5).

MLBs are classified into 5 major types: IMP, VIM, SPM, and GIM and SIM type enzymes. In Enterobacteriaceae only IMP and VIM enzymes have been found as yet. MBLs hydrolyze most beta-lactams (carbapenems and large spectrum cephalosporins), except aztreonam. This phenotype of multiple beta-lactam resistance and aztreonam susceptibility may be helpful for identification of these strains in the laboratory. If the strain is resistant to aztreonam it may be due to additional resistance mechanisms (efflux, other beta-lactamases, ESBL etc.). Their expression is not inducible.

The MBL enzymes are resistant to beta-lactamase inhibitors and susceptible to chelating agents like EDTA (2-MPA) and Dipicolinic acid (DPA).

Early detection of MBL-producing microorganisms is essential to prevent dissemination of these organisms. The enclosed tables, including strains of Enterobacteriaceae and Non-fermenters producing MBLs, show that MBL-producers (particularly in Enterobacteriaceae) may show low MIC values against carbapenems making it difficult for the laboratory to detect MBL-positive isolates.

Suspicious isolates (resistant to ceftazidime showing no synergy between clavulanate and third gen. cephalosporins and possibly showing reduced susceptibility to carbapenems) should be tested for carbapenemase activity using Imipenem, Meropenem and EDTA and Dipicolinic acid tests.

The first metallo-beta-lactamase producing strain of *E. coli* (in Spain) has been detected in Barcelona, using Imipenem+EDTA Neo-Sensitabs and E-test (3,8). The first metallo-beta-lactamase producing strain of *K. pneumoniae* was found in France (4).

MBL- producing gram-negatives have now emerged in Australia (15). The resistance gene bla-IMP4 appears highly mobile, this outbreak involved 5 different gram-negative genera. Diagnostic laboratories in Australia and other countries must be now in high alert, because early detection may limit the wide dispersal of MBL-genes.

Kyegong (27) and Miriagou (28) showed the efficiency of Dipicolinic acid (DPA) to detect metallo-β-lactamases in Enterobacteriaceae and non-fermenters. Miriagou found that the DPA/Imipenem synergy test was positive for all VIM-producing isolates of *Klebsiella*/Enterobacter and *P. mirabilis*, while EDTA based tests could not identify VIM-producing *P. mirabilis*

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Detection of beta lactamases

Page 10 / 13 - Document: 1.5.0

Carbapenemases. KPC+ MBL Confirm ID kit

### Acquired Metallo-beta-lactamases NON-FERMENTERS

MBL	3rd gen. Ceph MIC	AZT MIC µg/ml	IMP MIC µg/ml	MRP MIC µg/ml	Microorganisms	Genetic location
IMP 1-11	≥ 128	≤ 8/16	≥ 8	≥ 8	<i>Pseudomonas</i> spp. <i>Alcaligenes</i> spp. <i>Acinetobacter baumannii</i>	Chromosomal plasmid integron
IMP 12	≥ 128	32	32	128	<i>Pseudomonas putida</i>	
IMP 13-16	≥ 256	4-128	≥ 64	≥ 64	<i>Pseudomonas aeruginosa</i>	
VIM 1-3	R	S → R	2-128	1-128	<i>Achromobacter xylosoxidans</i> <i>Pseudomonas aeruginosa</i> <i>Pseudomonas putida</i> (VIM 2 and 4)	Chromosomal plasmid integron
VIM 4-11	> 256	S → R	32-256	32-256	<i>Acinetobacter baumannii</i> <i>Pseudomonas aeruginosa</i> , <i>A.baumannii</i>	
VIM 15-16	≥ 64	16-32	>128	≥128	<i>Pseudomonas aeruginosa</i>	
VIM-18	R	S → R	R	R	<i>Pseudomonas aeruginosa</i>	integron
SPM-1	≥ 256	4	R	R	<i>Pseudomonas aeruginosa</i>	Plasmid (not integron)
GIM-1	16 → 32	8-16	> 8	> 8	<i>Pseudomonas aeruginosa</i>	Integron
SIM-1	≥ 256	128	8-16	16	<i>Acinetobacter baumannii</i>	Integron
IND1-6	1-32-128	32-128	4-32-128	4-16-128	<i>Chryseobact indologenes</i>	Chromosomal (23)
AIM-1					<i>Pseudomonas aeruginosa</i>	

MBL are not inhibited by clavulanate, but are inhibited by EDTA or DPA

### Acquired Metallo-beta-lactamases ENTEROBACTERIACEAE

MBL	3rd gen. Ceph MIC	AZT MIC µg/ml	IMP MIC µg/ml	MRP MIC µg/ml	Microorganisms	Genetic location
IMP-1	≥ 32	< 0.5	2	0.5	<i>E. coli</i>	Integron  Plasmid
IMP-1	≥ 32	0.5 → R	4-128	4-128	<i>S. marcescens</i> , <i>K. pneumoniae</i> , <i>K. oxytoca</i> , <i>E. cloacae</i> / <i>E. aerogenes</i> , <i>Cit. freundii</i> , <i>P. rettgeri</i> , <i>M. morganii</i> , <i>Shigella flexneri</i>	
IMP-3	64	0.5	1	.	<i>Citrobacter youngae</i>	
IMP-4	256	.	3	6	<i>E. coli</i>	
IMP-6	> 128	0.25	2-8	64	<i>Serratia marcescens</i>	
IMP-6	> 128	128	32	> 128	<i>Enterobacter cloacae</i> , <i>Klebsiella pneumoniae</i> , <i>S.marcescens</i>	
IMP-8	R	S → R	0.5-8	0.25-4		
VIM-1	R	8-128	8-32	2-32	<i>E. coli</i> , <i>P. mirabilis</i> (integron) <i>C. koseri</i> , <i>K. oxytoca</i>	Plasmid
VIM-1	16-128	S → R	1-64	1-32	<i>Klebsiella pneumoniae</i> , <i>E. cloacae</i>	Plasmid (integron)
VIM-2	≥ 32	S → R	≥ 1	0.5 → > 2	<i>Citrobacter freundii</i> / <i>E.cloacae</i>	Plasmid
VIM-2	≥ 128	32	16-64	8-64	<i>Serratia marcescens</i> , <i>P. rettgerii</i>	Integron
VIM-2	8	16	4	0.1 (S)	<i>Klebsiella oxytoca</i>	Plasmid (integron)
VIM-4	≥ 32	4 → R	2-4	0.5-1	<i>K. pneumoniae</i> / <i>E. cloacae</i>	Plasmid
VIM-12	≥ 128	16	8	4	<i>K. pneumoniae</i>	Plasmid (16)
VIM-12	> 32	1	1	0.25	<i>E. coli</i>	Plasmid (22)
VIM-2 + GES7					<i>E. coli</i>	Integron
KHM-1	R	0.25	2	4	<i>C. freundii</i>	Plasmid
NDM-1	R				<i>K. pneumoniae</i> , <i>E.coli</i> , <i>C.freundii</i> (31,33)	Plasmid

MBL are not inhibited by clavulanate, but are inhibited by EDTA or DPA.

**Procedure for metallo-beta-lactamase (MBL) detection**

Some resistance profiles may suggest MBL production, for example:

a) *Pseudomonas aeruginosa*, *Pseudomonas* spp. and *Acinetobacter* spp.

All isolates non-susceptible to carbapenems and resistant to either ticarcillin, ticarcillin+clavulanate or ceftazidime should be tested for MBL production.

b) *Enterobacteriaceae*

For *E. coli*, *Klebsiella* spp., *P. mirabilis*, *Salmonella* spp. and *Shigella* spp.: All carbapenem S-I-R isolates that are resistant to cefoxitin and amoxicillin+clavulanate and are non-susceptible to ceftazidime (inhibition zone < 18 mm) should be tested for MBL production. In all other cases all isolates are non-susceptible to carbapenems (18).

**Enterobacteriaceae**

Apply one Dipicolinic Acid Diatabs (DPA) on an inoculated Mueller Hinton (MH) plate. Apply one Meropenem Neo-Sensitabs and one Ertapenem Neo-Sensitabs onto the plate on either side of the DPA, 5mm from the DPA (edge to edge). Apply Imipenem 10 µg + EDTA (IM10E) on an inoculated MH plate. Apply one Imipenem 10 µg Neo-Sensitabs.

**Non-fermenters**

Apply one DPA Diatabs on the MH plate. Apply one Imipenem Neo-Sensitabs and one Meropenem Neo-Sensitabs at either side of the DPA, 5 mm from the DPA (edge to edge). Apply Imipenem + EDTA (IM10E) on the inoculated MH plate. Apply one Imipenem Neo-Sensitabs.

**Interpretation (Double disk synergy test)**

The use of two chelating agents EDTA and DPA will enhance the detection of metallo-β-lactamases (MBL) in the clinical laboratory. A key hole or ghost zone between carbapenems (one or more) and Dipicolinic Acid indicates the presence of an MBL.

- An Imipenem + EDTA (10+750 µg) zone 7mm larger than Imipenem 10 µg indicates the presence of a Metallo-Beta-Lactamase.



A key hole or ghost zone between carbapenems (one or more) and Dipicolinic Acid indicates the presence of a MBL.

**Combined disk test. KPC and MBL Confirm ID kit**

Apply Meropenem, Meropenem+DPA on an inoculated MH plate. Interpretation: A Meropenem+DPA inhibition zone  $\geq 5$  mm than Meropenem alone indicates the presence of a metallo-beta-lactamase.



*P.aeruginosa* FN 8173 metallo- $\beta$ -lactamase positive

**Please note:**

Most MH agar brands contain physiological levels of Zn<sup>++</sup> ions and should be used for carbapenem testing. Iso-Sensitest agar has low levels of zinc ions and may give false susceptibility results for carbapenems in the presence of MBL. Strains of *Acinetobacter baumannii* producing certain oxacillinases may give a false positive metallo- $\beta$ -lactamase test result.

**References:**

- 1) Walsh T.R. et al: Evaluation of new E-test for detecting metallo-beta-lactamases in routine clinical testing. J. Clin. Microbiol., **40**, 2755-59, 2002.
- 2) Jong D. et al: Imipenem-EDTA disk method for differentiation of metallo-beta-lactamase producing clinical isolates of *Pseudomonas* spp. and *Acinetobacter* spp. J. Clin. Microbiol., **40**, 3798-3801, 2002.
- 3) Larrosa M. N. et al: *E. coli* multi-resistente productora de una metalo-beta-lactamasa (Spanish). SEIMC Congress, Bilbao, 16-19th May, 2004.
- 4) Lartigue MF et al: First detection of a carbapenem-hydrolyzing metalloenzyme in an Enterobacteriaceae isolate in France. Antimicrob. Agents Chemother., **48**, 4929-30, 2004.
- 5) Castanheira M et al: Molecular characterization of a beta-lactamase gene blaGIM-1, a new subclass of metallo-beta-lactamases. Antimicrob. Agents Chemother., **48**, 4654-61, 2004.
- 6) Jing-Jou Yan et al: Comparison of the double-disk, combined disk and E-test methods for detecting metallo-beta-lactamases in gram-negative bacilli. Diagn. Microbiol. Infect. Dis., **49**, 5-11, 2004.
- 7) Walsh T.R. et al: Metallo-beta-lactamases: The quite before the storm? Clin. Microbiol. Reviews, **18**, 306-325, 2005.

- 8) Tortola M.T. et al: First detection of a carbapenem-hydrolyzing metalloenzyme in two Enterobacteriaceae isolates in Spain. Antimicrob. Agents Chemother., **49**, 3492-4, 2005.
- 9) Chu Y.W. et al: EDTA susceptibility leading to false detection of metallo-beta-lactamase in *P. aeruginosa* by E-test and an Imipenem-EDTA disk method. Intl. J. Antimicrob. Agents, **26**, 340-341, 2005.
- 10) Lincopan N. et al: First isolation of metallo-β-lactamase producing multiresistant *K. pneumoniae* from a patient in Brazil. J. Clin. Microbiol., **43**, 516-9, 2005.
- 11) Pitout J.D.D. et al: Detection of *P. aeruginosa* producing metallo-β-lactamases in a large centralized laboratory. J. Clin. Microbiol., **43**, 3129-35, 2005.
- 12) Marqué S. et al: Regional occurrence of plasmid-mediated carbapenem-hydrolyzing oxacillinase OXA-58 in *Acinetobacter* spp. in Europe. J. Clin. Microbiol., **43**, 4885-8, 2005.
- 13) Choi Y.S. et al: Evaluation of Imipenem Disk Hodge Test and double disk synergy tests to detect SIM-1 type metallo-β-lactamase-producing *Acinetobacter* isolates. Presentation D-1705, 45th ICAAC, 2005.
- 14) Lee Kyungwon et al: Novel acquired metallo-β-lactamase gene blaSIM-1 in class 1 integron from *Ac. baumannii* clinical isolates from Korea. Antimicrob. Agents Chemother., **49**, 4485-91, 2005.
- 15) Peleg A.Y.: Dissemination of the metallo-β-lactamase gene blaIMP-4 among gram-negative pathogens in a clinical setting in Australia. CID, **41**, 1549-56, 2005.
- 16) Pournaras et al: VIM-12 a novel plasmid-mediated metallo-β-lactamase from *K. pneumoniae* that resembles a VIM-1/VIM-2 hybrid. Antimicrob. Agents Chemother., **49**, 5153-6, 2005.
- 17) Berges L. et al: Prospective evaluation of Imipenem/EDTA combined disk and E-test for detection of MBL producing *P. aeruginosa*. Clin. Microbiol. Infect., **12**, Suppl. 4, P1451, 2006.
- 18) Gornaglia G. et al: Metallo-β-lactamases as emerging resistance determinants in Gram-negative pathogens; open issues. Intl. J. Antimicrob. Agents, **29**, 380-88, 2007.
- 19) Franklin C. et al: Phenotypic detection of carbapenem-susceptible metallo-β-lactamase-producing Gram-negative bacilli in the clinical laboratory. J. Clin. Microbiol., **44**, 3139-3144, 2006.
- 20) Moland E.S. et al: Prevalence of newer β-lactamases in gram-negative clinical isolates collected in the U.S. from 2001 to 2002. J. Clin. Microbiol., **44**, 3318-3324, 2006.
- 21) Deshpande L.M. et al: Emergence of serine carbapenemases (KPC and SME) among clinical strains of Enterobacteriaceae in the U.S. Medical Centers: Report from the MYSTIC Program (1999-2005).
- 22) Iconomidis A. et al: First occurrence of an *E. coli* clinical isolate producing the VIM1/VIM2 hybrid metallo-β-lactamase VIM12. Antimicrob. Agents Chemother., **51**, 3038-39, 2007.
- 23) Perilli M. et al: Identification and characterisation of a new metallo-β-lactamase from a clinical isolate of *Chryseobacterium indologenes*. Antimicrob. Agents Chemother., **51**, 2988-2990, 2007.
- 24) Anderson K.F. et al: Evaluation of methods to identify the *K. pneumoniae* carbapenemase in Enterobacteriaceae. J. Clin. Microbiol., **45**, 2723-25, 2007.
- 25) Kim Soo-Young et al: Convenient test using a combination of chelating agents for the detection of metallo-beta-lactamases in the clinical laboratory. J. Clin. Microbiol., **45**, 2798-2801, 2007.
- 26) Giakkoupi P. et al: Supplementation of growth media with Zn<sup>2+</sup> facilitates detection of VIM-2 producing *P. aeruginosa*. JCM **46**, 1568-9, 2008.
- 27) Kyeong S.S. et al: Dipicolinic acid based disk methods for detection of metallo-β-lactamase producing *Pseudomonas* spp and *Acinetobacter* spp. Diagn. Microbiol. Infect Dis **62**, 102-105, 2008.
- 28) Miriagou V et al: Improved detection of VIM-1 production in *Proteus mirabilis* by an Imipenem-dipicolinic acid double disk synergy test. In press.
- 29) Multiresistant hospital bacteria linked to India and Pakistan. HPS Weekly Report News, July 8<sup>th</sup>, 2009.
- 30) Carvalho-Assef APD et al: First report of the carbapenem-resistant IMP-16 producing *Pseudomonas stutzeri*. Presentation C2-634, 49<sup>th</sup> ICAAC, 2009.
- 31) Pasteran F. et al: Spread of multiple clones of metallo-β-lactamases (MBL) producing *Pseudomonas putida*. Presentation C2-636, 49<sup>th</sup> ICAAC, 2009.
- 32) Treviño M. et al: Carbapenem-resistant *E. cloacae* and the emergence of metallo-β-lactamases producing strains in a third level hospital (Santiago de Compostela, Spain). Eur. J. Clin. Microbiol. Infect. Dis **28**, 1253-58, 2009.
- 33) Dongeun Yong et al: Characterization of a new metallo-β-lactamase gene NDM-1 and a novel erythromycin esterase gene carried on a unique genetic structure in *K. pneumoniae* sequence type 14 from India. Antimicrob. Agents Chemother **53**, 5046-5054, 2009.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 1 / 1 - Document: 1.6.0

## Detection of beta lactamases

## Detection of susceptible beta lactamases in one strain

### Detection of multiple beta-lactamases in one strain

Diagnostic problems posed by coexistence of different classes of beta-lactamases in a single bacterial isolate could be solved by the combined use of various phenotypic detection methods. See below example with multiresistant *K. pneumoniae* from Taiwan and USA.

	Neo-Sensitabs				
K. pneumoniae producing:	Cefoxitin	Cefepime	Ceftazidime+ Clavulanate or Cefepime+ Clavulanate synergy	D.P.A. + Meropenem or Imipenem+ EDTA synergy	Boronic acid Cefotaxime/ Ceftazidime or Cloxacillin Cefotaxime/Ceftazidime synergy
AmpC	R	S	negative	negative	POSITIVE
ESBL	S (V)	I / R	POSITIVE	negative	negative
Metallo-β-lactamase	R	I / R	negative	POSITIVE	negative (V)
KPC*	R	I / R	negative	negative	Negative (V)
AmpC + ESBL	R	I / R	POSITIVE	negative	POSITIVE
AmpC + metallo-β-lactamases	R	I / R	negative	POSITIVE	POSITIVE
AmpC + ESBL + metallo-β-lactamases	R	I / R	POSITIVE	POSITIVE	POSITIVE

\*KPC shows synergism between Boronic acid and Imipenem/Meropenem, but **no synergism** between Cloxacillin and Imipenem/Meropenem.

### References:

- 1) Jing-Jou Yan et al: Complexity of *Klebsiella pneumoniae* isolates resistant to both cephamycins and extended spectrum cephalosporins at a teaching hospital in Taiwan. J. Clin. Microbiol., **42**, 5337-40, 2004.
- 2) Smith Moland E. et al: *Klebsiella pneumoniae* isolate producing at least 8 different beta-lactamases, including AmpC and KPC beta lactamases. Antimicrob. Agents Chemother., **51**, 800-801, 2007.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 1 / 2 - Document: 1.7.0

## Detection of beta lactamases

Detection of beta-lactam resistance phenotypes in Enterobacteriaceae

### Detection of $\beta$ -lactam Resistance Phenotypes in Enterobacteriaceae

	AMP	AMC	CLOT	CAZ	FEP	FOX	CAZ+CLAV FEP+CLAV	IMI	IMI+ EDTA MRP+DPA	CLOXA Boronic	Comments
1) E. coli / P. mirabilis / Salmonella spp. / Shigella spp. / Klebsiella spp., /C. diversus											
Penicillinase low (E. coli)	I/R	S/I	S	S	S	S		S			
Penicillinase high	R	R	R	S	S	S		S			
Cephalosporinase low	S/I	S/I	I/R	S	S	S		S			
AmpC high / plasmid	R	R	R	I/R	S	I/R		S		synergy CAZ or CTX	FOX antagonism with CAZ or FEP indicates inducible
AmpC + ESBL	R	R	R	R	S/R	R	Synergy	S		synergy CAZ or CTX	
IRT	R	R	S/I	S	> CAZ (zone)			S			
Oxacillinase	R	R	S/I	S	≤ CAZ (zone)			S			
Penicillinase + Cephalosporinase	R	R	R	S/I	S	I/R		S			
Chromosomal K-1 high (Klebsiella oxytoca)	R	I/R	R	S	S/I	S	false synergy FEP+CLAV	S			
ESBL	R	V/R	I/R	V/R	V/R	S/R	synergy	S			
Metallo-β- lactamase	R	R	R	I/R	I/R	R	-	S/R	Synergy		AZT S <sup>R</sup>
Carbapenemase class A KPC	R	S/R	R	I/R	I/R	R		I/R (syn CLAV)		No synergy CLOXA with MRP	Synergy: IMI OR MRP with BORONIC
2) Enterobacter spp. / C. freundii											
Penicillinase	R	R	R	S	S	R		S			
AmpC derepressed	R	R	R	R	S	R		S		syn CTX syn CAZ	
AmpC + ESBL	R	R	R	R	S/R	R	Synergy	S		syn CTX syn CAZ	
ESBL, Metallo-β- lactamase and Carbapenemases	As for E. coli etc.										
3) M. morganii / Providencia spp. / Serratia spp.											
Penicillinase	R	R	R	S	S	R		S			
Amp C derepressed	R	R	R	R	S	I/R		S		syn CTX syn CAZ	
ESBL, Metallo-β- lactamase and Carbapenemases	As for E. coli etc.										
4) P. vulgaris / P. penneri											
Penicillinase	R	S	R	S	S	S		S			

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 2 / 2 - Document: 1.7.0

## Detection of beta lactamases

## Detection of beta-lactam resistance phenotypes in Enterobacteriaceae

Chromosomal β-lactamase derepressed	R	S	R	R	S	S		S			
-------------------------------------	---	---	---	---	---	---	--	---	--	--	--

AMP; Ampicillin, AMC; Amoxicillin – clavulanic acid, CLOT; Cephalothin, CAZ; Ceftazidime, FEP; Cefepime, FOX; Cefoxitin, CAZ+CLAV; Ceftazidime + clavulanic acid, FEP+CLAV; Cefepime + clavulanic acid, IMI; Imipenem, IMI+EDTA; Imipenem+EDTA, Cloxa ; Cloxacillin, AZT; Aztreonam; CTX=Cefotaxime;MRP=Meropenem.

ORIGINAL ROSCO DOCUMENT

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 1 / 2 - Document: 1.8.0

## Detection of beta lactamases

## Detection of beta-lactam resistance phenotypes in non-fermenters

### Detection of $\beta$ -lactam resistance phenotypes in non-fermenters

	TIC	TCC	PIP	CAZ	FEP	IMI	MRP	CAZ+CLAV FEP+CLAV	IMI+EDTA MRP+DPA	Remarks
<b>1) <i>Pseudomonas aeruginosa</i></b>										
Penicillinase (TEM 1-2, PSE 1-4)	R	S/I	I/R	S	S	S	S			
Amp C partial derep.	S/I	any	I	I	S/I	S	S			
Amp C derepressed	I/R	R	I/R	R	I/R	S	S			CAZ < AZT TIC > CTAX
Oxacillinase	I/R	I/R	any	S	S/R ( $\leq$ CAZ)	S	S			FEP $\leq$ CAZ TIC > CTAX
OXA-31, OXA 1, 4.	I/R	I/R	any	S	I/R	S	S			
Mex XY-OprM Efflux I/R (1)				S	I/R	S/R	S			CAZ $\geq$ FEP zone
PER 1-2, VEB-1 (ESBL)	R	S	S/I	R	R	S (PER-1) I/R (VEB-1)	S (PER-1) I/R (VEB-1)	Synergy		Phenotype PIP S, CAZ R indicates PER-1 enz.
ESBL	I/R	any	R	I/R	R	S	S	Synergy		
Metallo- $\beta$ -lactamase	R	R	R	R	R	I/R	I/R		Synergy	AZT, S <sup>R</sup>
Carbapenemase Class A KPC	R	S/I	R	I/R	I/R	I/R	I/R			AZT, S <sup>R</sup> Syn. MRP+BOR
Increased efflux	I/R	I/R	I/R	I/R	I/R	S <sup>R</sup>	I			CAZ > AZT
Loss of Opr porin	S	S	S	S	S	I/R	S/I			
<b>2) <i>Stenotrophomonas maltophilia</i></b>										
Beta-lactamase L-1	R	R	R	R	R	R	R	-	Synergy	AZT, S
Beta-lactamase L-2	-	S	I	R	R	S	-	-	-	AZT, R
Beta-lactamase L-1+L-2	R	R	R	R	R	R	R	-	Synergy	AZT, R
<b>3) <i>Acinetobacter baumannii</i></b>										
	TIC	PIP	PTZ	CAZ	FEP	IMI	MERO			
Penicillinase (TEM 1-2)	R	R	S	S	S	S	S			
Oxacillinase (OXA 21, 37)	R	R	S	S	S	S	S			FEP $\leq$ CAZ
Amp C partial derep.	I	I	S/I	I	S/I	S	S			
Amp C derepressed	R	R	I/R	R	I/R	S	S			
PER-1, VEB-1 (ESBL)	R	any	S	R	R	S	S	Synergy		Syn. TIC+CLAV Syn. PIP+TAZO
ESBL	R	R	any	R	R	S	S	Synergy		
Metallo- $\beta$ -lactamase	R	R	R	R	R	I/R	I/R		Synergy	AZT, S <sup>R</sup>

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 2 / 2 - Document: 1.8.0

## Detection of beta lactamases

## Detection of beta-lactam resistance phenotypes in non-fermenters

Carbapenemase Class D (Oxa 23-27, 40, 51, 58)	R	R	I/R	I/R	I/R	S/I/R	S/I/R			AZT, S <sup>R</sup>
Loss of porins	S	S	S	S	S	I/R	S/I			

TIC; Ticarcillin, TCC; Ticarcillin+clavulanic acid, PIP; Piperacillin, PTZ; Piperacillin+tazobactam, CAZ; Ceftazidime, FEP; Cefepime, IMI; Imipenem, MRP; Meropenem, CAZ+CLAV; Ceftazidime + clavulanic acid, FEP+CLAV; Cefepime + clavulanic acid, IMI+EDTA; Imipenem + EDTA, AZT; Aztreonam; BOR=Boronic acid.

### References:

1)Hocquet D. et al: Involvement of the MexXY-OprM Efflux system in emergence of Cefepime resistance in clinical strains of Ps. aeruginosa. Antimicrob. Ag. Chemother., 50, 1347-51, 2006.

### Screening of 16S rRNA Methylases (HLR to Aminoglycosides)

Unlike aminoglycoside-modifying enzymes that vary in their substrate profile, the acquired 16S rRNA methylases confer high level resistance (HLR) to almost all clinically important aminoglycosides. They have been identified in several nosocomial pathogens, including *P. aeruginosa*, *Serratia marcescens*, *E. coli*, *P. mirabilis*, *K. pneumoniae* and *Acinetobacter* spp., *Enterobacter cloacae*, *Citrobacter freundii* (8,9,11).

These enzymes (RmtA, RmtB, RmtC, ArmA) are capable of conferring very high levels of resistance (MIC > 512 µg/ml) against amikacin, gentamicin, isepamicin, netilmicin and tobramycin, while apramycin, neomycin and streptomycin are not affected. The responsible genes *armA*, *rmtA*, *rmtB*, *rmtC* and *rmtD* are located in self-transmissible plasmids (7).

#### Screening method

A high-level amikacin resistance (MIC > 512 µg/ml) corresponding to no-zone of inhibition around Amikacin 40 µg Neo-Sensitabs may be used as a marker for screening the 16S rRNA methylase producing strains.

The diffusion test is performed on MH-agar using a 0.5 McF inoculum and incubation at 35-37 °C overnight.

Strains of Enterobacteriaceae and non-fermenters (*P. aeruginosa* and *Acinetobacter* spp.) showing no-zone of inhibition around Amikacin 40 µg Neo-Sensitabs should be suspected of possessing 16S rRNA methylases.

#### a) Enterobacteriaceae

16S rRNA methylase **positive** strains will show:

Amikacin:	No zone of inhibition	Resistant
Gentamicin:	No zone of inhibition	Resistant
Netilmicin:	No zone of inhibition	Resistant
Tobramycin:	No zone of inhibition	Resistant
Neomycin 120 µg:	Zone of inhibition ≤ 20 mm or no zone	Resistant
Apramycin 100 µg:	Zone of inhibition ≥ 20 mm (6) (S)	Resistant

#### NpmA-enzyme

#### b) Non-fermenters

16S rRNA methylase **positive** strains will show:

Amikacin:	No zone of inhibition
Gentamicin:	No zone of inhibition
Netilmicin:	No zone of inhibition
Tobramycin:	No zone of inhibition
Neomycin 120 µg:	No zone or small zone
Streptomycin 100 µg:	Small zone in most cases

Galimand et al (5) found in 12 clinical isolates of Enterobacteriaceae the *armA* gene associated with ESBL beta-lactamase CTXM-3 (cefotaxime zone < ceftazidime zone) on a conjugative plasmid.

Bogaerts et al (9) investigated the presence of 16S rRNA methylase mediated high level resistance to aminoglycosides in clinical isolates of Enterobacteriaceae from 2 University Hospitals in Belgium. They screened for HLR to gentamicin, tobramycin and amikacin resistance and deleted by PCR, *armA* genes in 18 *K. pneumoniae*, *E. coli*, *E. aerogenes*, *E. cloacae*, and *C. amalonaticus*, whereas *rmtB* was detected in a single *E. coli* isolate. These strains were susceptible to Apramycin and Neomycin Neo-Sensitabs (except 2 strains). All 16S rRNA methylase positive strains produced ESBL's predominantly type CTX-M3 (13).

The concomitant presence of 16S rRNA methylase genes (*armA* or *rmtB*) and beta lactamase CTX-M among amikacin-resistant ESBL-producing *K. pneumoniae* isolates are widely spread in Taiwan (12).

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 2 / 2 - Document: 2.1.0

## Screening of 16S + rRNA methylases

## Screening of 16S rRNA methylases (HLR to Aminoglycosides)

The emergence of 16S rRNA methylases in Enterobacteriaceae and non-fermenters (*P. aeruginosa*, *Acinetobacter* spp.) in strains that already are ESBL positive, may result in the spread of multidrug-resistant isolates producing both ESBLs and 16S rRNA methylases becoming an important clinical problem.

Wachino et al (10) describes a new plasmid-mediated 16S rRNA methyltransferase NpmA isolated from *E. coli* and providing total aminoglycoside resistance (including apramycin, neomycin and streptomycin).

### References:

- 1) Galimand M et al: Plasmid-mediated high-level resistance to aminoglycosides in Enterobacteriaceae due to 16S rRNA methylation. *Antimicrob. Ag. Chemother.*, **47**, 2565-2571, 2003.
- 2) Yohei Doi et al: Plasmid-mediated 16S rRNA methylase in *Serratia marcescens* conferring HLR to aminoglycosides. *Antimicrob. Ag. Chemother.*, **48**, 491-6, 2004.
- 3) Jing-Jou Yan et al: Plasmid-mediated 16S rRNA methylases conferring HL aminoglycoside resistance in *E. coli* and *K. pneumoniae* isolates from 2 Taiwanese hospitals. *J. Antimicrob. Chemother.*, **54**, 1007-1012, 2004.
- 4) Kunkazu Yamane et al: Global spread of multiple aminoglycoside resistance genes. *Emerg. Infect. Dis.*, **11**, 951-953, 2005.
- 5) Galimand M. et al: Worldwide disseminated armA aminoglycoside resistance methylase gene is borne by composite transposon Tn1548. *Antimicrob. Ag. Chemother.*, **49**, 2949-53, 2005.
- 6) González-Zorn B. et al: Genetic basis of dissemination of armA. *J. Antimicrob. Chemother.*, **56**, 583-5, 2005.
- 7) Wachino Jun-idri et al: Novel plasmid mediated 16S rRNA methylase RmtC found in a *P. mirabilis* isolate demonstrating extraordinary high-level resistance against various aminoglycosides. *Antimicrob. Ag. Chemother.*, **50**, 178-184, 2006.
- 8) Lee S. et al: Dissemination of 16S rRNA methylase in AmpC producing *E. cloacae*, *C. freundii* and *S. marcescens* in Korea. *Clin Microbiol. Infect.*, **12**, Suppl. 4, P1274, 2006.
- 9) Bogaerts P. et al: Emergence of ArmA and RmtB aminoglycoside resistance rRNA methylases in Belgium. *JAC*, **59**, 459-464, 2007.
- 10) Wachino Jun-idri et al: *Antimicrob. Ag. Chemother* **51**, 4401-09, 2007.
- 11) Fritsche T. R. et al: Detection of methyltransferases conferring high-level resistance to aminoglycosides in Enterobacteriaceae from Europe, North America and Latin America. *Antimicrob. Agents Chemother* **52**, 1843-1845, 2008.
- 12) Ling Ma. et al: Widespread dissemination of aminoglycoside resistance genes armA and rmtB in *K. pneumoniae* isolates in Taiwan, producing CTX-M type ESBL-s. *Antimicrob. Ag. Chemother* **53**, 104-111, 2009.
- 13) Kim S. et al: 16S rRNA methylase genes in amikacin resistant ESBL producing *E. coli* and *K. pneumoniae* isolates at University Hospital in Korea. Abstract, Presentation C2-701, 49<sup>th</sup> ICAAC, 2009.

### Screening for plasmid-mediated quinolone resistance

The plasmid gene responsible for quinolone resistance (*qnrA*, *qnrB*, *qnrC*, *qnrD*, *qnrS*, *aac(6')*lb-cr and *QepA*) is carried on class 1 integrons of the In 4 family, an efficient mechanism for rapid horizontal and vertical dissemination of antibiotic resistance determinants among bacteria.

The plasmid mediated mechanisms have led to resistance to almost all clinical important antimicrobials, such as  $\beta$ -lactams, aminoglycosides, macrolides, phenicols, sulphonamides and trimethoprim.

The identification in the US of *qnr* in clinical strains of *K. pneumoniae* isolates besides producing plasmidic  $\beta$ -lactamases and ESBL's (7) and its discovery in strains of *E. coli* from Southeast Asia and *Salmonella* in Hong Kong indicates the emergence of this new mechanism of quinolones resistance in clinical strains.

It is important to indicate that a significant relation exists between quinolone resistance and resistance to 3. gen. cephalosporins (co-resistance): ESBL and/or plasmid mediated AmpC (6,14).

Poirel et al (10) have shown in vivo selection of fluoroquinolone-resistant *E. coli* isolates expressing plasmid mediated quinolone resistance and ESBL and physical linkage between ESBL and *qnrA*-encoding genes in the same integron.

Although *QnrA*, *QnrB*, *QnrS* produce low levels of quinolone resistance, it facilitates selection for a high level of quinolone resistance.

*QnrB*, another plasmid-mediated gene for quinolone resistance has been discovered in plasmids encoding the ESBL: CTA-15 from a *K. pneumoniae*. These strains show low-level resistance to quinolones and MIC of 16  $\mu$ g/ml towards nalidixic acid, and show similar multiresistance phenotypes as *qnrA* containing strains (11).

Lavigne et al (12) screened for *qnr* genes 112 clinical isolates of ESBL-producing *E. coli* from French hospitals in 2004. 7.7 % of CTX-M-producing *E. coli* presented a plasmid-mediated resistance to quinolones. All strains harboured a *qnrA* gene located on a class 1 integron.

Poirel et al (13) listed 186 ESBL positive Enterobacteriaceae. From them 2.2 % and 1.6 % carried a *QnrA1* and a *QnrS1* determinant respectively. The association of the *qnrA* gene with class 1 integrons was confirmed.

Hyunjoo Pai et al (14) screened *E. coli* and *K. pneumoniae* producing ESBLs or plasmid mediated AmpC beta-lactamases for the presence of *qnrA* and *qnrB* genes. *QnrB* was present in 54 of 54 DHA-1 producing *K. pneumoniae* isolates and 10 of 45 SHV-12 producing isolates. It is possible that *qnrB* contributes to the widespread distribution of DHA-1 (plasmid mediated AmpC) in areas, where 3rd generation cephalosporins and fluoroquinolones are widely used.

According to Lavilla (15) and Jones (16): the presence of *aac(6')*lb-cr were associated with quinolone resistance and aminoglycoside resistance (tobramycin is the best indicator).

Pitout et al (18) and Ruiz (20) found that isolates with *aac(6')*lb-cr were often associated with CTX-M-15.

**Screening procedure**

Perform antibiogram as usual (standard procedure): MH agar, inoculum McF 0.5, incubation at 35-37 °C for 18-24 hours.

Strains of Enterobacteriaceae should be suspected of plasmid-mediated quinolone resistance when showing unusual multiresistance phenotypes such as:

**Neo-Sensitabs**

Ampicillin	- no zone (HLR)
Sulphonamides	- no zone (HLR)
Trimethoprim	- no zone (HLR)
Trimethoprim+Sulfa	- no zone (HLR)
Streptomycins	- no zone (HLR)
Nalidixan 130 µg	- no zone or zone < 23 mm
Nalidixan 30 µg	- zone < 15 mm (MIC ≥16 µg/ml)
Norfloxacin 10 µg	- zone ≤ 23 mm
Ciprofloxacin 1 µg	- zone < 25 mm (MIC ≥0.125 µg/ml)
Ciprofloxacin 5 µg	- zone < 28 mm (MIC ≥0.125 µg/ml)
Ceftazidime	- zone < 20 mm
Chloramphenicol	- may show resistance
Tetracyclines	- may show resistance
aac(6')Ib-cr	- Tobra R (20), Kana R
QepA	- Genta R, Tobra R (20), Cipro R, Norflox R (15-17)

Suspected strains can be tested for the presence of the *qnr* gen by PCR.

It should be noted that strains showing the above-mentioned resistance phenotypes are most probably integron-carrying. Enterobacteriaceae and barrier precaution should be established to prevent further spread. A new chromosomal gen, *qnrM* from *S.maltophilia*, produce quinolone-resistance in *E.coli* (21).

In a selected group of ciprofloxacin and ceftazidime-resistant Enterobacteriaceae (mainly *K. pneumoniae* and *E. cloacae*), carriage of *qnrA* gene was 32 % (9). From those 73 % were ESBL-positive.

Cavaco et al (22) found that the combination of ciprofloxacin, norfloxacin and nalidixic acid is the best option for detecting *qnr* and *aac(6')Ib-cr* resistance mechanisms in Enterobacteriaceae.

#### **Detection of qnr and aac(6')Ib-cr and differentiation from mutations**

Screening with nalidixic acid is efficient for the detection of mutants, but it is not efficient for the detection of some isolates carrying qnr and aac(6')Ib-cr.

Transferable genes will be best detected using ciprofloxacin and norfloxacin, because these fluoroquinolones due to their chemical structures are attacked by strains carrying aac(6')Ib-cr, reducing their antimicrobial activity.

Combining the use of Nalidixic acid 30 µg, Norfloxacin 10 µg and Ciprofloxacin 1 µg Neo-Sensitabs it is possible to screen for the new resistance mechanisms qnrA, qnrS and aac(6')Ib-cr in strains of E.coli and Salmonella spp.

**Susceptible** strains E.coli and Salmonella spp. show zones of inhibition ≤ 25 mm with Nalidixic acid 30 µg Neo-Sensitabs (MIC ≤ 4-8 µg/ml).

With Ciprofloxacin 1 µg Neo-Sensitabs, susceptible strains show zones of inhibition ≥ 29 mm (MIC ≤ 0.016-0.03 µg/ml).

With Norfloxacin 10 µg Neo-Sensitabs, susceptible strains show zones of inhibition ≥ 33 mm (MIC ≤ 0.06 µg/ml).

Strains with **1 or 2 mutations** will show zones of inhibition ≤ 12 mm with Nalidixic acid 30 µg Neo-Sensitabs, while Ciprofloxacin will show zones ≤ 28 mm and Norfloxacin zones ≤ 32 mm.

Strains with **qnr or aac(6')Ib-cr** will show zones of inhibition ≥ 13 mm with Nalidixic acid 30 µg Neo-Sensitabs, zones of 16-25 mm (MIC 0.125-0.5 µg/ml) with Ciprofloxacin 1 µg Neo-Sensitabs and zones of ≤ 30 mm (MIC 0.25-2 µg/ml) with Norfloxacin 10 µg/ml Neo-Sensitabs.

#### **Conclusion:**

Strains with mutations are best detected with Nalidixic acid 30 µg Neo-Sensitabs. Strains with aac(6')Ib-cr are best detected using Ciprofloxacin 1 µg and Norfloxacin 10 µg Neo-Sensitabs, showing reduced zones of inhibition compared to fully susceptible strains. Strains possessing aac(6')Ib-cr, show also resistance to Tobramycin and Kanamycin, because these aminoglycosides are attacked by the mentioned enzyme. This is useful to differentiate them from the qnr enzymes.

#### **References:**

- 1) Minggui Wang et al: Emerging plasmid-mediated quinolone resistance associated with the *qnr* gene in *K. pneumoniae* clinical isolates in the U.S. Antimicrob. Ag. Chemother., **48**, 1295-9, 2004.
- 2) Rodriguez-Martinez J.M.: Mechanisms of plasmid-mediated resistance to quinolones. Enferm. Infecc. Microbiol. Clin., **23**, 25-31, 2005.
- 3) Hedi Mammeri et al: Emergence of plasmid-mediated quinolone resistance in *E. coli* in Europe. Antimicrob. Ag. Chemother., **49**, 71-6, 2005.
- 4) Xian-Zhi Li: Quinolone resistance in bacteria: emphasis on plasmid-mediated mechanisms. Intl. J. Antimicrob. Ag., **25**, 453-63, 2005.
- 5) Jin-Yong Jeong et al: Detection of *qnr* in clinical isolates of *E. coli* from Korea. Antimicrob. Ag. Chemother., **49**, 2522-4, 2005.
- 6) Poirel L. et al: Association of plasmid-mediated quinolone resistance with ESBL: VEB-1. Antimicrob. Ag. Chemother., **49**, 3091-4, 2005.
- 7) Robicsek A. et al: Broda distribution of plasmid-mediated quinolone resistance in the U.S. Antimicrob. Ag. Chemother., **49**, 3001-3, 2005.
- 8) Nordmann P, Poirel L: Emergence of plasmid-mediated resistance to quinolones in Enterobacteriaceae. J. Antimicrob. Chemother., **56**, 463-9.
- 9) Corkill J.E. et al: High prevalence of the plasmid mediated quinolone resistance determinant *anrA* in multidrug-resistant Enterobacter from blood cultures in Liverpool, UK. JAC, **56**, 1115-7, 2005.
- 10) Poirel L. et al: In vivo selection of fluoroquinolone-resistant *E. coli* isolates expressing plasmid-mediated quinolone resistance and ESBL. Antimicrob. Ag. Chemother., **50**, 1525-27, 2006.
- 11) Jacoby G.A. et al: QnrB, another plasmid-mediated gene for quinolone resistance. Antimicrob. Ag. Chemother., **50**, 1178-82, 2006.

- 12) Lavignes P. et al: qnrA in CTX-M producing E. coli isolates from France. Antimicrob. Agents Chemother., **50**, 4224-28, 2006.
- 13) Poirel L. et al: Prevalence and genetic analysis of plasmid mediated quinolone resistance determining 5 QnrA and QnrS in Enterobacteriaceae from a French University Hospital. Antimicrob. Agents Chemother., **50**, 3992-97, 2006.
- 14) Hyundoo Pai et al: Association of QnrB determinants and production of ESBL's or plasmid-mediated AmpC beta-lactamases in clinical isolates of K. pneumoniae. Antimicrob. Agents Chemother., **51**, 366-68, 2007.
- 15) Lavilla S. et al: Prevalence of the quinolone modifying enzyme aac(6')Ib-cr in ESBL-producing Enterobacterial isolates in Barcelona. Poster P-1523 ECCMID, april 2008.
- 16) Jones G. L. et al: The prevalence and distribution of plasmid-mediated quinolone resistance in clinical isolates of E.coli. Poster P-1530, ECCMID, Barcelona, april 2008.
- 17) Yamane K. et al: Plasmid-mediated qepA gene among E.coli clinical isolates from Japan. Antimicrob. Agents Chemother **52**, 1564-66, 2008.
- 18) Pitout J. D.D et al: Surveillance of plasmid-mediated quinolone resistance determinants in Enterobacteriaceae with the Calgary region Canada: the emergence of aac(6')Ib-cr. JAC **61**, 999-1002, 2008.
- 19) Fabregas A. et al: Quinolone resistance in the food chain. Int.J. Antimicrob. Agents, **31**, 307-315, 2008.
- 20) Ruiz E. et al: Detection of genes qnr and aac(6')Ib-cr in strains of E.coli and K. pneumoniae producing ESBL<sup>s</sup>. Poster 39P XIII Congress SEIMC, Madrid, may 2008.
- 21) Sanchez M. B. et al: A new chromosomal qnr gene, qnrM from S.maltophilia produces quinolone-resistance in E.coli. Presentation 3, XIII Congress SEIMC, may 2008, Madrid.
- 22) Cavaco L.M. et al: Evaluation of quinolones for use in detection of determinants of acquired quinolone resistance, including the new transmissible resistance mechanisms qnrA, qnrB, qnrS and aac(6')Ib-cr in E.coli and Salmonella enterica and determinations of wild type distributions. J. Clin. Microbiol. **47**, 2751-58, 2009.
- 23) Cavaco L.M. et al: qnrD, a novel gene conferring transferable quinolone resistance in Salmonella enterica serovar Kentucky and bovismorbificans strains of human origin. Antimicrob. Agents Chemother **53**, 603-608, 2009.

### **Detection of hypermutable strains**

Hypermutable or mutator strains are those that have an increased spontaneous mutation rate due to defects in genes involved in DNA repair or error avoidance systems.

Hypermutable strains have been described in natural populations of *E. coli*, *Shigella* spp, *Salmonella* spp. *P. aeruginosa*, *H. influenzae*, *Helicobacter pylori*, *N. meningitidis* and *S. aureus* among others. Hypermutable strains could have an important role in the emergence and spread of antibiotic resistance within bacterial populations, because they are substantially more resistant than nonmutators.

### **Detection of E.coli hypermutable strains**

Galan et al(1) and Conejo et al (2) described a diffusion method for detecting hypermutable strains in *E. coli*. A similar technique may be used with Neo-sensitabs.

MH-Blood agar plates are seeded with 100ul of an overnight Brain Heart Infusion broth culture. Fosfomycin Neo-sensitabs and Rifampicin 30 ug Neo-sensitabs are placed on the surface of the agar plate. Incubation takes place during 24 hours at 35-37 degrees Celsius and the inhibition zones are examined for the presence of colonies inside the inhibition zones.

Interpretation is as follows:

<u><b>Strong mutators:</b></u>	Strains with > 70 colonies in the Fosfomycin zone and > 10 colonies in the Rifampicin zone.
<u><b>Weak mutators:</b></u>	Strains with > 30 but <70 colonies in the Fosfomycin zone or > 70 colonies in the Fosfomycin zone but < 10 colonies in the Rifampicin zone.
<u><b>Nonmutators:</b></u>	Strains with < 30 colonies in the Fosfomycin zone and < 10 colonies in the rifampicin zone.

The same technique should be useful for *Shigella* and *salmonella* strains.

Denamur et al (3) found a high frequency of hypermutable strains in uropathogenic *E. coli* isolates, using nalidixic acid, fosfomycin, spectinomycin and streptomycin. The majority of hypermutable strains yielded colonies inside the inhibition zone (squatted colonies).

Baquero et al (4) found that hypermutable *E. coli* strains were found more frequently in a collection of ESBL-producing isolates, than in non-ESBL *E. coli*.

Ellington et al (5) found that hypermutator phenotypes were found among *E. coli* expressing ESBL CTX-M beta lactamases and they had an increased propensity to fosfomycin resistance.

### **Detection of P. aeruginosa hypermutable strains**

A technique similar to the one described for *E. coli* may be used for *P. aeruginosa* (6,7,8).

MH agar is inoculated with a 0.5 McFarland inoculum or McF 1.0 (mucoid colonies).

The following Neo-Sensitabs are added to the agar plate : Ceftazidime 30 ug, Ciprofloxacin 5 or 10 ug, Imipenem 10 or 15 ug, Meropenem 10 ug and Tobramycin 10 or 40 ug Neo-sensitabs.

Incubate for 24 hours (36 hours for slow growers) at 35-27 degrees.

The zones of inhibition are examined for the presence of resistant mutants (squatted colonies) inside the zones. When the diameter of the inhibition zones is reduced by > 5 mm (squatted colonies zone compared to the largest zone) the test is positive for hypermutability. With Tobramycin any resistant mutant inside the zone indicates a positive result.

Hypermutable strains will show the presence of resistant mutants inside the zone for 3 or more antibiotics.

Oliver et al(9) found that chronic infections caused by hypermutable *P. aeruginosa* of the sort frequently found in cystic fibrosis patients, will almost invariably be present in a high proportion at the onset of treatment, and consequently these strains should be considered resistant to all agents when they are used as monotherapy.

Henrichfreise et al (10) characterized a multiresistant hypermutable strain of *P. aeruginosa* from an ICU patient using Ciprofloxacin, Tobramycin, Ceftazidime, Imipenem and Meropenem by disk diffusion.

Wayne et al(11)found that the presence of hypermutator *P. aeruginosa* was associated with poorer lung function in cystic fibrosis patients.Mucoid isolates were more likely to be hypermutators.  
Macia et al (12)found that hypermutable *P. aeruginosa* strains from chronic respiratory infections had an increased susceptibility to colistin.

### **Hypermutable *Haemophilus influenzae***

Watson et al (13)found that hypermutable *H. influenzae* with mutations in *mutS* were prevalent in the cystic fibrosis lung environment.Rifampicin,Nalidixic acid and Spectinomycin were used to detect the hypermutable strains.Of 14 hypermutable *H influenzae* isolates,12 were from cystic fibrosis sputum.  
Perez-Vazquez et al (14)found that 39.3 % of 28 ciprofloxacin resistant *H.influenzae* were hypermutable and conclude that resistance to fluoroquinolones in *H. influenzae* is strongly associated to hypermutability.

### **References**

- 1)Galan JC et al : Fosfomycin and Rifampicin diffusion tests for detection of *E.coli* mutator strains : *J.Clin.Microbiol* **42**,4310-4312,2004.
- 2) Conejo MC et al: Detection of hypermutable *E.coli* strains in a collection of clinical isolates by the Fosfomycin-Rifampicin disk method.*J. Clin Microbiol.* **44**,2672, 2006.
- 3)Denamur E et al : "High frequency of mutator strains among human uropathogenic *E. coli* isolates. *J Bacteriol* **184**,605-609,2002
- 4)Baquero MR et al: Increased mutation frequencies in *E. coli* isolates harbouring ESBLs. *Antimicrob Ag Chemother.***49**,4754-4756,2005
- 5)Ellington MJ et al " Mutators among CTX-M beta-lactamase producing *E. coli* and risk for the emergence of fosfomycin resistance.*J.Antimicrob Chemother* **58**,848-852,2006
- 6)Macia MD et al : Detection and susceptibility testing of hypermutable *P. aeruginosa* strains with the E-test and disc diffusion.*Antimicrob Ag, Chemother.***48**,2665-2672,2004.
- 7)Macia MD et al : Hypermutation is a key factor in the development of multiple-antimicrobial resistance in *P. aeruginosa* strains causing chronic lung infections. *Antimicrob Ag. Chemother.***49**,3382-3386,2005
- 8) Manno G et al : " *P. aeruginosa* with hypermutable phenotype in cystic fibrosis patients.Poster P2100, ECCMID Barcelona,april 2008
- 9)Oliver A. et al " Hypermutation and the preexistence of antibiotic-resistant *P. aeruginosa* mutants: implications for susceptibility testing and treatment of chronic infections" *Antimicrob Agents Chemother* **48**,4226-4233,2004
- 10)Henrichfreise B et al:"Characterization of a multiresistant mutator strain of *P. aeruginosa* from an ICU patient" Poster p1267,ECCMID 2005.
- 11)Waine DJ et al : " Association between hypermutator phenotype,clinical variables,mucoid phenotype and antimicrobial resistance in *P. aeruginosa*".*J. Clin. Microbiol* **46**,3491-3493,2008
- 12)Macia MD et al." Increased susceptibility to colistin in hypermutable *P.aeruginosa* strains from chronic respiratory infections".*Antimicrob Agents Chemother.***51**,4531-4532,2008
- 13)Watson ME et al : " Hypermutable *Haemophilus influenzae* with mutations in *mutS* are found in cystic fibrosis sputum".*Microbiol.***150**,2947-2958,2004
- 14) Perez-Vazquez M et al : " Fluoroquinolone resistance in *H. influenzae* is associated to hypermutable strains. Presentation C2-1340,ICAAC 2006.

## **Detection of heteroresistant strains**

### **Introduction**

Heteroresistance refers to populations containing a majority of bacteria inhibited by concentrations below the susceptibility breakpoint, together with a small number of microorganisms (1 to 1000 to 1 10.000.000) that are resistant.

Heteroresistance has been described in *S. aureus* and coagulase negative staphylococci, particularly in MRSA, hVISA. Besides, in *Enterococcus faecium* (teicoplanin), pneumococci (penicillin), *Acinetobacter baumannii* (carbapenems and colistin), *Helicobacter pylori* (metronidazole, clarithromycin), *Klebsiella pneumoniae* (colistin), *Clostridium difficile* (metronidazole), *Ps aeruginosa* (carbapenems, piperacillin + tazobactam) *Cryptococcus neoformans* (fluconazole) among others.

The most common method used to detect heteroresistance is to see colonies inside the inhibition zones of disks (Neo-sensitabs), containing different antimicrobials. There is no doubt that the phenomenon of heteroresistance may be clinically relevant. In *S. aureus*, heteroresistance to oxacillin/methicillin is clinically relevant. The clinical importance of *S. aureus* isolates heteroresistant to glycopeptides, is under discussion, while the amount of hVISA strains is increasing worldwide. The use of antimicrobials in vivo may select heteroresistant mutants, that later on result in a population stable resistant.

Isolates of *A. baumannii* heteroresistant to colistin may show stable MICs of more than 8 – 16 µg/ml. These strains can be selected in vitro and probably can be selected in patients treated with colistin. *A. baumannii* isolates from patients previously treated with colistin, show a higher degree of heteroresistance.

### **Carbapenem heteroresistance in *A. baumannii***

Fernandez-Cuenca et al (1) indicates that *A. baumannii* heteroresistant to carbapenems should be detected in the laboratory, because it is probable that the resistant subpopulations, will be selected in the presence of imipenem or meropenem, resulting in a therapeutic failure. Microdilution and automatic methods (Wider, Vitek etc) cannot detect heteroresistance, because they are using too small inocula. To detect heteroresistance methods of diffusion on agar should be used and the inoculum should never be below to that corresponding to a 0.5 Mc Farland.

Pournaras et al (2) in 2005 describes the spread of isolates heteroresistant to carbapenems, they showed subcolonies present in the clear zone of inhibition. Resistant colonies were retested and again a subpopulation of resistant isolates was grown inside the zone of inhibition. Del Rosario Quintana et al (3) conclude that the automatic system Wider cannot detect the carbapenem heteroresistant strains of *A. baumannii*. Fernandez F et al (4) tested 30 clinical isolates of *A. baumannii* and heteroresistance was defined as presence of colonies inside the zones around imipenem disks. Heteroresistant strains showed MICs of 4-16 µg/ml towards imipenem and were associated to the presence of the beta-lactamase OXA-58 gen.

Gomez MC et al (5) tested 44 isolates of *A. baumannii* and concluded that 84 % of the isolates showing colonies inside the zone of imipenem, also possessed the OXA-58 gen, while none of the isolates without colonies inside the zone, possessed the gen OXA-58. Neou et al (6) tested 142 non-repetitive isolates of *A. baumannii*. Agar MICs for imipenem were 0.25 to 4 µg/ml. Colonies grown at 8 µg/ml did not show resistance stability when subcultured in drug-free medium. Agar MICs for meropenem were 0.25 to 4 µg/ml. Colonies grown at 8 to 32 µg/ml showed stability to meropenem resistance after 1 week subculture in drug-free medium, but they were susceptible to imipenem. The authors suggest that apparently carbapenem susceptible *A. baumannii* populations contain an amount of resistant meropenem subpopulations. The implementation of screening techniques to identify heteroresistant isolates is of significant

importance.

Observe the presence of subcolonies inside the zone of inhibition of Imipenem and Meropenem Neo-sensitabs

#### **Colistin heteroresistance in *A. baumannii***

Jian Li et al (7) tested 16 colistin susceptible clinical isolates of *A. baumannii* by population analysis profiles and by serial passaging with or without exposure to colistin. They demonstrated the presence of heterogeneous colistin-resistant *A. baumannii* in "colistin susceptible" isolates.

The authors conclude that colistin heteroresistant *A. baumannii* isolates cannot be discriminated from colistin susceptible by MIC measurements alone.

Colistin heteroresistant *A. baumannii* may be a preliminary stage that leads to the proliferation of resistant subpopulations upon exposure to colistin.

Hawley et al (8) conclude that the isolates exhibiting heteroresistance is significantly higher among isolates recovered from patients previously treated with colistin

Park et al (9) indicates that high colistin resistance rates in *Acinetobacter* have been reported from Korean hospitals and that was not due to clonal dissemination, but they arose independently.

Hawley et al (10) identified one colistin-dependent *A. baumannii* isolate. When plated on Mueller Hinton agar with a Colistin 10 ug disk, the isolate grew heavily immediately around the disk.

The Colistin 10 ug Neo-sensitabs prediffusion method ( 2 hours + 18/22 hours ) will detect both colistin heteroresistant and resistant strains.

#### **Carbapenem and PIP + TAZO Heteroresistance in *P. aeruginosa*.**

The presence of subcolonies inside the carbapenems disk zone indicates that heterogeneous subpopulations with reduced susceptibility to carbapenems may exist in a number of *P. aeruginosa* strains that appear to be carbapenem susceptible by conventional automated susceptibility methods.

Pournaras et al (11) tested 14 non-repetitive isolates of *P. aeruginosa* in which a few subcolonies appeared within the zone of inhibition of imipenem and meropenem disks. These isolates represented 27.5 % of the apparently carbapenem susceptible isolates.

Population analysis showed distinct subpopulations that grew in concentrations up to 18 ug/ml imipenem and 12 ug/ml meropenem. The heterogeneous subpopulations retained their resistance levels implying a rather stable expression of resistance. Conventional MIC dilution methods, using the standard 10,000 CFU per spot inoculum may miss carbapenem resistant mutants.

Pournaras et al (12) describe an isolate showing heteroresistance to Piperacillin + tazobactam. The isolate was reported as susceptible by automatic methods and by agar dilution. Nevertheless, the isolate exhibited distinct colonies within the inhibition zone around the piperacillin + tazobactam disk.

#### **Colistin resistant *P. aeruginosa***

Brannon et al (13) tested 19 colistin-resistant *P. aeruginosa* isolated from colistin-treated cystic fibrosis patients.

75 % of the colistin resistant CF strains were highly resistant to colistin ( MIC > 200 ug/ml) and the remaining were moderately resistant ( MIC > 2 ug/ml).

Montero et al (14) tested colistin-resistant *P. aeruginosa* (CORPA) isolated from 10 patients. In all cases, multidrug-resistant *P. aeruginosa* susceptible only to colistin and amikacin were isolated, before the emergence of CORPA. 9 of the patients had previously received prolonged courses of colistin ( mean 40 days).

Detection using Colistin 10 ug Neo-sensitabs and the prediffusion method ( 2 hours + 18/22 hours prediffusion)

#### **Colistin heteroresistant /resistant *Klebsiella pneumoniae* and carbapenem heteroresistant *E. aerogenes***

Poudyal et al (15) tested 22 multidrug-resistant clinical isolates of *K. pneumoniae*. 6 isolates were colistin-resistant with MICs  $\geq 32$  ug/ml. Colistin heteroresistance was observed in 15 of 16 isolates considered colistin-susceptible. Similar to our recent finding of colistin heteroresistance in *A. baumannii*, the MIC alone may not provide information to guide treatment, because heteroresistance is not detected by an MIC method.

Antoniadou et al (16) tested 18 colistin-resistant *K. pneumoniae* isolates from 13 patients over a 16 months period. Most of the isolates possessed ESBLs or metallo-beta-lactamases or both. Selective pressure due to extensive or inadequate colistin use may lead to the emergence of colistin resistance.

Papaioannou et al (17) in a study over 3 years (2005-2008) and almost 5000 isolates, found that 1.5 % of *A. baumannii*, 16 % of *Klebsiella pneumoniae* and 4.8 % of *P. aeruginosa* were colistin resistant.

Gordon et al (19) showed that the Microscan failed to detect heteroresistance to carbapenems in a patient with *E. aerogenes* bacteremia, while disk diffusion and E-test detected it.

Detection using Colistin 10 ug Neo-sensitabs and the prediffusion method ( 2 hours + 18/22 hours prediffusion) and Imipenem/Meropenem Neo-sensitabs by the diffusion method

#### **Metronidazole heteroresistance in *Clostridium difficile***

Pelaez et al (18) found that initially metronidazole-resistant *C. difficile* isolates became susceptible after thawing.;

However they presented slow-growing subpopulations within the inhibition zones of the metronidazole disk. The authors conclude that: resistance to metronidazole in toxigenic *C. difficile* is heterogeneous, and prolonged exposure to metronidazole can select for in vitro resistance. We recommend routine performance of the disk diffusion method, with primary fresh *C. difficile* isolates in order to ensure that metronidazole heteroresistant populations do not go undetected.

Similar heterogeneous metronidazole resistance has been observed in *Bacteroides* spp and *Helicobacter pylori*.

Detection using Metronidazole Neo-sensitabs. Observe the presence of colonies inside the zone of inhibition

#### **Vancomycin/Teicoplanin heteroresistance in *S. aureus* and coagulase negative staphylococci**

Is treated separately under the prediffusion method ( 2 hours + 18/22 hours ) prediffusion

#### **Vancomycin/Teicoplanin heteroresistance/resistance in *Enterococcus***

Is treated separately under the prediffusion method ( 2 hours + 18/22 hours) prediffusion

#### **References**

- 1) Fernandez-Cuenca F et al : Epidemiological and clinical features associated with the colonisation/infection by *A. baumannii* heteroresistant to carbapenems. Poster P-1655, ECCMID 2009.
- 2) Pournaras S et al : Heteroresistance to carbapenems in *A. baumannii*" JAC, **55**, 1055-1056, 2005.
- 3) Del Rosario Quintana C. et al " Detection of *A. baumannii* heteroresistant to carbapenems in a University Hospital ( Spanish)" Poster 231, XIII Meeting SEIMC, May 2009.
- 4) Fernandez F et al " Association between heteroresistance to carbapenems in *A. baumannii*, production of beta-lactamases (AmpC and oxacillinases) and sequences type ISABA" ( Spanish) Poster 12, XIII Meeting SEIMC, May 2009.
- 5) Gomez MC et al : Association between heterogeneous phenotypic resistance to carbapenems and production of oxacillinases, sequences type ISABA, porins and OPRD-like in clinical isolates of *A. baumannii* (Spanish). Poster 175. XIII Meeting SEIMC, May 2009.
- 6) Neou E et al "Investigation of carbapenem heteroresistance among susceptible *Acinetobacter baumannii*" Poster ECCMID Helsinki 2009.
- 7) Jian Li et al " Heteroresistance to colistin in multidrug-resistant *A. baumannii*" Antimicrob. Ag. Chemother. **50**, 2946-2950, 2006.
- 8) Hawley JS et al : Colistin heteroresistance in *Acinetobacter* and its association with previous colistin therapy" Antimicrob. Ag. Chemother **52**, 351-352, 2009

- 9) Park YK. Et al " Independent emergence of colistin-resistant *Acinetobacter* spp isolates from Korea" *Diagn Microbiol Infect. Dis* **64**,51-59,2009
- 10) Hawley JS et al: " Development of colistin-dependent *Acinetobacter baumannii-calcoaceticus* complex" *Antimicrob. Ag. Chemother* **51**,4529-4530,2007
- 11) Pournaras S et al : " Characterization of clinical isolates of *P. aeruginosa* heterogeneously resistant to carbapenems" *J. Med Microbiol* **56**,66-70,2007.
- 12) Pournaras S et al : Piperacillin/tazobactam-heteroresistant *P. aeruginosa* from urinary infection successfully treated by piperacillin/tazobactam" *JAC* 757-758,2008.
- 13) Brannon MK et al : " Polymyxin-resistant *P. aeruginosa* from colistin-treated Cystic fibrosis patients : diverse regulatory mutations associated with a specific lipid A phenotype" Presentation ICAAC 2006
- 14) Montero M et al " Emergence of colistin-resistant *P. aeruginosa* in a University Hospital. Presentation C1-1065, ICAAC 2008.
- 15) Poudyal A et al " In vitro pharmacodynamics of colistin against multidrug-resistant *Kl pneumoniae*" *J. Antimicrob. Chemother.* **62**,1311-1318,2008
- 16) Antoniadou A et al : " Colistin-resistant *Kl pneumoniae* emerging in intensive care unit patients: first report of a multiclonal cluster". *J Antimicrob Chemother* 59,786-790,2007
- 17) Papaioannou V et al : " Colistin resistance in multidrug resistant gram-negative bacteria". Presentation ECCMID 2009.
- 18) Pelaez T et al : " Metronidazole resistance in *Clostridium difficile* is heterogeneous". *J Clin Microbiol.* **46**,3028-3032,2008
- 19) Gordon NC et al: " Failure of the Microscan Walk Away system to detect heteroresistance to carbapenems in a patient with *E. aerogenes* bacteremia. *J Clin Microbiol* **47**,3024-3025,2009.

## Efflux pump inhibitors and NMP Diatabs

### Introduction

The accelerated evolution of antibiotic resistance to important human pathogens and the scarcity of new anti-infective drug families under development, makes that other ways are being tried.

Efflux is a general mechanism responsible for bacterial resistance to antibiotics. This active drug transport is involved in low intrinsic susceptibility, cross resistance to chemically unrelated classes of molecules and selection/acquisition of additional mechanisms of resistance.

As a consequence, inhibition of bacterial efflux mechanisms appears to be promising in order to a) increase the intracellular concentration of antibiotics, that are expelled by efflux pumps b) restore the drug susceptibility of resistant clinical strains and c) reduce the capability for acquired additional resistance(1).

The resultant efflux pump inhibitor/antibiotic combination drug should exhibit increased potency, enhanced spectrum of activity and reduced propensity for acquired resistance(2)

### Enterobacteriaceae

Kim et al(3) found that *Shigella flexneri* isolated in Korea and resistant to fluoroquinolones, showed an increased susceptibility

to ciprofloxacin, norfloxacin and ofloxacin in the presence of an efflux pump inhibitor (CCCP).

Freyre et al(4) studied the effect of an inhibitor of efflux pumps (PABN) on the MICs to different fluoroquinolones in clinical isolates of *E. coli*. Strains resistant to the fluoroquinolones showed an increased reduction of the MICs for levofloxacin and moxifloxacin.

Tran QT et al(5) studied the effect of the efflux pump inhibitor PABN on the resistance of *Enterobacter aerogenes* and *K. pneumoniae*.

Efflux was involved in resistance (chloramphenicol, sparfloxacin) in *E. aerogenes* isolates more frequently than in *K. pneumoniae*.

Kern WV et al (6) studied the efflux pump inhibitor 1-(1-naphthylmethyl)-piperazine (NMP) on clinical isolates of *E. coli*.

NMP was moderately active in reversing multidrug resistance in clinical isolates of *E. coli* and can partially restore fluoroquinolone susceptibility through inhibition of efflux pumps.

Schumacher A et al (7) studied the effect of NMP on drug susceptibility of Enterobacteriaceae (other than *E. coli*). NMP has shown to reverse multidrug resistance in *E. coli* overexpressing RND type efflux pumps. On other Enterobacteriaceae, NMP consistently reduced the MIC of linezolid in *C. freundii*, *Enterobacter aerogenes* and *K. pneumoniae* clinical isolates. Significant effects were also seen for levofloxacin, tetracycline and chloramphenicol in *E. aerogenes*, and for levofloxacin and tetracycline for *K. pneumoniae*.

Effect of NMP was more likely in isolates with decreased susceptibility to fluoroquinolones.

### Non-fermenters and Vibrio

Bina XR et al(8) studied the effect of NMP on antimicrobial susceptibility and virulence factor production in *Vibrio cholerae*.

NMP potentiated antimicrobial compounds that were substrates for the *V. cholerae* RND efflux systems. NMP inhibited the production of virulence factors cholera toxin and the toxin coregulated pilus.

Bean D et al(9) found that *A. baumannii* belonging to the multidrug resistant OXA-23 clone1 appeared to decrease in susceptibility to tigecycline in the presence of NMP. The converse was seen when NMP was combined with doxycycline, tetracycline or minocycline. The synergy seen between NMP and the tetracyclines must be due to NMP effect against the AdeABC resistance-nodulation division efflux pump.

Vila J et al (10) give a review of porins, efflux pumps and multidrug resistance in *Acinetobacter baumannii*. The efflux pump AdeABC of the family RND affects the following antimicrobials: aminoglycosides, beta-lactams, chloramphenicol, erythromycin, tetracyclines and reduced susceptibility to fluoroquinolones.

Lomovskaya O et al (11) identify and characterize the inhibitors of multidrug resistance efflux pumps in *P. aeruginosa*. They conclude that inhibition of efflux pumps in *P. aeruginosa* may significantly improve the

clinical performance of fluoroquinolones. It would suppress the emergence of *S. aureus* and pneumococci mutants resistant to ciprofloxacin.

#### NMP Diatabs

Rosco has developed a test (NMP Diatabs) for detecting efflux pumps, particularly of the RND family. The test is performed on Mueller Hinton Agar using McFarland 0.5 inoculum and the NMP Diatabs (9 mm) is placed near the corresponding Neo-sensitabs that should be tested for synergy. The distance between NMP Diatabs and Neo-sensitabs, will depend on the size of the inhibition zone of the corresponding Neo-sensitabs with the particular bacteria. If the zone is  $\leq 20$  mm, the NMP tablet should be placed 6-8 mm apart (edge to edge). If the zone is  $\geq 30$  mm the distance should be 10-12 mm apart. Synergy (enlargement of the zone of inhibition, keyhole zone or phantom zone) indicates the presence of an Efflux Pump, probably of the RND family.

#### References

- 1) Mahamoud A et al : " Antibiotic efflux pumps in gram-negative bacteria: the inhibitor response strategy. *J. Antimicrob Chemother* **59**, 1223-1229, 2007
- 2) Lomovskaya O et al : " Practical applications and feasibility of efflux pump inhibitors in the clinic. A vision for applied use". *Biochemical Pharmacol* **71**, 910-918, 2006
- 3) Kim JY et al : " Resistance to fluoroquinolones by the combination target site mutations and enhanced expression of genes for efflux pumps in *Shigella flexneri* and *Sh. Sonnei* strains isolated in Korea". *Clin Microbiol Infect* **14**, 760-765, 2008
- 4) Freyre C et al : " Effect of the inhibition of efflux pumps on Moxifloxacin and other fluoroquinolones on isolates of uropathogenic *E. coli*" (Spanish) Presentation 86. Meeting SEIMC June 2009
- 5) Tran QT et al : " Occurrence of efflux mechanism and cephalosporinase variant in a population of *Enterobacter aerogenes* and *K. pneumoniae* isolates producing ESBLs". *Antimicrob. Agents Chemother* **53**, 1652-1656, 2009
- 6) Kern WV et al : " Effect of 1-(1-naphthylmethyl)-piperazine, a novel putative efflux pump inhibitor, on antimicrobial drug susceptibility in clinical isolates of *E. coli*". *J. Antimicrob Chemother.* **57**, 339-343, 2006
- 7) Schumacher A et al : " Effect of 1-(1-naphthylmethyl)-piperazine, a novel putative efflux pump inhibitor, on antimicrobial drug susceptibility in clinical isolates of Enterobacteriaceae other than *E. coli*". *J. Antimicrob Chemother.* **57**, 344-348, 2006
- 8) Bina XR et al : " Effect of the efflux inhibitor 1-(1-phenylmethyl)-piperazine and PABN on antimicrobial susceptibility and virulence factor production in *Vibrio cholerae*". *J. Antimicrob Chemother.* **63**, 103-108, 2009
- 9) Bean D et al : " Paradoxical effect of 1-(1-naphthylmethyl)-piperazine on resistance to tetracyclines in multidrug resistant *A. baumannii*". *J. Antimicrob Chemother* **63**, 349-352, 2009
- 10) Vila J et al : " Porins, efflux pumps and multidrug resistance in *A. baumannii*". *J. Antimicrob Chemother.* **59**, 1210-1215, 2007
- 11) Lomovskaya O et al : " Identification and characterization of inhibitors of multidrug resistance efflux pumps in *P. aeruginosa*: novel agents for combination therapy". *Antimicrob Agents Chemother.* **45**, 105-116, 2001.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Page 1 / 1 - Document: 7.1.0

Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

### Introduction

High molecular weight antimicrobials (vancomycin, teicoplanin, daptomycin, colistin) diffuse poorly on agar media, resulting in difficulties in the interpretation of results when using the current disc diffusion method.

As a consequence, CLSI do not recommend the current diffusion method for vancomycin with staphylococci, daptomycin with staphylococci and enterococci and colistin with gram negative rods.

Teicoplanin has not yet been evaluated, but we expect it will perform like vancomycin. ROSCO Diagnostica has taken on the investigation of this problem and has developed a 2 + 18 hours (or 2 + 22 hours) prediffusion technique, permitting an easier differentiation between susceptible and resistant strains when testing against these antimicrobials.

### Principle

The principle of the prediffusion technique, was developed by a Danish microbiologist Frølund-Thomsen, several decades ago. The idea is to give the high molecular weight antimicrobial a longer period of time to diffuse into the agar before bacterial growth takes place.

In a current disc diffusion minute colonies of growth are visible after approximately 8 hours incubation at 35 degrees. As a consequence the antimicrobial has only approximately 8 hours to diffuse into the agar, because when minute colonies are formed, further antimicrobial diffusion will not affect the size of the inhibition zone. When using the 2 + 22 hours prediffusion the antimicrobial has 8 + 24 = 32 hours to diffuse into the agar, i. e. 4 times more than with the current diffusion method.

This results in a much larger zone size difference between 2 consecutive MIC values approximately 5 mm with the prediffusion method compared to 1.0-1.5 mm with the current disc diffusion method (see enclosed regression lines)

Another important point is that with the prediffusion method, the antimicrobial depot is eliminated after 2 hours. From this moment no further antimicrobial is added to the agar and the antimicrobial that has diffused during the 2 hours will continue further diffusion without any pressure from the depot.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Page 1 / 6 - Document: 7.2.0

Detection of VISA, GISA and hVISA using the Neo-Sensitabs prediffusion method

### DETECTION OF VISA, GISA AND hVISA USING THE NEO-SENSITABS PREDIFFUSION METHOD

1. One Neo-Sensitabs of each Vancomycin 30 µg and Teicoplanin 30 µg are placed on an uninoculated plate containing the susceptibility test medium Mueller-Hinton Agar.
2. After 2 hours at room temperature, the tablet (disc) is removed (by knocking the plate against the table), but prior to this a short name (VAN or TEI) is written on the back of the plate in order to make it possible to identify the antimicrobial.
3. Now the plate is maintained at room temperature for a further 18 to 22 hours (overnight).
4. The plate is now inoculated with the strain to be tested using a McFarland 0.5 inoculum. Additional antimicrobial discs (Neo-Sensitabs) may be added now, using a dispenser and thereafter the plate is incubated overnight at 35 degrees.
5. The zones of inhibition are then measured and compared with the corresponding zone breakpoints



*S. aureus* CB 182, Susceptible strain. V = Vancomycin, T = Teicoplanin, D = Daptomycin

#### Please note:

In the laboratory the prediffusion plate can be prepared the day before it is inoculated to avoid loss of time and results are available within 24 hours. Surplus of prediffused plates may be kept in the refrigerator for another 24 hours.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Page 2 / 6 - Document: 7.2.0

Detection of VISA, GISA and hVISA using the Neo-Sensitabs prediffusion method

### Interpretation

#### DETECTION of VISA, GISA, hVISA and hGISA

VISA, GISA and hVISA strains will show the following zones of inhibition using Mueller Hinton Agar and McFarland 0.5 inoculum and the prediffusion method:

hVISA, hGISA	VISA, GISA
Teicoplanin 30 ug inhibition zone < 20 mm <b>OR</b> Vancomycin 30 ug inhibition zone ≤ 22 mm	Teicoplanin 30 ug inhibition zone < 20 mm <b>AND</b> Vancomycin 30 ug inhibition zone ≤ 22 mm

Please notice that Teicoplanin in general is the most sensitive drug to detect these isolates. The current MIC methods and automatic systems are unable to detect hVISA strains, because they use too small inocula and consequently cannot detect heteroresistant isolates.

In the laboratory it may be useful to test vancomycin, teicoplanin and daptomycin together in order to be able to detect hVISA/VISA strains first and thereafter in order to find the best drug for treatment of MRSA and hVISA/VISA infections.

An 9 cm Mueller-Hinton Agar plate will be adequate for testing the 3 antimicrobials by the prediffusion method.



S. aureus ATCC 700698 Daptomycin, susceptible hVISA strain.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Page 3 / 6 - Document: 7.2.0

Detection of VISA, GISA and hVISA using the Neo-Sensitabs prediffusion method

## DETECTION OF COAGULASE NEGATIVE GLYCOPEPTIDE INTERMEDIATE STAPHYLOCOCCI

The same procedure as above is used for Teicoplanin 30 ug. Isolates showing zones of inhibition < 20 mm with Teicoplanin 30 ug (prediffusion method) should be reported as resistant to teicoplanin and possibly heteroresistant to vancomycin.




*S. aureus* AT 403, GISA strain, Daptomycin susceptible.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Page 4 / 6 - Document: 7.2.0

Detection of VISA, GISA and hVISA using the Neo-Sensitabs prediffusion method




# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Page 5 / 6 - Document: 7.2.0

Detection of VISA, GISA and hVISA using the Neo-Sensitabs prediffusion method



# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Page 6 / 6 - Document: 7.2.0

Detection of VISA, GISA and hVISA using the Neo-Sensitabs prediffusion method

### References

- 1) Nielsen SV, Casals JB : Detection of decreased susceptibility to glycopeptides in *S. aureus* using tablet(disc) prediffusion. 15th ECCMID, april 2005.
- 2) Ferreira Nunes AP et al : Heterogeneous resistance to vancomycin in *S. epidermidis*, *S. haemolyticus* and *S. warneri* clinical strains: characterisation of glycopeptide susceptibility profiles and cell wall thickening. *Int J Antimicrob Agents* **27**, 307-315, 2006.
- 3) Liu C et al : *Staphylococcus aureus* with heterogeneous resistance to vancomycin: epidemiology, clinical significance and critical assessment of diagnostic methods. *Antimicrob Agents Chemother*. **47**, 3040-3045, 2003
- 4) Sandberg A et al : First report of heterogeneous vancomycin intermediate *S. aureus* (hVISA, VISA) in Denmark. Poster AR-07. ICAAC 2004.
- 5) Jae-Hoon Song et al : Emergence in Asian countries of *S. aureus* with reduced susceptibility to vancomycin. *Antimicrob Agents Chemother* **48**, 4926-28, 2004.
- 6) Denis O et al : Emergence of vancomycin-intermediate *S. aureus* in a Belgian hospital: microbiological and clinical features. *J Antimicrob Chemother*. **50**, 383-391, 2002.
- 7) Lodise : Relationship between vancomycin MIC and failure among patients with MRSA bacteremia treated with vancomycin. *Antimicrob Ag. Chemother* **52**, 3315-3320, 2008.
- 8) In-Gyu Bae et al : hVISA phenotype in bloodstream MRSA isolates from an international cohort of patients with infective endocarditis : prevalence, genotype and clinical significance. *J. Infect Dis* **200**, 1355-1366, 2009.
- 9) Musta AC et al : Vancomycin MIC plus heteroresistance and outcome of MRSA bacteremia : trends over 11 years. *J. Clin Microbiol* **47**, 1640-44, 2009.
- 10) Lulitanond A et al : The first VISA strains isolated from patients in Thailand. *J. Clin. Microbiol* **47**, 2311-16, 2009.
- 11) Sader HS et al : Occurrence of vancomycin-tolerant and hVISA strains among *S. aureus* causing bloodstream infections in 9 USA hospitals. *J. Antimicrob Chemother* **64**, 1024-28, 2009.
- 12) Pillai SK et al : Development of reduced vancomycin susceptibility in MSSA. *Clin Infect Dis* **49**, 1169-1174, 2009.
- 13) Katayama Y. et al: Selection of heterogeneous intermediate *S. aureus* by imipenem. *Antimicrob. Ag. Chemother*, **53**, 3190-96, 2009.
- 14) Fusco D. N. et al: Clinical failure of vancomycin in a dialysis patient with methicillin susceptible, vancomycin-heteroresistant *S. aureus*. *Diagn. Microbiol. Infect. Dis.* **65**, 180-183, 2009.
- 15) Perichon B. et al: Van-A type vancomycin resistant *S. aureus*. *Antimicrob. Ag. Chemother* **53**, 4580-87, 2009.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Page 1 / 5 - Document: 7.3.0

### Detection of vancomycin resistant enterococci (VRE)

#### Detection of vancomycin resistant enterococci ( VRE )

Although clinical laboratories can reliably detect high level resistance to vancomycin, there are reports of poor proficiency in the detection of low-level inducible resistance to vancomycin.

Pendle et al (8) found that the Vitek 2 correctly identified only 2 of 16 isolates of *E. faecium* vancomycin-resistant isolates (VRE) with low-level vancomycin resistance. The remaining were not identified or identified as *aerococcus viridans*. Disk testing using the CLSI method, was unreliable for detection of vancomycin resistance upon primary isolation.

Heteroresistance of *E. faecium* to vancomycin, was reported by Alam et al (1) in 2001. The heterogeneity of the VRE isolate was observed only by Etest (colonies inside the zone), automated Microscan or microdilution were not adequate for detecting heteroresistance.

In some European countries up to 5 % of non-hospitalised persons could be asymptotically colonised with VRE (2). Hospital outbreaks of VRE are almost exclusively caused by a specific genogroup of vancomycin-resistant *E. faecium* characterised by co-resistance to ampicillin (and ciprofloxacin) and the presence of the variant *esp* gene (3, 7).

VanB phenotype and *vanA* genotype *Enterococcus faecium* with heterogeneous expression of glycopeptide resistance, are reported from Korea (5) and China (14). Outbreaks of vancomycin resistant *E. faecium* expressing VanD-like phenotype associated with a *vanA* genotype are reported from France (4) and Korea (9). Interpretation of MICs for vancomycin was difficult since most of the isolates presented heterogeneous expression of resistance, with colony growth in the E-test zone. Was better visible after 48 hours incubation. Colonies growing inside the zone of inhibition once retested, showed a homogeneous phenotype of resistance to vancomycin. Teicoplanin was not effective against vanD-*vanA* VRE strains both in vitro and in vivo.

Neves (11) in Brazil reports the emergence of *vanA* genotype among *E. gallinarum* from patients in a University Hospital in Rio de Janeiro. Resistance to teicoplanin, when testing motile enterococci, may be a clue to the presence of the *vanA* gene.

Dobbs et al (6) report the nosocomial spread of *Enterococcus faecium* resistant to both vancomycin and linezolid in a tertiary care Medical Center. Exposure to fluoroquinolones may be important for the emergence of linezolid resistance.

Souli et al (13) in a prevalence survey of faecal carriage in a tertiary care university hospital in Athens (Greece) found 6 isolates of vancomycin and linezolid resistant *vanA* positive *E. faecium*.

#### Procedure

- 1) One Teicoplanin 30 ug and one Vancomycin 30 ug Neo-sensitabs are placed on an uninoculated plate containing plain MH agar.
- 2) After 2 hours at room temperature, the tablets are removed (knocking the plate against the table) and the shortnames (VAN and TEI) are written on the back of the plate for identification purposes.
- 3) The plate is maintained at room temperature for further 18-22 hours (overnight)
- 4) The plate is now inoculated with the strain to be tested using McFarland 0,5 inoculum and additional antimicrobial disks (Neo-sensitabs) may be added using a dispenser. Thereafter, incubate at 35 degrees overnight.
- 5) The zones of inhibition are measured and compared with the corresponding zone breakpoints.

#### Interpretation

**Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™**  
**Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar**

Page 2 / 5 - Document: 7.3.0

**Detection of vancomycin resistant enterococci (VRE)**

VanA: Vancomycin R and Teicoplanin R : no zone of inhibition  
Van B: Vancomycin zone < 16 mm(hazy edge) and Teicoplanin zone >20 mm.  
VanC: Vancomycin zone < 12 mm (sharp edge). Teicoplanin zone >20 mm.

VanB Phenotype and vanA genotype : Vancomycin no zone, Teicoplanin zone < 6 mm( R ).  
VanC Phenotype and vanA genotype : Vancomycin no zone, Teicoplanin zone <12 mm ( R )  
VanD Phenotype and vanA genotype : Vancomycin no zone, Teicoplanin zone < 16mm ( R )


ORIGINAL ROSCO DOCUMENT

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Page 3 / 5 - Document: 7.3.0


Detection of vancomycin resistant enterococci (VRE)



**Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™**  
**Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar**

Page 4 / 5 - Document: 7.3.0

**Detection of vancomycin resistant enterococci (VRE)**



E.faecalis ATCC 51299 (van B) compared with susceptible E.faecalis

V=Vancomycin, T=Teicoplanin

**References**

- 1) Alam MR et al : Heteroresistance to vancomycin in Enterococcus faecium. J Clin Microbiol. **39**,3379-3381,2001.
- 2) Mascini EM et al : Vancomycin-resistant enterococci: consequences for therapy and infection control. Clin Microbiol Infect **11**, suppl 4, 43-56,2005.
- 3) Willems RJL et al : Global spread of vancomycin-resistant E. faecium from distinct nosocomial genetic complex. Emerg Infect Dis **11**,821-827,2005.
- 4) Naas T et al : First nosocomial outbreak of vancomycin-resistant E. faecium expressing a VanD-like phenotype associated with a vanA genotype. J. Clin Microbiol, **43**,3642-49,2005.
- 5) Jae-Hoon Song : High frequency of vancomycin-resistant E. faecium isolates with VanB phenotype and vanA genotype in Korean hospitals. Diagn Microbiol Infect Dis , **56**,401-406,2006.
- 6) Dobbs TE et al Nosocomial spread of E. faecium resistant to vancomycin and linezolid in a tertiary care medical center. J. Clin Microbiol **44**,3368-3370,2006.
- 7) Werner G et al : Emergence and spread of vancomycin resistance among enterococci in Europe. Eurosurveillance, **13**,1-16 november 2008.
- 8) Pendle S et al.: Difficulties in detection and identification of E. faecium with low-level inducible resistance to vancomycin, during a hospital outbreak. Clin Microbiol Infect **14**,853-857,2008.
- 9) Jae-Hoon Song et al : Clinical implications of vancomycin-resistant E. faecium(VRE) with VanD phenotype and vanA genotype. J. Antimicrob Chemother. **61**,838-844,2008.
- 10) Khan SA et al : Heteroresistance to vancomycin and novel point mutations in Tn1546 of E. faecium ATCC 51559. Int J Antimicrob Agents **31**,27-36,2008.
- 11) Neves FPG et al : Emergence of the vanA genotype among E. gallinarum isolates colonising the intestinal tract of patients in a university hospital in Rio de Janeiro, Brazil. Int J Antimicrob Agents **33**,211-215,2009.

**Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™**  
**Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar**

Page 5 / 5 - Document: 7.3.0

**Detection of vancomycin resistant enterococci (VRE)**

- 12) San Millan A et al : Van-B type E. faecium clinical isolate succesively inducibly resistant to,dependent on, and constitutively resistant to vancomycin.Antimicr Ag Chemother. **53**,1974-1982,2009.
- 13) Souli M et al : Colonisation with vancomycin and linezolid-resistant E. faecium in a university hospital:molecular epidemiology and risk factor analysis. Int J Antimicr Agents **33**,137-142,2009.
- 14) Li Gu et al : A new Tn1546 type vanB phenotype-vanA genotype vancomycin-resistant E. faecium isolates in mainland China.Diagn Microbiol Infect Dis,**63**,70-75,2009.

ORIGINAL ROSCO DOCUMENT

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 1 / 4 - Document: 7.4.0

## Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

## Detection of staphylococci and enterococci non-susceptible to daptomycin using the Neo-Sensitabs prediffusion method

### Detection of staphylococci and enterococci non-susceptible to daptomycin using the Neo-Sensitabs prediffusion method

Daptomycin, due to its large molecular weight, diffuses slowly into agar, resulting in small differences in the size of inhibition zones between susceptible and non-susceptible isolates.

Extending the diffusion time by a prediffusion method and eliminating the daptomycin depot after 2 hours, allow for better discrimination of non-susceptible isolates showing MICs near the susceptible breakpoint.

Data from a recent study (Katz et al 1) suggests that the 2+18 hours prediffusion method using Daptomycin (DAPCa) Neo-Sensitabs is a promising method for discriminating between daptomycin susceptible and non-susceptible isolates.

### Principle

When using the 2+18/22 hours prediffusion, daptomycin has 4 times longer time to diffuse into the agar, before bacterial colonies are formed, than with the current diffusion method. This results in a much larger zone size difference between 2 consecutive MIC, as it can be seen in the enclosed regression line.

Another important feature is that with the prediffusion method the antimicrobial depot is eliminated after 2 hours, avoiding further antibiotic diffusion from the depot into the agar.

### Procedure

- 1) One Daptomycin 30 µg Neo-Sensitabs (DAPCa) is placed on an uninoculated plate, containing plain MH agar.
- 2) After 2 hours at room temperature, the tablet is removed (knocking the plate against the table) and a short name (DAP) is written on the back of the plate for identification purposes.
- 3) The plate is maintained at room temperature for further 18-22 hours (overnight)
- 4) The plate is now inoculated with the strains to be tested using a McFarland 0.5 inoculum. Additional antimicrobial disks (Neo-Sensitabs) may be added now using a dispenser. Thereafter incubate at 35 degrees, overnight.
- 5) The zones of inhibition are measured and compared with the corresponding zone breakpoints.

### Interpretation

#### a) Staphylococci

Susceptible:	zone of inhibition $\geq 22$ mm (corresponding to MIC $\leq 1$ µg/ml)
Non-susceptible:	zone of inhibition $< 20$ mm (corresponding to MIC $\geq 2$ µg/ml)
Indeterminate:	zone 20-21 mm.

#### b) Enterococci

Susceptible:	zone of inhibition $\geq 12$ mm (corresponding to MIC $\leq 4$ µg/ml)
Resistant:	no zone of inhibition


**Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™**  
**Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar**

Page 2 / 4 - Document: 7.4.0

**Detection of staphylococci and enterococci non-susceptible to daptomycin using the Neo-Sensitabs prediffusion method**



*S. aureus* ATC 403, GISA strain, Daptomycin susceptible




*S. aureus* ATC 408 GISA strain Daptomycin non-susceptible.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 3 / 4 - Document: 7.4.0

Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Detection of staphylococci and enterococci non-susceptible to daptomycin using the Neo-Sensitabs prediffusion method



## References

Revision date: 22.03.10  
Rev./Appr. by: KM / KM  
Replaces date: 11.02.10

Establish date: 03.12.09  
Establish by: km  
Revision no.: 3

**ROSCO**  
DIAGNOSTICA

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 4 / 4 - Document: 7.4.0

**Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar**

**Detection of staphylococci and enterococci non-susceptible to daptomycin using the Neo-Sensitabs prediffusion method**

- 1) Katz BD et al : A new prediffusion method for the detection of Daptomycin (DAP) non-susceptible strains using Neo-sensitabs. Presentation D-226, ICAAC september 2007, Chicago, USA
- 2) Katz BD et al : Detection of daptomycin non-susceptible strains using Neo-sensitabs prediffusion method. Diagn Microbiol. Infect Dis. **61**, 315-320, 2008.
- 3) Jevitt LA et al : Evaluation of FDA approved disk diffusion breakpoint for daptomycin . Presentation D-48, 44<sup>th</sup> ICAAC, 2004.
- 4) Koeth LM et al : Comparison of Daptomycin E-test MICs on European Mueller-Hinton and Iso-sensitest Agars against 20 S. aureus. Presentation D-2221, 48<sup>th</sup> ICAAC, 2008.
- 5) Morrissey I : Susceptibility of European staphylococci to daptomycin. Presentation C1-152, 48<sup>th</sup> ICAAC, 2008.
- 6) Marco F et al : Daptomycin is effective in treatment of experimental endocarditis due to MRSA and GISA strains. Antimicrob. Ag. Chemotherapy **52**, 2538-2543, 2008.
- 7) Castanheira M et al : Update of the in vitro activity of daptomycin tested against 6710 gram-positive cocci isolated in North America (2006). Diagn Microbiol Infect Dis **61**, 235-239, 2008.
- 8) Campanile F et al : In vitro activity of daptomycin against methicillin and multiresistant S. haemolyticus invasive isolates carrying different mec complexes. Diagn Microbiol Infect Dis **61**, 227-231, 2008.
- 9) Leonard SN et al : Evaluation of vancomycin and daptomycin against MRSA and hVISA in an in vitro PK/PD model with simulated endocardial vegetations. Presentation P\_1057, ECCMID 2008.
- 10) Koeth L et al : Comparison of daptomycin E-test MICs on Mueller-Hinton, Iso-sensitest and BHI agars with and without calcium supplementation to broth microdilution MICs against 20 S. aureus isolates. Presentation P-593, ECCMID 2008.
- 11) Jones RN et al : Update on daptomycin activity and spectrum when tested against gram-positive strains collected in European Medical Centers. Presentation 1638, ECCMID 2009.
- 12) Sader HS et al : Potent antimicrobial activity of daptomycin tested against S. aureus with vancomycin MIC of 2 µg/ml isolated in the USA and European hospitals. Presentation 1639, ECCMID 2009.
- 13) Traczewski MM et al : Inhibitory and bactericidal activities of daptomycin, vancomycin and teicoplanin against MRSA isolates collected from 1985 to 2007. Antimicrob. Ag. Chemother **53**, 1735-38, 2009.
- 14) Friedrich L et al : Evidence for daptomycin E-test lot-related MIC elevations for S. aureus. Diagn Microbiol Infect Dis **65**, 306-311, 2009.
- 15) Tenover FC et al : Characterisation of S. aureus strain with progressive loss of susceptibility to vancomycin and daptomycin during therapy. Int J Antimicrob Agents **33**, 564-68, 2009
- 16) Sader HS et al : Nine-hospital study comparing broth microdilution and Etest method results for vancomycin and daptomycin against MRSA. Antimicrob Ag Chemother, **53**, 3162-65, 2009.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Page 1 / 5 - Document: 7.5.0

### Detection of colistin resistant *P. aeruginosa*, *A. baumannii* and Enterobacteriaceae

#### Detection of Colistin resistant *P. aeruginosa*, *A. baumannii* and Enterobacteriaceae

It is well-known the detection of *P. aeruginosa* showing colistin MICs of 4-16 ug/ml (R) cannot be accurately detected using the current disk diffusion method. Only high level resistance (MIC  $\geq$  128 ug/ml), will show no zone of inhibition around a 10 ug disk, and detected.

Testing of *Acinetobacter*, *B. cepacia* and *Stenotrophomonas* is unreliable by the current disk diffusion method.

Testing of one compound (colistin or polymyxin B) predict the results for the other. Colistin is preferred for testing, because it is likely to be used clinically and it is also slightly more sensitive in detecting resistance.

Siqueira et al (3) studied the importance of medium, calcium concentration and inoculum when testing polymyxin B against *Acinetobacter* spp and *P. aeruginosa*. Higher calcium content resulted in higher MICs

Lo-Ten-Foe (17), van der Heijden (5) compared different methods : disc diffusion, E-test, broth microdilution, Vitek 2 and agar dilution for colistin susceptibility testing against *P. aeruginosa*, *A. baumannii*, and heteroresistant *E. cloacae*.

The Vitek 2 displayed low sensitivity in the detection of heteroresistant subpopulations of *E. cloacae*. The current disk diffusion was an unreliable method for testing colistin.

Hiramatsu et al (6) compared disk diffusion, E-test and agar dilution for testing the susceptibility of metallo-beta-lactamases producing *P. aeruginosa* to polymyxins. The authors found poor concordance between disk diffusion and E-test methods compared to agar dilution

Goi-Barisic (7) and Reis (1) found that the current disk diffusion test with colistin was not adequate for routine testing of *A. baumannii*.

Ko et al (13) report the emergence of colistin-resistant *K. pneumoniae* isolates from blood in Korean hospitals.

Landman et al (2) report the citywide emergence of *P. aeruginosa* strains with reduced susceptibility to polymyxin B. 61 % of the isolates were susceptible to imipenem and 5% of the isolates had reduced susceptibility to polymyxin B (MIC 4 – 8 ug/ml).

Kitzis et al (8) report very high frequency of mutation to colistin resistance in gram negative bacteria, and resistance could not be detected by the current disc-diffusion for most resistant mutants.

Savini et al (15) report the isolation of a colistin-resistant *Hafnia alvei*, representing a serious clinical and microbiological concern.

Halaby et al (14) report the emergence of colistin resistance during use of selective decontamination of the digestive tract in an ICU, masked by the use of inappropriate laboratory methods. Colistin resistance among ESBL producing *K. pneumoniae* isolates emerged after the introduction of selective decontamination and the routinely used disc diffusion method failed to detect colistin resistance.

Rosco Diagnostica has developed the 2+18 (22) hours prediffusion method, permitting the detection of colistin resistant strains.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Page 2 / 5 - Document: 7.5.0

Detection of colistin resistant *P.aeruginosa*,  
*A. baumannii* and Enterobacteriaceae

### Procedure


- 1) One Colistin 10 ug Neo-sensitabs is placed on an uninoculated plate, containing plain MH agar.
- 2) After 2 hours at room temperature, the tablet is removed ( by knocking the plate against the table) and the shortname COL is written on the back of the plate, for identification purposes.
- 3) The plate is maintained at room temperature for further 18-22 hours (overnight).
- 4) The plate is now inoculated by the strain to be tested, using Mc Farland 0.5 inoculum and additional antimicrobial disks (Neo-sensitabs) may be added using a dispenser. Thereafter incubate at 35 degrees overnight.
- 5) The zones of inhibition are measured and compared with the corresponding zone breakpoints.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Page 3 / 5 - Document: 7.5.0

Detection of colistin resistant *P.aeruginosa*, *A. baumannii* and Enterobacteriaceae




# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Page 4 / 5 - Document: 7.5.0

Detection of colistin resistant *P.aeruginosa*,  
*A. baumannii* and Enterobacteriaceae



Colistin susceptible E.coli ATCC 25922

C=Colistin

## Interpretation

Heteroresistant strains will currently show colonies inside the inhibition zone.

Susceptible strains :	inhibition zone $\geq 15$ mm ( MIC $\leq 2$ ug/ml).
Intermediate	(indeterminate) : 14- 11 mm
Resistant :	$\leq 10$ mm ( MIC $> 4$ ug/ml)

## References

- 1) ReisAO et al : Polymyxin-resistant Acinetobacter spp isolates : What is next ? Emerg Infect Dis **9**,1025-26,2003.
- 2) Landman D et al : Citywide emergence of P. aeruginosa strains with reduced susceptibility to polymyxin B. J. Antimicrob Chemother. **55**, 954-57,2005.
- 3) Siqueira I et al : Importance of medium,calcium concentration,time of incubation and inoculum for testing polymyxin B against Acinetobacter spp and P. aeruginosa.Presentation D-1702,45<sup>th</sup> ICAAC 2005.
- 4) Jian Li et al : Heteroresistance to colistin in multidrug-resistant A. baumannii.Antimicrob Ag Chemother. **50**,2946-2950,2006
- 5) Van der Heijden IM et al : Comparison of disc diffusion,Etest and broth microdilution for testing susceptibility for carbapenem-resistant P. aeruginosa to polymyxins. Annals Clin Microbiol and Antimicrob. **6**, 1-8,2007.
- 6) Hiramatsu K et al : Comparison of disk diffusion,Etest and agar dilution for testing the susceptibility of metallo-beta-lactamase producing P. aeruginosa to polymyxins. Presentation D-2236, 48<sup>th</sup> ICAAC,2008.
- 7) Goic-Barisic I et al : In vitro testing of carbapenem-resistant A. baumannii to colistin. ECCMID 2008.
- 8) Kitzis MD et al : Very high frequency mutation to colistin resistance in gram negative bacteria.Presentation E-2038, 48<sup>th</sup> ICAAC,2008.Poudyal A et al : In vitro pharmacodynamics of colistin against multidrug-resistant K. pneumoniae.J. Antimicrob Chemother **62**, 1311-18,2008.
- 9) Montero M et al : Emergence of colistin-resistant P. aeruginosa in a University Hospital. Presentation C1-1065,48<sup>th</sup> ICAAC,2008.
- 10) Antoniadou A et al : Colistin-resistant K. pneumoniae emerging in intensive care unit patients:first report of a multiclonal cluster. J Antimicrob Chemother **59**,786-790,2007.
- 11) Park YK et al : Independent emergence of colistin-resistant Acinetobacter spp isolates from Korea. Diagn Microbiol Infect Dis **64**, 51-59,2009.

**Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™**  
**Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar**

Page 5 / 5 - Document: 7.5.0

**Detection of colistin resistant *P.aeruginosa*,  
*A. baumannii* and Enterobacteriaceae**

- 12) Papaioannou V et al : Colistin resistance in multidrug resistant gram-negative bacteria. Presentation ECCMID 2009.
- 13) Ko SK et al : Emergence of colistin-resistant *K. pneumoniae* isolates from blood in Korean Hospitals and their non-clonality. Presentation C2-099. 49<sup>th</sup> ICAAC 2009.
- 14) Halaby T et al : Emergence of colistin resistance during use of selective decontamination of the digestive tract in an ICU, masked by the use of inappropriate laboratory methods. Presentation C2-100, 49<sup>th</sup> ICAAC, 2009.
- 15) Savini V et al : Isolation of colistin-resistant *Hafnia alvei*. J. Med Microbiol **58**, 278-280, 2009.
- 16) Borda N et al : Comparison of methods : Diffusion, Prediffusion and E-test on isolates of *Acinetobacter baumannii-calcoaceticus* complex, against colistin. (In press).
- 17) Lo-Ten-Foe JR et al: Comparative evaluation of the Vitek2, disk diffusion, Etest, broth microdilution and agar dilution susceptibility testing methods for colistin in clinical isolates including heteroresistant *Enterobacter cloacae* and *A. baumannii* strains. Antimicrob Ag Chemother. **51**, 3726-2730, 2007.

### Resistance mechanisms and automated methods

It is well-known that automated methods, due to their use of small inocula have difficulties in detecting heteroresistant populations, leading to false susceptible results and consequently very major errors. Also well-known is the difficulties of automated methods in testing beta-lactam antimicrobials and detection of beta-lactamases.

We have screened papers and presentations during 2005 to 2009, including only well recognised authors and magazines, such as Clinical Microbiology & Infection, J. Clinical Microbiology, J. Medical Microbiology, Int. J. Antimicrobial Agents and presentations at the ICAAC and ECCMID.

#### MicroScan

Jones et al (1) report from CAP: false automated system results for VRSA, VISA and piperacillin/tazobactam with *P. aeruginosa*. Suboptimal MRSA and ESBL detection in commercial automated systems.

Sader et al (2) report an unacceptable rate ( 21 – 32 %) of Very Major errors for piperacillin/tazobactam with *P. aeruginosa*.

Unda et al ( 4) report false resistance to amikacin with the MicroScan WalkAway 96 system as well as major errors with gentamicin and tobramycin; concluding that the system was not reliable for susceptibility testing of *P. aeruginosa* to amikacin.

Sader et al ( 7) studied the accuracy of MicroScan , Vitek and Vitek 2 for susceptibility testing of *P. aeruginosa* to 5 broad spectrum beta-lactams. All systems tested exhibited a high, unacceptable level of very major (false susceptible) errors for piperacillin/tazobactam ( 19-27%) and minor error rates were elevated for cefepime and aztreonam ( 8-32 %), leading to consistent trends towards false resistance.

Jones et al (9) detected errors by automated systems in CAP surveys. Numerous false susceptible and false resistant were noted for each antimicrobial. False positive ESBL results for a CMY-2 AmpC producing strain. Besides, false susceptibility to linezolid and false resistance to cefepime (2005), false susceptibility to amox+ clav and piperacillin/tazobactam ( 2006).

Juretschko et al (10) reported unacceptable levels of errors (minor, major and very major) with false susceptibility to piperacillin-tazobactam and piperacillin and false resistance to aztreonam, cefepime and ceftazidime in *P. aeruginosa*.

Kulah et al (11) studied imipenem resistance in *A. baumannii* and MicroScan showed the worst performance with 25 % very major errors and 44.6 % minor errors.

Juretschko et al (15) found that 25 % of all cystic fibrosis isolates identified as MSSA by automated methods were *mecA* positive and therefore MRSA (false susceptibility, very major error).

Gordon et al (20) describe the failure of MicroScan to detect heteroresistance to carbapenems in a patient with *E. aerogenes* bacteremia and conclude : reliance on automated susceptibility testing, while is more rapid than disc diffusion methods, may result in a delay in detecting or inability to detect the development of resistance.

#### Phoenix

Juretschko et al (10) reported false resistance to aztreonam and cefepime for *P. aeruginosa*

Kehl et al ( 12) found that the Phoenix automated system made false positive identification of vancomycin resistant enterococci (VRE).

Arslan U et al(14) studied automated systems to detect HLR to gentamicin in enterococci.4 of 5 isolates with gentamicin MICs 256 ug/ml,were found susceptible by Phoenix.Strains with aac-aph gene were reported false susceptible .

Juretschko et al ( 15) found that 25 % of all cystic fibrosis isolates identified as MSSA by Phoenix were mecA positive and therefore MRSA (very major error).

Boyd et al(16) report misidentification of *K. pneumoniae* carbapenemase (KPC) as ESBL by the Phoenix instrument and conclude that KPC can remain undetected or misinterpreted as ESBL by automated systems.Routine phenotypic testing should be performed.

Boyd et al (18) mention the erroneous reporting of ESBL producers by Phoenix and Vitek 2.

9 of 22 non-ESBL were reported as ESBL by Phoenix (specificity 59 %).The authors conclude that false positive ESBL may be reported by Phoenix and that implementation of phenotypic testing should be considered in ESBL "positive" isolates by automated methods.

Fisher et al ( 19) compared Phoenix with disc diffusion for identification of ESBL,AmpC and KPC producers.The Phoenix system misclassified nearly 50% of the isolates as ESBL positive.Besides a high false positive detection rate for AmpC. Of 8 KPC producing isolates the Phoenix system reported 50% as meropenem and imipenem susceptible (very major error).The authors conclude that we need reliable methods to detect these important resistance mechanisms in the clinical laboratory.

### Vitek 2

Jones et al (1) report from CAP 30 % false susceptibility with piperacillin/tazobactam and *P. aeruginosa*.

Sader et al (2) report an unacceptable rate (21-32 %) of Very major errors for piperacillin/tazobactam with *P. aeruginosa*..

Ceyssens et al (3) report that VITEK-2 overestimates the number of ESBLs among Enterobacteriaceae,that may lead to unnecessary use of carbapenems.

Karatuna et al ( 5)studied the reliability of carbapenem resistance in *P. aeruginosa*,by automated systems (Vitek2)and concluded that laboratories using automated systems should consider using at least a second method to validate intermediate or resistant results for carbapenems tested against *Pseudomonas* spp.

Tokatlidou D et al ( 6) observed the overdetection of imipenem resistance by Vitek 2 in VIM 12 producing *K. pneumoniae* with Etest low level carbapenem MICs .

Sader et al (7) found an unacceptable level of very major(false susceptible) errors for piperacillin/tazobactam( 19-27 %)as well elevated minor error rates for cefepime and aztreonam,when testing against *P. aeruginosa*.

Navon-Venezia et al (8) found that when testing cefepime susceptibility of ESBL producing *E. coli* and *K. pneumoniae* the Vitek 2 results, were major errors.The high proportion errors and major errors negate the use of Vitek 2 for cefepime MIC determination when ESBL positive strains are considered.

Jones et al (9) studied errors of automated methods according to CAP surveys.False susceptibility to imipenem when testing a SME-1 producing *S. marcescens*.Besides false susceptibility to meropenem and false resistance to ceftazidime ( 2003).False susceptibility to piper/tazo , amikacin imipenem and piperacillin (2004). False susceptibility to linezolid (2005).False resistance to ceftazidime (2006).

Juretschko et al ( 10) reported unacceptable levels of error with false susceptibility to piperacillin/tazobactam and imipenem and false resistance to aztreonam,cefepime and ceftazidime.

Pitout et al (13) report that Vitek Legacy and Vitek 2 failed to detect piperacillin/tazobactam resistance in 74 % of 101 piper/tazo resistant,ESBL producing strains,especially CTX-M-15 producing isolates that co-produced OXA-1.

Arslan et al (14) found that Vitek 2 failed to detect HLR gentamicin resistance in *E. faecium* isolates with aac-aph gene (false susceptible,very major error).

Juretschko et al ( 15) conclude that the VITEK 2 misidentified 25 % of MRSA isolated from cystic fibrosis patients as MSSA (very major error).

Boyd et al (16) report misidentification of *K. pneumoniae* carbapenemases (KPC) as ESBL by the Vitek 2 and conclude that KPC isolates can evade detection with automated methods and be misinterpreted as ESBL.

Bulik et al (17) compared meropenem MICs for KPC producing *K. pneumoniae* by Etest,broth microdilution and Vitek2.When comparing Vitek2 MICs to broth microdilution,only 22 % of isolates were in agreement.

Boyd et al ( 18) indicates the erroneous reporting of ESBL producers by Vitek 2. 13 of 22 non-ESBLs were reported as ESBL positive by the Vitek2 (specificity 40.9%).Implementation of phenotypic testing should be considered.

### References

- 1) Jones RN et al : Accuracy for antibacterial and antifungal susceptibility tests : report from the College of American Pathologists (CAP) microbiology proficiency survey program for 2004-2005. Presentation D-1654,45<sup>th</sup> ICAAC 2005.
- 2) Sader HS et al : Influence of errors caused by automated systems in the antimicrobial susceptibility rates of *P. aeruginosa* against selected beta-lactam antimicrobials.Presentation C-007,45<sup>th</sup> ICAAC 2005.
- 3) Ceyssens C et al : Overestimation of ESBLs by VITEK-2 and possible consequences for antibiotic policy.Clin Microbiol Infect **12**,suppl 4,Abstract P452,2006.
- 4) Unda F et al : False resistance to amikacin in *P. aeruginosa* with the MicroScan WalkAway 96 system.Clin Microbiol Infect **12**,suppl 4,Abstract P444,2006
- 5) Karatuna O et al : Reliability of carbapenem resistance in *Pseudomonas* spp detected by automated systems.Clin Microbiol Infect **12**,suppl 4 ,Abstract P1453,2006.
- 6) Tkatlidou D et al : Overdetection of imipenem resistance by VITEK 2 in VIM 12 producing *K. pneumoniae* with Etest low level carbapenem MICs.Clin Microbiol Infect **12**,suppl 4, abstract P1270,2006.
- 7) Sader HJ et al : Accuracy of 3 automated systems (MicroScan WalkAway,VITEK and VITEK 2) for susceptibility testing of *P. aeruginosa* against 5 broad-spectrum beta-lactam agents.J. Clin Microbiol. **44**,1101-1104,2006.
- 8) Navon-Venezia S et al : Cefepime susceptibility testing of ESBL producing *E. coli* and *K. pneumoniae*,by Vitek 2 results in major errors.Presentation D-0456,46<sup>th</sup> ICAAC,2006.
- 9) Jones RN et al : Detection of antimicrobial susceptibility testing errors by automated systems:College of American Pathologists microbiology surveys 2003-2006.Presentation D-223, 47<sup>th</sup> ICAAC 2007.
- 10) Juretschko S et al : Accuracies of beta-lactam susceptibility test results for *P. aeruginosa* with 4 automated systems (BDPhoenix,Microscan WalkAway,Viek and Vitek 2).J. Clin Microbiol **45**, 1339-1342,2007.
- 11) Kulah C et al : Comparison of 3 automated systems (BD Phoenix,Microscan WalkAway and Vitek 2) for detecting imipenem resistance in *A. baumannii*.ECCMID 2008.
- 12) Kehl SC et al : False-positive identification of vancomycin-resistant *Enterococcus* from the BD Phoenix instrument.Presentation C-022,ASM Meeting,june 2008.
- 13) Pitout JDD et al : Antimicrobial susceptibility of well-characterised multiresistant CTX-M producing *E. coli* : failure of automated systems to detect resistance to piperacillin/tazobactam.Int J Antimicrob Agents,**32**,333-338,2008.
- 14) Arslan U et al : Vitek 2 and Phoenix fail to detect high level gentamicin resistant *Enterococcus faecium* isolates with aac-aph gene. Presentation D-2241, 48<sup>th</sup> ICAAC 2008.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 4 / 4 - Document: 8.1.0

## Resistance mechanisms and automated methods

## Resistance mechanisms and automated methods

- 15) Juretschko S et al : Failure of automated identification systems to detect methicillin resistance of *S. aureus* in patients with cystic fibrosis. Presentation D-789, 49<sup>th</sup> ICAAC 2009.
- 16) Boyd NK et al : Misidentification of *K. pneumoniae* carbapenemases (KPC) as ESBL producers by automated testing instruments. Presentation D-731, 49<sup>th</sup> ICAAC 2009.
- 17) Bulik CC et al : Comparison of meropenem MICs for carbapenemase producing *K. pneumoniae* (KPC) by Etest, broth microdilution and VITEK 2 methods. Presentation D-721, 49<sup>th</sup> ICAAC 2009.
- 18) Boyd NK et al : Erroneous reporting of ESBL producers from Vitek2 and Phoenix automated systems. Presentation D-1441, 49<sup>th</sup> ICAAC 2009.
- 19) Fisher MA et al : Performance of the Phoenix bacterial identification system compared with disc diffusion methods for identifying ESBL, AmpC and KPC producers. *J. Med Microbiol* **58**, 774-778, 2009.
- 20) Gordon NC et al : Failure of MicroScan WalkAway system to detect heteroresistance to carbapenems in a patient with *E. aerogenes* bacteremia. *J. Clin. Microbiol* **47**, 3024-3025, 2009.

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Detection of resistance mechanisms

## Detection of resistance mechanisms (general)

### Detection of Resistance Mechanisms (General)

Mechanisms of resistance include production of inactivating enzymes, alteration of drug targets, and altered drug uptake or efflux. Find enclosed in the table below the antibiotics recommended to detect certain resistance mechanisms (1).

Antibiotic (Neo-Sensitabs)	Phenotype	Mechanism of resistance	Bacteria
<b>1) Beta-lactams</b>			
Penicillin + pH indicator (Beta-lactamase - D.T.)	Penicillin resistance	Penicillinase	Staphylococci, Haemophilus, Gonococci
Oxacillin 1 µg res.	Resistance to all beta-lactams	Additional PBP	Staphylococci
Cefoxitin res.	Resistance to all beta-lactams	mecA	Staphylococcus aureus
Ampicillin res.	a) Resistance to penicillins and beta-lactams inhibitor comb.	Altered PBPs	Enterococci
Oxacillin 1 µg (zone < 20 mm) (zone < 14 mm) (zone < 12 mm)	Penicillin resistance	PBP alteration	Pneumococci, Streptococci, Gonococci
Ceftizoxime res.	Resistance to third generation cephalosporins	PBP alteration	Pneumococci
Ampicillin 2.5 µg (zone < 20 mm)	Resistance to AMP, AMX, AM+CL, CCLOR, CEFUR (BLNAR strains)	PBP alteration	Haemophilus
Amoxycillin and Amoxycillin+Clavulanate	b) Penicillin resistance AM+CL synergy	Beta-lactamase (BRO-1, BRO-2)	Moraxella catarrhalis
Ceftazidime, Ceftriaxone (zone < 24 mm)		Screening ESBL	Klebsiella spp, E. coli, Salmonella
Cefpodoxime I/R		Screening ESBL	E. coli, Klebsiella Salmonella
Cefotaxime, Ceftriaxone Ceftazidime/Cefepime and Amoxycillin+Clavulanate	c) Synergy between CTX, CTR, CAZ and AMC, CP+CL, (double disk synergy)	Extended spectrum beta-lactamase (ESBL)	Enterobacteriaceae d)
Ceftazidime+Clavulanate	CAZ+CL zone ≥ 5 mm than CAZ alone	ESBL	Enterobacteriaceae
Cefepime and Amoxycillin + Clavulanate	Synergy between FEP and AMC	ESBL	Enterobacter, Serratia, Citrobacter freundii
Cefepime+Clavulanate	FEP+CL zone ≥ 5 mm than FEP alone	ESBL (confirmatory)	Enterobacteriaceae
Cefotaxime+Clavulanate	CTX+CL zone ≥ 5 mm than CTX alone	ESBL	Enterobacteriaceae
Cefoxitin + Cephalosporins Imipenem + Cephalosporins	Antagonism, indicates cephalosporin resistance	Inducible cephalosporinase AmpC	Enterobacteriaceae
Cefepime+Clavulanate Cefoxitin/3rd gen. cepha.	No synergy No antagonism Cefoxitin R, Ceftazidime R	Plasmid mediated AmpC	Enterobacteriaceae
Amoxycillin+Clavulanate I/R Cefazolin	Amoxycillin+Clavulanate R Cefazolin S	Inhibitor resistant TEM β-lactamase	E. coli Klebsiella
Cefotaxime+Boronic Ceftazidime+Boronic	zone ≥ 5 mm than CTX or CAZ alone	AmpC beta lactamase	Enterobacteriaceae

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Detection of resistance mechanisms

## Detection of resistance mechanisms (general)

Antibiotic (Neo-Sensitabs)	Phenotype	Mechanism of resistance	Bacteria
Cefoxitin res.	Antibiotic resistance	Porin alteration	E. coli Klebsiella
Aztreonam, Ceftazidime, Cefepime and Ticarcillin + Clavulanate	Synergy between TC+Cl and AZT, FEP, CAZ.	ESBL	Ps. aeruginosa
Imipenem+EDTA Meropenem+DPA	Synergy between Imipenem and EDTA Synergy MRP and DPA	Metallo-β-lactamase	Ps. aeruginosa Acinetobacter E. coli
Cloxacillin	Synergy between cefotaxime, ceftazidime and cloxacillin	AmpC beta-lactamase	Enterobacteriaceae
Dipicolinic acid	Synergism with Meropenem and/or Imipenem	Metallo-β-lactamase	Enterobacteriaceae Non-fermenters
Boronic acid	Synergism with Cefotaxime and/or Ceftazidime	AmpC beta-lactamase	Enterobacteriaceae
Meropenem+Boronic Meropenem+Cloxa	Synergy (≥5mm) No synergy	KPC beta lactamase	K.pneumoniae P.aeruginosa
Meropenem+Boronic Meropenem+Cloxa	Synergy synergy	AmpC+porin loss	Klebsiella Enterobacter
<b>2) Aminoglycosides</b>			
Kanamycin res.	Amikacin and Isepamicin resistance	APH(3'), ANT(4')	Staphylococci
Gentamicin res.	Resistance to aminoglycosides except streptomycin	APH(2'')-AAC(6')	Staphylococci
Kanamycin 500 µg (zone < 14 mm)	HLR to amikacin (no synergy with penicillins)	APH(3'), ANT(4')	Enterococci (HLR)
Gentamicin 250 µg (zone < 14 mm)	HLR to all aminoglycosides	APH(2'')-AAC(6')	Enterococci (HLR)
Streptomycin 500 µg (zone < 14 mm)	Streptomycin resistance		Enterococci (HLR)
Amikacin + Tobramycin res.	Resistance to aminoglycosides	APH(3')-VI	Acinetobacter
Netilmicin + Tobramycin res.	Resistance to aminoglycosides	AAC(3)	Pseudomonas
<b>3) Others</b>			
Erythromycin + Clindamycin	Inducible MLS resistance (antagonism)	Ribosomal methylation	Staphylococci, Streptococci,
Nalidixan I/R	Reduced sensitivity to quinolones	DNA gyrase	Enterobacteriaceae Vibrio cholerae Haemophilus/Moraxella
Nalidixan I/R	Reduced sensitivity to quinolones	DNA gyrase	Gonococci Meningococci
Ciprofloxacin 0.5 µg (zone < 20 mm)	Quinolone resistance	DNA-gyrase	Gonococci, Haemophilus
Vancomycin Teicoplanin	Vancomycin resistance 2+18 hours' prediffusion	Van A, Van B VISA, hVISA	Enterococci, Staphylococci
Metronidazole	Imidazole resistance	Reductase	Anaerobes

- for non beta-lactamase producing enterococci.
- Synergy: AM+CL zone > 5 mm larger than AMOXY (resistance to penicillin, amoxycillin and ampicillin, susceptible to amoxycillin+clavulanate).

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Page 3 / 3 - Document: 9.1.0

## Detection of resistance mechanisms

## Detection of resistance mechanisms (general)

- c) Beta-lactam resistance (except cephamycins and carbapenems).
- d) Except *Proteus penneri* and *P. vulgaris*.

**Note:** The mentioned zone sizes are valid for McFarland 0.5 inoculum.

### References:

- 1) Clinical Microbiology and Infection. Vol 2. Suppl. 1. December 1996.
- 2) Hakonen A. et al: Detection of decreased fluoroquinolone susceptibility in Salmonellas and validation of Nalidixic acid screening test. J. Clin. Microbiol., 37, 3572-77, 1999.

ORIGINAL ROSCO DOCUMENT

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Intrinsic (natural) resistance

## Intrinsic (natural) resistance

### Intrinsic (Natural) Resistance

Antimicrobial resistance can be classified as either intrinsic or acquired. Intrinsic resistance may be related to inherent or natural characteristics in a bacteria and may be used for recognition of a bacterial species, and results of in vitro susceptibility testing is not relevant to report as treatment options. The most relevant drug related natural resistance in a group or species is listed below:

BACTERIA	NATURAL RESISTANCE
Enterobacteriaceae	
Enterobacteriaceae in general	Penicillinase stable penicillins, Macrolides, Fucidin, Rifampicin, Glycopeptides
Enterobacteria, group 2 K. pneumoniae, K. oxytoca, C. diversus, Esch. hermannii	Aminopenicillins, Carboxypenicillins
Enterobacteria, group 3 E. cloacae, E. aerogenes, C. freundii, S. marcescens, M. morganii, Prov. rettgeri, P. stuartii, Hafnia alvei (except P. vulgaris AMC S)	Aminopenicillins, Amoxicillin+Clavulanate, 1 <sup>st</sup> gen Cephalosporins
Citrobacter freundii	Aminopenicillins, Amoxicillin+Clavulanate, 1 <sup>st</sup> gen Cephalosporins, Cefoxitin,
Citrobacter koseri (diversus)	Aminopenicillins, Carboxypenicillins
Enterobacter aerogenes, E. cloacae	Aminopenicillins, Amoxicillin+Clavulanate, Cefoxitin, 1 <sup>st</sup> gen Cephalosporins, Nitrofurantoin
Klebsiella pneumoniae, K. oxytoca	Aminopenicillins, Carboxypenicillins
Morganella morganii	Aminopenicillins, Amoxicillin+Clavulanate, 1 <sup>st</sup> and 2 <sup>nd</sup> gen Cephalosporins, Cefoxitin, Polymyxins, Tetracyclines, Nitrofurantoin, Fosfomycin
Proteus mirabilis	Polymyxins, Tetracyclines, Nitrofurantoin
Proteus vulgaris, P. penneri	Aminopenicillins, Carboxypenicillins Cefuroxime, Polymyxins, Tetracyclines, Nitrofurantoin
Providencia rettgeri	Aminopenicillins, Polymyxins, Tetracyclines, Nitrofurantoin, Amoxycillin+Clavulanate.
Providencia stuartii	Aminopenicillins, Amoxicillin+Clavulanate, Polymyxins, Tetracyclines, Nitrofurantoin, Gentamicin, Tobramycin, Netilmicin
Salmonella spp.	1 <sup>st</sup> and 2 <sup>nd</sup> gen Cephalosporins, Cefuroxime (active in vitro, not active in vivo) Aminoglycosides (in vivo)
Serratia marcescens	Aminopenicillins, Amoxicillin+Clavulanate, 1 <sup>st</sup> and 2 <sup>nd</sup> gen Cephalosporins Polymyxins
Shigella spp.	1 <sup>st</sup> and 2 <sup>nd</sup> gen Cephalosporins, Aminoglycosides (in vivo)
Yersinia enterocolitica	Aminopenicillins, Carboxypenicillins, Amoxicillin+Clavulanate, 1 <sup>st</sup> and 2 <sup>nd</sup> gen Cephalosporins, Cefoxitin
Gram positive cocci in general	Aztreonam, Nalidixic acid, Polymyxins
Staphylococci	
Staphylococcus spp. in general	Nalidixic acid, Polymyxins
S. saprophyticus	Novobiocin, Fosfomycin
Penicillin resistant staphylococci (Oxa S)	Penicillin, Aminopenicillins, Ureidopenicillins, Carboxypenicillins
Methicillin resistant staphylococci	All beta-lactams
Micrococcus spp.	Nitrofurantoin, Mupirocin
Streptococci/enterococci	
Streptococcus spp.	Polymyxins, Nalidixic acid, Aminoglycosides (low level)
Enterococcus faecalis	Cephalosporins, Clindamycin, Mupirocin, Aminoglycosides (low level - HLR test), Novobiocin, Trim+Sulfa (in vivo)
E. faecium	Cephalosporins, Aminoglycosides (low level - HLR test), Nitrofurantoin, Trim+Sulfa (in vivo)
E. gallinarum/casseliflavus	Vancomycin (MIC 4-16 µg/ml)
Arcanobacterium spp.	Bacitracin, Mupirocin, Optochin
Pediococcus/Leuconostoc	Glycopeptides
Lactobacillus/Erysipelothrix	

# Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

## Intrinsic (natural) resistance

## Intrinsic (natural) resistance

BACTERIA	NATURAL RESISTANCE
Non fermenters	
Acinetobacter baumannii/calcoaceticus	Aminopenicillins, 1 <sup>st</sup> and 2 <sup>nd</sup> gen Cephalosporins, Chloramphenicol, Trimethoprim, Fosfomycin, Nitrofurantoin
Achromobacter xylosoxidans	Aminopenicillins, 1 <sup>st</sup> , 2 <sup>nd</sup> and 3 <sup>rd</sup> gen Cephalosporins, Aminoglycosides, Aztreonam
Alc denitrificans	Cefotaxime
Burkholderia cepacia	Aminopenicillins, Ureidopenicillins, Carboxypenicillins, Amoxicillin+Clavulanate, 1 <sup>st</sup> and 2 <sup>nd</sup> gen Cephalosporins, Quinolones, Aminoglycosides, Polymyxins, Nitrofurantoin, Fosfomycin, Chloramphenicol, Imipenem, Trimethoprim
Elisabethkingia meningoseptica	Aminoglycosides, Carboxypenicillins, 1 <sup>st</sup> , 2 <sup>nd</sup> and 3 <sup>rd</sup> gen. Cephalosporins, Polymyxins, Tetracyclines, Chloramphenicol, Ticarcillin+Clavulanate, Quinolones, Imipenem
Ochrobactrum anthropi	Ureidopenicillins, Carboxypenicillins, Ticarcillin+Clavulanate, 3 <sup>rd</sup> gen. Cephalosporins, Aztreonam, Ertapenem
Pseudomonas aeruginosa	Aminopenicillins, Amoxicillin+Clavulanate, 1 <sup>st</sup> and 2 <sup>nd</sup> gen. Cephalosporins, Cefotaxime, Ceftriaxone, Chloramphenicol, Nalidixic acid, Trim+Sulfa, Tetracyclines, Nitrofurantoin
Stenotrophomonas maltophilia	Ureidopenicillins, Carboxypenicillins, 1 <sup>st</sup> and 2 <sup>nd</sup> gen Cephalosporins, Imipenem, Cefotaxime, Aztreonam, Aminoglycosides, Tetracyclines (except Minocycline) Fosfomycin
Listeria	Oxacillin, Cephalosporins, Aztreonam, Polymyxins, Nalidixic acid, Clindamycin, Fosfomycin
Neisseria/Branhamella	
Branhamella catarrhalis	Lincomycin, Clindamycin, Trimethoprim
Gonococci, meningococci	Lincomycin, Clindamycin, Polymyxins, Trimethoprim, Vancomycin
Campylobacter/Helicobacter	
Campylobacter spp.	Vancomycin, Trimethoprim, Polymyxins, Lincomycin, Novobiocin, Aztreonam
Helicobacter pylori	Vancomycin, Polymyxins, Nalidixic acid, Trimetoprim, Sulfonamides
Corynebacterium in general	Fosfomycin, Mupirocin, Polymyxins, Nalidixic acid
C. jeikeium/urealyticum	All Penicillins, 1 <sup>st</sup> , 2 <sup>nd</sup> and 3 <sup>rd</sup> gen. Cephalosporins, Amoxicillin+Clavulanate., Imipenem, Meropenem, Aminoglycosides, Chloramphenicol, Nalidixic acid, Trim+Sulfa, Polymyxins, Fosfomycin, Mupirocin, Macrolides
<b>Anaerobes</b> in general	Aminoglycosides, Aztreonam (except Fusobacteria), Trimethoprim, Nalidixic acid
Bacteroides fragilis group	Aminoglycosides, Vancomycin, Aminopenicillins, 1 <sup>st</sup> , and 2 <sup>nd</sup> gen. Cephalosporins, Polymyxins, Glycopeptides, Fosfomycin, Aztreonam, Oxgall
Clostridium spp.	Kanamycin, Trimethoprim, Aztreonam, Polymyxins, Fosfomycin
Fusobacteria spp.	Nalidixic acid, Vancomycin, Macrolides (low level)
Porphyromonas spp	Polymyxins, Fosfomycin, Aminoglycosides
Prevotella spp.	Glycopeptides, Fosfomycin, Aminoglycosides
Actinomyces/Propionibacterium	1 <sup>st</sup> , and 2 <sup>nd</sup> gen. Cephalosporins, Polymyxins, Metronidazole
Mobiluncus spp.	Metronidazole
Peptostreptococcus/Eubacterium	Polymyxins, Fosfomycin
Veillonella spp.	Macrolides (low level), Glycopeptides

### References:

- 1) Livermore DM et al. Interpretative reading: recognizing the unusual and inferring resistance mechanisms from resistance phenotypes. J Antimicrob Chemother. 2001; 48 Suppl 1:87-102.
- 2) Communiqué January 2008: Société Française de Microbiologie (CA-SFM)