


No	Subject	Rev. date	Document
1	Performance of Susceptibility Testing		
	Front page and preface	29-11-2012	1.1.0
	Characteristics of Neo-Sensitabs/Inoculum standardisation/Incubation and reading of plates	01-09-2009	1.2.0
	Storage of NEO-SENSITABS	16-09-2010	1.3.0
	Procedure according to CLSI for bacteria and yeasts/Inoculum standardisation/Incubation and reading of plates	18-03-2009	1.4.0
	Prediffusion method	31-01-2013	1.5.0
	Interpretation of results	16-09-2010	1.6.0
2	NEO-SENSITABS Range (Cartridges with Spring)		
	Antibacterials	06-02-2013	2.1.0
	Antifungals	11-04-2013	2.2.0
3	Interpretation Zones and MIC Breakpoints according to CLSI		
	Enterobacteriaceae	16-01-2013	3.1.0
	<i>P. aeruginosa</i> , <i>Acinetobacter</i> spp., <i>B. cepacia</i> , <i>S. maltophilia</i>	24-01-2013	3.2.0
	Staphylococci	24-01-2013	3.3.0
	Enterococci	24-01-2013	3.4.0
	Pneumococci	15-01-2013	3.5.0
	Streptococci (other than <i>S. pneumoniae</i>)	24-01-2013	3.6.0
	<i>Haemophilus</i> spp.	29-01-2013	3.7.0
	<i>Moraxella catarrhalis</i>	04-12-2012	3.8.0
	Meningococci	15-01-2013	3.9.0
	Gonococci	25-01-2013	3.10.0
	<i>Campylobacter</i>	24-01-2013	3.11.0
	<i>Helicobacter pylori</i>	15-01-2013	3.12.0
	<i>Vibrio cholerae</i>	04-12-2012	3.13.0
	Anaerobes	24-01-2013	3.14.0
	Yeasts (CLSI's M44-A2 Method of Diffusion on Agar)	03-12-2013	3.15.0
	Quality control and Control Limits on Mueller-Hinton Agar for Nonfastidious Organisms	15-01-2013	3.16.0
	Control Limits on Mueller-Hinton with Blood and $\pm$ Supplements for Fastidious Organisms	15-01-2013	3.16.1
	<i>Corynebacteria</i>	18-01-2013	3.17.0
4	EUCAST and ROSCO Interpretation		
	Susceptibility testing of Enterobacteriaceae using EUCAST potency Neo-Sensitabs and EUCAST breakpoints	18-04-2011	4.1.0
	<i>P. aeruginosa</i> , <i>Acinetobacter</i> spp, <i>S.maltophilia</i>	18-04-2011	4.2.0
	Staphylococci	18-04-2011	4.3.0
	Enterococci	18-04-2011	4.4.0
	Pneumococci	18-04-2011	4.5.0
	Streptococci (other than <i>S.pneumoniae</i>)	18-04-2011	4.6.0
	<i>Haemophilus</i> spp.	18-04-2011	4.7.0
	<i>Moraxella catarrhalis</i>	18-04-2011	4.8.0
	Meningococci	18-04-2011	4.9.0
	Gonococci	18-04-2011	4.10.0
	Yeast (tentative)	18-04-2011	4.12.0
	EUCAST recommended strains for internal Quality Control and control limits on MH agar and McF 0.5 inoculum	18-04-2011	4.13.0
	Quality control and control limits on MH agar+5% horse blood and 20mg/l β NAD according to EUCAST	18-04-2011	4.13.1

No	Subject	Rev. date	Document
5	EUCAST interpretation using MH agar and McF 1.0 inoculum		
	Anaerobes	20-02-2013	5.1.0
6	Interpretation according to MIC Breakpoints of SFM (France) Using CLSI Potency Neo-Sensitabs		
	Rapidly Growing Bacteria	18-01-2013	6.1.0
	Haemophilus spp., S. pneumoniae, Streptococcus spp., N. gonorrhoeae, N. meningitidis, Campylobacter spp. and Anaerobes	18-01-2013	6.2.0
7	Interpretation according to MIC Breakpoints of BSAC (UK and Ireland) using BSAC and CLSI potency Neo-Sensitabs		
	Rapidly growing bacteria	24-01-2013	7.1.0
	Haemophilus spp. S. pneumoniae, Streptococcus spp., N. gonorrhoeae, N. meningitidis, Moraxella catarrhalis, Coryneforms, Campylobacter spp., Pasteurella multocida and anaerobes	24-01-2013	7.2.0
8	Interpretation according to MIC breakpoints of BSAC (UK and Ireland) using BSAC potency Neo-Sensitabs		
	BSAC urines	24-01-2013	8.1.0
9	Veterinary practice according to CLSI breakpoints		
	Veterinary practice	11-04-2013	9.1.0

NEO-SENSITABS™

User's Guide

Susceptibility testing



EUCAST and CLSI potency Neo-Sensitabs™


2013

EUCAST and CLSI-potency Neo-Sensitabs™ User's Guide 2013

EUCAST and CLSI-potency NEO-SENSITABS™ User's Guide 2013 contains updated text, tables and references, all necessary information when using Neo-Sensitabs tablets for susceptibility testing.

Totally new interpretation tables using Mueller Hinton Agar and McFarland 0.5 inoculum according to the MIC breakpoints recommended by EUCAST are included.

The different interpretation tables following the CLSI (formerly NCCLS) recommendations have been updated according to the latest information of the CLSI described in "Performance Standards for Antimicrobial Disk Susceptibility Testing", 23rd Informational Suppl., **M100-S23**, 2013.



Furthermore, the User's Guide includes updated Zone Diameter Interpretative Standards according to national recommendation groups including interpretation according to MIC breakpoints recommended by the SFM and BSAC has been updated.

Other new aspects that have been revised in the updated User's Guide of the EUCAST and CLSI-potency NEO-SENSITABS™ are:

- 2 + 18 hours (or 2+22) hours prediffusion method for high molecular weight antimicrobials: Colistin, Daptomycin, Vancomycin and Teicoplanin has been included.
- Interpretative zone breakpoints for 2+18 hours' prediffusion method for Daptomycin, Vancomycin and Teicoplanin in the different tables.
- Prediffusion method with Colistin is recommended for both *A. baumannii* and *P. aeruginosa* strains.
- Kits for detection of ESBLs using Neo-Sensitabs have been developed.
- Kit for detecting Plasmid-mediated AmpC beta-lactamases, using Cefotaxime, Ceftazidime and Cloxacillin has been developed.
- Kits for detecting carbapenemases class A (KPC and GES enzymes), class D (metallo-betalactamases) and OXA-48 have been developed.
- Kit for detecting metallo-β-lactamases using 2 chelating agents EDTA and Dipicolinic acid has been developed.
- New drugs that will be introduced in the near future are: Ceftazidime+Avibactam, Ceftaroline+Avibactam, Ceftolozane+Tazobactam, CXA-101, Tedizolid (TR-700), Monosulfactam (BAL30072), Fidaxomicin (OPT-80), Solithromycin (CEM-101), Deformylase inhibitor and new Fluoroquinolones.

The User's Guide for EUCAST and CLSI-potency Neo-Sensitabs™ 2013, is available on our website www.rosco.dk and updated information is continuously included.

ROSCO Diagnostica A/S is welcoming any feedback and questions on susceptibility testing from users directly (info@rosco.dk) or through our representatives.

ROSCO Diagnostica A/S

Characteristics of Neo-Sensitabs

Neo-Sensitabs are 9 mm in diameter and each is print coded for safe identification. The tablets are manufactured with the aid of microbial inert auxiliary substances by a dry process using crystalline antimicrobials. This procedure results in very uniform tablets that are homogenous in their content of active ingredients and have an extraordinary stability, usually not less than 4 years shelf-life at room temperature (chapter 1).

Neo-Sensitabs and the MIC-breakpoints (chapter 3) are standardized according to potencies recommended by CLSI- Clinical and Laboratory Standards Institute (formerly NCCLS, EUCAST, SFM France (chapter 5) and BSAC (chapter 6). The size of the zone therefore is equivalent to the size of a paper disk, if the potency is the same.

Recently EUCAST (chapter 4) informed that they will develop a disk diffusion method on Mueller-Hinton agar. Neo-Sensitabs will follow the EUCAST recommendations.

All antimicrobials have received letter codes in order to achieve optimal recognition and zone measurements with automatic instruments a 5-digit code for each Neo-Sensitabs type has been chosen. The new Neo-Sensitabs with potency according to CLSI have of course different code than the old Neo-Sensitabs e.g. Erythromycin 78 µg (code: ERYTR) and Erythromycin 15 µg (code: ERY15).

The new ranges of Neo-Sensitabs are produced in cartridges with spring and are used together with the new dispenser from Rosco. The interpretation tables are available in this supplement; however, some additional information concerning Neo-Sensitabs is still only available in the User's Guide for Neo-Sensitabs Ed. 19, 2007/2008 (soon Ed. 20, 2009).

Dispenser

The dispensers are to be used with Rosco Neo-Sensitabs **cartridges with a spring**, where the potencies of the tablets are according to recommendations of CLSI (formerly NCCLS).

Models available from Rosco:

- 1) Adaptable to 8-10 cm petri dishes delivering up to 7 Neo-Sensitabs at a time (Dispenser 101).
- 2) Adaptable for square petri dishes (12 x 12 cm) delivering up to 16 Neo-Sensitabs at a time (Dispenser 104).
- 3) A single cartridge dispenser.

The Dispensers are made of hard plastic and are operated by pushing the handle down, and the Neo-Sensitabs will be transferred to the agar surface.

When using several dispensers, the colour code on the top of the handles can be used for differentiation.

The holes in the bottom plate ensure that the tablets are placed onto the agar surface in a pre-determined pattern.

The dispensers are easy to disassemble for inside cleaning. They must be cleaned occasionally (wipe with ethanol and hot water to remove dust from the Neo-Sensitabs tablets).

The tablets come packed in cartridges (tubes), matching the dispenser top-plate holes. Each cartridge accommodates 50 Neo-Sensitabs and a red block that prevents dispensing when one cartridge is empty. Insert the cartridges gently and carefully one by one through the top-plates holes. For easy identification the bottom of each tube is labeled with a short name (5-digit alpha-numeric code, also marked on each tablet) of the antimicrobial contained in the cartridge.

Further information is available on www.rosco.dk

Storage instructions

1. On receipt check the temperature symbol on the outer container.
Neo-Sensitabs with a 2 °C to 8 °C symbol should be stored in a refrigerator and Neo-Sensitabs with a 25 °C as an upper temperature symbol on the outer container should be stored at room temperature ($\leq 25^{\circ}\text{C}$).
2. If Neo-Sensitabs are stored in the refrigerator, allow the cartridge to reach room temperature ($\leq 25^{\circ}\text{C}$) before being opened, i.e. 30 – 60 minutes, in order to avoid water condensation on the tablets.
3. Neo-Sensitabs with the temperature symbol 2 °C to 8 °C may be left at room temperature ($\leq 25^{\circ}\text{C}$) for up to 2 months, without essential loss of activity.
4. **Opened cartridges** placed in a dispenser must be used within 2 months for Neo-Sensitabs with the temperature symbol 2 °C to 8 °C, and within 12 months for Neo-Sensitabs with the temperature symbol below 25 °C. **The dispenser must be kept at room temperature.**

The stability of antimicrobials in paper disks is decreased compared to Neo-Sensitabs. The CLSI (1) recommends frozen storage of paper disks containing beta-lactam antibiotics. In case of Imipenem, Cefaclor and Clavulanic acid combinations, paper disks should be stored frozen until the day of use. In a comparative stability study between Neo-Sensitabs and Oxoid paper disks (2), it was observed that disks containing Ticarcillin 75 µg, lost activity after 15 days at 4-6 °C, while Ampicillin 10 µg and Amoxycillin + Clavulanate 20+10 µg disks lost activity after one month at 4-6 °C. The corresponding Neo-Sensitabs were stable at room temperature (and at 4-6 °C) during the study period of six months.

Steward et al (3) noticed overdetection of imipenem/meropenem resistance in the project ICARE, most probably due to inactivation of the reagents used (Vitek, disk diffusion etc.) and recommended the use of a second independent antimicrobial susceptibility testing method to validate carbapenem intermediate and resistant strains.

References:

- 1) NCCLS: Performance Standards for Antimicrobial Disk Susceptibility Testing, 10th Ed. **M2-A10**, January 2009
- 2) del Cuerpo M. et al: Stability of beta-lactam antibiotics in paper disks and tablets used in the diffusion test. Rev. Esp. Quimioter., **10**, nr. 3, 1997 (Spanish).
- 3) Steward C.D. et al: Antimicrobial susceptibility testing of carbapenems: multicenter validity testing and accuracy levels of 5 antimicrobial test methods for detecting resistance in Enterobacteriaceae and Pseudomonas aeruginosa isolates. J. Clin. Microbiol., **41**, 351-358, 2003.

Performance of Susceptibility Testing

Procedure according to CLSI for bacteria and yeasts/Inoculum standardisation/Incubation and reading of plates

Inoculum

When using the technique of Kirby-Bauer, the inoculum is standardized according to the method described by the CLSI (chapter 2 and 3), which results in confluent growth for bacteriae and semi-confluent growth for most of the *Candida* species isolates.

Inoculum standardisation

Direct colony suspension method:

Suspend several morphologically similar colonies from an 18-24 h agar plate (non selective) into 4-5 ml sterile saline solution, and then immediately adjusting the turbidity to match that of the BaSO₄ standard (0.5 McFarland). For *Candida* spp. that are subcultured onto blood agar or Sabouraud dextrose agar five distinct colonies from a 24-hour-old culture are suspended in 5 ml sterile saline (0.145 mol/L; 8.5 g/L NaCl; 0.85% saline).

- Within 15 minutes, dip a sterile cotton swab into the adjusted suspension and remove inoculum from the swab by exerting firm pressure on the inside of the tube.
- Within 15 minutes swabs are used to inoculate the test plates.
- Inoculate the dried surface of the appropriate agar plate by streaking the swab over the entire surface. Allow the surface to dry 3-5 minutes or maximum 15 minutes before applying Neo-Sensitabs to the media.
- Select appropriate tablets e.g. such as recommended by CLSI (NCCLS).² Use no more than nine Neo-Sensitabs per 150 mm plate or four Neo-Sensitabs per 100 mm plate when testing *Candida* spp., *H. influenzae*, *N. gonorrhoeae*, and *Streptococcus* spp.
- Neo-Sensitabs is dispensed onto the surface of the inoculated agar plate.

Incubation and reading of plates

- Within 15 minutes, place the agar side up in a 35 °C (± 2 °C) incubator. For some strains special recommendations are noted in the tables e.g. *Haemophilus* spp., *N. gonorrhoeae*, *S. pneumoniae* and other streptococci should be incubated in an atmosphere enriched with 5 % CO₂. Incubate no more than 5 plates in a stack. Plates in the middle of the stack will take longer to reach the desired incubation temperature than plates at the top and the bottom.
- Examine the plates after 16-18 hours incubation (20-24 h for *Candida* spp. *N. gonorrhoeae*, *S. pneumoniae* and other streptococci). Full 24-hour incubation is recommended for the detection of Methicillin resistant *Staphylococcus aureus* (MRSA) and *Enterococcus* spp. for vancomycin resistance. Hold the plate up to transmitted light and examine the oxacillin, linezolid and vancomycin zones for light growth (minute colonies) of methicillin, linezolid or vancomycin resistant colonies, respectively, within apparent zones of inhibition. Any discernible growth within the zone of inhibition is indicative of methicillin, linezolid or vancomycin-resistance. The edges of the zones of inhibition contain a large number of small colonies when using Trimethoprim, Sulphonamides and Trimethoprim + Sulfamethoxazole tablets. In this case zones of inhibition are measured up to colonies of normal size (disregard slight growth and measure the more obvious margin).
- For some antimicrobial agents (chloramphenicol, clindamycin, erythromycin and tetracycline) the zones of inhibition will contain a gradient of growth. In this case zone diameters should be read half-way between the start of inhibition and complete inhibition.
- With *Proteus* spp. ignore the thin veil of swarming growth in an otherwise obvious zone of inhibition. The diameters of the zones of complete inhibition are measured as determined by gross visual inspection. Zones are measured to the nearest whole millimeter.
- The measured zone diameters of inhibition are compared with the zone interpretative tables for the individual antibiotics in order to determine the agent(s) most suitable for use in antimicrobial therapy.

Pre-diffusion Method (2 + 18 or 2+22 hours) for Antimicrobials Diffusing Poorly on Agar

High molecular weight antimicrobials (Vancomycin, Teicoplanin, Daptomycin, Colistin) diffuse poorly on agar media, resulting in difficult to interpret results when using the disc diffusion method. Rosco Diagnostica has developed a 2 + 18 or 2 + 22 hour prediffusion technique, permitting an easier differentiation between susceptible and resistant strains when testing against these antimicrobials.

Procedure

One Neo-Sensitabs of the antimicrobial to be tested is placed on an uninoculated plate containing the susceptibility test medium (Mueller-Hinton plain or BHI Agar + 5 % blood). The plate is then placed in the incubator at 37°C for 2 hours. After 2 hours the tablet (disc) is removed (by knocking the plate against the table) and the plate is maintained at room temperature for further 18 or 22 hours (overnight).

The plate is now inoculated with the strain to be tested using a McFarland 0.5 inoculum.

Additional antibiotic discs (Neo-Sensitabs) may be added now using a dispenser (if MH agar is used) and thereafter the plate is incubated overnight at 35-37 °C.

The zones of inhibition are then measured. Zone breakpoints are tentative and for research use only.

Notice: In the laboratory, the prediffusion plate can be prepared the day before it is inoculated, in which case there is no loss of time and results are obtained within 24 hours.

Interpretation

IA) Detection of Visa/GISA/hVISA strains (medium BHI + 5 % blood or MH plain)

VISA/GISA/hVISA strains will show the following zones of inhibition:

BHI+5% blood

Teicoplanin 30 µg: inhibition zone < 20 mm and/or
Vancomycin 30 µg: inhibition zone < 20 mm

MH plain

Teicoplanin R < 20 mm
Vancomycin R ≤ 20 mm

IB) Detection of GISCN/hGISCH strains (medium BHI + 5 % blod)

Teicoplanin 30 µg: inhibition zone < 20 mm

GISCN = Glycopeptide intermediate staphylococci, coagulase negative (2).

Strains showing zones of inhibition < 20 mm around Teicoplanin 30 µg should be reported as heteroresistant to both Teicoplanin and Vancomycin.

IC) Detection of VanA, VanB and VanC in enterococci

Use Mueller-Hinton Agar (without blood add), McFarland 0.5 inoculum.

Vancomycin 30 µg (2+18 hours' prediffusion method):

Susceptible: zone > 16 mm (sharp edge)
VanB: zone < 16 mm (hazy edge)
VanC: zone < 12mm (sharp edge)

The VanA genotype will show **no zone** of inhibition in the current diffusion test with Vancomycin 30 µg Neo-Sensitabs.

ID) Detection of VRE, vanB phenotype and vanA genotype

Use Mueller-Hinton Agar, McFarland 0.5 inoculum

Vancomycin 30 µg : no zone

Teicoplanin 30 µg : zone < 16 mm (MIC 4-12 µg/ml). Report as R

II) Daptomycin testing (medium used Mueller-Hinton plain)

a) Staphylococci:

Daptomycin 30 µg Neo-Sensitabs: Susceptible zone ≥ 22 mm (corresponding to an MIC of
≤ 1 µg/ml). Resistant zone < 20 mm (MIC ≥ 2 µg/ml)

b) Enterococci:

Daptomycin 30 µg Neo-Sensitabs: Susceptible zone ≥ 12 mm (corresponding to an MIC of
≤ 4 µg/ml).

III) Colistin testing (medium used Mueller-Hinton plain)

Colistin 10 µg Neo-Sensitabs: Susceptible zone ≥ 15 mm (corresponding to an MIC of ≤ 2 µg/ml).

References:

- 1) Nielsen S.V., Casals J.B.: Detection of decreased susceptibility to glycopeptides in *S. aureus* using tablet (disc) prediffusion. 15th Eur. Cong. Clin. Microbiol. Inf. Dis. (ECCMID), April 2005.
- 2) Ferreira Nunes AP et al: Heterogeneous resistance to vancomycin in *S. epidermidis*, *S. haemolyticus* and *S. warneri* clinical strains: characterisation of glycopeptide susceptibility profiles and cell wall thickening. Intl. J. Antimicrob. Ag., **27**, 307-315, 2006.
- 3) Katz B.D. et al: A new pre-diffusion method for the detection of Daptomycin (DAP) non-susceptible strains using Neo-Sensitabs. Presentation D-226, ICAAC September 2007, Chicago, USA.
- 4) Borda N. et al: Comparison of methods: diffusion (DF), prediffusion (PDF) and E-test on isolates of *Ac. baumannii-calcoaceticus* complex (Abc) against colistin. 2007 (in press).
- 5) Katz B.D. et al: Detection of daptomycin- non-susceptible strains using Neo-Sensitabs trade prediffusion method. Diagn. Microbiol. Infect. Dis. **61**. 315-320, 2008.
- 6) Koeth L. et al: Multisite evaluation of the Daptomycin Neo-Sensitab prediffusion method against 20 *S. aureus*. ECCMID, Milan 2011.
- 7) Boyen F. et al: Disk prediffusion is a reliable method for testing Colistin susceptibility in porcine *E. coli* strains. Vet. Microbiol. **144**, 359-362, 2010.

Note:

Detailed description of the prediffusion methods for Vancomycin, Teicoplanin, Daptomycin and Colistin see: "Detection of resistance mechanisms using Neo-Sensitabs and Diatabs" documents 7.1.0; 7.2.0; 7.3.0; 7.4.0 and 7.5.0.

Susceptible (S):

The infection due to the strain tested may be expected to respond to a normal dosage of this antimicrobial.

Intermediate (I)

The intermediate category implies clinical applicability in body sites where the drugs are concentrated (e.g. urine) or when high dosage of an antimicrobial can be used (e.g. betalactams).

The intermediate category also comprises a "buffer zone" which should prevent small uncontrolled technical factors from causing major discrepancies in interpretations; thus, when a zone falls within the intermediate range, the results may be considered equivocal, and if alternative drugs are not available MIC testing may be indicated.

Resistant (R):

The antimicrobial cannot be recommended for treatment in this case.

If only "S" criteria are specified:

For some organism/antimicrobial combinations, the absence of resistant strains precludes defining any category other than susceptible. For strains yielding results suggestive of "non susceptible", organism identification and antimicrobial susceptibility test results should be confirmed. Subsequently the isolates should be submitted to a Reference Laboratory for further testing.

References:

- 1) Gylling Pedersen O.: Standardizing, manufacture, and control of Neo-Sensitabs. Acta Clin. Belg., **28**, 139-149, 1973.
- 2) CLSI: Performance Standards for Antimicrobial Susceptibility Testing 20th Inf. Suppl. **M100-S20**, 2010.
- 3) NCCLS: Performance Standards for Antimicrobial Disk Susceptibility Testing. 8th Ed. **M2-A9**, 2006.
- 4) NCCLS: Method for Antifungal Disk Diffusion Susceptibility Testing. **M44-A2**, 2008.
- 5) Kauppila J. et al: Comparison of Neo-Sensitabs with paper disks in the routine disk diffusion antimicrobial susceptibility testing. ECCMID poster P-865, Barcelona 2008.
- 6) Kahlmeter G: Implementation of European breakpoints and the future of EUCAST. ECCMID abstract 512, april 2008, Barcelona.

EUCAST-and CLSI potency NEO-SENSITABS™ NEO-SENSITABS Range (Cartridges with Spring)

Page 1 / 4 - Document: 2.1.0

Antibacterials

Antibacterials

Neo-Sensitabs	Identification Code Neo-Sensitabs	Diffusible Amount of Antimicrobial
A. Penicillins (Penams)		
PENICILLIN 1 UNIT EUCAST	PENG1	1 U
PENICILLIN	PEN10	10 U
AMPICILLIN	AMP10	10 µg
AMOXYCILLIN 25 µg BSAC (U)	AMX25	25 µg
AMOXYCILLIN	AMOXY	30 µg
AMPICILLIN 2 µg EUCAST	AMP.2	2 µg
OXACILLIN 1 µg	OXA.1	1 µg
MECILLINAM (Amdinocillin)	MEC10	10 µg
TICARCILLIN	TIC75	75 µg
PIPERACILLIN	PIPER	100 µg
PIPERACILLIN EUCAST	PIP30	30 µg
TEMOCILLIN	TEMOC	30 µg
B. Beta-lactam / Beta-lactamase Inhibitor Combinations		
AMOXYCILLIN+CLAVULANATE 20+10 µg	AMC30	20+10 µg
AMOXYCILLIN+CLAVULANATE 2+1 µg (Augmentin)	AMC.3	1+2 µg
AMPICILLIN+SULBACTAM	SAM20	10+10 µg
TICARCILLIN+CLAVULANATE (Timentin)	TIM85	75+10 µg
PIPERACILLIN+TAZOBACTAM	PI+TZ	100+10 µg
PIPERACILLIN+TAZOBACTAM BSAC	PTZ85	75+10 µg
CEFOTAXIME+CLAVULANATE	CTX+C	30+10 µg
CEFTAZIDIME+CLAVULANATE	CAZ+C	30+10 µg
CEFEPIME+CLAVULANATE	FEP+C	30+10 µg
CEFPODOXIME+CLAVULANATE*	CPD+C	10+1 µg
PIPERACILLIN +TAZOBACTAM EUCAST	PTZ36	30+6 µg
C. (1) Cephalosporins (Cephems)		
CEPHALOTHIN	CEP30	30 µg
CEFACTOR	CCL30	30 µg
CEFAZOLIN	CFZ30	30 µg
CEPHALEXIN BSAC (U)	CFLEX	30 µg
CEFADROXIL	CDX30	30 µg
CEFUROXIME	CXM30	30 µg
CEFIXIME	CFM.5	5 µg
CEFPODOXIME	CPD10	10 µg
CEFOTAXIME	CTX30	30 µg
CEFTAZIDIME	CAZ30	30 µg
CEFTRIAXONE	CTR30	30 µg
CEFTIZOXIME	ZOX30	30 µg
CEFEPIME	FEP30	30 µg
CEFOTAXIME EUCAST	CTX.5	5 µg
CEFTAZIDIME EUCAST	CAZ10	10 µg
CEFTIBUTEN EUCAST	CTB30	30 µg
C. (2) Cephamycins and Oxacephems		
CEFOXITIN	CFO30	30 µg
CEFOXITIN**	CFO10	10 µg

EUCAST-and CLSI potency NEO-SENSITABS™ NEO-SENSITABS Range (Cartridges with Spring)

Neo-Sensitabs	Identification Code Neo-Sensitabs	Diffusible Amount of Antimicrobial
C. (3) Cephalosporins active against MRSA		
CEFTOBIPROLE (investigational drug)	CFBIP	30 µg
CEFTAROLINE (investigational drug)	CPT30	30 µg
D. Penems and Carbapenems		
IMIPENEM	IMI10	10 µg
MEROPENEM	MRP10	10 µg
ERTAPENEM	ETP10	10 µg
DORIPENEM	DOR10	10 µg
E. Monobactams		
AZTREONAM	AZT30	30 µg
F. Aminoglycosides		
STREPTOMYCINS 10 µg	STR10	10 µg
STREPTOMYCIN 500 µg (HLR)**	ST500	500 µg
KANAMYCIN 30 µg	KAN30	30 µg
KANAMYCIN 500 µg (HLR)**	KA500	500 µg
NEOMYCIN (Framycetin)*	NEOMY	120 µg
AMIKACIN	AMI30	30 µg
GENTAMICIN 10 µg	GEN10	10 µg
GENTAMICIN 30 µg EUCAST	GEN30	30 µg
GENTAMICIN 250 µg (HLR)**	GN250	250 µg
NETILMICIN 10 µg EUCAST	NET10	10 µg
NETILMICIN 30 µg	NET30	30 µg
TOBRAMYCIN	TOB10	10 µg
SPECTINOMYCIN	SPECT	200 µg
G. Tetracyclines		
TETRACYCLINES 30 µg (Oxytetracycline)	TET30	30 µg
DOXYCYCLINE	DOX30	30 µg
MINOCYCLINE	MIN30	30 µg
TIGECYCLINE	TIG15	15 µg
H. Chloramphenicol and Derivatives		
CHLORAMPHENICOL 30 µg	CLR30	30 µg
I. Macrolides, Lincosamides, Streptogramins, Ketolides and Oxazolidinones		
ERYTHROMYCIN	ERY15	15 µg
AZITHROMYCIN	AZI15	15 µg
CLARITHROMYCIN	CLA15	15 µg
CLINDAMYCIN	CLI.2	2 µg
QUINUPRISTIN/DALFOPRISTIN	SYN15	15 µg
TELITHROMYCIN	TEL15	15 µg
LINEZOLID	LINEZ	30 µg
LINEZOLID EUCAST	LIZ10	10 µg
J. (1) Glycopeptides		
VANCOMYCIN 5 µg EUCAST	VAN.5	5 µg
VANCOMYCIN 30 µg	VAN30	30 µg
J. (2) Lipoglycopeptides		
TELAVANCIN (investigational drug)	TLV30	30 µg
TEICoplanin	TPN30	30 µg

EUCAST-and CLSI potency NEO-SENSITABS™ NEO-SENSITABS Range (Cartridges with Spring)

Page 3 / 4 - Document: 2.1.0

Antibacterials

Neo-Sensitabs	Identification Code Neo-Sensitabs	Diffusible Amount of Antimicrobial
DALBAVANCIN (investigational drug)	DAL60	60 µg
J. (3) Cyclic Lipopeptides		
a) Gram positive spectrum: DAPTOMYCIN (+ Ca)	DAPCa	30 µg
b) Gram negative spectrum: COLISTIN 10 µg	CO.10	10 µg
POLYMYXINS (not for susceptibility testing)	CO150	150 µg
K. Sulphonamides and Similars		
SULPHONAMIDES*	SULFA	240 µg
TRIMETHOPRIM	TRIM5	5 µg
TRIMETHOPRIM BSAC (U)	TP2.5	2.5 µg
TRIMETHOPRIM+SULFA	SxT25	1.25+23.75 µg
L. Nitrofurans		
NITROFURANTOIN	NI300	300 µg
NITROFURANTOIN EUCAST	NI100	100 µg
NITROFURANTOIN BSAC (U)	NI200	200 µg
FURAZOLIDONE*	FURAZ	50 µg
M. Quinolone Derivatives		
NALIDIXAN	NAL30	30 µg
CIPROFLOXACIN 5 µg	CIPR5	5 µg
CIPROFLOXACIN 1 µg	CIPR1	1 µg
MOXIFLOXACIN	MOXIF	5 µg
GATIFLOXACIN	GATIF	5 µg
LEVOFLOXACIN	LEVOF	5 µg
NORFLOXACIN	NORFX	10 µg
NORFLOXACIN BSAC (U)	NOR.2	2 µg
OFLOXACIN	OFL.5	5 µg
N. Others		
BACITRACIN (not for susceptibility testing)	BACIT	40 µg
FOSFOMYCIN (+G6-P)	FO200	200 µg
FUCIDIN 10 µg EUCAST	FUC10	10 µg
METRONIDAZOLE 5 µg EUCAST	MTR.5	5 µg
METRONIDAZOLE 16 µg*	MTR16	16 µg
MUPIROCIN 10 µg*	MUPIR	10 µg
MUPIROCIN 200 µg EUCAST	MP200	200 µg
NOVOBIOCIN 5 µg*	NOV05	5 µg
RETAPAMULIN	RETA2	2 µg
RIFAMPICIN	RIF.5	5 µg
Special tests		
CLOXACILLIN (AmpC test) Diatabs**	CLOXA	-
PHENYLBORONIC ACID (AmpC and KPC test) Diatabs**	BORON	-
DIPICOLINIC ACID (metallo-β- lactamases) Diatabs**	D.P.A	-
IMIPENEM+EDTA **	IM10E	10+750 µg

Note:

Revision date: 06.02.13
Rev./Appr. by: Maria / TP
Replaces date: 16.09.10

Establish date: 01.08.07
Establish by: JBC
Revision no.: 9

ROSCO
DIAGNOSTICA

**EUCAST-and CLSI potency
NEO-SENSITABS™
NEO-SENSITABS Range (Cartridges
with Spring)**

Page 4 / 4 - Document: 2.1.0

Antibacterials

- * There are no potency recommendations from CLSI so far.
- ** Special potency tablets (discs) for detection of resistance mechanisms

ORIGINAL ROSCO DOCUMENT

Revision date: 06.02.13
Rev./Appr. by: Maria / TP
Replaces date: 16.09.10

Establish date: 01.08.07
Establish by: JBC
Revision no.: 9

ROSCO
DIAGNOSTICA

**EUCAST-and CLSI potency
NEO-SENSITABS™
NEO-SENSITABS Range (Cartridges
with Spring)**

Page 1 / 1 - Document: 2.2.0

Antifungals

Neo-Sensitabs	Identification Code Neo-Sensitabs	Diffusible Amount of Antimicrobial
AMPHOTERICIN B*	AMPHO	10 µg
CICLOPIROX*	CICLO	50 µg
CLOTRIMAZOLE*	CLOTR	10 µg
ECONAZOLE*	ECONZ	10 µg
FLUCONAZOLE	FLUCZ	25 µg
5-FLUOROCYTOSINE 1 µg*	FLU.1	1 µg
ISOCONAZOLE*	ISOCN	10 µg
ITRACONAZOLE*	ITRAC	10 µg
KETOCONAZOLE*	KETOC	15 µg
MICONAZOLE*	MICON	10 µg
NATAMYCIN*	NATAM	50 µg
NYSTATIN*	NYSTA	50 µg
TERBINAFINE*	TERBI	30 µg
VORICONAZOLE	VOR.1	1 µg
CASPOFUNGIN	CASP5	5 µg
POSACONAZOLE	POSAC	5 µg

* There are no potency recommendations from CLSI so far.

EUCAST-and CLSI potency NEO-SENSITABS™ Interpretation Zones and MIC Breakpoints according to CLSI

ENTEROBACTERIACEAE

Zone diameter interpretative criteria and MIC breakpoints according to CLSI (formerly NCCLS) (chapter 1) when testing Enterobacteriaceae are listed in the table below. (1)

Table 3.1-1 Interpretation for Enterobacteriaceae

Mueller-Hinton agar. Inoculum: McFarland 0.5. Incubation for 16-18 hours ambient air at 35°C ± 2 degrees.

				Zone diameter in mm			Break-points MIC µg/ml	
NEO-SENSITABS		POTENCY	CODE	S	I	R	S	R
k)	Amikacin	30 µg	AMI30	≥ 17	16-15	≤ 14	≤ 16	≥ 32
	Amoxycillin+Clav.	20+10 µg	AMC30	≥ 18	17-14	≤ 13	≤ 8/4	≥ 32/16
e)	Ampicillin	10 µg	AMP10	≥ 17	16-14	≤ 13	≤ 8	≥ 32
	Ampicillin+Sulbactam	10+10 µg	SAM20	≥ 15	14-12	≤ 11	≤ 8/4	≥ 32/16
m)	Aztreonam	30 µg	AZT30	≥ 21	20-18	≤ 17	≤ 4	≥ 16
b)	ESBL screening			-	-	≤ 27	Screening ESBL	
d)	Cefazolin	30 µg	CFZ30	≥ 23	22-20	≤ 19	≤ 2	≥ 8
	Cefepime	30 µg	FEP30	≥ 18	17-15	≤ 14	≤ 8	≥ 32
	Cefepime+Clavulanate	30+10 µg	FEP+C	Detection of ESBL				
	Cefixime	5 µg	CFM.5	≥ 19	18-16	≤ 15	≤ 1	≤ 1
d)	Cefotaxime	30 µg	CTX30	≥ 26	25-23	≤ 22	≥ 1	≥ 4
b)	ESBL screening			-	-	≤ 27	Screening ESBL	
	Cefotaxime+Clav.	30+10 µg	CTX+C	Detection of ESBL				
	Cefoxitin	30 µg	CFO30	≥ 18	17-15	≤ 14	≤ 8	≥ 32
	Cefpodoxime	10 µg	CPD10	≥ 21	20-18	≤ 17	≤ 2	≥ 8
	ESBL screening			-	-	≤ 17	Screening ESBL	
	Ceftaroline	30 µg	CPT30	≥ 23	22-20	≤ 19	≤ 0.5	≥ 2
m) *)	Ceftazidime	30 µg	CAZ30	≥ 21	20-18	≤ 17	≤ 4	≥ 16
b)	ESBL screening			-	-	≤ 22	Screening ESBL	
	Ceftazidime+Clav.	30+10 µg	CAZ+C	Detection of ESBL				
	Ceftizoxime	30 µg	ZOX30	≥ 25	24-22	≤ 21	≤ 1	≥ 4
d)	Ceftriaxone	30 µg	CTR30	≥ 23	22-20	≤ 19	≤ 1	≥ 4
b)	ESBL screening			-	-	≤ 25	Screening ESBL	
	Cefuroxime (parenteral)	30 µg	CXM30	≥ 18	17-15	≤ 14	≤ 8	≥ 32
	Cefuroxime (oral)	30 µg	CXM30	≥ 23	22-15	≤ 14	≤ 4	≥ 32
d)	Cephalothin	30 µg	CEP30	≥ 18	17-15	≤ 14	≤ 8	≥ 32
	Chloramphenicol	30 µg	CLR30	≥ 18	17-13	≤ 12	≤ 8	≥ 32
	Ciprofloxacin	5 µg	CIPR5	≥ 21	20-16	≤ 15	≤ 1	≥ 4
c)	Salmonella spp.			≥ 31	30-21	≤ 20	≤ 0.06	≥ 1
	Ciprofloxacin	1 µg	CIPR1					
c)	Salmonella spp.			≥ 26	-	< 26	≤ 0.06	≥ 0.12
h)	Cloxacillin		CLOXA	Detection of plasmid mediated AmpC			beta-lactamases	
	Colistin	10 µg	Co. 10					
	2+18 hour pre-diffusion			≥ 15	-	<15	≤ 2	>4
i)	Doxycycline	30 µg	DOX30	≥ 14	13-11	≤ 10	≤ 4	≥ 16
	Doripenem	10 µg	DOR10	≥ 23	22-20	≤ 19	≤ 0.25	≥ 1
j)	Ertapenem	10 µg	ETP10	≥ 23	21-19	≤ 18	≤ 0.5	≥ 1

EUCAST-and CLSI potency NEO-SENSITABS™

Interpretation Zones and MIC Breakpoints according to CLSI

Enterobacteriaceae

NEO-SENSITABS				Zone diameter in mm			Break-points MIC µg/ml	
POTENCY				S	I	R	S	R
CODE								
Carbapenemase screen				-	-	≤ 22	≤ 0.5	-
							(possible carbapenemase)	
a)	Fosfomycin (inj)	200 µg	FO200	≥ 22	21-18	< 18	≤ 16	≥ 32
	Fosfomycin (U)	200 µg	FO200	≥ 16	15-13	≤ 12	≤ 64	≥ 256
	Gatifloxacin	5 µg	GATIF	≥ 18	17-15	≤ 14	≤ 2	≥ 8
k)	Gentamicin	10 µg	GEN10	≥ 15	14-13	≤ 12	≤ 4	≥ 8
g) j)	Imipenem	10 µg	IMI10	≥ 23	22-20	≤ 19	≤ 1	≥ 4
Carbapenemase screen				-	-	<23	≤ 1	-
							(possible carbapenemase)	
g)	Imipenem+EDTA 10+750 µg		IM10E	Detection of metallo-β-lactamases				
	Kanamycin	30 µg	KAN30	≥ 18	17-14	≤ 13	≤ 6	≥ 25
	Levofloxacin	5 µg	LEVOF	≥ 17	16-14	≤ 13	≤ 2	≥ 8
	Salmonella spp.			≥ 29	28-19	≤ 18	≤ 0.125	≥ 2
a)	Mecillinam (U)	10 µg	MEC10	≥ 15	14-12	≤ 11	≤ 8	≥ 32
j)	Meropenem	10 µg	MRP10	≥ 23	22-20	≤ 19	≤ 1	≥ 4
Carbapenemase screen				-	-	<25	≤ 0.5	-
							(possible carbapenemase)	
i)	Minocycline	30 µg	MIN30	≥ 16	15-13	≤ 12	≤ 4	≥ 16
	Moxifloxacin	5 µg	MOXIF	≥ 19	18-16	≤ 15	≤ 2	≥ 8
a)	Nalidixan (U)	30 µg	NAL30	≥ 19	18-14	≤ 13	≤ 16	≥ 32
c)	Screening quinolones			-	-	< 15	-	Reduced suscept. to quinolones
	Netilmicin	30 µg	NET30	≥ 15	14-13	≤ 12	≤ 12	≥ 32
a) f)	Nitrofurantoin (U)	300 µg	NI300	≥ 17	16-15	≤ 14	≤ 32	≥ 128
a)	Norfloxacin (U)	10 µg	NORFX	≥ 17	16-13	≤ 12	≤ 4	≥ 16
	Ofloxacin	5 µg	OFL.5	≥ 16	15-13	≤ 12	≤ 2	≥ 8
	Salmonella spp.						≤ 0.12	≥ 2
	Piperacillin	100 µg	PIPA	≥ 21	20-18	≤ 17	≤ 16	≥ 128
	Piperacillin+Tazobactam	100+10µg	PI+TZ	≥ 21	20-18	≤ 17	≤ 16/4	≥ 128/4
	Streptomycin	10 µg	STR10	≥ 15	14-12	≤ 11	-	-
a)	Sulphonamides (U)	240 µg	SULFA	≥ 17	16-13	≤ 12	≤ 256	≥ 512
i)	Temocillin	30 µg	TEMOC	≥ 18	17-15	≤ 14	≤ 16	≥ 32
i)	Tetracyclines	30 µg	TET30	≥ 15	14-12	≤ 11	≤ 4	≥ 16
	Ticarcillin	75 µg	TIC75	≥ 20	19-15	≤ 14	≤ 16	≥ 128
	Ticarcillin+Clavulanate	75+10 µg	TIM85	≥ 20	19-15	≤ 14	≤ 16/2	≥ 128/2
n)	Tigecycline	15 µg	TIG15	≥ 19	18-15	≤ 14	≤ 2	≥ 8
k)	Tobramycin	10 µg	TOB10	≥ 15	14-13	≤ 12	≤ 4	≥ 8
a)	Trimethoprim (U)	5 µg	TRIM5	≥ 16	15-11	≤ 10	≤ 4	≥ 16
	Trimethoprim+Sulfa	1.25+23.75 µg	SxT25	≥ 16	15-11	≤ 10	≤ 2/38	≥ 8/152
Boronic acid				BORON			Detection of AmpC and KPC Beta.lactamases	
Dipicolinic acid				D.P.A			Detection of metallo-β-lactamases	

a) (U) For urinary tract infections only.

- b) Strains of *Klebsiella*, *E. coli*, *Salmonella* that **produce ESBL**, may be clinically resistant to therapy with penicillins, cephalosporins or aztreonam, despite apparent in vitro susceptibility, see "Detection of resistance mechanisms using Neo-Sensitabs™ and DiaTabs™".

Strains showing zones ≤ 27 mm with Aztreonam and/or Cefotaxime, ≤ 17 mm with Cefpodoxime, ≤ 22 mm with Ceftazidime and/or ≤ 25 mm with Ceftriaxone should be suspected of producing ESBL. For confirmatory tests, use Ceftazidime $\pm$ Clavulanate, Cefotaxime $\pm$ Clavulanate, Cefepime $\pm$ Clavulanate and/or Cefpodoxime $\pm$ Clavulanate.

An increase in zone diameter of ≥ 5 mm for the combination Cefpodoxime+Clavulanate, Cefotaxime+Clavulanate, Ceftazidime+Clavulanate or Cefepime+Clavulanate compared to Cefpodoxime/Cefotaxime/ Ceftazidime/Cefepime alone is confirmatory of the presence of an ESBL. If resistant to cefotaxime, ceftazidime and ceftriaxone, negative for ESBL and susceptible for cefepime, report as found.

Using the new CLSI breakpoints for cephalosporins (9, 10) a substantial number of *E. coli* and *P. mirabilis* containing ESBLs would be reported as susceptible to CAZ, FEP and AZT.

For *K. oxytoca* and *C. koseri* showing **R** to Cefuroxime, Ceftriaxone and Piperacillin/Tazobactam, but susceptible to ceftazidime, are not ESBL producers (false positives may be obtained using 3rd generation Cephalosporins). Test for ESBL production with Ceftazidime and Clavulanate.

- c) Strains resistant to nalidixic acid show a decreased susceptibility to quinolones (MIC CIPRO ≥ 0.125 µg/ml). Nalidixic makes a good screening for detection of decreased fluoroquinolone susceptibility in ***Salmonella* spp.** (4).

According to Hakanen et al. (2) and Parry et al. (7) *Salmonella enterica* isolates from Southeast Asia may show a new quinolone resistance pattern: NAL susceptible and CIPRO reduced susceptibility (MIC ≥ 0.12 µg/ml), therefore test both Nalidixan and Ciprofloxacin 1 µg Neo-Sensitabs.

If resistant to Ciprofloxacin, report as **R** to all fluoroquinolones.

- d) For Enterobacteriaceae isolated from the **CSF test cefotaxime** (or ceftriaxone) instead of cephalothin (or cefazolin).
- e) *Klebsiella* and *Enterobacter* spp. should always be reported as **R** to ampicillin.
- f) *Klebsiella*, *Enterobacter* and *Proteus* spp. should always be reported as **R** to nitrofurantoin.
- g) Isolates with MIC > 0.5 µg/ml against Imipenem and highly resistant to ceftazidime (KPC and GES enzymes) should be suspected of possessing carbapenemases. For detecting carbapenemases Ambler classes A and D. Synergy between BORON and carbapenems indicates KPC. For detecting metallo-beta-lactamases test, Dipicolinic acid with Meropenem and Ertapenem and Imipenem+EDTA against Imipenem Neo-Sensitabs or use the KPC+MBL Confirm ID kit. Further information is giving in User's Guide "Detection of resistance mechanisms using Neo-Sensitabs™ and DiaTabs™" (www.rosco.dk).
- h) Enterobacteriaceae suspicious of possessing plasmid mediated Amp C beta lactamases (Cefoxitin R, Ceftazidime R, Cefepime S and Carbapenems S) should be tested for synergy between Boronic acid and Ceftazidime+Clavulanate and Cefotaxime+Clavulanate and/or between Cefotaxime/Ceftazidime and Cloxacillin 500 µg (distance 5-10 mm from edge to edge) by using combined Neo-Sensitabs disc tests (AmpC Confirm ID kit). See further details in "Detection of resistance mechanisms using Neo-Sensitabs™ and DiaTabs™" (www.rosco.dk).
- i) Sader et al. recommends smaller inhibition zones for tetracyclines. (3)

- j) Results with one carbapenem cannot be extrapolated to the other. If a strain is resistant to Imipenem or Meropenem, report as resistant to ertapenem. Screen values for possible carbapenemases (metallo-beta-lactamases) are ≤ 23 mm for Ertapenem, Meropenem and Imipenem. (6)
- k) - If I/R to Tobramycin and susceptible to Gentamicin, report Amikacin as **I** (Intermediate).
- If **I** to Gentamicin and susceptible to other aminoglycosides, report as Gentamicin **R**.
- If **I** to Tobramycin, Gentamicin **R** and Amikacin **S**, report as Tobramycin **R**.
- If Gentamicin **R** and Tobramycin **R**, report as resistant to netilmicin.
- l) Rodriguez-Villalobos et al found that Temocillin Neo-Sensitabs showed better discrimination between susceptible and resistant and lower number of discrepancies than paper disks (BD). (5)
- m) New CLSI proposed breakpoints for Cephalosporins, Aztreonam and Enterobacteriaceae (Teleconference, August 2008).
- n) Tigecycline 15 µg disc has poor resolution (1-2 mm) between disc modal zones for adjacent MICs. (8)

References:

- 1) CLSI: Performance Standards for Antimicrobial Susceptibility Testing 23rd Inf. Suppl. **M100-S23**, 2013.
- 2) Hakanen et al J. Clin. Microbiol., **43**, 5775-8, 2005.
- 3) Sader H.S. et al: Reevaluation of CLSI Disk diffusion breakpoints for Tetracyclines for testing Enterobacteriaceae. J. Clin. Microbiol., **45**, 1640-1643, 2007.
- 4) Aznar E. et al: Detection of decreased susceptibility to fluoroquinolones in Salmonella spp. by 5 different methods including real-time PCR. Int. J. Antimicrob. Agents, **30**, 67-71, 2007.
- 5) Rodriguez-Villalobos H. et al: Comparison of 4 commercial methods for antimicrobial susceptibility testing of Temocillin. ECCMID, abstract R2411, Barcelona, april 2008.
- 6) Rhomberg PR et al: Regression analysis of MIC versus disk diffusion zone diameters for 3 carbapenems tested against Enterobacteriaceae isolates harboring serine carbapenemases with matched control strains. ASM Gen. Meet. June 2008, presentation C-031.
- 7) Parry C.M et al: Suitable disc antimicrobial susceptibility breakpoints defining *S. enterica* serovar *Typhi* isolates with reduced susceptibility to fluoroquinolones. Antimicrob. Ag. Chemother. **54**, 5201-5208, 2010.
- 8) Thean Yen Tan et al: Influence of different Mueller-Hinton agars and media ages on Etest susceptibility testing of tigecycline. Diagn. Microbiol. Infect. Dis. **68**, 93-95, 2010.
- 9) Wang et al: Susceptibility of ESBL-producing Enterobacteriaceae with the new CLSI breakpoints. Presentation D-1533, 50th ICAAC, Sept. 2010.
- 10) Watz N. et al: Impact of the new CLSI cephalosporin breakpoints on results reporting in an adult (Stanford) and pediatric (LPCH) hospital. Presentation D-1531, 50th ICAAC, Sept. 2010.

PSEUDOMONAS AERUGINOSA

Zone diameter interpretative criteria and MIC breakpoints according to CLSI (formerly NCCLS) (1) when testing *P. aeruginosa*, *Acinetobacter* spp., *S. maltophilia*, *B. cepacia*, are listed in the tables below.

Table 3.2-1 Interpretation for *Pseudomonas aeruginosa*

Mueller-Hinton Agar. Inoculum McFarland 0.5. Incubation at 35 °C ± 2 degrees ambient air for 16-18 hours.

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
a)	Amikacin 30 µg	AMI30	≥ 17	16-15	≤ 14	≤ 8	≥ 32
	Aztreonam 30 µg	AZT30	≥ 22	21-16	≤ 15	≤ 8	≥ 32
	Cefepime 30 µg	FEP30	≥ 18	17-15	≤ 14	≤ 8	≥ 32
d)	Cefepime+Clavulanate FEP+C	30+10 µg	detection of ESBL				
	Ceftazidime 30 µg	CAZ30	≥ 18	17-15	≤ 14	≤ 8	≥ 32
d)	Ceftazidime+Clavulanate CAZ+C	30+10µg	detection of ESBL				
	Ciprofloxacin 5 µg	CIPR5	≥ 21	20-16	≤ 15	≤ 1	≥ 4
c)	Colistin 10 µg	Co.10					
c)	2+18 hours' prediffusion method		≥ 15	14-11	≤ 10	≤ 2	≥ 8
f)	Doripenem 10 µg	DOR10	≥ 19	18-17	≤ 16	≤ 2	≥ 8
	Fosfomycin 200 µg	FO200	≥ 22	21-19	≤ 18	≤ 16	≥ 32
	Gatifloxacin 5 µg	GATIF	≥ 18	17-15	≤ 14	≤ 2	≥ 8
a)	Gentamicin 10 µg	GEN10	≥ 15	14-13	≤ 12	≤ 4	≥ 8
e,f)	Imipenem 10 µg	IMI10	≥ 19	18-16	≤ 15	≤ 2	≥ 8
b)	Imipenem+EDTA 10+750µg	IM10E	detection of metallo-β-lactamases				
	Levofloxacin 5 µg	LEVOF	≥ 17	16-14	≤ 13	≤ 2	≥ 8
e,f)	Meropenem 10 µg	MRP10	≥ 19	18-16	≤ 15	≤ 2	≥ 8
	Minocycline 30 µg	MIN30	≥ 19	18-15	≤ 14	≤ 4	≥ 16
	Netilmicin 30 µg	NET30	≥ 15	14-13	≤ 12	≤ 4	≥ 8
	Ofloxacin (U) 5 µg	OFL.5	≥ 16	15-13	≤ 12	≤ 2	≥ 8
	Piperacillin 100 µg	PIPRA	≥ 21	20-15	≤ 14	≤ 16	≥ 128
	Piperacillin+Tazobactam 100+10µg PI+TZ		≥ 21	20-15	≤ 14	≤ 16/4	≥ 128/4
	Tetracyclines 30 µg	TET30	≥ 19	18-15	≤ 14	≤ 4	≥ 16
	Ticarcillin 75 µg	TIC75	≥ 20	19-15	≤ 14	≤ 16	≥ 128
	Ticarcillin+Clavulanate 75+10 µg TIM85		≥ 24	23-17	≤ 16	≤ 16/2	≥ 128/2
a)	Tobramycin 10 µg	TOB10	≥ 15	14-13	≤ 12	≤ 2	≥ 8
	Dipicolinic acid	D.P.A	Detection of metallo-β-beta lacta- mases				

- *Pseudomonas aeruginosa* isolated from patients with cystic fibrosis requires incubation up to 24 hours.

Clinical isolates of *P. aeruginosa* heterogeneously resistant to carbapenems have been isolated (1). Subcolonies appearing within the zone of inhibition of Imipenem/Meropenem show higher MIC values. Automated systems and conventional agar dilution MICs using the standard 10⁴ CFU per spot inoculum may miss carbapenem-resistant mutants.

Detection of *P. aeruginosa* resistant to Colistin in University Hospital in Barcelona (3).

Unacceptable levels of error in beta-lactam susceptibility were detected using 4 automated methods (2). The authors suggest that clinical laboratories using automated systems should consider accurate alternative methods (agar diffusion methods) for routine use.

Cabot et al (4) conclude that isolates non-susceptible to Ceftazidime or Piperacillin-Tazobactam were found to hyperproduce AmpC, MexAB-OprM, while particularly MexX-OprM overexpression was found among Cefepime-non-susceptible isolates.

Overexpression of MexXY-OprM is known to confer low-level resistance to aminoglycosides.

All Imipenem resistant isolates showed inactivating mutations in oprD.

References:

- 1) Pournaras S. et al: Characterization of clinical isolates of *P. aeruginosa* heterogeneously resistant to carbapenems. J. Med. Microbiol, **56**, 66-70, 2007.
- 2) Juretschko S. et al: Accuracies of beta-lactam susceptibility test results for *P. aeruginosa* with 4 automated systems (BD Phoenix, Microscan WalkAway, Vitek and Vitek2). J. Clin. Microbiol, **45**, 1339-42, 2007.
- 3) Montero M. et al: Detection of *P. aeruginosa* resistant to Colistin in a University Hospital. Poster 65-P, XIII Congress, SEIMC, Madrid, May 2008.
- 4) Gabot et al: Overexpression of AmpC and efflux pumps in *P. aeruginosa* isolates from bloodstream infections: prevalence and impact on resistance in a Spanish multicenter study. Antimicrob. Agents Chemother. **55**, 1906-1911, 2011.

ACINETOBACTER spp.

Table 3.2-1

Interpretation for *Acinetobacter* spp.

Mueller-Hinton Agar. Inoculum McFarland 0.5. Incubation at 35 °C ± 2 degrees ambient air for 20-24 hours.

NEO-SENSITABS CODE				Zone diameter in mm			Break-points MIC µg/ml	
POTENCY				S	I	R	S	R
a)	Amikacin	30 µg	AMI30	≥ 17	16-15	≤ 14	≤ 16	≥ 32
	Ampicillin	10+10 µg	SAM20	≥ 15	14-12	≤ 11	≤ 8/4	≥ 32/16
	+Sulbactam							
	Aztreonam	30 µg	AZT30	≥ 22	21-16	≤ 15	≤ 8	≥ 32
	Cefepime	30 µg	FEP30	≥ 18	17-15	≤ 14	≤ 8	≥ 32
d)	Cefepime	30+10 µg	FEP+C	detection of ESBL				
	+Clavulanate							
d)	Ceftazidime	30 µg	CAZ30	≥ 18	17-15	≤ 14	≤ 8	≥ 32
	Ceftazidime	30+10µg	CAZ+C	detection of ESBL				
	+Clavulanate							
	Cefotaxime	30 µg	CTX30	≥ 23	22-15	≤ 14	≤ 8	≥ 64
	Ceftriaxone	30 µg	CTR30	≥ 21	17-14	≤ 13	≤ 8	≥ 64
	Chloramphenicol	30 µg		≥ 18	17-13	≤ 12	≤ 8	≥ 32
	CLR30							
c)	Ciprofloxacin	5 µg	CIPR5	≥ 21	20-16	≤ 15	≤ 1	≥ 8
c)	Colistin (7)	10 µg	Co.10					
	2+18 hours' prediffusion			≥ 15	14-11	≤ 10	≤ 2	≥ 8
	Doripenem	10 µg	DOR10	≥ 18	17-15	≤ 14	≤ 1	≥ 4
	Doxycycline	30 µg	DOX30	≥ 14	13-11	≤ 10	≤ 4	≥ 16
	Gatifloxacin	5 µg	GATIF	≥ 18	17-15	≤ 14	≤ 2	≥ 8
a)	Gentamicin	10 µg	GEN10	≥ 15	14-13	≤ 12	≤ 4	≥ 8
*)	Imipenem	10 µg	IMI10	≥ 16	15-14	≤ 13	≤ 4	≥ 16
b)	Imipenem	10+750 µg	IM10E	detection of metallo-β-lactamases				
	+EDTA							
	Levofloxacin	5 µg	LEVOF	≥ 17	16-14	≤ 13	≤ 2	≥ 8
*)	Meropenem	10 µg	MRP10	≥ 16	15-14	≤ 13	≤ 4	≥ 16
	Minocycline	30 µg	MIN30	≥ 16	15-13	≤ 12	≤ 4	≥ 16
	Piperacillin	100 µg	PIPR	≥ 21	20-18	≤ 17	≤ 16	≥ 128
	Piperacillin	100+10 µg	PI+TZ	≥ 21	20-18	≤ 17	≤ 16/4	≥ 128/4
	+Tazobactam							
	Rifampicin	5 µg	RIF.5	≥ 16	15-12	≤ 11	≤ 4	≥ 16
	Tetracyclines	30 µg	TET30	≥ 15	14-12	≤ 11	≤ 4	≥ 16
	Ticarcillin	75 µg	TIC75	≥ 20	19-15	≤ 14	≤ 16	≥ 128
	Ticarcillin	75+10 µg	TIM85	≥ 20	19-15	≤ 14	≤ 16/2	≥ 128/2
	+Clavulanate							
**)	Tigecycline	15 µg	TIG15	≥ 22	21-20	≤ 19	≤ 1	>2
a)	Tobramycin	10 µg	TOB10	≥ 15	14-13	≤ 12	≤ 4	≥ 8
	Trimethoprim (U) 5 µg		TRIM5	≥ 16	15-11	≤ 10		
	Trimethoprim	1.25+23.75 µg	SxT25	≥ 16	15-11	≤ 10	≤ 2/38	≥ 8/152
	+Sulfa							
	Dipicolinic acid		D.P.A	Detect. of	metallo - beta	lactama-ses		

Multiresistant *Acinetobacter baumannii*, resistant to all beta-lactams including imipenem/meropenem, should be suspected of producing carbapenemases; either Class D enzymes (oxacillinases OXA-23, OXA-24, OXA-58) or Class A enzymes (metallo-β-lactamases).

There are differences in activity between tigecycline and minocycline, consequently, they should be tested separately (11).

EUCAST-and CLSI potency NEO-SENSITABS™

Page 4 / 5 - Document: 3.2.0

Interpretation Zones and MIC Breakpoints according to CLSI

P. aeruginosa, *Acinetobacter* spp., *B. cepacia*, *S. maltophilia*

Clinical isolates of Colistin hetero-resistant *E. cloacae* and *A. baumannii* were not detected by Vitek-2 (9). Colistin-dependent *A. baumannii* were not detected by broth microdilution. The strain grows around a 10 µg colistin disk.

*Screen values for possible carbapenemases (OXA and MBL) are <21mm for Ertapenem and <22mm for Meropenem and Imipenem.

**Tentative breakpoints. Oxoid MHA gives average 3.5 mm smaller zones than BD-MHA and should not be used for testing tigecycline (12).

BURKHOLDERIA CEPACIA AND S. MALTOPHILIA

Table 3.2-2

Interpretation for *Burkholderia cepacia*

Mueller-Hinton Agar. Inoculum McFarland 0.5. Incubation at 33-35°C ambient air for 20-24 hours.

NEO-SENSITABS CODE	POTENCY		Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
c) Ceftazidime 30 µg	CAZ30		≥ 21	20-18	≤ 17	≤ 8	≥ 32
c) Colistin 10 µg	Co.10		≥ 15	14-11	≤ 10	≤ 2	≥ 4
c) 2+18 hours' prediffusion			≥ 20	19-16	≤ 15	≤ 4	≥ 16
Doripenem 10 µg	DOR10		≥ 20	19-16	≤ 15	≤ 4	≥ 16
Meropenem 10 µg	MRP10		≥ 20	19-16	≤ 15	≤ 4	≥ 16
Temocillin 30 µg	TEMOC		≥ 18	17-15	≤ 14	≤ 16	≥ 32
Tigecycline 15 µg	TIG15		≥ 20	19-17	≤ 16	≤ 2	≥ 8
Trimethoprim+ Sulfa 1.25+23.75µg	SxT25		≥ 16	15-11	≤ 10	≤ 2/38	≥ 8/152

Table 3.2-3

Interpretation for *Stenotrophomonas maltophilia*

Mueller-Hinton Agar. Inoculum McFarland 0.5. Incubation at 33-35°C ambient air for 20-24 hours.

NEO-SENSITABS CODE	POTENCY		Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
Fosfomycin (+ G6P) 200 µg	FO200		≥ 20	-	-	≤ 32	-
Levofloxacin 5 µg	LEVOF		≥ 17	16-14	≤ 13	≤ 2	≥ 8
Minocycline 30 µg	MIN30		≥ 19	18-15	≤ 14	≤ 4	≥ 16
Moxifloxacin 5 µg	MOXIF		≥ 18	17-15	≤ 14	≤ 1	≥ 4
Tigecycline 15 µg	TIG15		≥ 20	19-17	≤ 16	≤ 2	≥ 8
Trimethoprim+ Sulfa 1.25+23.75µg	SxT25		≥ 16	15-11	≤ 10	≤ 2/38	≥ 8/152

- a) *S. maltophilia* is intrinsic resistant towards aminoglycosides (report R).
For *S. maltophilia* and *B. cepacia* incubation at 30 °C for 24 hours might be more appropriate for detecting beta lactam resistance (2, 5, 8).
For *Pseudomonas* if **I/R** to Tobramycin and susceptible to Gentamicin, report Amikacin as resistant. Resistance to Trimethoprim + Sulfamethazole is increasing in Brazil (27%). (13)
- b) For detection of carbapenemases (including metallo-beta-lactamases), see "Detection of resistance mechanisms using Neo-Sensitabs™ and DiaTabs™".
- c) For accurate detection of colistin resistance in multidrug resistant strains, use the 2+18 hours' prediffusion method. This technique permits a good separation between susceptible and resistant strains: Place Colistin 10 µg Neo-Sensitabs on a non-inoculated Mueller-Hinton plate and incubate for 2 hours at room temperature ($\leq 25^{\circ}\text{C}$). Thereafter eliminate the tablet by knocking the plate against a table and leave the plate for 18 hours at room temperature of $\leq 25^{\circ}\text{C}$ (next day). Inoculate the plate with the strain to be tested. Incubate for 18-24 hours at 33-35 °C and read the inhibition zones (see Ch. 1.5).
Colistin 10 µg disk testing without prediffusion can be used as screening test for high level resistance with *P. aeruginosa* (MIC ≥ 128 µg/ml corresponds to no zone of inhibition).
- d) For detection of strains producing ESBLs. Further information in User's Guide "Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™".
- e) For *Acinetobacter*/*Pseudomonas*: results for one carbapenem (Imipenem, Meropenem) cannot be extrapolated to others.

References:

- 1) CLSI: Performance Standards for Antimicrobial Susceptibility Testing 23rd Inf. Suppl. **M100-S23**, 2013.
- 2) Wheat P.F. et al: Effect of temperature on antimicrobial susceptibilities of *Ps. maltophilia*. J. Clin. Pathol., **38**, 1055-58, 1985.
- 3) Bonfiglio G, Livermore D.M.: Effect of media composition on the susceptibility of *X. maltophilia* to beta-lactam antibiotics. J. Antimicrob. Chemother., **28**, 837-842, 1991.
- 4) Hsueh Po-Ren et al: Pandrug resistant *Acinetobacter baumannii* causing nosocomial infections in a University Hospital, Taiwan. Emerg. Infect. Dis., **8**, August 2002.
- 5) King A.: Susceptibility testing of *S. maltophilia*: Effect of temperature and medium on results. Clin. Microbiol. Infect., **40**, Suppl. 3, 335, 2004.
- 6) Pournaras S. et al: Heteroresistance to carbapenems in *Acinetobacter baumannii*. J. Antimicrob. Chemother., **55**, 1055-6, 2005.
- 7) Borda N. et al: Comparison of methods: Diffusion (DF) Prediffusion (PDF) and Epsilometric on isolates of *Acinetobacter baumannii*-calcoaceticus complex (Abc) against Colistin. 2007 (in press) Spanish.
- 8) Nicodemo A.C. et al: Antimicrobial therapy for *S. maltophilia* infections. Eur. J. Clin. Microbiol. Inf. Dis., **26**, 229-237, 2007.
- 9) Lo-Ten-Foe JR. et al: Comparative evaluation of the Vitek 2, disk diffusion, E-test, broth micro- and agar dilution susceptibility testing methods for Colistin in clinical isolates, including hetero resistant *E. cloacae* and *A. baumannii* strains. Antimicrobial Ag Chemother **51**, 3726-3730, 2007.
- 10) Hawley J.S. et al: Development of Colistin- dependent *Acinetobacter baumannii*-Ac calcoaceticus complex. Antimicrob. Ag. Chemother **51**, 4529-30, 2007.
- 11) Moland E.S et al: In vitro activity of Tigecycline against MDR *Acinetobacter baumannii* and selection of Tigecycline-Amikacin synergy. Antimicrob. Ag. Chemother **52**, 2940-42, 2008.
- 12) Thamlikitkul V. et al: Effect of different MH-agars on tigecycline disk diffusion susceptibility for *Acinetobacter* spp. J. Antimicrob. Chemother **61**, 847-8, 2008.
- 13) Paez J. et al: Trimethoprim + Sulfamethazole resistant and class 1 integron detection in clinical isolates of *S. maltophilia* in Brazil. Presentation C2-114, 50th ICAAC, Sept 2010.

EUCAST-and CLSI potency NEO-SENSITABS™ Interpretation Zones and MIC Breakpoints according to CLSI

STAPHYLOCOCCI

Zone diameter interpretative criteria and MIC breakpoints according to CLSI (formerly NCCLS) (1) when testing staphylococci are listed in the table below.

Table 3.3-1 Interpretation for *Staphylococcus* spp.

Mueller-Hinton agar. Inoculum: McFarland 0.5. Incubation for 16-18 hours ambient air at 35 °C ± 2 degrees, however, incubate for Oxacillin and Vancomycin for 24 hours at 33-35 °C.

				Zone diameter in mm			Break-points MIC µg/ml	
NEO-SENSITABS		POTENCY	CODE	S	I	R	S	R
d)	Amikacin	30 µg	AMI30	≥ 17	16-15	≤ 14	≤ 16	≥ 32
	Amoxycillin+Clav.	20+10 µg	AMC30		Test Cefoxitin			
	Ampicillin	10 µg	AMP10		Test Penicillin			
i)	Azithromycin	15 µg	AZI15	≥ 18	17-14	≤ 13	≤ 2	≥ 8
	Cefoxitin (10,11,13,16)	10 µg	CFO10	≥ 17	-	≤ 16	-	mecA pos.
	<i>S. aureus</i> / <i>S. lugdunensis</i>			≥ 21	-	≤ 20	-	mecA pos.
a) b)	Cefoxitin	30 µg	CFO30	≥ 22	-	≤ 21	≤ 4	mecA pos.
	<i>S. aureus</i> / <i>S. lugdunensis</i>			≥ 25	-	≤ 24	≤ 0.25	≥ 0.5
	Coagulase neg. staph.						(oxa)	
b)	Cefepime	30 µg	FEP30		Test Cefoxitin		≤ 8	≥ 32
b)	Cefotaxime	30 µg	CTX30		Test Cefoxitin		≤ 8	≥ 64
	Ceftaroline	30 µg	CPT30	≥ 24	23-21	≤ 20	≤ 1	≥ 4
b)	Ceftazidime	30 µg	CAZ30		Test Cefoxitin		≤ 8	≥ 32
a)	Ceftizoxime	30 µg	ZOX30					
	<i>S. aureus</i>			≥ 20	-	≤ 15	Oxa S	mecA pos. (Oxa R)
b)	Ceftriaxone	30 µg	CTR30		Test Cefoxitin		≤ 8	≥ 64
b)	Cefuroxime parenteral	30 µg	CXM30		Test Cefoxitin		≤ 8	≥ 32
b)	Cefuroxime oral	30 µg	CXM30		Test Cefoxitin		≤ 4	≥ 32
b)	Cephalothin	30 µg	CEP30		Test Cefoxitin		≤ 8	≥ 32
	Chloramphenicol	30 µg	CLR30	≥ 18	17-13	≤ 12	≤ 8	≥ 32
	Ciprofloxacin	5 µg	CIPR5	≥ 21	20-16	≤ 15	≤ 1	≥ 4
	Clarithromycin	15 µg	CLA15	≥ 18	17-14	≤ 13	≤ 2	≥ 8
i)	Clindamycin	2 µg	CLI.2	≥ 21	20-15	≤ 14	≤ 0.5	≥ 4
h)	Daptomycin	30 µg	DAPCa					
	2+18 hours' prediffusion			≥ 22	-	< 20	≤ 1	≥ 2
	Doripenem	10 µg	DOR10	≥ 30	-	-	≤ 0.5	-
	Doxycycline	30 µg	DOX30	≥ 16	15-13	≤ 12	≤ 4	≥ 16
i)	Erythromycin	15 µg	ERY15	≥ 23	22-14	≤ 13	≤ 0.5	≥ 8
	Fucidin	10 µg	FUC10	≥ 21	-	≤ 20	≤ 1	> 1
	Fucidin	100 µg	FUCID	≥ 28	27-24	≤ 23	≤ 1	≥ 4
c)	Gatifloxacin	5 µg	GATIF	≥ 23	22-20	≤ 19	≤ 0.5	≥ 2
d)	Gentamicin	10 µg	GEN10	≥ 15	14-13	≤ 12	≤ 4	≥ 8
b)	Imipenem	10 µg	IMI10		Test Cefoxitin		≤ 4	≥ 16
e) d)	Kanamycin	30 µg	KAN30	≥ 18	17-14	≤ 13	≤ 6	≥ 25
c)	Levofloxacin	5 µg	LEVOF	≥ 19	18-16	≤ 15	≤ 1	≥ 4
j)	Linezolid	30 µg	LINEZ	≥ 21	-	-	≤ 4	-
b)	Meropenem	10 µg	MRP10		Test Cefoxitin		≤ 4	≥ 16
	Minocycline	30 µg	MIN30	≥ 19	18-15	≤ 14	≤ 4	≥ 16
c)	Moxifloxacin	5 µg	MOXIF	≥ 24	23-21	≤ 20	≤ 0.5	≥ 2
	Mupirocin	10 µg	MUPIR	≥ 16	15-10	-	≤ 4	128-LLR
	Staphylococci			-	-	no zone		>256-HLR
d) f)	Netilmicin	30 µg	NET30	≥ 15	14-13	≤ 12	≤ 12	≥ 32
	Nitrofurantoin (U)	300 µg	NI300	≥ 17	16-15	≤ 14	≤ 32	≥ 128
	Novobiocin	5 µg	NOVO5	≥ 14	-	≤ 13	-	-

Interpretation Zones and MIC Breakpoints according to CLSI

Staphylococci

				Zone diameter in mm			Break-points MIC µg/ml	
NEO-SENSITABS		POTENCY	CODE	S	I	R	S	R
c)	Norfloxacin (U)	10 µg	NORFX	≥ 17	16-13	≤ 12	≤ 4	≥ 16
c)	Ofloxacin	5 µg	OFL.5	≥ 18	17-15	≤ 14	≤ 1	≥ 4
a)	Oxacillin	1 µg	OXA.1					
	<i>S. aureus/ S. lugdunensis</i>			≥ 13	12-11	≤ 10	≤ 2	≥ 4
	Coagulase neg. staph.			≥ 18	-	≤ 17	≤ 0.25	≥ 0.5
a)	Penicillin	10 U	PEN10	≥ 29	-	≤ 28	≤ 0.12	beta lactamase
	Quinupristin-Dalfoprist.	15 µg	SYN15	≥ 19	18-16	≤ 15	≤ 1	≥ 4
	Retapamulin	2 µg	RETA2	≥ 20	19-17	≤ 16	≤ 0.5	≥ 2
	Rifampicin	5 µg	RIF.5	≥ 20	19-17	≤ 16	≤ 1	≥ 4
	Sulphonamides (U)	240 µg	SULFA	≥ 17	16-13	≤ 12	≤ 100	≥ 350
	Teicoplanin	30 µg	TPN30					
g)	2+18 hours' prediffusion (BHI blood)			≥ 18	-	<18	≤ 4	GISA/TRCNS
	2+18 hours' prediffusion, MH plain			≥ 22	-	< 20	≤ 4	GISA
	Telithromycin	15 µg	TEL15	≥ 22	21-19	≤ 18	≤ 1	≥ 4
	Tetracyclines	30 µg	TET30	≥ 19	18-15	≤ 14	≤ 4	≥ 16
	Tigecycline	15 µg	TIG15	≥ 19	-	-	≤ 0.5	-
d)	Tobramycin	10 µg	TOB10	≥ 15	14-13	≤ 12	≤ 4	≥ 8
	Trimethoprim (U)	5 µg	TRIM5	≥ 16	15-11	≤ 10	≤ 4	≥ 16
	Trimethoprim+Sulfa	1.25+23.75 µg	SxT25	≥ 16	15-11	≤ 10	≤ 2/38	≥ 8/152
	Vancomycin	30 µg	VAN30	-	-	No zone	-	>32 VRSA
g)	2+18 hours' prediffusion (BHI blood)			≥ 22	-	< 20	≤ 2	VISA/GISA
	2+18 hours' prediffusion, MH plain			≥ 21	-	≤ 20	≤ 2	VISA/GISA

- a) Penicillin zone edge test: Robert Skov presented data showing, that the Penicillin zone edge test was more sensitive than the Nitrocefin test, using Penicillin 10 unit disc. Sharp zone edges indicate beta-lactamase production, while fussy zone edges indicate no beta-lactamase production. Staphylococci and beta-lactams: test only Penicillin, Cefoxitin (Ceftizoxime) and Oxacillin. Isolates that test as mecA negative (Cefoxitin S), but which show Oxacillin total resistance (no zone), should be reported as resistant. Isolates that test mecA positive and are Cefoxitin S and Oxacillin S, should be reported as susceptible to Oxacillin (PBP not functioning properly). If both Cefoxitin and Oxacillin are tested against *S. aureus* or *S. lugdunensis* and either result is resistant, the organism is reported as Oxacillin resistant. mecA-independent beta-lactam resistance (17) (Cefoxitin S, Oxacillin R, Ceftobiprole R) is likely to become more important as Ceftobiprole/Ceftaroline are developed and used clinically. Therefore test both Cefoxitin and Oxacillin. If Oxacillin resistant, report the isolate as resistant regardless of the result given by Cefoxitin.
- b) Cefoxitin and *S. aureus* incubation 18-24 h at 35 °C. For coagulase neg. staph. incubate for 24 h. Results may be reported after 18 h, if resistant (10, 13). Cefoxitin (Oxacillin) resistant staphylococci should be considered resistant to all available beta-lactam antibiotics and combinations with beta-lactamase inhibitors.
- c) There is cross resistance between quinolones against staphylococci.
- d) Staphylococci resistant to Gentamicin should be reported as resistant to Netilmicin, Tobramycin and Amikacin (APH (2'') + AAC (6'')). Tobramycin sensitive strains will currently be susceptible to Gentamicin and Netilmicin. Strains resistant to Tobramycin, are reported as resistant to Kanamycin and Amikacin.
- e) Interpretation valid for Amikacin and Isepamicin. Isolates resistant to Kanamycin should be reported as resistant to Amikacin (no matter zone size around Amikacin)
- f) Netilmicin may give false sensitive results with coagulase negative staphylococci. Use Gentamicin. Isolates sensitive to Gentamicin will currently be susceptible to Netilmicin.

- g) The CLSI Staphylococci Working Group has agreed to remove the disk diffusion zone breakpoint for Staphylococci, because the current diffusion method does not distinguish between Vanco S and hVISA or VISA strain. VRSA can be detected with the current diffusion method (without prediffusion) because they produce no zone of inhibition around Vancomycin 30 µg Neo-Sensitabs™.

For detection of hVISA/GISA/VISA strains use the 2+18 hours' prediffusion method with Teicoplanin 30 µg and Vancomycin 30 µg (see "Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™"). TRCNS = Teicoplanin resistant coagulase negative staphylococci (ex. *S. haemolyticus* (8)). Fitzgibbon et al: (12) showed that 5 µg Teicoplanin on BHI agar and MCF 0.5 inoculum, gave the best results as screening method for h-GISA strains.

- h) For detection of Daptomycin resistance, use the 2+18 hours' prediffusion method (see "Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™").
- i) Results with Erythromycin are also valid for Clarithromycin and Azithromycin. Macrolide-resistant isolates may have inducible clindamycin resistance that can be detected by the **D-zone test**, using Erythromycin 15 µg and Clindamycin 2 µg placed approximately 15 mm apart (edge to edge).
- j) For reading the inhibition zone, use transmitted light and measure the inner edge of the zone of inhibition. Staphylococci showing **R** to Linezolid, Chloramphenicol, Clindamycin and Quinu/Dalfopristin may possess cfr methyltransferases (15).

Vancomycin susceptibility among S.aureus (14)

Vancomycin	Name	MIC µg/ml
Susceptible	VSSA	≤ 2 µg/ml
Heteroresistant*	hVISA/hGISA	1-2 µg/ml
Intermediate	VISA/GISA	4-8 µg/ml
Resistant	VRSA	≥ 16 µg/ml

*Consist of sub-populations ($\leq 10^{-6}$) that may grow in media containing >2 µg/ml Vancomycin. For detecting hVISA/hGISA and VISA/GISA use the 2+18 hour pre-diffusion method with Vancomycin 30 µg/ml and Teicoplanin 30 µg/ml Neo-Sensitabs.

Infections involving hVISA are a problem, since these strains are reported by laboratories to be susceptible to vancomycin on the basis of current recommended MIC breakpoints. Therefore, the true prevalence of hVISA is unknown.

Stegger et al (18) found that isolates containing mecA(LGA251), now called mecC, are phenotypically resistant to beta-lactams, but have not been recognized as classical MRSA with conventional PCRs for mecA, owing to the different nucleotide composition. The mecC producing isolates were resistant to Cefoxitin, indicating that the diffusion method with Cefoxitin is reliable for any type of *S. aureus*.

Kotsakis et al (19) describe the emergence of a mecA negative *S. lugdunensis* with a modified PBP 1A/1B that expresses resistance to all beta-lactams. The strain was resistant to Oxacillin and Cefoxitin, but was not detected by the mecA test.

Lindquist et al (20) demonstrates the presence of pseudo-SCC elements resembling SCCmec type II among Multiresistant-MSSA. All 54 isolates were PCR positive for the nuc gene and negative by the mecA gene. Methods testing the nuc gene (GeneOhm MRSA) gave false resistant results.

References:

- 1) CLSI: Performance Standards for Antimicrobial Susceptibility Testing 23rd Inf. Suppl. **M100-S23**, 2013.

- 2) Felten A. et al: Evaluation of 3 techniques for detection of low-level methicillin-resistant *S. aureus* (MRSA): a disk diffusion method with Cefoxitin and Moxalactam, the Vitek 2 system and the MRSA-screen latex agglutination test. *J. Clin. Microbiol.*, **40**, 2766-2771, 2002.
- 3) Tenover F.C. et al: Vancomycin-resistant *S. aureus* isolate from a patient in Pennsylvania. *Antimicrob. Ag. Chemother.*, **48**, 275-280, 2004.
- 4) Howden B.P. et al: Identification of a new agar dilution screening method for the accurate detection of hVISA, 44th Annual ICAAC, presentation D-59, 2004.
- 5) Sandberg A et al: First report of heterogeneous vancomycin intermediate *S. aureus* (hVISA/VISA) in Denmark. 44th Annual ICAAC, presentation AR-07, 2004.
- 6) Nielsen S.V., Casals, J.B.: Detection of decreased susceptibility to glycopeptides in *S. aureus* using tablet (disc) prediffusion. 15th Eur. Cong. Clin. Microbiol. Infect. Dis. (ECCMID), April 2005.
- 7) Steward C.D. et al.: Testing for induction of Clindamycin resistance in Erythromycin-resistant isolates of *S. aureus*. *J. Clin Microbiol*, **43**, 171-21, 2005.
- 8) Falcone M. et al: Teicoplanin use and emergence of *S. haemolyticus*: is there a link? *C.M.I.* **12**, 96-97, 2005.
- 9) Wagenvoort J.H.T. et al: Superior discriminative performance of ceftizoxime disk diffusion test for detecting MRSA. *Eur. J. Clin. Microbiol. Infect. Dis.*, **25**, 405-6, 2006.
- 10) Report for CLSI, September 2006.
- 11) Skov R. et al: Phenotypic detection of methicillin resistance in *S. aureus* by disk diffusion, testing and E-test on MHA. *J. Clin. Microbiol.*, **44**, 4395-99, 2006.
- 12) Fitzgibbon M.H. et al: Investigation of reduced susceptibility to glycopeptides among MRSA isolates from patients in Ireland and evaluation of agar screening methods for detection of hGISA. *J.Clin. Microbiol.* **45**, 3263-69, 2007.
- 13) Skov R et al: Evaluation of 10, 30 and 60 µg/ml Cefoxitin Neo-Sensitabs for prediction of Methicillin resistance in *S.aureus* ESCMID. Abstract nr. 1134-03-150, Copenhagen 2005.
- 14) Sakoulas G. Et al: Increasing antibiotic resistance among MRSA strains. *CID*, **46**, S-360-5367, 2008.
- 15) Fritsche T.R. et al: Detection of cfr rRNA methyltransferases among *S.aureus* and coag. neg. *Staphylococci* recovered from human infections. ASM meeting, June 2008, presentation A-047.
- 16) Roisin S. et al: Evaluation of new Vitek-2 card and disk diffusion method for determining susceptibility of *S.aureus* to Oxacillin. *JCM*, **46** 2525-2528, 2008.
- 17) Banerjee R. et al: A mecA-negative strain of methicillin-resistant *Staphylococcus aureus* with high-level beta-lactam resistance contains mutations in three genes. *Antimicrob. Agents Chemother.* **54**, 4900-4902, 2010.
- 18) Stegger M et al: Rapid detection, differentiation and typing of MRSA harbouring either mecA or the new mecA homologue mecA(LGA251). *Clin. Microbiol. Infect.* **18**, 395-400, 2011.
- 19) Kotsakis SD et al: *S. lugdunensis* strain with a modified PBP 1A/1B expressing resistance to beta-lactams. *Eur. J. Clin. Microbiol. Infect. Dis.* **31**, 169-172, 2012.
- 20) Lindquist M et al: Detection and characterization of SCCmec remnants in multiresistant MSSA causing a clonal outbreak in a Swedish county. *Eur. J. Clin. Microbiol. Infect. Dis.* **31**, 141-147, 2012.

SCREENING TESTS FOR STAPHYLOCOCCI

Inoculum: Mc Farland 0.5 and incubation at 35°C for 18-24 hours

Test	Organism group	Neo-S potency	Test method	Medium	Interpretation
Beta-lactamase	Staphylococci	Penicillin 10 units	Disk diffusion	MH	R: ≤28 mm
Oxacillin res. mecA	S.aureus S.lugdu-nensis	Cefoxitin 30 µg and Oxacillin 1 µg	Disk diffusion	MH	S: ≥22 mm R: ≤21 mm (cefoxitin) S: ≥13 mm R: ≤10 mm (oxacillin)
	Coagulase neg. staph.				S: ≥25 mm R: ≤24 mm (cefoxitin) S: ≥18 mm R: ≤17 mm (oxacillin)
Vancomycin/Teicoplanin reduced susceptibility	S.aureus Coagulase neg. staph.	Vancomycin 30 µg Teicoplanin 30 µg	Pre-diffusion 2+18 hours	BHI+5% blood or MH plain	Vanco: R: <22 mm (hVISA/VISA/GISA) Teico: R: <20 mm R: <18 mm (hVISA/VISA/GISA) (BHI blood) TRCNS
Daptomycin reduced susceptibility	S.aureus	Daptomycin 30 µg	Pre-diffusion 2+18 hours	MH	S: ≥22 mm R: <20 mm
Inducible Clindamycin resistance (D-zone)	Staphylococci	Erythromycin 15 µg Clindamycin 2 µg	Disk diffusion	MH	Flattering of the disk adjacent to the erythromycin disk = inducible clindamycin resistance. Report clindamycin R
Aminoglycosides	Staphylococci	Kanamycin 30 µg Gentamicin 10 µg Tobramycin 10 µg	Disk diffusion	MH	Kana 30: S: ≥18 R: ≤13 Genta 10: S: ≥15 R: ≤12 Tobra 10: S: ≤15 R: ≤12 Results with kanamycin are valid for amikacin and isepamicin. Strains resistant to genta are reported as resistant to tobra, netilmicin and amikacin (APH (2 ⁿ) + AAC (6 ¹)). Isolates susceptible to gentamicin will currently be susceptible to netilmicin. Strains resistant to tobramycin are reported as resistant to kanamycin and amikacin.

ENTEROCOCCI

Zone diameter interpretative criteria and MIC breakpoints according to CLSI (formerly NCCLS) (1) when testing enterococci are listed in the table below.

Ligu et al (9) isolated 23 vanA genotype VREF and found that 12 of them (52%) had vanB phenotype (Teicoplanin MICs 2-12 µg/ml).

Table 3.4-1 Interpretation for *Enterococcus* spp.

Mueller-Hinton agar. Inoculum: McFarland 0.5. Incubation for 18-24 hours ambient air at 35 °C ± 2 degrees.

				Zone diameter in mm			Break-points MIC µg/ml	
NEO-SENSITABS	POTENCY	CODE		S	I	R	S	R
e) Ampicillin	10 µg	AMP10		≥ 17	-	≤ 16	≤ 8	≥ 16
Chloramphenicol	30 µg	CLR30		≥ 18	17-13	≤ 12	≤ 8	≥ 32
Ciprofloxacin (U)	5 µg	CIPR5		≥ 21	20-16	≤ 15	≤ 1	≥ 4
g) Daptomycin	30 µg	DAPCa		≥ 12	-	-	≤ 4	-
2+18 hours' prediffusion				≥ 12	-	-	≤ 4	-
Doxycycline	30 µg	DOX30		≥ 16	15-13	≤ 12	≤ 4	≥ 16
f) Erythromycin	15 µg	ERY15		≥ 23	22-14	≤ 13	≤ 0.5	≥ 8
Fosfomycin (U)	200 µg	FO200		≥ 16	15-13	≤ 12	≤ 64	≥ 256
Gatifloxacin (U)	5 µg	GATIF		≥ 18	17-15	≤ 14	≤ 2	≥ 8
h) Imipenem	10 µg	IMI10		≥ 16	15-14	≤ 13	≤ 4	≥ 16
Levofloxacin (U)	5 µg	LEVOF		≥ 17	16-14	≤ 13	≤ 2	≥ 8
Linezolid	30 µg	LINEZ		≥ 23	22-21	≤ 20	≤ 2	≥ 8
Minocycline	30 µg	MIN30		≥ 19	18-15	≤ 14	≤ 4	≥ 16
Nitrofurantoin (U)	300 µg	NI300		≥ 17	16-15	≤ 14	≤ 32	≥ 128
Norfloxacin (U)	10 µg	NORFX		≥ 17	16-13	≤ 12	≤ 4	≥ 16
Penicillin	10 U	PEN10		≥ 15	-	≤ 14	≤ 8	≥ 16
a) Quinupristin/Dalfopristin	15 µg	SYN15		≥ 19	18-16	≤ 15	≤ 1	≥ 4
Rifampicin	5 µg	RIF.5		≥ 20	19-17	≤ 16	≤ 1	≥ 4
c) Teicoplanin	30 µg	TPN30		≥ 14	13-11	≤ 10	≤ 8	≥ 32
2+18 hour pre-diffusion				≥ 16	-	<16	≤ 8	TRE
Tetracyclines	30 µg	TET30		≥ 19	18-15	≤ 14	≤ 4	≥ 16
Tigecycline	15 µg	TIG15		≥ 19	-	-	≤ 0.25	-
c) Vancomycin	30 µg	VAN30		>16	15-14	< 14	≤ 4	≥ 32
(2+18 hours' prediffusion)				≥16	15-12	< 12	≤ 4	VRE
b) d) Gentamicin	250 µg	GN250		-	-	< 14	-	HLR
b) Streptomycin	500 µg	ST500		-	-	< 14	-	HLR

- E. faecalis* is intrinsically resistant to Quinupristin/Dalfopristin.
- For detection of High Level Aminoglycoside Resistance (HLR) use high level Neo-Sensitabs (Gentamicin 250 µg and Streptomycin 500 µg).
- With Vancomycin and Teicoplanin plates should be incubated for full 24 hours and examined carefully. Vancomycin resistant strains show a hazy zone edge while sensitive strains show a sharp zone edge.

- d) Gentamicin HLR in enterococci indicates resistance to all aminoglycosides, except Streptomycin. *E. faecium* shows intrinsic resistance towards Kanamycin, Tobramycin and Netilmicin due to the production of the enzyme AAC (6').
- e) Results with Ampicillin and *E. faecalis* are valid for Amoxicillin and Piperacillin. Reduced susceptibility to Ampicillin is common in *E. faecium*. If resistant to Ampicillin, report as resistant to ureidopenicillins. Martin et al. (10) report Ampicillin S/Imipenem R and Ampicillin R/Imipenem S in 5% of clinical isolates.
- f) Results with Erythromycin are valid also for Azithromycin and Clarithromycin.
- g) For detection of Daptomycin resistance, use the 2+18 hours' prediffusion method (valid only for *E. faecalis* susceptible to Vancomycin).
- h) Ampicillin susceptible but Penicillin (MIC ≥ 16 µg/ml) and Imipenem (MIC ≥ 8 µg/ml) resistant strains of *E. faecalis* were widely disseminated in a Greek hospital (31.4% with this phenotype)(7,8). Ampicillin and Imipenem should be tested separately. Guardabassi et al. (11) found the same pattern in Denmark.

For more information on glycopeptide resistance see: "Detection of resistance mechanisms using Neo-Sensitabs™ and DiaTabs™".

Gagnon et al (12) found during routine screening for VRE 6 Enterococcus faecium isolates positive for vanA by PCR, but susceptible to Vancomycin and Teicoplanin by phenotypic testing. The 6 isolates were missing the vanR and vanS genes that are responsible for the activation of resistance genes. False resistance results with PCR testing.

Sarti et al (13) describe the isolation of 8 beta-lactamase-producing multiresistant *E. faecium* isolates from 2010. The beta-lactamase activity was demonstrated at a high inoculum and in the presence of Methicillin after overnight incubation.

References:

- 1) CLSI: Performance Standards for Antimicrobial Susceptibility Testing 23rd Inf. Suppl. **M100-S23**, 2013.
- 2) Jenkins R.D. et al: False susceptibility of enterococci to aminoglycosides with blood-enriched Mueller-Hinton Agar for disk susceptibility testing. J. Clin. Microbiol., **22**, 369-374, 1985.
- 3) Torres C. et al: Detection of Aminoglycoside-penicillin synergy against *Enterococcus faecium* using high content aminoglycoside disks. Eur. J. Clin. Microbiol. Infect. Dis., **14**, 878-882, 1995.
- 4) Centinkaya Y. et al: Vancomycin resistant enterococci. Clin. Microbiol. Reviews, **13**, 686-707, 2000.
- 5) Corso A. et al: First report of Van A *E. gallinarum* dissemination within an intensive care unit in Argentina. Intern. J. Antimicrob. Ag., **25**, 51-56, 2005.
- 6) Willems R.J.L. et al: Global spread of vancomycin-resistant *E. faecium* from distinct nosocomial genetic complex. Emerg. Infect. Dis., **11**, 821-828, 2005.
- 7) Lemming L. et al: Characterization of atypical *Enterococcus faecalis* resistant to penicillin but susceptible to ampicillin from bacteremic episodes. 45th ICAAC presentation D1651, 2005.
- 8) Metzidie E. et al: Spread of an unusual penicillin and imipenem resistant but ampicillin susceptible phenotype among *E. faecalis* clinical isolates. JAC, **57**, 158-59, 2006.
- 9) Ligu et al: A new Tn 1546 type of vanB phenotype- vanA genotype vancomycin-resistant *E. faecium* isolates in mainland China. Diagn. Microbiol. & Infect. Dis. **63**, 70-75, 2009.
- 10) Martin et al: Carbapenem and Ampicillin resistance in Enterococci. Presentation C2-109, 50th ICAAC, Sept. 2010.
- 11) Guardabassi L. et al: Gentamicin-resistant Enterococcus faecalis sequence type 6 with reduced Penicillin susceptibility: Diagnostic and Therapeutic implications. J. Clin. Microbiol. **48**, 3820-3821, 2010.
- 12) Gagnon S et al: vanA containing *E. faecium* susceptible to Vancomycin and Teicoplanin because of major nucleotide deletions in Tn1546. J. Antimicrob. Chemother. **66**, 2758-2762, 2011.

- 13) Sarti M et al: Polyclonal diffusion of beta-lactamase-producing *Enterococcus faecium*. J. Clin. Microbiol. **50**, 169-172, 2012.

ORIGINAL ROSCO DOCUMENT

Interpretation Zones and MIC Breakpoints according to CLSI

Pneumococci

PNEUMOCOCCI

Zone diameter interpretative criteria and MIC breakpoints according to CLSI (formerly NCCLS) (1) when testing *Streptococcus pneumoniae* are listed in the table below.

The current tablet diffusion procedure described for Neo-Sensitabs should be followed, but no more than 9 Neo-Sensitabs should be placed on a large (150 mm) plate or 4 Neo-Sensitabs on a 90-100 mm plate. The recommended medium is Mueller-Hinton agar supplemented with 5% defibrinated sheep blood.

Soriano et al (9) found *S. pneumoniae* isolates with penicillin/cefotaxime, MICs of 16 µg/ml, in strains from Romania.

Table 3.5-1 Interpretation for *Streptococcus pneumoniae*

Mueller-Hinton + 5 % blood. Inoculum: McFarland 0.5. Incubation at 35 °C ± 2 degrees in 5-7 % CO₂ for 20-24 hours.

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
a)	Azithromycin 15 µg	AZI15	≥ 18	17-14	≤ 13	≤ 0.5	≥ 2
	Ceftaroline 30 µg	CPT30	≥ 26	-	-	≤ 0.5	-
	Ceftizoxime 30 µg	ZOX30	≥ 26	-	≤ 25	≤ 0.5	-
	(Cefotaxime, Ceftriaxone)					(3rd gen. cepha.)	
	Chloramphenicol 30 µg	CLR30	≥ 21	-	≤ 20	≤ 4	≥ 8
	Clarithromycin 15 µg	CLA15	≥ 21	20-17	≤ 16	≤ 0.25	≥ 1
	Clindamycin 2 µg	CLI.2	≥ 19	18-16	≤ 15	≤ 0.25	≥ 1
	Doripenem 10 µg	DOR10	-	-	-	≤ 1	-
	Doxycycline 30 µg	DOX30	≥ 28	27-25	≤ 24	≤ 0.25	≥ 1
	Erythromycin 15 µg	ERY15	≥ 21	20-16	≤ 15	≤ 0.25	≥ 1
e)	Gatifloxacin 5 µg	GATIF	≥ 21	20-18	≤ 17	≤ 1	≥ 4
e)	Imipenem 10 µg	IMI10	≥ 28	27-25	≤ 24	≤ 0.12	≥ 1
e)	Levofloxacin 5 µg	LEVOF	≥ 17	16-14	≤ 13	≤ 2	≥ 8
	Linezolid 30 µg	LINEZ	≥ 21	-	-	≤ 2	-
e)	Meropenem 10 µg	MER10	≥ 26	25-23	≤ 22	≤ 0.25	≥ 1
e)	Moxifloxacin 5 µg	MOXIF	≥ 18	17-15	≤ 14	≤ 1	≥ 4
e)	Norfloxacin 10 µg	LEVOF			< 12	Reduced Susceptibility to quinolones	>8
b)	Ofloxacin 5 µg	OFL.5	≥ 16	15-13	≤ 12	≤ 2	≥ 8
	Oxacillin 1 µg	OXA.1	≥ 20	≤ 19	≤ 19	≤ 0.06	≥ 0.12
	(penicillin)					(pen)	(I/R) (pen)
	Quinupristin/Dalfopristin 15 µg	SYN15	≥ 19	18-16	≤ 15	≤ 1	≥ 4
	Rifampicin 5 µg	RIF.05	≥ 19	18-17	≤ 16	≤ 1	≥ 4
	Teicoplanin 30 µg	TPN30	≥ 16	-	-	≤ 1	-
	Telithromycin 15 µg	TEL15	≥ 19	18-16	≤ 15	≤ 1	≥ 4
	Tetracyclines 30 µg	TET30	≥ 24	23-21	≤ 20	≤ 1	≥ 4
	Tigecycline 15 µg	TIG15	≥ 19	-	-	≤ 0.25	-
	Trimethoprim +Sulfa 1.25+23.75 µg	SxT25	≥ 19	18-16	≤ 15	≤ 0.5/9.5	≥ 4/76
c)	Vancomycin 30 µg	VAN30	≥ 17	-	-	≤ 1	-
	non-meningeal criteria						
	Amoxycillin+Clav. 20+10 µg	AM+CL	≥ 23	22-20	≤ 19	≤ 2/1	≥ 8/4
	Ampicillin 10 µg	AMP10	≥ 20	19-17	≤ 16	≤ 2	≥ 8
	(valid for Amoxycillin and Penicillin G)						
	Oxacillin 1 µg		≥ 20	19-10	no zone	≤ 0.06	≥ 2
	(valid for Penicillin-Voral)						
	Cefepime 30 µg	FEP30	≥ 24	23-22	≤ 21	≤ 1	≥ 4

Interpretation Zones and MIC Breakpoints according to CLSI

Pneumococci

				Zone diameter in mm			Break-points MIC µg/ml	
NEO-SENSITABS	POTENCY	CODE	S	I	R	S	R	
d)	Cefotaxime	30 µg	CTX30	≥ 24	23-22	≤ 21	≤ 1	≥ 4
	Cefpodoxime	10 µg	CPD10	≥ 23	22-19	≤ 18	≤ 0.5	≥ 2
	Ceftriaxone	30 µg	CTR30	≥ 24	23-22	≤ 21	≤ 1	≥ 4
	Cefuroxime (oral)	30 µg	CXM30	≥ 24	23-22	≤ 21	≤ 1	≥ 4

- a) Cefotaxime and Ceftriaxone must not be tested against pneumococci (CSF) by the diffusion method. A surrogate test: Ceftizoxime is used instead. Ceftizoxime detects reduced sensitivity to 3rd and 4th generation cephalosporins.
- b) Oxacillin 1 µg is used for screening susceptibility towards Penicillin is reported as S or I/R to Penicillin (not Oxacillin).
Non-meningeal pneumococcal isolates with Penicillin MIC ≤ 2 µg/ml (Ampicillin 10 µg zone ≥ 20 mm) can be considered susceptible to Amoxycillin, Amoxycillin+ Clavulanate, Cefepime, Cefotaxime, Ceftriaxone and Ertapenem (1).
Non-meningeal pneumococci with Penicillin MIC ≤ 0.06 µg/ml (Oxacillin zone ≥ 20 mm) can also be considered susceptible to Ampicillin (parenteral), Ampicillin+ Sulbactam, Imipenem and Meropenem.
 Penicillin and Cefotaxime/Ceftriaxone or Meropenem MIC's should be determined for isolates with Oxacillin zones ≤ 19mm (1).
- c) For use in acute otitis media, acute sinusitis and community acquired pneumonia (non-meningitis).
- d) Interpretation criteria for Cefotaxime requires doses appropriate for serious pneumococcal infections; e.g. at least 1 g (adults) or 50 mg/kg (children) every 8 hours or more frequently.
- Penicillin-resistant strains from the CSF should be considered resistant to Ampicillin, Amoxycillin, Amoxycillin+Clavulanate, and first and second generation cephalosporins.
 - Emergence in USA of pneumococci with very high level resistance to Penicillin (MIC ≥ 8 µg/ml), associated with high Amoxicillin MIC's (≥ 8 µg/ml) and high Cefotaxime MIC's.
 - Chiu et al (7) observed an increasing Ceftriaxone resistance in *S. pneumoniae* from Taiwan.
- e) Currently used CLSI MIC breakpoints for fluoroquinolones define many pneumococci isolates as susceptible even though they harbour QRDR mutations (8). Microbiological resistance breakpoints are as follows: Levofloxacin (MIC > 1 µg/ml ~ zone ≤ 20 mm), Moxifloxacin (MIC > 0.12 µg/ml ~ zone ≤ 25 mm), Gatifloxacin (MIC > 0.25 µg/ml ~ zone ≤ 24 mm) and Norfloxacin (MIC > 8 µg/ml ~ zone < 12mm). If resistant to Moxifloxacin, reports as resistant to all fluoroquinolones. Isolates susceptible to Levofloxacin are susceptible to Moxifloxacin, but not vice versa.
- CLSI recommends Penicillin, Cefotaxime (or Ceftriaxone) and Meropenem to be tested by a reliable MIC method and reported routinely with **CSF isolates** of *S. pneumoniae*.
- f) Macrolide-resistant isolates may have inducible Clindamycin resistance, which can be detected by the D-zone test (See 3.3.0 page 5).
- g) According to Rosco's Regression line analysis the following interpretation corresponds to the recommended MIC's: S ≥ 24 mm, I: 23-21 mm, R ≤ 20 mm

References:

- 1) CLSI: Performance Standards for Antimicrobial Susceptibility Testing 23rd Inf. Suppl. **M100-S23**, 2013.

- 2) Williams-Bouyer N. et al: Predicting susceptibility of *Streptococcus pneumoniae* to Ceftriaxone and Cefotaxime by Cefuroxime and Ceftizoxime disk diffusion testing. J. Clin. Microbiol., **37**, 3707-3710, 1999.
- 3) Pérez-Trallero E. et al: Fluoroquinolone and macrolide treatment failure in pneumococcal pneumonia and selection of multidrug resistant isolates. Emerg. Infect. Dis., **9**, 1159-1161, 2003.
- 4) Schrag S.J. et al: Emergence of *S. pneumoniae* with very high level resistance to Penicillin. Antimicrob. Ag. Chemother., **48**, 3016-23, 2004.
- 5) Lim S. et al: Antimicrobial susceptibility breakpoints and first-step parC mutations in *S. pneumoniae*: redefining fluoroquinolone resistance. Emerg. Infect. Dis., **9**, 833-7, 2003.
- 6) Faccone D. et al: Emergence of *S. pneumoniae* clinical isolate highly resistant to Telithromycin and fluoroquinolones. J. Clin. Microbiol., **43**, 5800-5803, 2005.
- 7) Chiu C.H. et al: Increasing Ceftriaxone resistance and multiple alterations of PBPs among Penicillin-resistant *S. pneumoniae* isolates in Taiwan. Antimicrob. Ag. Chemother., **51**, 3404-06, 2007.
- 8) de Cueto M. et al: Fatal Levofloxacin failure in treatment of a bacteremic patient infected with *S. pneumoniae*, with a pre-existing parC mutation. J. Clin. Microbiol. **46**, 1558-1560, 2008.
- 9) Soriano F. et al: Breakthrough in Penicillin resistance? *S. pneumoniae* isolates with Penicillin/Cefotaxime MICs of 16 µg/ml and their genotypic and geographical relatedness. J. Antimicrob. Chemother **62**. 1234-40, 2008.
- 10) Dallas S. D et al: Development of Doxycycline MIC and disk diffusion interpretative breakpoints for *S. pneumoniae*. Presentation 091, 111th General Meet. ASM 2011.

Zone diameter interpretative criteria and MIC breakpoints according to CLSI (formerly NCCLS) (1) when testing *Streptococci* other than *S. pneumoniae* are listed in the table below

The current tablet diffusion procedure described for Neo-Sensitabs should be followed, but no more than 9 Neo-Sensitabs should be placed on a large (150 mm) plate or 4 Neo-Sensitabs on a 90-100 mm plate. The recommended medium is Mueller-Hinton agar supplemented with 5% defibrinated sheep blood.

Small colony-forming beta-haemolytic strains with group A, C, F or G antigens (*S. anginosus*, previously termed *S. milleri*) are considered part of the viridans group. The viridans group also includes: *S. mitis*, *S. oralis*, *S. sanguinis*, *S. salivarius*, *S. intermedius*, *S. constellatus*, *S. mutans* and *S. bovis*.

Table 3.6-1

Interpretation for Streptococci (except *S. pneumoniae*)

Mueller-Hinton + 5 % blood. Inoculum: McFarland 0.5. Incubation at 35 °C ± 2 degrees in 5-7 % CO₂ for 20-24 hours.

				Zone diameter in mm			Break-points MIC µg/ml	
NEO-SENSITABS	POTENCY	CODE		S	I	R	S	R
a) Ampicillin	10 µg	AMP10						
Beta-haemolytic				≥ 24	-	-	≤ 0.25	-
Viridans				≥ 24	23-17	≤ 16	≤ 0.25	≥ 8
c) e) Azithromycin	15 µg	AZI15		≥ 18	17-14	≤ 13	≤ 0.5	≥ 2
Cefepime	30 µg	FEP30						
Beta-haemolytic				≥ 24	-	-	≤ 0.5	-
Viridans				≥ 24	23-22	≤ 21	≤ 1	≥ 4
Cefotaxime	30 µg	CTX30						
Beta-haemolytic				≥ 24	-	-	≤ 0.5	-
Viridans				≥ 28	27-24	≤ 23	≤ 1	≥ 4
Ceftaroline	30 µg	CPT30						
<i>S. pyogenes</i>				≥ 26	-	-	≤ 0.5	-
<i>S. agalactiae</i>				≥ 26	-	-	≤ 0.5	-
Ceftibuten (screen penicillin group B)	30 µg	CTB30		≥ 20	-	<20	-	≥ 0.25
Ceftizoxime (screen penicillin group B)	30 µg	ZOX30		-	-	<30	-	≥ 0.25
Ceftriaxone	30 µg	CTR30						
Beta-haemolytic				≥ 24	-	-	≤ 0.5	-
Viridans				≥ 27	26-25	≤ 24	≤ 1	≥ 4
Chloramphenicol	30 µg	CLR30		≥ 21	20-18	≤ 17	≤ 4	≥ 16
c) e) Clarithromycin	15 µg	CLA15		≥ 21	20-17	≤ 16	≤ 0.25	≥ 1
c) Clindamycin	2 µg	CLI.2		≥ 19	18-16	≤ 15	≤ 0.25	≥ 1
Daptomycin 2+18 h prediffusion	30 µg	DAPCa						
				≥ 22	-	-	≤ 1	-
Doripenem	10 µg	DOR10						
Beta-haemolytic				≥ 24	-	-	≤ 0.12	-
Doxycycline (viridans)	30 µg	DOX30		≥ 23	22-19	≤ 18	≤ 2	≥ 8
c) Erythromycin	15 µg	ERY15		≥ 21	20-16	≤ 15	≤ 0.25	≥ 1
Gatifloxacin	5 µg	GATIF		≥ 21	20-18	≤ 17	≤ 1	≥ 4

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
Imipenem	10 µg	IMI10	≥ 26	25-23	≤ 22	≤ 0.25	≥ 1
Levofloxacin	5 µg	LEVOF	≥ 17	16-14	≤ 13	≤ 2	≥ 8
Linezolid	30 µg	LINEZ	≥ 21	-	-	≤ 2	-
Meropenem	10 µg	MER10	≥ 26	25-23	≤ 22	≤ 0.5	-
Minocycline	30 µg	MIN30	≥ 23	22-19	≤ 18	≤ 2	≥ 8
Moxifloxacin	5 µg	MOXIF	≥ 21	20-18	≤ 17	≤ 1	≥ 4
Ofloxacin	5 µg	OFL.5	≥ 16	15-13	≤ 12	≤ 2	≥ 8
b) Oxacillin (penicillin screening)	1 µg	OXA.1	≥ 14	≤ 13	≤ 13	≤ 0.12 (pen) (I/R)	≥ 0.25
a) Penicillin	10 U	PEN10					
Beta-haemolytic			≥ 24	-	-	≤ 0.12	-
Viridans			≥ 26	25-13	≤ 12	≤ 0.12	≥ 4
Quinupristin/Dalfopristin	15 µg	SYN15	≥ 19	18-16	≤ 15	≤ 1	≥ 4
Retapamulin	2 µg	RETA2	≥ 15	-	-	≤ 0.25	-
Rifampicin	5 µg	RIF.5	≥ 20	19-17	≤ 16	≤ 1	≥ 4
Teicoplanin	30 µg	TPN30	≥ 16	-	-	≤ 1	-
Telithromycin	15 µg	TEL15	≥ 19	18-16	≤ 15	≤ 1	≥ 4
Tetracyclines	30 µg	TET30	≥ 23	22-19	≤ 18	≤ 2	≥ 8
Tigecycline	15 µg	TIG15	≥ 19	-	-	≤ 0.25	-
Trimethoprim	5 µg	TRIM5	≥ 17	16-15	≤ 14	≤ 2	> 2
f) Trimethoprim+Sulfa	1.25+23.75 µg	SxT25	≥ 19	18-16	≤ 15	≤ 0.5/9.5	≥ 4/76
Vancomycin	30 µg	VAN30	≥ 17	-	-	≤ 1	-

- Viridans streptococci isolated from blood or CSF should be tested for Penicillin or Ampicillin susceptibility using an MIC method.
- Oxacillin 1 µg Neo-Sensitabs is useful for screening for Penicillin susceptibility in streptococci.
- Macrolide resistant isolates may have inducible Clindamycin resistance, which can be detected by the D-zone test (See 3.3.0 page 5). Group B streptococci are susceptible to Ampicillin, Penicillin and Cefazolin, but may be resistant to Clindamycin and/or Erythromycin. In case of severe Penicillin allergy they should be tested and reported (D-zone test). Bonfiglio et al (10) found *S.agalactiae*, Erythromycin S and Clindamycin R.
- If resistant to Levofloxacin, report as resistant to all fluoroquinolones.
- Result for Erythromycin is also valid for Clarithromycin and Azithromycin.

The first report of group B streptococci with confirmed reduced susceptibility to Penicillin (MIC 0.12-0.25 µg/ml) in blood isolates from Hong Kong (7).

Kimura et al (8,11) found that Ceftizoxime and Ceftibuten are the best surrogates for detecting Group B streptococci with reduced susceptibility to Penicillin, using the diffusion method.

Marvaud et al (9) reported the emergence of high level resistance to Gentamicin in *S.pyogenes*.

- Please note that *S. pyogenes* (group A) and *S. agalactiae* (group B) are currently resistant to Trimethoprim+Sulfa.

Nagano et al (12) isolated multiple group B Streptococci with reduced Penicillin susceptibility, suspecting a possible nosocomial spread. They showed Penicillin MICs 0.25 to 1 µg/ml and were resistant to Ceftizoxime and Ceftibuten (good screening for detecting them in routine). In Canada the development of Penicillin non-susceptibility in group B streptococci has been described in vivo in adult patients.

References:

- 1) CLSI: Performance Standards for Antimicrobial Susceptibility Testing 23rd Inf. Suppl. **M100-S23**, 2013.
- 2) NCCLS: Performance Standards for Antimicrobial Susceptibility Tests 10th Ed, **M2-A10**, 2009.
- 3) Alcaide F. et al: In vitro activities of 22 beta-lactam antibiotics against Penicillin resistant and Penicillin susceptible viridans group streptococci isolated from blood. *Antimicrob. Agents Chemother.* **39**, 2243-2247, 1995.
- 4) Jiunn-Jong Wu et al: High incidence of Erythromycin-resistant streptococci in Taiwan. *Antimicrob. Agents Chemother.* **41**, 844-846, 1997.
- 5) Descheemaeker P. et al: Macrolide resistance and Erythromycin resistance determinants among Belgian *S. pyogenes* and *S. pneumoniae* isolates. *J. Antimicrob. Chemother.*, **45**, 167-173, 2000.
- 6) Seppälä H. et al: Three different phenotypes of Erythromycin resistant *Streptococcus pyogenes* in Finland. *J. Antimicrob. Chemother.*, **32**, 885-891, 1993.
- 7) Yiu Wai Chu et al: Invasive group B streptococcus isolates showing reduced susceptibility to Penicillin in Hong Kong. *J. Antimicrob. Chemother* **60**, 1407-1409, 2007.
- 8) Kimura K. et al: First molecular characterization of Group B streptococci with reduced Penicillin susceptibility. *Antimicrob. Ag, Chemother* **52**, 2890-2897, 2008.
- 9) Marvaud J. C. et al: Emergence of high level resistance to Gentamicin in *S. pyogenes*. Abstract Presentation C1-619, 49th ICAAC, 2009.
- 10) Bonfiglio L. et al: Emergence of unrelated isolates of *S. agalactiae* harboring the *inuB* gene. Abstract Presentation C1-616, 49th ICAAC, 2009.
- 11) Kimura K. et al: Practical disk diffusion for detecting group B streptococcus with reduced Penicillin susceptibility. *J. Clin. Microbiol.* **47**, 4154-4157, 2009.
- 12) Nagano N et al: Nosocomial spread of multidrug-resistant group B Streptococci with reduced Penicillin susceptibility belonging to clonal complex 1. *J. Antimicrob. Chemother.* **67**, 849-856, 2012.

EUCAST-and CLSI potency NEO-SENSITABS™

Interpretation Zones and MIC Breakpoints according to CLSI

Haemophilus spp.

Zone diameter interpretative criteria and MIC breakpoints according to CLSI (formerly NCCLS) (1) when testing *Haemophilus* spp. are listed in the table below.

Note that too high inoculum may lead to false-resistant results with some beta-lactam antibiotics. Apply in general no more than 9 Neo-Sensitabs to the surface of a 140-150 mm plate and no more than 4 Neo-Sensitabs on a 90-100 mm plate.

The zone margin should be considered as the area showing no obvious growth visible with the unaided eye. Faint growth, or tiny colonies that may appear to fade from the more obvious zone, should be ignored in the measurement (1).

Table 3.7-1 Interpretation for *Haemophilus influenzae/parainfluenzae*

HTM-agar. Inoculum: McFarland 0.5. Incubation at 35-37 °C in 5-7 % CO₂ for 16-18 hours.

NEO-SENSITABS			Zone diameter in mm			Break-points MIC µg/ml	
	POTENCY	CODE	S	I	R	S	R
	Amoxycillin+Clav	2+1 µg AMC.3	≥ 16	-	< 16	-	-
	Amoxycillin+Clav.	20+10 µg AMC30	≥ 20	-	≤ 19	≤ 4/2	≥ 8/4
b)	Ampicillin	2 µg AMP.2	≥ 18	17-16	< 16	≤ 1	≥ 4
	Ampicillin	10 µg AMP10	≥ 22	21-19	≤ 18	≤ 1	≥ 4
	Ampicillin+Sulbactam	10+10 µg SAM20	≥ 20	-	≤ 19	≤ 2/1	≥ 4/2
c)	Azithromycin	15 µg AZI15	≥ 12	-	-	≤ 4	-
	Aztreonam	30 µg AZT30	≥ 26	-	-	≤ 2	-
	Cefaclor	30 µg CCL30	≥ 22	21-17	≤ 16	≤ 8	≥ 32 (amp R)
c)	Cefepime	30 µg FEP30	≥ 26	-	-	≤ 2	-
c)	Cefixime	5 µg CFM30	≥ 21	-	-	≤ 1	-
c)	Cefotaxime	30 µg	≥ 26	-	-	≤ 2	-
	Cefonicid	30 µg CFCID	≥ 20	19-17	≤ 16	≤ 4	≥ 16
c)	Cefpodoxime	10 µg CPD10	≥ 21	-	-	≤ 2	-
	Ceftaroline	30 µg CPT30	≥ 30	-	-	≤ 0.5	-
c)	Ceftazidime	30 µg CAZ30	≥ 26	-	-	≤ 2	-
c)	Ceftizoxime	30 µg ZOX30	≥ 26	-	-	≤ 2	-
	Ceftriaxone	30 µg CTR30	≥ 26	-	-	≤ 2	-
	Cefuroxime	30 µg CXM30	≥ 20	19-17	≤ 16	≤ 4	≥ 16
	Cephalothin	30 µg CEP30	≥ 22	-	≤ 21	≤ 8	(amp R)
	Chloramphenicol	30 µg CLR30	≥ 29	28-26	≤ 25	≤ 2	≥ 8
a) c)	Ciprofloxacin	5 µg CIPR5	≥ 21	-	-	≤ 1	-
	Clarithromycin	15 µg CLA15	≥ 13	12-11	≤ 10	≤ 8	≥ 32
	Doripenem	10 µg DOR10	≥ 16	-	-	≤ 1	-
	Doxycycline	30 µg DOX30	≥ 27	26-24	≤ 23	≤ 2	≥ 8
	Ertapenem	10 µg ETP10	≥ 19	-	-	≤ 0.5	-
a) c)	Gatifloxacin	5 µg GATIF	≥ 18	-	-	≤ 1	-
c)	Imipenem	10 µg IMI10	≥ 16	-	-	≤ 4	-
a)	Levofloxacin	5 µg LEVOF	≥ 17	-	-	≤ 2	-
c)	Meropenem	10 µg MRP10	≥ 20	-	-	≤ 0.5	-
a)	Moxifloxacin	5 µg MOXIF	≥ 18	-	-	≤ 1	-
a)	Nalidixan (screen quinolones)	30 µg NAL30	-	-	< 19	-	decreased suscept. to quinolones
a) c)	Ofloxacin	5 µg OFL.5	≥ 16	-	-	≤ 2	-
	Penicillin	10 U PEN10	≥ 20	19-17	≤ 16	≤ 1	≥ 4
	Piperacillin/Tazobactam		≥ 21	-	-	≤ 1	-
	Rifampicin	5 µg RIF.5	≥ 20	19-17	≤ 16	≤ 1	≥ 4

Interpretation Zones and MIC Breakpoints according to CLSI

Haemophilus spp.

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
Telithromycin	15 µg	TEL15	≥ 15	14-12	≤ 11	≤ 4	≥ 16
Tetracyclines	30 µg	TET30	≥ 29	28-26	≤ 25	≤ 2	≥ 8
Tigecycline (8)	15 µg	TIG15	≥ 19	-	-	≤ 1	-
Trimethoprim	5 µg	TRIM5	≥ 23	22-20	≤ 19	≤ 1	≥ 4
Trimethoprim+Sulfa	1.25+23.75µg	SxT25	≥ 16	15-11	≤ 10	≤ 0.5/9.5	≥ 4/76

- a) Strains resistant to nalidixic acid as well as strains showing zones < 30 mm with the quinolones should be suspected of having reduced susceptibility to quinolones. The current CLSI/NCCLS MIC breakpoints seem too high for detecting low level resistance.
- b) Ampicillin 2 µg Neo-Sensitabs best detect the beta-lactamase-negative, Ampicillin-resistant (BLNAR) strains. BLNAR should be reported as resistant to Amoxycillin +Clav., Ampicillin+Sulbactam as well as to 1st and 2nd generation cephalosporins (Cefaclor, Cefuroxime). Garcia-Cobos et al (10) in Spain (2007) found an increase to 12.8% of BLNAR strains. Garcia-Cobos et al (11) suggests new Amoxicillin breakpoints in order to detect BLNAR isolates: MIC ≤ 0.5 µg/ml (**S**, no resistance mechanisms), MIC 1-2 µg/ml (**I**, treatable) and MIC ≥ 4 µg/ml (**R**)
Ampicillin 2 µg will show zones of inhibition ≥18 mm (MIC ≤0.5 µg/ml). Cefaclor and Cephalexin (12) are also effective in detecting BLNAR strains.

The majority of *H. influenzae* that are resistant to Ampicillin and Amoxycillin are beta-lactamase producing strains (plasmidic TEM or ROB-1).

Beta-lactamase producing strains are easily detected

- 1) with a rapid beta-lactamase test e.g. Beta-lactamase (acido) Diagnostic Tablets (ROSCO Ref. 45521), using several colonies, because producers and non-producers may coexist in the sample,
- 2) using Amoxycillin+Clavulanate Neo-Sensitabs compared to Amoxycillin alone (same potency). Synergism (larger zone with Amoxycillin+Clavulanate) will be seen in the presence of a beta-lactamase.

ROB-1 resistance is widespread in Australia, including Cefaclor resistant strains (13). PBP3 mutations that conferred reduced susceptibility to Amoxicillin and 3rd generation cephalosporins were found in clinical *H. parainfluenzae* isolates (15). Inhibitor-resistant TEM-34 was first detected in this pathogen.

Differentiation of strains showing beta-lactam resistance:

<i>H. influenzae</i>	AMP 2	AM+CL 20+10 µg	Beta- lactamase	Cefpodoxime+Clav.	Cefaclor
Susceptible	Zone > 18 mm	Zone S	Neg.	No synergism	≥ 23 mm
BLNAR	Zone < 18 mm	-	Neg.	No synergism	≥ 22 mm
Beta-lactamase pos. TEM, ROB-1	No zone	Zone > 20 mm (S)	Pos.	No synergism	≥ 23 mm Variable (ROB-1)
AM+CL resistant BLPACR	No zone	Zone < 20 mm	Pos.	No synergism	zone < 20 mm
ESBL positive	No zone	S or R	Pos/neg.	Synergism	V


BLPACR isolates were defined as AMC, MIC ≥ 4 µg/ml. Both BLNAR and BLPACR are susceptible to Ceftobiprole (14) and Ceftriaxone.

- c) The CLSI has not yet defined other categories than "S" due to the current absence of resistant strains.
- d) If Gentamicin susceptible, report susceptible to Amikacin, Tobramycin and Netilmicin.
- e) *H. influenzae* tested against Pip/Tazo demonstrated uniformly high susceptibility ($MIC \leq 0.5 \mu\text{g/ml}$) including BLNAR strains worldwide (9) and BLPACR.

Søndergaard et al (16) conclude that: Ampicillin 2 μg and Penicillin 1 Unit show an excellent performance in detecting BLNAR isolates. Discs containing Amox + Clav 2 + 1 μg could be used for the detection of N526K in beta-lactamase-producing isolates (BLPACR).

Algorithm: Penicillin 1 unit zones $\geq 12 \text{ mm}$ are considered susceptible to penicillins/cephalosporins. Smaller zones are tested for beta-lactamase. Beta-lactamase positive showing zones $< 23 \text{ mm}$ with Cefaclor 30 μg are considered positive for mutational resistance and considered resistant to Cefuroxime.

Garcia-Cobos et al (17) have studied how to detect the different beta-lactam resistance mechanisms in *H. influenzae*. The algorithm could be:



BLNAR = Beta-lactamase negative Ampicillin resistant

BLPAR = Beta-lactamase positive Ampicillin resistant

BLPACR = Beta-lactamase positive Amoxycillin + Clavulanate resistant

References:

- 1) NCCLS: Performance Standards for Antimicrobial Susceptibility Tests 10th Ed, **M2-A10**, 2009.
- 2) CLSI: Performance Standards for Antimicrobial Susceptibility Testing 23rd Inf. Suppl. **M100-S23**, 2013.
- 3) Zerva L. et al: Reevaluation of interpretative criteria for *H.influenzae* by using Meropenem (10 μg), Imipenem (10 μg) and Ampicillin (2 and 10 μg) disks. J. Clin. Microbiol., **34**, 1970-1974, 1996.
- 4) Campos J. et al: Long-term persistence of Ciprofloxacin-resistant *H.influenzae* in patients with cystic fibrosis. J.I.D., **174**, 1345-1347, 1996.
- 5) Pitout M.J. et al: Characterization of ESBL activity in *Haemophilus parainfluenzae*. 42 ICAAC. Presentation C2-645 (2002).

- 6) Perez-Vázquez M. et al: Laboratory detection of *H. influenzae* with decreased susceptibility to Nalidixic acid, Ciprofloxacin, Levofloxacin and Moxifloxacin due to *gyrA* and *ParC* mutations. J. Clin. Microbiol., **42**, 1185-91, 2004.
- 7) Kärpänöja P. et al: Disc diffusion susceptibility testing of *Haemophilus influenzae* by NCCLS methodology, using low-strength Ampicillin and co-amoxycavulanate discs, J. Antimicrob. Chemother. **53**, 660-663, 2004.
- 8) Brown S.D. et al: Comparative in vitro antimicrobial activity of Tigecycline a new glycycline compound, in freshly prepared medium and Quality Control. J. Clin. Microbiol., **45**, 2173-79, 2007.
- 9) Hirakata Y. et al: Proposed susceptibility testing criteria for *H. influenzae* tested against Piperacillin+Tazobactam. Report using results from the SENTRY program and Japan surveillance isolates. Presentation C-040. ASM meeting, June 2008.
- 10) Garcia-Cobos S et al: Antibiotic resistance in *H. influenzae* decreased except for BLNAR isolates, in parallel with community antibiotic consumption in Spain from 1997 to 2007. Antimicrob. Ag Chemother **52**, 2760-2766, 2008.
- 11) Garcia-Cobos S et al: Low BLNAR *H. influenzae* strains are best detected by testing Amoxicillin susceptibility by the broth microdilution method. Antimicrob. Ag. Chemother **52**, 2407-2414, 2008.
- 12) Endo R. et al: A simple disk diffusion test using Cephalexin to identify beta-lactamase negative Ampicillin resistant *H. influenzae*. Abstract, Presentation D-1432, 49th ICAAC, 2009.
- 13) Tristram S.G et al: bla_{ROB-1} presence on pB1000 in *H. influenzae* is widespread, and variable Cefaclor resistance is associated with altered penicillin-binding proteins. Antimicrob. Agents Chemother. **54**, 4945-47, 2010.
- 14) Davies T.A et al: Characterization of beta-lactamase negative Ampicillin-resistant (BLNAR) and beta-lactamase positive Amoxicillin+Clavulanate resistant (BLPACR) *H. influenzae*. Presentation C2-117, 50th ICAAC, Sept. 2010.
- 15) Garcia-Cobos S. et al: Mechanisms of resistance to beta-lactam antibiotics in *H. parainfluenzae*: first report of an inhibitor-resistant beta-lactamase TEM-34 in this pathogen. Presentation C1-1424, 50th ICAAC, Sept. 2010.
- 16) Søndergaard A et al: Detection of N526K-substituted PBP 3 conferring low-level mutational resistance to beta-lactams in *Haemophilus influenzae* by disc diffusion testing on MH agar according to EUCAST guidelines. J. Antimicrob. Chemother. **67**. 1401-1404, 2012.
- 17) Garcia-Cobos et al: Evaluation of the EUCAST disc diffusion susceptibility testing method for *H. influenzae* based on the resistance mechanism to beta-lactam antibiotics. J. Antimicrob. Chemother. **68**, 159-163, 2013.

Zone diameter interpretative criteria and breakpoints when testing *Moraxella catarrhalis* are listed in the table below. The CLSI (formerly NCCLS) has now established MIC breakpoints or zone-diameter interpretative standards for *Moraxella catarrhalis*.

Table 3.8-1 Interpretation for *Moraxella catarrhalis*

Mueller-Hinton + 5 % blood+20 mg/l β-NAD. Inoculum: McFarland 0.5. Incubation for 16-18 hours at 35 °C ± 2 degrees. CO₂ may be used.

			Zone diameter in mm			Break-points MIC µg/ml	
NEO-SENSITABS	POTENCY	CODE	S	I	R	S	R
a)	Amoxycillin+Clav.	20+10 µg AMC30	≥ 20	-	≤ 19	≤ 4/2	≥ 8/4
	Ampicillin	10 µg AMP10	≥ 30	29-27	≤ 26	≤ 0.12	penase
	Azithromycin	15 µg AZI15	≥ 18	17-13	≤ 12	≤ 2	≥ 8
	Cefuroxime (oral)	30 µg CXM30	≥ 28	27-24	≤ 23	≤ 1	≥ 4
	Chloramphenicol	30 µg CLR30	≥ 26	25-23	≤ 22	≤ 2	≥ 8
	Ciprofloxacin	5 µg CIPR5	≥ 24	-	≤ 23	≤ 1	-
	Clarithromycin	15 µg CLA15	≥ 18	17-13	≤ 12	≤ 2	≥ 8
	Clindamycin	2 µg CLI.2	≥ 20	19-17	≤ 16	≤ 0.5	≥ 4
	Doxycycline	30 µg DOX30	≥ 24	23-21	≤ 20	≤ 2	≥ 8
	Erythromycin	15 µg ERY15	≥ 22	21-17	≤ 16	≤ 0.5	≥ 8
	Gatifloxacin	5 µg GATIF	≥ 30	29-27	≤ 26	≤ 0.25	≥ 0.5
	Levofloxacin	5 µg LEVOF	≥ 22	-	≤ 21	≤ 2	-
	Moxifloxacin	5 µg MOXIF	≥ 30	29-27	≤ 26	≤ 0.25	≥ 0.5
	Nalidixan	30 µg NAL30	-	-	< 18	-	Reduced susceptibility to quinolones
	screen quinolones		-	-	< 18	-	-
a)	Ofloxacin	5 µg OFL.5	≥ 24	-	≤ 23	≤ 1	-
	Penicillin	10 U PEN10	≥ 32	31-27	≤ 26	≤ 0.06	penase
	Tetracyclines	30 µg TET30	≥ 24	23-21	≤ 20	≤ 2	≥ 8
	Trimethoprim+Sulfa	1.25+23.75µg SxT25	≥ 16	15-11	≤ 10	≤ 0.5/9.5	≥ 2/38

a) Beta-lactamase-positive strains may be detected using Amoxycillin+Clavulanate Neo-Sensitabs and Ampicillin Neo-Sensitabs. Synergism (>8 mm larger zone with Amoxycillin+Clavulanate) is seen in the presence of beta-lactamases (BRO-1 or BRO-2). They should be reported as resistant to Penicillin, Ampicillin, Amoxycillin, Ticarcillin and Piperacillin.

b) Report resistant to Penicillin, Ampicillin and Amoxicillin (rare strains do not produce penicillinase).

References:

- Westley Catlin B.: *Branhamella catarrhalis*: an organism gaining respect as a pathogen. Clinical Microbiology Reviews, **Vol. 3**, No. 4, 293-320, 1990.
- Chaibi E.B. et al: Beta lactamases de *Branhamella catarrhalis* et leurs implications phenotypiques. Res. Microbiol., **146**, 761-771, 1995.
- Verduin C.M. et al: *Moraxella catarrhalis*: from emerging to established pathogen. Clin. Microbiol. Rev., **15**, 125-144, 2002.

Zone diameter interpretative criteria and MIC breakpoints when testing meningococci are listed in the table below. The CLSI (formerly NCCLS) has established MIC breakpoints for meningococci, and inhibition zone interpretative standards have now been approved by CLSI (1).

Table 3.9-1 Interpretation for *Neisseria meningitidis*

Mueller-Hinton agar + 5 % blood. McFarland 0.5. Incubation at 35 °C ± 2 degrees with 5-7 % CO₂ for 20 to 24 hours.

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
	Cefotaxime	30 µg	CTX30	≥ 34	-	≤ 0.12	-
	Ceftriaxone	30 µg	CTR30	≥ 34	-	≤ 0.12	-
	Chloramphenicol	30 µg	CLR30	≥ 26	25-20	≤ 2	≥ 8
	Meropenem	10 µg	MRP10	≥ 30	-	≤ 0.25	-
c)	Oxacillin (screen penicillin)	1 µg	OXA.1	≥ 10	-	≤ 0.06 (pen)	MIC
e)	Mecillinam (Screen pen./amp.)	10 µg	MEC10	≥ 23	-	≤ 0.06 (pen)	-
				≥ 23	-	≤ 0.12 (amp)	-
a)	Azithromycin	15 µg	AZI15	≥ 20	-	≤ 2	-
a)	Ciprofloxacin	5 µg	CIPR5	≥ 35	34-33	≤ 0.03	≥ 0.12
	Ciprofloxacin	1 µg	CIPR1	≥ 30	29-26	≤ 0.03	≥ 0.12
a)	Doxycycline	30 µg	DOX30	≥ 26	-	≤ 2	-
a)	Levofloxacin	5 µg	LEVOF	≥ 35	34-33	≤ 0.03	≥ 0.12
a)	Minocycline	30 µg	MIN30	≥ 26	-	≤ 2	-
d)	Nalidixan (screen quinolones)	30 µg	NAL30	-	-	-	≥ 8
a)	Ofloxacin	5 µg	OFL.5	≥ 35	34-33	≤ 0.03	≥ 0.12
a)	Rifampicin	5 µg	RIF.5	≥ 25	24-20	≤ 0.5	≥ 2
a) b)	Trimethoprim+Sulfa	1.25/23.75	SxT25	≥ 30	29-26	≤ 0.12/2.3	≥ 0.5/9.5

a) Used for prophylaxis only (not treatment).

b) Trimethoprim+Sulfa predict susceptibility and resistance to Trimethoprim/Sulfamethoxazole and sulphonamides.

c) Oxacillin 1 µg is useful to screen for beta-lactamase negative meningococci with decreased susceptibility to Penicillin (chromosomal resistance). If the zone is < 10 mm (Oxacillin 1 µg), perform an MIC test for Penicillin.

d) Nalidixic acid is useful to screen for strains with reduced susceptibility to quinolones. If the zone is < 25 mm, control MIC's of Ciprofloxacin, Ofloxacin etc. Strains showing decreased susceptibility to Ciprofloxacin have been reported in Spain (9), Hong Kong (5, 7), USA (8) and Italy (10).

e) Mecillinam is used as a surrogate disk to screen for diminished Penicillin and Ampicillin susceptibility (6).

References:

- 1) CLSI: Performance Standards for Antimicrobial Susceptibility Testing 23rd Inf. Suppl. **M100-S13**, 2013.
- 2) Campos J. et al: Detection of relatively Penicillin G resistant *N. meningitidis* by disk susceptibility testing . Antimicrob. Ag. Chemother., **31**, 1478-1482, (1987).
- 3) Shultz T.R. et al: An invasive isolate of *N. meningitidis* showing decreased susceptibility to quinolones. Antimicro. Ag. Chemother., **44**, 1116, 2000.
- 4) Richter S.S. et al: Neisseria meningitidis with decreased susceptibility to Penicillin: report from the SENTRY antimicrobial surveillance program. North America, 1998-99. Diagn. Microbiol.Infect. Dis., **41**, 83 - 88, 2001.
- 5) Acala B. et al: *N. meningitidis* showing decreased susceptibility to ciprofloxacin: first report in Spain. J. Antimicro. Chemoter., **53**, 409, 2004.
- 6) Jorgensen J.H. et al: Multilaboratory evaluation of disk diffusion antimicrobial susceptibility testing of *N. meningitidis* isolates. J. Clin. Microbiol, **44**, 1744-54, 2006.
- 7) Chu Y. W. et al: A blood isolate of *N. meningitidis* showing reduced susceptibility to quinolones in Hong Kong. Int. J. Antimicrob. Ag., **30**, 94-95, 2007.
- 8) CDC. MMWR: Emergence of fluoroquinolone resistant *N. meningitidis*. Minnesota and North-Dakota, 2007-2008. MMWR **57** (7) 173-175, 2008.
- 9) Enriquez R. et al: Nalidixic acid disk for laboratory detection of Ciprofloxacin resistance in *N. meningitidis*. Antimicro. Ag. Chemoter **53**, 796-797, 2009.
- 10) Lapadula G. et al: Imported Ciprofloxacin resistant *N. meningitidis*. Emerg. Infect. Dis. **15**, 1852-54, 2009.

Zone diameter interpretative criteria and breakpoints according to CLSI (formerly NCCLS) (1,2) when testing gonococci are listed in the table below.

Table 3.10-1 Interpretation for *Neisseria gonorrhoeae*

GC Agar base and 1 % defined growth supplement. McFarland 0.5. Incubation at 35 °C ± 2 degrees in 5 % CO₂ for 20-24 hours.

				Zone diameter in mm			Break-points MIC µg/ml	
NEO-SENSITABS	POTENCY	CODE		S	I	R	S	R
e) Azithromycin	15 µg	AZI15		≥ 24	-	≤ 23	≤ 1	≥ 2
Cefepime	30 µg	FEP30		≥ 31	-	-	≤ 0.5	-
d) Cefixime	5 µg	CFM.5		≥ 31	-	-	≤ 0.25	-
Cefotaxime	30 µg	CTX30		≥ 31	-	-	≤ 0.5	-
Cefotetan	30 µg	CFTTN		≥ 26	25-20	≤ 19	≤ 2	≥ 8
Cefoxitin	30 µg	CFO30		≥ 28	27-24	≤ 23	≤ 2	≥ 8
Cefpodoxime	10 µg	CPD10		≥ 29	-	-	≤ 0.5	-
d) Ceftriaxone	30 µg	CTR30		≥ 35	-	-	≤ 0.25	-
Cefuroxime	30 µg	CXM30		≥ 31	30-26	≤ 25	≤ 1	≥ 4
Ciprofloxacin	1 µg	CIPR1		≥ 32	31-21	≤ 20	≤ 0.06	≥ 1
b) Ciprofloxacin	5 µg	CIPR5		≥ 41	40-28	≤ 27	≤ 0.06	≥ 1
Doxycycline	30 µg	DOX30		≥ 38	37-31	≤ 30	≤ 0.25	≥ 2
a) Nalidixan	30 µg	NAL30		-	-	< 20	decreased susceptibility to quinolones	
(screen quinolones)								
Ofloxacin	5 µg	OFL.5		≥ 31	30-25	≤ 24	≤ 0.25	≥ 2
c) Oxacillin	1 µg	OXA.1		≥ 12	-	-	≤ 0.06	-
(screen penicillin)								
d) Penicillin	10 U	PEN10		≥ 47	46-27	≤ 26	≤ 0.06	≥ 2
Spectinomycin	200 µg	SPECT		≥ 23	22-20	≤ 19	≤ 32	≥ 128
Tetracyclines	30 µg	TET30		≥ 38	37-31	≤ 30	≤ 0.25	≥ 2

- Nalidixic acid is useful to detect strains with reduced susceptibility to quinolones, however, strains of gonococci susceptible to Nalidixan acid and resistant to Ciprofloxacin were detected in the UK (5). Test both Nalidixan and Ciprofloxacin.
- Ciprofloxacin resistant gonococci are presumably resistant to all fluoroquinolones.
- Oxacillin 1 µg Neo-Sensitabs is useful to detect beta-lactamase negative gonococci, with decreased susceptibility to Penicillin (chromosomal resistance).
- A positive beta-lactamase test predicts resistance to Penicillin, Ampicillin, Amoxycillin, Ticarcillin and Piperacillin. Emergence in Japan of beta-lactamase negative strains, resistant to Penicillin (MIC 2-8 µg/ml) and decreased susceptibility to Cefixime (MIC ≥ 0.5 µg/ml) and Ceftriaxone (MIC ≥ 0.125 µg/ml). They contain fragments of PBP-2 from *N. cinerea* and *N. perflava* (5).
- Patients treated with Azithromycin for Chlamydia cannot be assumed to have been adequately treated for gonococcal infection. According to UK Department of Health Azithromycin should not be used to treat gonorrhoea (6).
- Multidrug-resistant gonococci with reduced susceptibility to 3rd generation cephalosporins found in Greece (Tzelepi 7).

Unemo et al (8) indicates that the first *N. gonorrhoeae* strain highly resistant to Ceftriaxone and Cefixime was isolated in Japan. A similar strain was isolated in France. F89 has a Cefixime MIC of 4

µg/ml and Ceftriaxone MIC of 1 to 2 µg/ml. The strain was susceptible to Spectinomycin, but *N. gonorrhoeae* appears to be emerging as a superbug.

Katz et al (9) reports the first isolate of *N. gonorrhoeae* with high level resistance to Azithromycin, identified in the United States. The emergence of both cephalosporin and macrolide resistance would severely limit our treatment options for gonorrhoea.

Unemo et al (10) describe a high level Spectinomycin-resistant (MIC > 1.024 µg/ml) *N. gonorrhoeae* strain from Norway with a novel resistance mechanism. The isolate was susceptible to Ceftriaxone and Cefixime.

References:

- 1) NCCLS: Performance Standards for Antimicrobial Susceptibility Tests 10th Ed, **M2-A10**, 2009.
- 2) CLSI: Performance Standards for Antimicrobial Susceptibility Testing. 23rd Inf. Suppl. **M100-S23**, 2013.
- 3) Gransden W.R. et al: Decreased susceptibility of *Neisseria gonorrhoeae* to Ciprofloxacin. The Lancet **335**, 51, 1990.
- 4) Anon.: Dramatic increase in Ciprofloxacin-resistant gonorrhoea in England and Wales. CDR Weekly, 10th April 2003.
- 5) Regunathan P.L. et al: Nalidixic acid susceptible, Ciprofloxacin-resistant *N. gonorrhoeae* strain in the UK. JAC, **56**, 437, 2005.
- 6) Health Protection Report: High level of Azithromycin resistance in *N. gonorrhoeae*, 4. April 2008.
- 7) Tzelepi E. et al: Cluster of multidrug-resistant *N. gonorrhoeae* with reduced susceptibility to the newer cephalosporins in Northern Greece. JAC **62**, 637-639, 2008.
- 8) Unemo M et al: High-level Cefixime and Ceftriaxone-resistant *N. gonorrhoeae* in France: Novel penA mosaic allele in a successful international clone causes treatment failure. Antimicrob. Agents Chemother. **56**, 1273-1280, 2012.
- 9) Katz AR et al: *N. gonorrhoeae* with high-level resistance to Azithromycin: Case report of the first isolate identified in the United States. Clin. Infect. Dis. **54**, 841-844, 2012.
- 10) Unemo M. et al: *Neisseria gonorrhoeae* strain with high-level resistance to Spectinomycin due to a novel resistance mechanism (mutated ribosomal protein S5) verified in Norway. Antimicrob. Agents Chemother. **57**, 1057-1061, 2013.

Zone diameter interpretative criteria and MIC breakpoints when testing campylobacter are listed in the table below. The CLSI (formerly NCCLS) has established MIC breakpoints for Erythromycin, Ciprofloxacin, Tetracycline and Doxycycline (1). No zone diameter interpretative standards from CLSI are available for campylobacter so far.

Not more than 9 Neo-Sensitabs should be placed on a 150 mm agar plate or 4 Neo-Sensitabs on a 90-100 mm plate.

The plates are incubated at 36-37°C in Campylobacter atmosphere + catalysator for 42 to 48 hours, or at 42 °C for 20-24 hours, and the inhibition zones measured to the nearest mm.

Table 3.11-1 Interpretation for *Campylobacter* spp.

Mueller-Hinton agar with 5 % blood added. McFarland 1.0. Incubation in Campylobacter atmosphere + catalysator at 36 °C for 42-48 hours or at 42 °C for 20-24 hours.

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
Azithromycin	15 µg	AZI15	≥ 22	21-18	≤ 17	≤ 1	≥ 4
Amoxycillin	25 µg	AMX25	≥ 28	27-24	≤ 23	≤ 2	≥ 8
Amoxycillin+Clavulanate	(20+10) µg		≥ 23	22-20	≤ 19	≤ 2	≥ 8
Ampicillin	10 µg	AMP10	≥ 22	21-19	≤ 18	≤ 2	≥ 8
Cephalothin (ID)	30 µg	CEP30	-	-	< 14 (ID)	-	> 32
Chloramphenicol	30 µg	CLR30	≥ 24	23-21	≤ 20	≤ 4	≥ 16
Ciprofloxacin (6,7)	1 µg	CIPR1	≥ 18	14-17	≤ 13	≤ 1	≥ 4
b) Ciprofloxacin (6,7)	5 µg	CIPR5	≥ 25	24-21	≤ 20	≤ 1	≥ 4
Doxycycline	30 µg	DOX30	≥ 25	24-21	≤ 20	≤ 2	≥ 8
b) Erythromycin (6,7)	15 µg	ERY15	≥ 16	15-13	≤ 12	≤ 8	≥ 32
Erythromycin (10)	15 µg	ERY15	≥ 23	-	< 23	≤ 1	> 1
Gentamicin	10 µg	GEN10	≥ 16	15	≤ 14	≤ 4	≥ 16
Imipenem	10 µg	IMI10	≥ 26	25-23	≤ 22	≤ 1	≥ 2
Meropenem	10 µg	MRP10	≥ 26	25-23	≤ 22	≤ 1	≥ 2
a) Nalidixan	30 µg	NAL30					
identification (ID)			-	-	< 14 (ID)	-	> 128
screen quinolones			-	-	< 18	-	≥ 32
Tetracyclines	30 µg	TET30	≥ 25	24-21	≤ 20	≤ 4	≥ 16

ID = for identification purposes only.

- a) Strains resistant to Nalidixic acid show a decreased susceptibility to fluoroquinolones.
- For *Campylobacter* spp. the absence of zone of inhibition around beta-lactams, aminoglycosides, macrolides or quinolones indicates high level resistance.
- b) According to CLSI, strains showing no zone of inhibition with Ciprofloxacin and Erythromycin should be reported as resistant to these antibiotics.

QC ranges for *Campylobacter jejuni* ATCC 33560 (9)

Neo-Sensitabs™	Potency	Zone diameter in mm.
Ciprofloxacin	5 µg	32-45 mm
Erythromycin	15 µg	26-38 mm
Nalidixic acid	30 µg	25-34 mm
Tetracycline	30 µg	29-39 mm

Note: The proposed zone diameters are for 24-hour readings only

Lehtopolku et al. (10) recommend a lower MIC breakpoint and larger zones of inhibition for Erythromycin and *Campylobacter*. Erythromycin resistant strains (zone < 23mm) were uniformly resistant.

Perez Pomata et al. (11) conclude that the diffusion method is adequate for testing Erythromycin and Ciprofloxacin against isolates of *Campylobacter jejuni*.

Gaudreau (12) has published studies about the disk diffusion method for Erythromycin and Ciprofloxacin susceptibility testing of *campylobacter* spp. Here results match the CLSI Erythromycin and Ciprofloxacin disk diffusion method standardization: all isolates without zone of inhibition around Erythromycin 15 µg and Ciprofloxacin 5 µg disks were resistant to macrolides and Ciprofloxacin.

References:

- 1) CLSI: Methods for Antimicrobial Dilution and Disk Susceptibility Testing of infrequently isolated fastidious bacteria. Approved Guideline. Suppl **M45-A**, 2006.
- 2) Vanhoof R. et al: Disk sensitivity testing for *Campylobacter jejuni*. Eur. J. Clin. Microbiol., **3**, 160-162, 1984.
- 3) Engberg J. et al: Comparison of 2 agar dilution methods and 3 agar diffusion methods, including E-test, for antibiotic susceptibility testing of thermophilic *Campylobacter* species. Clin. Microbiol. Infect., **5**, 580-584, 1999.
- 4) Lucero C. et al: *Campylobacter* spp.: Comparison of disc diffusion and agar dilution methods for susceptibility testing. ICAAC 2003, Presentation D-238.
- 5) McDermott P.F. et al: Development of a standardized susceptibility test for *campylobacter* with Q.C. ranges for Ciprofloxacin, Doxycycline, Erythromycin, Gentamicin and Meropenem. Microb. Drug Resist., **10**, 124-131, 2004.
- 6) Schönberg-Norio D. et al: Activities of Telithromycin, Erythromycin, fluoroquinolones, and Doxycycline against *Campylobacter* strains isolated from Finnish subjects. Antimicrob. Ag. Chemother., **50**, 1086-8, 2006.
- 7) Gaudreau C. et al: Comparison of disk diffusion and agar dilution methods for Erythromycin and Ciprofloxacin susceptibility testing of *C. jejuni* subsp. *jejuni*. Antimicrob. Ag. Chemother., **51**, 1524-26, 2007.
- 8) Gaudreau C. et al: Comparison of disk diffusion and agar dilution methods for Erythromycin, Ciprofloxacin and Tetracycline susceptibility testing of *Camp. jejuni* subsp *jejuni*. Antimicrob. Ag. Chemother **52**, 4475-77, 2008.
- 9) Disk Diffusion QC ranges for Ciprofloxacin, Nalidixic acid, Erythromycin and Tetracycline against *C. jejuni*. Summary Subcom. Vet. Antimicrob. Suscept. Testing 26-27, January 2010, page 7.
- 10) Lehtopolku et al.: Antimicrobial Susceptibilities of Multidrug Resistant *C. jejuni* and *C. coli* strains: in vitro activity of 20 antimicrobial agents. Antimicrob. Ag. Chemother., **54**, 1232-1236, 2010
- 11) Perez Pomata MT et al.: Susceptibility of *C. jejuni* Against Erythromycin and Ciprofloxacin by the diffusion method (Spanish). Presentation 033, SEIMC 2010, Barcelona.
- 12) Gaudreau C: Disk diffusion method for Erythromycin and Ciprofloxacin susceptibility testing of *C. jejuni* and *C. coli*. J. Clin. Microbiol. **51**, 380, 2013.

Zone diameter interpretative criteria and MIC breakpoints when testing *Helicobacter pylori* are listed in the table below. Only breakpoints for Clarithromycin are available from CLSI (formerly NCCLS) (1). Other MIC breakpoints or zone diameter interpretative standards for *Helicobacter pylori* have not yet been established by CLSI.

Table 3.12-1 Interpretation for *Helicobacter pylori*

Mueller-Hinton agar with 5-10 % blood added. McFarland 3 or 4. Incubation at 35-37 °C in a microaerophilic atmosphere (Campylobacter) for 72 hours.

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
Ampicillin	10 µg	AMP10	≥ 26	25-21	≤ 20	≤ 0.5	>1
Azithromycin	15 µg	AZI15	≥ 26	25-21	≤ 20	≤ 0.25	≥ 1
Amoxycillin (use Ampicillin 10 µg)							
Ciprofloxacin	5 µg	CIPR5	≥ 26	25-23	≤ 22	≤ 0.5	≥ 1
Clarithromycin	15 µg	CLA15	≥ 26	25-21	≤ 20	≤ 0.25	≥ 1
Doxycycline	30 µg	DOX30	≥ 30	29-26	< 26	≤ 1	>1
Erythromycin	15 µg	ERY15	≥ 20	19-17	≤ 16	≤ 1	≥ 4
Levofloxacin	5 µg	LEVOF	≥ 26	-	< 26	≤ 0.5	≥ 1
Metronidazole	16 µg	MTR16	≥ 26	26-23	≤ 22	≤ 8	≥ 16
Moxifloxacin	5 µg	MOXIF	≥ 26	-	< 26	≤ 0.5	≥ 1
Rifampicin	5 µg	RIF.5	≥ 21	-	< 21	-	-
Tetracyclines	30 µg	TET30	≥ 30	29-26	< 26	≤ 1	>1

Because of the increasing resistance of *H. pylori* against Metronidazole and Clarithromycin, alternative regimens including newer fluoroquinolones have been suggested. Bogaerts et al. (8) showed that *H. pylori* resistance to fluoroquinolones is occurring at high frequency in Belgium and consequently the susceptibility of the infecting isolates to fluoroquinolones should be determined before administering these agents.

Rasmussen et al (10) found an excellent correlation between E-test and Neo-Sensitabs for Metronidazole and Ciprofloxacin and acceptable correlation for Ampicillin, Tetracycline and Clarithromycin.

Wüppenhorst et al. (11) isolated one strain of *H. pylori* being resistant to Metronidazole (MIC ≥ 256 µg/ml), Clarithromycin (MIC 16 µg/ml), Levofloxacin (MIC ≥ 32 µg/ml) and Rifampicin (MIC ≥ 32 µg/ml), but susceptible to Amoxycillin (MIC 0.05 µg/ml). Multiresistant clinical *H. pylori* isolates exist in Germany and will probably increase in the future.

References:

- 1) CLSI: Performance Standards for Antimicrobial Susceptibility Testing. 23rd Inf. Suppl. **M100-S23**, 2013.
- 2) Huaxiang X. et al: Standardisation of Disk Diffusion Test and clinical significance for susceptibility testing of Metronidazole against *Helicobacter pylori*. Antimicrob. Ag. Chemother., **38**, 2357-2361, 1994.
- 3) Sørberg M. et al: Risk of development of in vitro resistance to Amoxycillin, Clarithromycin and Metronidazole in *Helicobacter pylori*. Antimicrob. Ag. Chemother., **42**, 1222-1228, 1998.
- 4) Kalach N. et al: High levels of resistance to Metronidazole and Clarithromycin in *Helicobacter pylori* strains in children. J. Clin. Microbiol., **39**, 394-397, 2001.
- 5) Glupczynski Y. et al: European multicenter survey of in vitro antimicrobial resistance in *Helicobacter pylori*. EJCMID **20**, 820-823, 2001.

- 6) McNulty C. et al: Helicobacter pylori susceptibility testing by disc diffusion. J. Antimicrob. Chemother., **49**, 601-609, 2002.
- 7) Tankovic J. et al: Single and double mutations in gyrA but not in gyrB are associated with low and high-level fluoroquinolone resistance in *H. pylori*. Antimicrob. Ag. Chemother., **47**, 3942-44, 2003.
- 8) Bogaerts P. et al: Prevalence and mechanisms of resistance to fluoroquinolones in *H. pylori* strains from patients living in Belgium. Helicobacter, **11**, 441-5, 2006.
- 9) Henriksen T.H. et al: A simple method for determining Metronidazole resistance of *H. pylori*. J. Clin. Microbiol., **35**, 1424-26, 1997.
- 10) Rasmussen L. et al: Comparison of the E-test and the Neo-Sensitabs method for antimicrobial susceptibility testing of *H. pylori*. ECCMID poster P-1179, Barcelona 2008.
- 11) Wüppenhorst N. et al: Isolation and eradication of clinical isolate of Helicobacter *pylori* resistant to five antimicrobials in Germany. J. Antimicrob. Chemother. **66**, 222-223, 2011.

EUCAST-and CLSI potency NEO-SENSITABS™ Interpretation Zones and MIC Breakpoints according to CLSI

Zone diameter interpretative criteria and MIC breakpoints according to CLSI (formerly NCCLS) (1) when testing *Vibrio cholerae* are listed in the table below.

Table 3.13-1 Interpretation for *Vibrio cholerae*

Mueller-Hinton agar. Inoculum McFarland 0.5. Incubation 35 °C ± 2 degrees ambient air for 16-18 hours.

				Zone diameter in mm			Break-points MIC µg/ml	
NEO-SENSITABS	POTENCY	CODE		S	I	R	S	R
a) Ampicillin	10 µg	AMP10		≥ 17	16-14	≤ 13	≤ 8	≥ 32
Tetracyclines	30 µg	TET30		≥ 19	18-15	≤ 14	≤ 4	≥ 16
Sulphonamides	240 µg	SULFA		≥ 23	22-20	≤ 19	≤ 100	≥ 350
Trimethoprim+Sulfa	1.25+23.75µg	SxT25		≥ 16	15-11	≤ 10	≤ 2/38	≥ 8/152
Chloramphenicol	30 µg	CLR30		≥ 18	17-13	≤ 12	≤ 8	≥ 32
Doxycycline	30 µg	DOX30		≥ 16	15-13	≤ 12	≤ 4	≥ 16
b) Nalidixan	30 µg	NAL30		-	-	< 15	-	reduced susceptibility to quinolones
- screening quinolones								
Azithromycin	15µg	AZI15						

a) Results for Ampicillin are used to predict susceptibility to Amoxycillin.

b) Strains resistant to Nalidixic acid show a decreased susceptibility to fluoroquinolones.

References:

- 1) CLSI: Performance Standards for Antimicrobial Susceptibility Testing. 23rd Inf. Suppl. **M100-S23**, 2013.
- 2) Mukhopadhyay A.K et al: Emergence of Fluoroquinolone resistance in strains of *V. cholerae* isolated from hospitalized patients with acute diarrhoea in Calcutta, India. Antim. Ag. Chemother., **42**, 206-207, 1998.
- 3) Garg et al: Emergence of fluoroquinolone-resistant strains of *Vibrio cholerae* 01 Biotype El Tor among hospitalized patients with cholera in Calcutta, India. Antim. Ag. Chemother., **45**, 1605-1606, 2001.
- 4) Lai-King Ng et al: Can E-test be used to determine *Vibrio cholerae* susceptibility to Erythromycin? Antimicrob. Ag. Chemother., **47**, 1479-80, 2003.

Interpretation Zones and MIC Breakpoints according to CLSI

Anaerobes

Organisms recognized as virulent or commonly resistant, should be considered for testing. These include the *Bacteroides fragilis* group, *Prevotella* and *Porphyromonas* spp., *Clostridium perfringens*, *Cl. ramosum*, *Cl. septicum* and *Cl. difficile*.

The recommended procedure for susceptibility testing of anaerobes by the diffusion method is the following:

- Use supplemented Brucella blood agar, it supports good growth for essentially all anaerobes. Brucella agar base is supplemented with 5 µg/ml haemin, 5% defibrinated sheep blood and 1 µg/ml vitamin K1 (haemin and vitamin K1 may be added before sterilisation).
- Direct suspension of colonies in broth to achieve a turbidity equivalent to a 1.0 McFarland standard (3×10^8 CFU/ml). Streak the surface of the agar with a cotton swab.
- Allow the inoculated plate to remain at room temperature of $\leq 25^\circ\text{C}$ (5-10 min.) until the surface of the agar looks dry. For some fastidious isolates that do not grow on control plates, pre-reduction of plates in an anaerobic environment may be necessary. Apply Neo-Sensitabs tablets.
- Invert the inoculated plate and incubate at 35°C in an anaerobic jar or alternative anaerobic environment, for 24-48 hours.

Zone diameter interpretative standards are correlated to the MIC breakpoints recommended by the CLSI for anaerobic bacteria (1), but zone diameter interpretative standards have not yet been established by CLSI.

Table 3.14-1 Interpretation for Anaerobes

Supplemented Brucella Blood agar. Inoculum: McFarland 1.0. Incubation at 35°C in anaerobic environment for 24-48 hours.

NEO-SENSITABS				Zone diameter in mm			Break-points MIC µg/ml	
				S	I	R	S	R
b)	Amoxycillin+Clav.	20+10 µg	AMC30	≥ 24	23-21	≤ 20	$\leq 4 / 2$	$\geq 16 / 8$
	Cefoxitin	30 µg	CFO30	≥ 16	15-12	≤ 11	≤ 16	≥ 64
	Chloramphenicol	30 µg	CLR30	≥ 23	22-20	≤ 19	≤ 8	≥ 32
	Clindamycin	2 µg	CLI.2	≥ 15	14-10	No zone	≤ 2	≥ 8
	Doxycycline	30 µg	DOX30	≥ 22	21-19	≤ 18	≤ 4	≥ 16
	Doripenem	10 µg	DOR10				≤ 2	≥ 8
a)	Ertapenem	10 µg	ERTAP	≥ 23	22-20	≤ 19	≤ 4	≥ 16
	Fucidin	10 µg	FUC10	≥ 18	-	≤ 18	≤ 2	≥ 4
a)	Imipenem	10 µg	IMI10	≥ 22	21-18	≤ 17	≤ 4	≥ 16
c)	Imipenem+EDTA	10+750 µg	IMI10E	Detection of metallo-β-lactamases				
	Linezolid	30 µg	LINEZ	≥ 21	-	< 21	≤ 4	> 4
c	Meropenem	10 µg	MRP10	≥ 22	21-18	≤ 17	≤ 4	≥ 16
b)	Metronidazole	16 µg	MTR16	≥ 22	21-19	≤ 18	≤ 8	≥ 32
	Clostridium difficile			≥ 28	-	-	≤ 2	-
	Moxifloxacin	5 µg	MOXIF	≥ 23	22-20	≤ 19	≤ 2	≥ 8
	Penicillin (Ampicillin)	10 U	PEN10	≥ 25	24-21	≤ 20	≤ 0.5	≥ 2
	Piperacillin+Tazobactam	100+10µg	PI+TZ	≥ 24	23-21	≤ 20	$\leq 32 / 4$	$\geq 128 / 4$
	Teicoplanin	30 µg	TPN30	≥ 18	-	-	≤ 4	-
	Ticarillin+Clavulanate	75+10 µg	TIM85	≥ 24	23-21	≤ 20	$\leq 32 / 2$	$\geq 128 / 2$
	Tigecycline	15 µg	TIG15	≥ 22	21-16	≤ 14	≤ 2	≥ 8
	Vancomycin	30 µg	VAN30	≥ 18	-	-	≤ 4	-

- The carbapenems show typically MIC's of ≤ 0.12 µg/ml against *B. fragilis* (inhibition zones ≥ 30 mm). Strains showing MIC's of 2-4 µg/ml (zones: 22-26 mm) may sometimes conduct to therapeutic failure (increased expression of *cfiA* (6)).
- The SFM recommends Metronidazole Neo-Sensitabs for sensitivity testing of anaerobes, due to stability problems with paper disks.
Baines et al (8) found that *Cl. difficile* ribotype 001 had reduced susceptibility to Metronidazole (MIC 3.5 µg/ml). To detect these strains, the susceptible MIC breakpoint should be lowered to $S \leq 2$ µg/ml (zone ≥ 28 mm (S)).
- See description of technique for the detection of metallo-β-lactamases in "Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™".
- Peptococcus/Bacteroides* spp: If resistant to Erythromycin, but susceptible to Clindamycin, report resistant to Clindamycin (MLS_B phenotype).

Wybo et al (7) in a Belgian study found that the prevalence of *cfiA* resistance gene was 6.7% in *B. fragilis*. Using the CLSI breakpoints (4/16 µg/ml), 3 of the 9 strains would be reported as susceptible to Meropenem, but with the EUCAST breakpoints (2/8 µg/ml) no false susceptible were obtained.

Quality Control Limits for anaerobes

Suppl. Brucella Blood Agar
Inoc. McF 1.0 Anaerobic incubation

Neo-Sensitabs	CODE	Bact. Fragilis ATCC 25285	MIC µg/ml	B.thetaiotaomicron ATCC29741	MIC µg/ml
Amoxycillin+Clav. 20+10 µg	AMC30	30-38 mm	0.5	33-38 mm	0.5
Cefoxitin 30 µg	CFO30	21-28 mm	8	17-24 mm	16
Cloramphenicol 30 µg	CLR30	26-32 mm	8	26-32 mm	16
Clindamycin 2 µg	CLI2	14-20 mm	1	14-20 mm	4
Imipenem 10 µg	IMI10	32-40 mm	0.06	32-40 mm	0.12
Meropenem 10 µg	MRP10	29-35 mm	0.12	29-35 mm	0.12/0.25
Metronidazole 16 µg	MTR16	29-36 mm	0.5	30-36 mm	1

Goldstein et al (9) indicates that susceptibility testing of anaerobes is performed in a minority of laboratories. Beta-lactamase is the most common resistance mechanism. Particularly in *Bacteroides fragilis* group are common but also found in *Clostridia*, *Porphyromonas* and *Fusobacterium*. Carbapenemases (metallo-beta-lactamases) are found in the *B. fragilis* group and some *Prevotella* spp. Metronidazole resistance has been reported worldwide in the *B. fragilis* group. Fluoroquinolone resistance (moxifloxacin) is found in the *B. fragilis* group and in *C. difficile*.

References:

- NCCLS Methods for Antimicrobial Susceptibility Testing of Anaerobic Bacteria, 8th Ed. Approved Standard, **M11-A8**, 2008.
- Edwards R. et al: Mechanism responsible for reduced susceptibility to Imipenem in *Bacteroides fragilis*. J. Antimicrob. Chemother., **38**, 941-951, 1996.
- Cormican M.G. et al: False resistance to Metronidazole by E-test, among anaerobic bacteria. Diagn. Microbiol. Infect. Dis., **24**, 117-119, 1996.
- Fang H. et al: Detection of Imipenem-resistant and Metronidazole-resistant *Bact. fragilis* group strains in fecal samples. Clin. Microbiol. Infect., **5**, 753-8, 1999.
- Aldrige K.E.: Detection of carbapenemase production in the *B. fragilis* group using the MBL E-test strip: correlation with broth microdilution testing. Clin. Microbio. Inf., **11**, Suppl. 2, 625, 2005.

EUCAST-and CLSI potency NEO-SENSITABS™

Page 3 / 3 - Document: 3.14.0

Interpretation Zones and MIC Breakpoints according to CLSI

Anaerobes

- 6) Lei Ang et al: Carbapenem resistance in *Bacteroides fragilis*. JAC, **59**, 1042-43, 2007.
- 7) Wybo et al: Carbapenem (cfiA)-resistance gene in *B. fragilis* group strains in Belgium. Poster P-912. ECCMID, Barcelona, april 2008.
- 8) Baines S. et al: Emergence of reduced susceptibility to Metronidazole in *Clostridium difficile*. ECCMID, abstract O463, april, 2008, Barcelona.
- 9) Goldstein EJC et al: Resistance trends in antimicrobial susceptibility of anaerobic bacteria. Clin. Microbiol. Newsletter. **33**, 1-15, 2011.

ORIGINAL ROSCO DOCUMENT

Only Fluconazole and Voriconazole are included in the M44-A guideline from CLSI (formerly NCCLS)(1) and interpretation zones are only available for Fluconazole so far. When MIC breakpoints are available from CLSI (2) they are used in the table below.

Table 13.5-1 Interpretation for Yeasts

Mueller-Hinton Glucose Methylene Blue agar. Inoculum McFarland 0.5 undiluted. Incubation at 35 °C ± 2 °C ambient air for 20-24 hours. MICs according to M27-S3 (2007).

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
Amphotericin B	10 µg	AMPHO	≥ 15	14-10	< 10	≤ 1	≥ 2
Caspofungin *) (24)	5 µg	CASP5					
C. albicans, C. tropicalis			≥ 17	16-15	≤ 14	≤ 0.25	≥ 1
C. parapsilosis,* C. guilliermondii			≥ 13	12-11	≤ 10	≤ 2	> 2
C. glabrata, C. krusei			≥ 16	-	-	≤ 0.5	-
Fluconazole **) (27)	25 µg	FLUCZ					
C. albicans, C. tropicalis, C. parapsilosis			≥ 17	16-14	≤ 13	≤ 2	≥ 8
C. glabrata			-	≥ 15 (SDD)	≤ 14 (SDD)	≤ 32 (SDD)	≥ 64
Itraconazole (11)	10 µg	ITRAC	≥ 23	22-14 (SDD)	≤ 13	≤ 0.12	> 0.5
Ketoconazole	15 µg	KETOC	≥ 28	27-21	≤ 20	≤ 0.12	≥ 0.5
Posaconazole (9,12)	5 µg	POSAC	≥ 17	16-14 (SDD)	≤ 13	≤ 1	> 2
Voriconazole (23)	1 µg	VOR.1					
C. albicans, C. tropicalis, C. parapsilosis			≥ 17	16-15	≤ 14	≤ 0.12	≥ 1
C. krusei			≥ 15	14-13	≤ 12	≤ 0.5	≥ 2
C. glabrata			≥ 16	-	≤ 15	≤ 0.5	> 0.5

*) Tentative. There is cross-resistance between Caspofungin and the other echinocandins: Anidulafungin and Micafungin (19).

**) *C. krusei* should be reported as resistant to Fluconazole (no matter the zone).

For further information on Susceptibility Testing of Yeasts, see Neo-Sensitabs User's Guide.

Table 13.5-2 Mold Disk Diffusion Testing (14, 17,21)

Mueller-Hinton Plain. Incubation for 16-24 hours (zygomycetes), 24 hours (Aspergillus), 48 hours (other spp.). Temperature: 35-37°C

Inoculum:

Prepare a suspension of sporulating colonies in 0.85% saline, add 1 drop of Tween 20. Allow heavy particles to settle for 3-5 minutes and the upper suspension is treated for 15 seconds in a vortex mixer. The density of the suspension is read on a spectrophotometer at 530 nm wave length and the optical density adjusted at 0.09 to 0.13 for Aspergillus.

EUCAST-and CLSI potency NEO-SENSITABS™

Page 2 / 4 - Document: 3.15.0

Interpretation Zones and MIC Breakpoints according to CLSI

Yeasts (CLSI's M44-A2 Method of Diffusion
on Agar)

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
Amphotericin B (zygomycetes only)	10 µg	AMPHO	≥ 15	-	-	≤ 1	-
Caspofungin	5 µg	CASP5	≥ 17	-	-	≤ 1	-
Itraconazole	10 µg	ITRAC	≥ 17	-	-	≤ 1	≥ 2
Posaconazole	5 µg	POSAC	≥ 23	-	-	≤ 0.25	≥ 0.5
Voriconazole	1 µg	VOR.1	≥ 17	-	-	≤ 1	≥ 2

Note: The base medium should not be supplemented with neither 2% glucose nor 0.5% methylene blue dye.

Itraconazole susceptible strains should be reported as susceptible to both, Posaconazole and Voriconazole.

Interpretation table for Local treatment

In local treatment of fungal infections, a high concentration of antifungal is placed at site of the infection. Consequently other MIC breakpoints and zone interpretations should be used in those cases.

Local Treatment MH Glucose Methylene Blue Agar or Shadomy McFarland 0.5 inoculum			
Susceptible	≥ 20 mm	≥ 15 mm	≥ 10 mm
Intermediate	12-19 mm	10-14 mm	-
Resistant	≤ 11 mm	no zone	no zone
	Ciclopirox Clotrimazole Econazole Fluconazole Isoconazole Ketoconazole Miconazole Tioconazole Terbinafine	Natamycin Nystatin Itraconazole	Griseofulvin

Fluorocytosine cannot be tested on MH-agar (antagonists), but has to be tested on Shadomy agar or similar.

Table 13.5-3 Quality Control Zone Diameters (mm) Ranges

Mueller-Hinton Glucose Methylene Blue agar. Inoculum McFarland 0.5 undiluted. Incubation at 35 °C ± 2 °C for 20-24 hours.

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm		
			<i>C. albicans</i> ATCC 90028	<i>C. parapsilosis</i> ATCC 22019	<i>C. krusei</i> ATCC 6258
Amphotericin B	10 µg	AMPHO	20-27	22-29	18-25
Fluconazole	25 µg	FLUCZ	28-39	22-33	-
Itraconazole	10 µg	ITRAC	21-30	19-26	16-22
Ketoconazole	15 µg	KETOC	31-42	35-45	22-29
Voriconazole	1 µg	VOR.1	31-42	28-37	23-31
Caspofungin	5 µg	CASP5	15-22	13-23	16-22
Posaconazole	5 µg	POSAC	24-34	25-36	23-31

Carrillo-Munoz et al (30) determined Posaconazole susceptibility of clinical yeast isolates with Neo-Sensitabs and a microdilution method. Complete agreement between Posaconazole Neo-Sensitabs and the microdilution was 92.3 % after 24 hours incubation. The authors conclude that the Agar Diffusion with Posaconazole Neo-sensitabs can improve Posaconazole susceptibility testing due to its excellent correlation and reduced percentage of disagreements in comparison with microdilution testing.

References:

- 1) NCCLS. Methods for Antifungal Disk Diffusion Susceptibility Testing on Yeasts. Approved Guideline M44-A2, 2008.
- 2) NCCLS. Reference Methods for Broth Dilution Antifungal Susceptibility Testing of yeasts. Approved Standard, 2nd Ed. M27-A2, 2002.
- 3) Vandenbossche I. et al: Susceptibility testing of Fluconazole by the NCCLS broth microdilution method. E-Test and Disk diffusion for application in the routine laboratory. J. Clin. Microbiol, **40**, 918-921, 2002.
- 4) Pfaller M.A. et al: Evaluation of the NCCLS M44-P Disk diffusion method for determining Fluconazole susceptibility of 276 clinical isolates of *Cryptococcus neoformans*. ICAAC 2003, presentation M-1204.
- 5) Barry A. et al: Quality Control for Fluconazole disk susceptibility tests on Mueller-Hinton Agar with glucose and methylene blue. J.Clin.Microbiol., **41**, 3410-12, 2003.
- 6) Rementeria A. et al: Utility of Neo-S. tablets of Fluconazole and Voriconazole for in vitro susceptibility testing of *Candida* spp. with the NCCLS M44-P method of diffusion on agar. VII Congress Mycology. Salamanca, July 2004 (Spanish).
- 7) Espinel-Ingroff A. et al: Correlation between Neo-S. tablets on 3 media, NCCLS reference disk diffusion and broth microdilution methods for testing *Candida* spp. and *Cryptococcus neoformans* with Fluconazole and Voriconazole. Abstract, ECCMID 2005.
- 8) Espinel-Ingroff A, Canton E.: Comparison of 3 antifungal susceptibility methods for testing *Candida* spp. & *Cryptococcus neoformans* with Caspofungin (CAS) and Amphotericin B (AMB). Abstract, presentation M-1604, 45th ICAAC 2005.
- 9) Sims C.R. et al: Correlation between microdilution, E-test and disk diffusion methods for antifungal susceptibility testing of Posaconazole against *Candida* spp. J. Clin. Microbiol., **44**, 2105-08, 2006.
- 10) Brown S. et al: Q.C. limits for Posaconazole disk susceptibility tests on MH agar with glucose and methylene blue. J. Clin. Microbiol., **45**, 222-23, 2007.
- 11) Espinel-Ingroff A. et al: Correlation of Neo-Sensitabs tablet diffusion assay results on 3 different agar media with CLSI broth microdilution M27-A2 and disk diffusion M44-A results for testing

- susceptibilities of *Candida* spp. and *Cryptococcus neoformans* to Amphotericin B, Caspofungin, Fluconazole, Itraconazole and Voriconazole. *J. Clin. Microbiol.* **45**, 858-64, 2007.
- 12) Diekema D.J. et al: Evaluation of E-test and disk diffusion methods with broth microdilution antifungal susceptibility testing of Posaconazole against clinical isolates of *Candida*. *J. Clin. Microbiol.*, 2007 (in press).
 - 13) Carrillo-Muñoz A.J. et al: Activity of Caspofungin and Voriconazole against clinical isolates of *Candida* and other medically important yeasts by the CLSI M44A disk diffusion method with Neo-Sensitabs tablets. *Chemotherapy*, **54** (1) 38-42, 2008.
 - 14) Espinel-Ingróff A. et al: Multicenter evaluation of a new disk agar diffusion method for susceptibility testing of filamentous fungi with Voriconazole, Posaconazole, Itraconazole, Amphotericin B and Caspofungin. *J. Clin. Microbiol.*, **45**, 1811-20, 2007.
 - 15) Espinel-Ingróff A.: Standardized disk diffusion method for yeasts. *Clin. Microbiol. Newsl.*, **29**, 97-100, 2007.
 - 16) Carrillo-Muñoz A. J. et al: Antifungal activity in vitro of Posaconazole. Method of diffusion on agar (Neo-Sensitabs) based on document M44-A, against yeast. III. Forum Micológico, Dec 15th 2007, Alicante, Spain (spanish).
 - 17) Espinel Ingróff A. et al: Comparison of Neo-S tablet diffusion assay with CLSI broth microdilution M38-A and disk diffusion methods for testing susceptibility filamentous fungi with Amphotericin B, Caspofungin, Itraconazole, Posaconazole and Voriconazole. *J. Clin. Microbiol.* **46**, 1793-1803, 2008.
 - 18) Duran-Valle M T et al: Resistance to Voriconazole and/or Fluconazole in yeast isolated in a Madrid hospital (2005-2007) etc. XII Congress SEIMC, presentation 33, Madrid, May 2008.
 - 19) Pfaller M. A. et al: Correlation of MIC with outcome for *Candida* spp. tested against Caspofungin, Amphotericin B and Micafungin: analysis and proposal for interpretative MIC breakpoints. *JCM* **46**, 2620-2629, 2008.
 - 20) Pfaller M. A. et al: Wild-type MIC distribution and epidemiological cutoff values for *Aspergillus fumigatus* and 3 triazoles as determined by the CLSI both microdilution methods. *J. Clin. Microbiol.* **47**, 3142-46, 2009.
 - 21) M51-A Method for antifungal disk diffusion susceptibility testing of filamentous fungi approved guideline. CLSI 2009.
 - 22) Pfaller M. et al: Clinical breakpoints for the Echinocandins and *Candida*, revisited. Presentation M-369, 50th ICAAC, Sept. 2010.
 - 23) Pfaller M. et al: Clinical breakpoints for Voriconazole and *Candida* spp. revisited: Review of microbiologic, molecular, pharmacodynamic and clinical data as they pertain to the development of species-specific interpretative criteria. *Diagn. Microbiol. Infect. Dis.* **70**, 330-343, 2011.
 - 24) Pfaller M. et al: Clinical breakpoints for Echinocandins and *Candida* revisited: Integration of molecular, clinical and microbiological data to arrive at species-specific interpretative criteria. *Drug Resistance Updates*. **14**, 164-176, 2011.
 - 25) Ochiuzzi M. E et al: Correlation of E-test and Neo-Sensitabs diffusion assays on Mueller-Hinton-methylene blue agar, with broth microdilution reference method (CLSI-M27-A2) for testing susceptibilities of *Cryptococcus neoformans* to Amphotericin B and Fluconazole. *Med. Mycology*. **48**, 893-896, 2010.
 - 26) Carrillo-Muñoz A. J. et al: Antifungal activity of Posaconazole against *Candida* spp. and non-*Candida* clinical yeast isolates. *Rev. Esp. Quimioter.* **23**, 122-125, 2010.
 - 27) Nweze E. I et al: Agar-based disk diffusion assay for susceptibility testing of Dermatophytes. *J. Clin. Microbiol.* **48**, 3750-52, 2010.
 - 28) Pfaller M. A et al: Wild-type MIC distributions, epidemiological cutoff values and species-specific clinical breakpoints for Fluconazole and *Candida*: Time for harmonization of CLSI and EUCAST broth microdilution methods. *Drug Resistance Updates*. **13**, 180-195, 2010.
 - 29) Pfaller M. A et al: Echinocandin and triazole antifungal susceptibility profiles for *Candida* spp., *Cryptococcus neoformans* and *Asp. fumigatus*; application of new CLSI clinical breakpoints and epidemiologic cutoff values to characterize resistance in SENTRY Antimicrobial Surveillance Program (2009). *Diag. Microbiol. Infect. Dis.* **69**, 45-50, 2011.
 - 30) Carrillo-Munoz AJ et al: Posaconazole susceptibility of clinical yeast isolates determined by an agar diffusion and microdilution method. *Int. J. Antimicrob. Agents*. **37**, 271-273, 2011.

EUCAST-and CLSI potency NEO-SENSITABS™

Page 1 / 3 - Document: 3.16.0

Interpretation Zones and MIC Breakpoints according to CLSI

Quality control and Control Limits on Mueller-Hinton Agar for Nonfastidious Organisms

The acceptable limits for quality control strains, using agars and inoculum according to CLSI (Kirby-Bauer) are seen in the tables below (1,2).

Inoculum according to CLSI (Kirby-Bauer)

			Zone diameter in mm				
NEO-SENSITABS		CODE	<i>E. coli</i> ATCC 25922	<i>S. aureus</i> ATCC 25923	<i>Ps.</i> <i>aeruginosa</i> ATCC 27853	<i>E. coli</i> ATCC 35218	<i>E. faecalis</i> ATCC 29212
AMIKACIN	30 µg	AMI30	19-26	20-26	18-26	-	
AMOXYCILLIN	25 µg	AMX25	21-26	-	-	-	
AMOXYCILLIN+ CLAV.	20+10 µg	AMC30	18-24	28-36	-	17-22	
AMPICILLIN	10 µg	AMP10	16-22	27-35	-	9	
AMPICILLIN+ SULBACTAM	10+10 µg	SAM20	19-24	29-37		13-19	
AZITHROMYCIN	15 µg	AZI15	-	21-26	-	-	
AZTREONAM	30 µg	AZT30	28-36	-	23-29	-	
CEFTOBIPROLE	30 µg	CFBIP	30-36	26-34	23-31	-	
CEFAZOLIN	30 µg	CFZ30	21-27	29-35	-	-	
CEFEPIME	30 µg	FEP30	31-37	23-29	24-30	-	
CEFIXIME	5 µg	CFM.5	23-27	-	-	-	
CEFOTAXIME	30 µg	CTX30	29-35	25-31	18-22	-	
CEFOTAXIME	5 µg	CTX5	25-31	-	-	-	
CEFOXITIN	30 µg	CFO30	23-29	23-29	-	-	
CEFPODOXIME	10 µg	CPD10	23-28	19-25	-	-	
CEFTAROLINE	30 µg	CPT30	26-34	26-35	-	-	
CEFTAZIDIME	30 µg	CAZ30	25-32	16-20	22-29	-	
CEFTAZIDIME	10 µg	CAZ10	23-29	-	21-27	-	
CEFTIZOXIME	30 µg	ZOX30	30-36	27-35	12-17	-	
CEFTRIAXONE	30 µg	CTR30	29-35	22-28	17-23	-	
CEFUROXIME	30 µg	CXM30	20-26	27-35	-	-	
CEPHALOTHIN	30 µg	CEP30	15-21	29-37	-	-	
CHLORAMPHENICOL	30 µg	CLR30	21-27	19-26	-	-	
CIPROFLOXACIN	5 µg	CIPR5	30-40	22-30	25-33	-	
CLARITHROMYCIN	15 µg	CLA15	-	26-32	-	-	
CLINDAMYCIN	2 µg	CLI.2	-	24-30	-	-	
COLISTIN	10 µg	CO.10	-	-	-	-	
2+18 hours prediffusion			22-28	-	17-23	-	
DALBAVANCIN	60 µg	DAL60	-				
4+20 hours prediffusion			-	20-26	-	-	19-25
DAPTOMYCIN	30 µg	DAPCa	-	-	-	-	
2+18 hours prediffusion			-	22-28	-	-	
DORIPENEM	10 µg	DOR10	27-35	33-42	28-35	-	
DOXYCYCLINE	30 µg	DOX30	18-24	23-29	-	-	
ERTAPENEM	10 µg	ETP10	29-36	24-31	13-21	-	
ERYTHROMYCIN	15 µg	ERY15	-	22-30	-	-	
FOSFOMYCIN	200 µg	FO200	22-30	25-33	-	-	
FUCIDIN	10 µg	FUC10	-	24-32	-	-	
GATIFLOXACIN	5 µg	GATIF	30-37	27-33	20-28	-	

EUCAST-and CLSI potency NEO-SENSITABS™

Page 2 / 3 - Document: 3.16.0

Interpretation Zones and MIC Breakpoints according to CLSI

Quality control and Control Limits on Mueller-Hinton Agar for Nonfastidious Organisms

			Zone diameter in mm				
NEO-SENSITABS		CODE	<i>E. coli</i> ATCC 25922	<i>S. aureus</i> ATCC 25923	<i>Ps. aeruginosa</i> ATCC 27853	<i>E. coli</i> ATCC 35218	<i>E. faecalis</i> ATCC 29212
GENTAMICIN	10 µg	GEN10	19-26	19-27	17-23	-	
GENTAMICIN	250 µg		-	-	-	-	
ICLAPRIM	5 µg		14-22	25-33	-	-	
IMIPENEM	10 µg	IMI10	26-32	-	20-28	-	
KANAMYCIN	30 µg	KAN30	17-25	19-26	-	-	
LEVOFLOXACIN	5 µg	LEVOF	29-37	25-30	19-26	-	
LINEZOLID	30 µg	LINEZ	-	25-32	-	-	
MECILLINAM	10 µg	MEC10	24-30	-	-	-	
MEROPENEM	10 µg	MRP10	28-34	29-37	27-33	-	
MINOCYCLINE	30 µg	MIN30	19-25	25-30	-	-	
MOXIFLOXACIN	5 µg	MOXIF	28-35	28-35	17-25	-	
MUPIROCIN	10 µg	MUPIR	-	21-26	-	-	
MUPIROCIN	200 µg	MP200	-	29-38	-	-	
NALIDIXAN	30 µg	NAL30	22-28	-	-	-	
NETILMICIN	30 µg	NET30	22-30	22-31	17-23	-	
NITROFURANTOIN	300 µg	NI300	20-25	18-22	-	-	
NORFLOXACIN	10 µg	NORFX	28-35	17-28	22-29	-	
NOVOBIOCIN	5 µg	NOV.5	-	18-25	-	-	
OFLOXACIN	5 µg	OFL.5	29-33	24-28	17-21	-	
OXACILLIN	1 µg	OXA.1	-	18-24	-	-	
PIPERACILLIN+TAZOBACTAM	30+6 µg	PTZ36	21-27	-	23-29	-	
PENICILLIN	10U	PEN10	-	26-37	-	-	
PIPERACILLIN	100 µg	PIPA	24-30	-	25-33	12-18	
PIPERACILLIN+ TAZOBACTAM	100+10 µg	PI+TZ	24-30	27-36	25-33	24-30	
POLYMYXINS	150 µg	CO150	19-24	-	20-25	-	
QUINUPRISTIN/ DALFOPRISTIN	15 µg	SYN15	-	21-28	-	-	
RIFAMPICIN	5 µg	RIF.5	-	26-34	-	-	
RETAPAMULIN	2 µg	RETA2	-	23-30	-	-	
STREPTOMYCIN	10 µg	STR10	12-20	14-22	-	-	
SULPHONAMIDES	240 µg	SULFA	18-25	23-33	-	-	
TEICoplanin	30 µg	TPN30	-	15-21	-	-	
2+18 hours prediffusion			-	22-28	-	-	
TELITHROMYCIN	15 µg	TEL15	-	24-30	-	-	
TEMOCILLIN	30 µg	TEMOC	20-26	-	-	22-28	
TETRACYCLINES	30 µg	TET30	18-25	24-30	-	-	
TICARCILLIN	75 µg	TIC75	24-30	-	21-27	9	
TICARCILLIN+ CLAVULANATE	75+10 µg	TIM85	24-30	29-37	20-28	21-25	
TIGECYCLINE	15 µg	TIG15	20-27	20-25	9-13	-	
TOBRAMYCIN	10 µg	TOB10	18-26	19-29	20-26	-	

EUCAST-and CLSI potency NEO-SENSITABS™

Page 3 / 3 - Document: 3.16.0

Interpretation Zones and MIC
Breakpoints according to CLSI

Quality control and Control Limits on
Mueller-Hinton Agar for Nonfastidious
Organisms

Zone diameter in mm

NEO-SENSITABS		CODE	<i>E. coli</i> ATCC 25922	<i>S. aureus</i> ATCC 25923	<i>Ps.</i> <i>aeruginosa</i> ATCC 27853	<i>E. coli</i> ATCC 35218	<i>E. faecalis</i> ATCC 29212
TRIMETHOPRIM	5 µg	TRIM5	21-28	19-26	-	-	
TRIMETHOPRIM+ SULFA	1.25+23.75 µg	TR+SU	23-29	24-32	-	-	
VANCOMYCIN	30 µg	VAN30	-	17-21	-	-	
2+18 hours prediffusion			-	20-27	-	-	

K.pneumoniae
ATCC 700603

AMOXYCILLIN+CLAV	20+10 µg	AMC30	15-19				
AZTREONAM	30 µg	AZT30	9-17				
CEFOTAXIME	30 µg	CTX30	17-25				
CEFPODOXIME	10 µg	CPD10	9-16				
CEFTAZIDIME	30 µg	CAZ30	10-18				
CEFTRIAXONE	30 µg	CTR30	16-24				

EUCAST-and CLSI potency NEO-SENSITABS™

Page 1 / 2 - Document: 3.16.1

Interpretation Zones and MIC Breakpoints according to CLSI

Control Limits on Mueller-Hinton with Blood and ± Supplements for Fastidious Organisms

Inoculum according to CLSI (Kirby-Bauer)

			Zone diameter in mm			
NEO-SENSITABS		CODE	<i>Haemoph. influenzae</i> ATCC 49247	<i>Haemoph. influenzae</i> ATCC 49766	<i>Neisseria gonorrhoeae</i> ATCC 49226	<i>Streptoc. pneumoniae</i> ATCC 49619
AMOXYCILLIN	30 µg	AMOX				36-42
AMOXYCILLIN+	20+10 µg	AMC30	15-23	-	-	-
CLAVULANATE						
AMPICILLIN	10 µg	AMP10	13-21	-	-	30-36
AMPICILLIN	2.5 µg	AMP.L	9-14	-	-	25-31
AMPICILLIN+	10+10 µg	SAM20	14-22	-	-	-
SULBACTAM						
AZITHROMYCIN	15 µg	AZI15	13-21	-	27-36	19-25
AZTREONAM	30 µg	AZT30	30-38	-	-	-
CEFTOBIPROLE	30 µg	CFBIP	28-36	-	-	32-39
CEFEPIME	30 µg	FEP30	25-31	-	37-46	31-37
CEFIXIME	5 µg	CFM.5	25-33	-	37-45	16-23
CEFOTAXIME	30 µg	CTX30	31-39	-	38-48	31-39
CEFOTAXIME	5 µg	CTX5				28-34
CEFOXITIN	30 µg	CFO30	-	-	33-41	-
CEFPODOXIME	10 µg	CPD10	25-31	-	35-43	28-34
CEFTAROLINE	30 µg	CPT30	29-39	-	-	31-41
CEFTAZIDIME	30 µg	CAZ30	27-35	-	35-43	-
CEFTIZOXIME	30 µg	ZOX30	29-39	-	42-51	28-34
CEFTRIAZONE	30 µg	CTR30	31-39	-	39-51	30-35
CEFUROXIME	30 µg	CXM30	-	28-36	33-41	28-34
CEPHALOTHIN	30 µg	CEP30	-	-	-	26-32
CHLORAMPHENICOL	30 µg	CLR30	31-40	-	-	23-27
CIPROFLOXACIN	5 µg	CIPR5	34-42	-	48-58	-
CIPROFLOXACIN	1 µg	CIPR1	--	-	38-48	-
CLARITHROMYCIN	15 µg	CLA15	11-17	-	-	25-31
CLINDAMYCIN	2 µg	CLI.2	-	-	-	19-25
DOXYCYCLINE	30 µg	DOX30	-	-	-	25-34
DORIPENEM	10 µg	DOR10	21-31	-	-	30-38
ERTAPENEM	10 µg	ETP10	20-28	27-33	-	28-35
ERYTHROMYCIN	15 µg	ERY15				25-30
GATIFLOXACIN	5 µg	GATIF	33-41	-	45-56	24-31
IMIPENEM	10 µg	IMI10	21-29	-	-	34-42
LEVOFLOXACIN	5 µg	LEVOF	32-40	-	-	20-25
LINEZOLID	30 µg	LINEZ	-	-	-	25-34
MEROPENEM	10 µg	MRP10	20-28	-	-	28-35
MOXIFLOXACIN	5 µg	MOXIF	31-39	-	-	25-31
NITROFURANTOIN	300 µg	NI300	-	-	-	23-29
NORFLOXACIN	10 µg	NORFX	-	-	-	15-21
OFLOXACIN	5 µg	OFL.5	31-40	-	43-51	16-21
OXACILLIN	1 µg	OXA.1	-	-	-	9-14
PIPERACILLIN+TAZOBACTAM	30+6 µg	PTZ36	-	-	-	26-32
PENICILLIN	10U	PEN10	-	-	26-34	24-30
PIPERACILLIN+	100+10 µg	PI+TZ	33-38	-	-	-
TAZOBACTAM						

EUCAST-and CLSI potency NEO-SENSITABS™

Page 2 / 2 - Document: 3.16.1

Interpretation Zones and MIC Breakpoints according to CLSI

Control Limits on Mueller-Hinton with Blood and ± Supplements for Fastidious Organisms

NEO-SENSITABS			Zone diameter in mm			
			Haemoph. influenzae ATCC 49247	Haemoph. influenzae ATCC 49766	Neisseria gonorrhoeae ATCC 49226	Streptoc. pneumoniae ATCC 49619
QUINUPRISTIN/ DALFOPRISTIN	15 µg	SYN15	15-21	-	-	19-24
RETAPAMULIN	2 µg	RETA2	-	-	-	13-19
RIFAMPICIN	30 µg	RIF.5	22-30	-	-	25-30
SPECTINOMYCIN	200 µg	SPECT	-	-	26-35	-
TELITHROMYCIN	15 µg	TEL15	17-23	-	-	27-33
TETRACYCLINES	30 µg	TET30	14-22	-	30-42	27-31
TIGECYCLINE	15 µg	TIG15	23-31	-	30-40	23-29
TRIMETHOPRIM+ SULFA	1.25+23.75 µg	TR+SU	24-32	-	-	20-28
VANCOMYCIN	30 µg	VAN30	-	-	-	20-27

Q.C. ranges for antimicrobial agents against QC strains for disk diffusion testing are set using data from predefined structured multi-laboratory studies (1). It involves at least 7 independent laboratories, testing on 3 separate lots of medium from 2 different manufactures at least 30 times, using 2 separate disk lots from 2 manufacturers.

REFERENCES

- 1) Turnidge J. et al.: Statistical Methods for Establishing Q.C. Ranges for Antimicrobial Agents CLSI Susceptibility Testing. Antimicrob. Ag. Chemother., 51, 2483-2488, 2007.

Zone diameter interpretative criteria ROSCO and breakpoints when testing *Corynebacteria* are listed in the table below. No zone breakpoints are available from CLSI.

Table 3.17-1 Interpretation for *Corynebacteria*

Mueller-Hinton + 5 % blood+20 mg/l β-NAD. Inoculum: McFarland 0.5. Incubation for 24-48 hours at 35 °C ± 2 degrees.

				Zone diameter in mm			Break-points MIC µg/ml	
NEO-SENSITABS		POTENCY	CODE	S	I	R	S	R
a)	Cefepime	30 µg	FEP30	≥ 26	25-22	< 22	≤ 1	> 4
a)	Cefotaxime	30 µg	CTX30	≥ 26	25-22	< 22	≤ 1	> 4
a)	Ceftriaxone	30 µg	CTR30	≥ 26	25-22	< 22	≤ 1	> 4
	Ciprofloxacin	5 µg	CIPR5	≥ 26	25-19	≤ 18	≤ 1	≥ 4
	Clindamycin	2 µg	CLI.2	≥ 21	20-15	≤ 14	≤ 0.5	≥ 4
	Daptomycin	30 µg	DAPCa					
	2+18 hours' prediffusion			≥ 22	-	-	≤ 1	-
	Doxycycline	30 µg	DOX30	≥ 20	19-17	≤ 16	≤ 4	≥ 16
	Erythromycin	15 µg	ERY15	≥ 23	22-14	≤ 13	≤ 0.5	≥ 8
	Gentamicin	10 µg	GEN10	≥ 18	17-15	≤ 14	≤ 4	≥ 16
a) b)	Imipenem	10 µg	IMI10	≥ 26	25-23	< 22	≤ 0.25	> 1
	Linezolid	30 µg	LINEZ	≥ 23	22-21	≤ 20	≤ 2	≥ 8
a) b)	Meropenem	10 µg	MRP10	≥ 26	25-23	≤ 22	≤ 0.25	> 1
a)	Penicillin	10 U	PEN10	≥ 22	21-19	≤ 18	≤ 1	≥ 4
	Quinupristin-Dalfopristin	15 µg	SYN15	≥ 19	18-16	≤ 15	≤ 1	≥ 4
	Rifampicin	5 µg	RIF.5	≥ 26	25-21	≤ 20	≤ 1	≥ 4
	Teicoplanin	30 µg	TPN30	≥ 16	-	-	≤ 4	-
	Tetracycline	30 µg	TET30	≥ 24	23-19	≤ 18	≤ 4	≥ 16
	Trimethoprim	1.25 µg	SxT25	≥ 16	15-11	≤ 10	≤ 2/38	≥ 8/152
	+Sulfa	+23.75 µg						
	Vancomycin	30 µg	VAN30	≥ 17	-	-	≤ 2	-

a) Interpretative criteria may not apply to meningitis

b) For detection of possible resistance, ROSCO MIC breakpoints ($S \leq 0.25$ ug/ml, $R > 1$ ug/ml) are recommended for Imipenem and Meropenem. M45-A2 recommends $S \leq 4$ ug/ml and $R \geq 16$ ug/ml.

MIC breakpoints for carbapenems have been placed at $S \leq 0.25$ ug/ml in order to be able to detect isolates with increased resistance to carbapenems.

Testing of isolates from normally sterile sources (blood cultures, deep tissue, prosthetic devices) may be warranted, especially in immunodeficient patients.

Cephalosporins MIC interpretative criteria are adopted from those of *Streptococcus spp.*, Linezolid adopted from *Enterococcus spp.*; remaining from those of *Staphylococcus spp.* (CLSI, M45-A2)

Resistant results can be reported after 24 hours. Beta-lactam susceptible results should be reincubated and reported after 48 hours.

Penicillin results are predictive of susceptibility to cephalosporins and carbapenems.

Zone diameters tentative for 1 year.

***For further information
see:
[http://www.eucast.org/
eucast_disk_diffusion_
test/breakpoints/](http://www.eucast.org/eucast_disk_diffusion_test/breakpoints/)***

***For further information
see:***

***[http://www.eucast.org/
eucast_disk_diffusion_
test/breakpoints/](http://www.eucast.org/eucast_disk_diffusion_test/breakpoints/)***

***For further information
see:
[http://www.eucast.org/
eucast_disk_diffusion_
test/breakpoints/](http://www.eucast.org/eucast_disk_diffusion_test/breakpoints/)***

***For further information
see:
[http://www.eucast.org/
eucast_disk_diffusion_
test/breakpoints/](http://www.eucast.org/eucast_disk_diffusion_test/breakpoints/)***

***For further information
see:
[http://www.eucast.org/
eucast_disk_diffusion_
test/breakpoints/](http://www.eucast.org/eucast_disk_diffusion_test/breakpoints/)***

***For further
information see:
[http://www.eucast.org
/eucast_disk_diffusio
n_test/breakpoints/](http://www.eucast.org/eucast_disk_diffusion_test/breakpoints/)***

***For further information
see:
[http://www.eucast.org/
eucast_disk_diffusion_t
est/breakpoints/](http://www.eucast.org/eucast_disk_diffusion_test/breakpoints/)***

***For further information
see:
[http://www.eucast.org/
eucast_disk_diffusion_
test/breakpoints/](http://www.eucast.org/eucast_disk_diffusion_test/breakpoints/)***

***For further
information see:
[http://www.eucast.org/
eucast_disk_diffusion
_test/breakpoints/](http://www.eucast.org/eucast_disk_diffusion_test/breakpoints/)***

***For further information
see:
[http://www.eucast.org/
eucast_disk_diffusion
_test/breakpoints/](http://www.eucast.org/eucast_disk_diffusion_test/breakpoints/)***

***For further information
see:
[http://www.eucast.org/
eucast_disk_diffusion_
test/breakpoints/](http://www.eucast.org/eucast_disk_diffusion_test/breakpoints/)***

***For further information
see:***

***[http://www.eucast.org/
eucast_disk_diffusion
_test/breakpoints/](http://www.eucast.org/eucast_disk_diffusion_test/breakpoints/)***

***For further
information see:
[http://www.eucast.org
/eucast_disk_diffusio
n_test/breakpoints/](http://www.eucast.org/eucast_disk_diffusion_test/breakpoints/)***

Zone diameter interpretative standards are correlated to the MIC breakpoints recommended by EUCAST for anaerobic bacteria (1), but zone diameter interpretative standards have not yet been established by EUCAST. The table will be updated with when these are established.

Table 3.14-1 Interpretation for *Anaerobes*

Supplemented Brucella Blood agar. Inoculum: McFarland 1.0. Incubation at 35 °C in anaerobic environment for 24-48 hours.

NEO-SENSITABS				Zone diameter in mm			Break-points MIC µg/ml	
				S	I	R	S	R
Amoxycillin+Clav.	20+10 µg	AMC30		≥ 24	23-21	< 21	≤ 4	≥ 8
Chloramphenicol	30 µg	CLR30		≥ 23	22-21	< 21	≤ 8	> 8
Clindamycin	2 µg	CLI.2		≥ 13	-	< 13	≤ 4	> 4
Doripenem	10 µg	DOR10		≥ 26	-	< 26	≤ 1	> 1
Ertapenem	10 µg	ERTAP		≥ 26	-	< 26	≤ 1	> 1
Ertapenem Screen MBL						< 18		possible
Imipenem	10 µg	IMI10		≥ 24	23-19	< 19	≤ 2	≥ 8
Meropenem	10 µg	MRP10		≥ 24	23-19	< 19	≤ 2	≥ 8
Metronidazole	16 µg	MTR16		≥ 24	23-21	< 21	≤ 4	> 4
Metronidazole	5 µg	MTR.5		≥ 20	-	< 20	≤ 4	> 4
Penicillin (Ampicillin)	1 U	PENG1		≥ 21	20-18	< 18	≤ 0.25	≥ 0.5
Piperacillin+Tazobactam	30+6µg	PIZ36		≥ 20	19-17	< 17	≤ 8	≥ 16
Ticarcillin+Clavulanate	75+10 µg	TIM85		≥ 26	25-22	< 22	≤ 8	≥ 16
Vancomycin	5 µg	VAN.5		≥ 15	-	< 15	≤ 2	> 2

Zone diameter interpretative standards and MIC breakpoints recommended by EUCAST for *Clostridium difficile* (1)

Table 3.14-2 Interpretation for *Clostridium difficile*

Supplemented Brucella Blood agar. Inoculum: McFarland 1.0. Incubation at 35 °C in anaerobic environment for 24-48 hours.

NEO-SENSITABS				Zone diameter in mm			Break-points MIC µg/ml	
				S	I	R	S	R
b)	Metronidazole	16 µg	MTR16	≥ 26	-	< 26	≤ 2	> 2
	Metronidazole	5 µg	MTR.5	≥ 23	-	< 23	≤ 2	> 2
	Moxifloxacin	5 µg	MOXIF	≥ 20	-	< 20	≤ 4	> 4
	Vancomycin	5 µg	VAN.5	≥ 19	-	< 19	≤ 2	> 2

References:

- 1) EUCAST January 2013.

EUCAST-and CLSI potency NEO-SENSITABS™

Interpretation according to MIC Breakpoints of SFM (France) Using CLSI Potency Neo-Sensitabs

Rapidly growing bacteria

Interpretation zones and MIC breakpoints adapted to recommendations by the "Comité de l'Antibiogramme de la Société Française de Microbiologie" (2012)

Inoculum, media and incubation conditions acc. to SFM (2012)

				Zone diameter in mm			Concentrations critiques	
NEO-SENSITABS	POTENCY	CODE		S	I	R	S	R
	Amikacin	30 µg	AMI30	≥ 17	16-15	< 15	≤ 8	> 16
c)	Amoxycillin+Clav.	20+10 µg	AMC30	≥ 21	20-16	< 16	≤ 4/2	> 8/2
a)	Ampicillin	10 µg	AMP10					
	Enterobacteriaceae, <i>E. faecalis</i>			≥ 19	18-16	< 16	≤ 4	> 8
b)	Ampicillin+Sulbactam	10+10 µg	SAM20	≥ 19	18-16	< 16	≤ 4/8	> 8/8
	Azithromycin	15 µg	AZI15	≥ 22	21-17	< 17	≤ 0.5	> 4
	<i>Salmonella/Shigella</i>						≤ 16	
o) g)	Aztreonam	30 µg	AZT30	≥ 27	26-21	< 21	≤ 1	> 8
	<i>Ps. aeruginosa</i>			≥ 27	26-19	< 19	≤ 1	> 16
o) g)	Cefepime	30 µg*	FEP30	≥ 24	23-21	< 21	≤ 1	> 4
	<i>Ps. aeruginosa</i>			≥ 19	-	< 19	≤ 8	> 8
	<i>Acinetobacter</i> spp.			≥ 21	20-19	< 19	≤ 4	> 8
	Cefepime+Clavulanate	30+10µg**	CP+CL	Detection of ESBL				
g)	Cefixime	5 µg	CFM.5	≥ 22	21-19	< 19	≤ 1	> 2
g)	Cefotaxime	30 µg	CTX30	≥ 26	25-23	< 23	≤ 1	> 2
	Cefotaxime+Clav.	30+10µg**	CZ+CL	Detection of ESBL				
	Cefoxitin	30 µg	CFO30	≥ 22	21-15	< 15	≤ 8	> 32
	<i>Staphylococcus</i> spp			≥ 27	-	< 25	Oxa S	Mec A pos.
	Ceftobiprole	30 µg		≥ 24	-	-	≤ 2	-
g) m)	Cefpodoxime	10 µg	CPD10	≥ 24	23-21	< 21	≤ 1	> 2
	(Screen ESBL)							
o) g)	Ceftazidime	30 µg	CAZ30	≥ 26	25-21	< 21	≤ 1	> 4
	<i>Ps. aeruginosa</i>			≥ 19	-	< 19	≤ 8	> 8
	<i>Acinetobacter</i> spp.			≥ 21	20-19	< 19	≤ 4	> 8
	Ceftazidime+Clav.	30+10µg**	CZ+CL	Detection of ESBL				
g)	Ceftizoxime	30 µg	ZOX30	≥ 26	25-23	< 23	≤ 1	> 2
	Ceftriaxone	30 µg	CTR30	≥ 26	25-23	< 23	≤ 1	> 2
	Cefuroxime (parenteral)	30 µg	CXM30	≥ 22	-	< 22	≤ 8	> 8
	Cefuroxime (oral)	30 µg*	CXM30	≥ 30	29-23	< 23	≤ 1	> 4
	Cephalothin	30 µg	CEP30	≥ 18	17-12	< 12	≤ 8	> 32
	Chloramphenicol	30 µg	CLR30	≥ 23	-	< 23	≤ 8	> 8
	Ciprofloxacin	5 µg	CIPR5	≥ 25	24-22	< 22	≤ 0.5	> 1
	<i>P. aeruginosa</i>			≥ 25	24-22	< 22	≤ 0.5	> 1
	<i>Staphylococcus</i> spp; <i>Acinetobacter</i> spp.			≥ 22	-	< 22	≤ 1	> 1
	Clindamycin	2 µg*	CLI.2	≥ 23	22-20	< 20	≤ 0.25	> 0.5
u)	Cloxacillin		CLOXA	Detection of plasmid med. Amp C				
t)	Colistin	10 µg**	CO.10					
t) x)	2 + 18 h. pre-diffusion			≥ 15	14-11	≤ 10	≤ 2	> 2
	<i>Ps. aeruginosa</i> (prediffusion)			≥ 15	-	< 15	≤ 4	> 4
	<i>Acinetobacter</i> (prediffusion)			≥ 20	-	-	≤ 2	-
	Doxycycline	30 µg*	DOX30	≥ 19	18-17	< 17	≤ 4	> 8
x)	Daptomycin	30 µg	DAPCa					

EUCAST-and CLSI potency NEO-SENSITABS™

Page 2 / 5 - Document: 6.1.0

Interpretation according to MIC
Breakpoints of SFM (France) Using
CLSI Potency Neo-Sensitabs

Rapidly Growing Bacteria

NEO-SENSITABS			POTENCY	CODE	Zone diameter in mm			Concentrations critiques		
					S	I	R	S	R	
t)	2 + 18 h. prediffusion				≥ 22	-	< 20	≤ 1	> 1	
	<i>Staphylococcus</i> spp.				≥ 12	-	≤ 11	≤ 4	> 8	
	<i>Enterococcus</i> spp.				≥ 24	23-19	< 19	≤ 1	> 4	
	Doripenem	10 µg	DOR10	≥ 28	27-26	< 26	≤ 0.5	> 1		
	Ertapenem	10 µg	ETP10	≥ 22	21-19	< 19	≤ 1	> 2		
r)	Fosfomycin	200 µg [#]	FO200	≥ 18	-	< 18	≤ 32	> 32		
	Fucidin	10 µg [#]	FUC10	≥ 24	-	< 24	≤ 1	> 1		
n)	Gatifloxacin	5 µg	GATIF	≥ 21	20-18	< 18	≤ 1	≥ 2		
	Gentamicin	10 µg [#]	GEN10	≥ 18	17-16	< 16	≤ 2	> 4		
	<i>Staphylococcus</i> spp.				≥ 20	-	< 20	≤ 1	> 1	
	<i>P. aeruginosa</i> , <i>Acinetobacter</i> spp.				≥ 16	-	< 16	≤ 4	> 4	
h)	Gentamicin (HNR)	250 µg ^{**}	GN250	≥ 14	-	< 12	≤ 250	> 500		
	Imipenem	10 µg	IMI10	≥ 24	23-17	< 17	≤ 2	> 8		
	Pseudomonas			≥ 22	21-17	< 17	≤ 4	> 8		
	Imipenem+EDTA	10+750 µg ^{**}	IM10E	Detection of metallo-β-lactamases						
q)	Kanamycin	30 µg [#]	KAN30	≥ 17	16-15	< 15	≤ 8	> 16		
h)	Kanamycin (HNR)	500 µg ^{**}	KA500	≥ 14	-	< 10	≤ 250	> 500		
	Levofloxacin	5 µg	LEVOF	≥ 20	19-17	< 17	≤ 1	> 2		
	<i>Acinetobacter</i> spp; <i>P. aeruginosa</i>				≥ 20	19-17	< 17	≤ 1	> 2	
	<i>Staphylococcus</i> spp., <i>Enterococcus</i> spp.				≥ 20	19-17	< 17	≤ 1	> 2	
	Linezolid	30 µg	LINEZ	≥ 24	-	< 24	≤ 4	> 4		
	<i>Staphylococcus</i> spp., <i>Enterococcus</i> spp.				≥ 24	-	< 24	≤ 4	> 4	
	Moxifloxacin	5 µg	MOXIF	≥ 24	23-21	< 21	≤ 0.5	> 1		
	<i>Acinetobacter</i> spp., <i>Enterococcus</i> spp.				≥ 21	20-18	< 18	≤ 1	> 2	
	Mecillinam (U)	10 µg	MEC10	≥ 24	23-22	< 22	≤ 8	> 8		
	Meropenem	10 µg	MRP10	≥ 22	21-15	< 15	≤ 2	> 8		
	Minocycline	30 µg [#]	MIN30	≥ 19	18-17	< 17	≤ 4	> 8		
	<i>Staphylococci</i>				≥ 23	22-21	< 21	≤ 0.5	> 1	
	Mupirocin (staph)	10 µg ^{**}	MUPIR	≥ 19	18-10 (LLR)	no zone	≤ 2	> 256 (HLR)		
	p)	Nalidixan (U)	30 µg	NAL30	≥ 20	19-15	< 15	≤ 8	> 16	
		Enterobacteriaceae				-	-	< 19	Reduced susceptibility to quinolones	
		Netilmicin	30 µg	NET30	≥ 21	20-19	< 19	≤ 2	> 4	
<i>P. aeruginosa</i> , <i>Acinetobacter</i> spp.				≥ 19	-	< 19	≤ 4	> 4		
Nitrofurantoin (U)		300 µg	NI300	≥ 15	-	< 15	≤ 64	> 64		
Norfloxacin (U)		10 µg [#]	NORFX	≥ 25	24-22	< 22	≤ 0.5	> 1		
Ofloxacin		5 µg	OFL.5	≥ 25	24-22	< 22	≤ 0.5	> 1		
f)	Oxacillin 1 µg	1 µg [#]	OXA.1							
	<i>Staph. aureus</i>				≥ 16	-	< 16	≤ 2	> 2	
	Coag. neg. staph.				≥ 20	-	< 20	≤ 0.25	> 0.5	

EUCAST-and CLSI potency NEO-SENSITABS™

Interpretation according to MIC Breakpoints of SFM (France) Using CLSI Potency Neo-Sensitabs

Rapidly Growing Bacteria

			Zone diameter in mm			Concentrations critiques	
NEO-SENSITABS	POTENCY	CODE	S	I	R	S	R
Penicillin	10 U	PEN10	≥ 29	-	< 29	≤ 0.12	beta-lactamase
<i>Staphylococcus</i> spp.							
Piperacillin	100 µg [#]	PIPRA					
a) Enterobacteriaceae			≥ 20	19-16	< 16	≤ 8	> 16
<i>P. aeruginosa</i> , <i>Acinetobacter</i> spp.			≥ 18	-	< 18	≤ 16	> 16
Piperacillin+Tazobactam	100+10µg [#] PI+TZ						
e) Enterobacteriaceae			≥ 21	20-17	< 17	≤ 8/4	> 16/4
e) <i>P. aeruginosa</i> , <i>Acinetobacter</i> spp.			≥ 19	-	< 19	≤ 16/4	> 16/4
Quinu/Dalfopristin	15 µg	SYN15	≥ 22	21-19	< 19	≤ 1	> 2
Enterococcus			≥ 22	21-16	< 16	≤ 1	> 4
Rifampicin	5 µg	RIF.5					
<i>Staphylococcus</i> spp.			≥ 26	25-22	< 22	≤ 0.06	> 0.5
<i>Ps. aeruginosa</i> ; <i>Acinetobacter</i>			≥ 16	15-12	< 12	≤ 4	> 16
Streptomycin	10 µg [#]	STR10	≥ 15	14-13	< 13	≤ 8	> 16
h) Streptomycin (HNR)	500 µg ^{**}	ST500	≥ 14	-	< 12	≤ 250	> 500
Sulphonamides (U)	240 µg ^{**}	SULFA	≥ 17	16-12	< 12	≤ 64	> 256
Teicoplanin	30 µg	TPN30	≥ 17	-	-	≤ 2	> 8
k) s) 2 + 18 h. prediffusion			≥ 22	-	< 20	≤ 2	> 8
t) x) <i>Staphylococcus</i> spp.							(VISA/GISA)
(MH agar plain)							
Telithromycin	15 µg	TEL15	≥ 21	20-17	< 17	≤ 0.5	> 2
Tetracyclines	30 µg [#]	TET30	≥ 19	18-17	< 17	≤ 4	> 8
<i>Staphylococci</i>			≥ 23	22-21	< 21	≤ 1	> 2
a) Ticarcillin	75 µg	TIC75	≥ 24	23-22	< 22	≤ 8	> 16
<i>Pseudomonas</i>			≥ 22	-	< 22	≤ 16	> 16
d) o) Ticarcillin+Clavulanate	75+10 µg	TI+CL	≥ 24	23-22	< 22	≤ 8/2	> 16/2
<i>Pseudomonas</i>			≥ 22	-	< 22	≤ 16/2	> 16/2
Tigecycline (enterob)	15 µg	TIG15	≥ 21	20-19	< 19	≤ 1	> 2
<i>Staphylococcus</i>			≥ 22	-	< 22	≤ 0.5	> 0.5
Enterococcus			≥ 22	-	< 22	≤ 0.5	> 0.5
<i>Acinetobacter</i>			≥ 22	-	-	≤ 1	-
Tobramycin	10 µg	TOB10	≥ 18	17-16	< 16	≤ 2	> 4
<i>Staphylococcus</i> spp.			≥ 20	-	< 20	≤ 1	> 1
<i>P. aeruginosa</i> , <i>Acinetobacter</i> spp.			≥ 16	-	< 16	≤ 4	≥ 4
Trimethoprim (U)	5 µg	TRIM5	≥ 20	19-16	< 16	≤ 2	> 4
i) Trimethoprim+Sulfa	1.25+23.75µg	SxT25	≥ 16	15-13	< 13	≤ 2/38	> 4/76

EUCAST-and CLSI potency NEO-SENSITABS™

Page 4 / 5 - Document: 6.1.0

Interpretation according to MIC Breakpoints of SFM (France) Using CLSI Potency Neo-Sensitabs

Rapidly Growing Bacteria

NEO-SENSITABS				Zone diameter in mm			Concentrations critiques	
POTENCY		CODE	S	I	R	S	R	
k) s)	Vancomycin	30 µg	VAN30	≥ 17	-	-	≤ 2	> 8
t) x)	2 + 18 h. prediffusion <i>Staphylococcus</i> spp.(MH agar plain)			≥ 22	-	≤ 20	≤ 2	> 2
	<i>Enterococcus</i> spp.			≥ 16	-	< 16	-	(VISA/GISA), VRE
	Phenylboronic acid	BORON		Detection of AmpC and KPC beta lactamases				
	Dipicolinic acid	D.P.A		Detection of metallo beta lactamases				
	ESBL Confirm Kits		Confirmation of ESBL's					
	ESBL+ AmpC Screen Kit		Sreening of ESBL and AmpC					
	Total ESBL + AmpC Confirm Kit		Confirmation of ESBL, AmpC and ESBL+AmpC					
	AmpC Confirm Kit		Confirmation of AmpC					
	KPC, MBL Confirm Kit		Confirmation of KPC, MBL					
	KPC, MBL, OXA-48 Confirm Kit		Confirmation of KPC, MBL and OXA-48					

Rapidly growing bacteria includes Enterobacteriaceae, *Pseudomonas* spp., *Acinetobacter* spp., *Staphylococcus* spp., *Enterococcus* spp.

#) Differences of potency recommended by SFM and CLSI (if any available).

*) There are no potency recommendations from CLSI so far.

**) Special potency Neo-Sensitabs for detection of resistance mechanisms.

- Klebsiella* spp. produces a natural low level beta-lactamase that inactivates amino-, carboxy- and ureido-penicillins. They may appear susceptible in vitro, but they should be reported as Intermediate to carboxy- and ureido-penicillins.
- Critical concentrations of Ampicillin with a fixed concentration of Sulbactam (8 µg/ml).
- Critical concentrations of Amoxycillin with a fixed concentration of Clavulanic acid (2 µg/ml).
- Critical concentrations of Ticarcillin with a fixed concentration of Clavulanic acid (2 µg/ml).
- Critical concentrations of Piperacillin with a fixed concentration of Tazobactam (4 µg/ml).
- For detecting Methicillin/Oxacillin resistance in staphylococci follow the instructions enclosed with NeoSensitabs™. Test Cefoxitin Neo-Sensitabs and Oxacillin 1 µg. Strains resistant to Oxacillin should be reported as resistant to all beta-lactams, even if they appear susceptible in vitro.
- Strains of *Klebsiella* spp. and *E. coli* may be clinically resistant to cephalosporins and Aztreonam therapy by producing ESBL (extended spectrum beta-lactamase). Read more in "Detection of resistant mechanisms using Neo-Sensitabs™ and Diatabs™.
- These high content Neo-Sensitabs are used to detect high level resistance to the aminoglycosides.

- i) The interpretation of results is valid for other combinations of Trimethoprim+ Sulphonamide.
- k) Strains showing inhibition zones smaller than the limit for susceptible should be tested by an MIC method.
- l) Used for the detection of metallo- β -lactamases in gram-negatives. See further information in user guide "Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™".
- m) Strains showing zone < 21 mm with Cefpodoxime Neo-Sensitabs, should be suspected of producing ESBL (*E. coli*, *Klebsiella*, *Salmonella*).

For ESBL confirmatory tests see user guide "Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™" ESBL Confirm Kits.

- n) Results with Gentamicin and staphylococci are also valid for Netilmicin. Staphylococci that are resistant to Gentamicin, should be reported as resistant against both Netilmicin and Tobramycin (enzymes APH(2'') + AAC(6'')).
- o) Synergism between Ticarcillin+Clavulanate and Aztreonam/Ceftazidime/Cefepime permits the detection of ESBL producing strains in *Ps. aeruginosa*.
- p) Nalidixic acid is useful to detect strains with reduced susceptibility to quinolones in 1) Enterobacteriaceae (Nali zone < 19 mm) 2) *Haemophilus influenzae* (Nali zone < 21 mm) 3) Gonococci (Nali zone < 25 mm) and 4) *Vibrio cholerae* (Nali zone < 19 mm).
- q) With staphylococci, the interpretation is valid for Amikacin and Isepamicin.
- r) Do not take into account the presence of colonies inside the zone of inhibition. Resistant strains show homogeneous resistance.
- s) For detection of VISA, GISA, hVISA strains, see user guide "Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™", on detection of staphylococci with decreased susceptibility to Vancomycin. hVISA, GISA, VRE and Daptomycin Kit.
- t) Dipicolinic acid is used for detection of metallo- β -lactamases. For further information about methodology, see user guide "Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™". Total MBL Confirm Kit; KPC, MBL Confirm Kit.
- u) Cloxacillin and phenylboronic Acid are used for the detection of plasmid-mediated AmpC beta-lactamases. See description of the procedure in user guide "Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™" AmpC Confirm Kit.
- x) Special technique for susceptibility testing of high molecular weight antimicrobials (Colistin, Daptomycin, Teicoplanin, Vancomycin): 2+18 hours' prediffusion method permitting a good separation between susceptible and resistant strains. Description of the procedure can be found in document **1.5.0**.

EUCAST-and CLSI potency NEO-SENSITABS™

Page 1 / 5 - Document: 6.2.0

Interpretation according to MIC Breakpoints of SFM (France) Using CLSI Potency Neo-Sensitabs

**Haemophilus spp., S. pneumoniae,
Streptococcus spp., N. gonorrhoeae, N.
meningitidis, Campylobacter spp. and
Anaerobes**

Fastidious organisms

Interpretations zones and MIC breakpoints adapted to recommendations by the "Comité de l'Antibiogramme de la Société Française de Microbiologie" (2012)

Inoculum, media and incubation conditions acc. to SFM (2012)

				Zone diameter in mm			Concentrations critiques	
NEO-SENSITABS	POTENCY	CODE		S	I	R	S	R
e) Ampicillin	10 µg	AMP10						
			<i>Haemophilus</i> spp.	Use Amp 2.5 µg			-	> 1
			<i>Campylobacter</i> spp.	≥ 19	18-14	< 14	≤ 4	> 16
e) Ampicillin	2 µg**	AMP.2						
			<i>Haemophilus</i> spp.	≥ 20	-	< 20	≤ 1	> 1
	10 U	PEN10						
			<i>N. gonorrhoeae</i>	≥ 34	33-23	< 23	≤ 0.06	> 1
	1 µg#	OXA.1						
a) Oxacillin			<i>S. pneumoniae</i>	≥ 20	-	< 20	≤ 0.06 (pen)	> 1 MIC test
d) Streptococcus spp.				≥ 15	-	< 15	≤ 0.25 (pen)	MIC test
l) N. gonorrhoeae				≥ 14	-	< 14	≤ 0.06 (pen)	> 1 MIC test
g) N. meningitidis				≥ 12	-	< 12	≤ 0.06 (pen)	> 1 MIC test
Amoxycillin+Clav.	20+10 µg	AMC30						
			<i>Haemophilus</i> spp./Moraxella	≥ 24	-	< 24	≤ 1	> 1
			<i>Campylobacter</i> spp.	≥ 21	20-14	< 14	≤ 4/2	> 16/2
			Anaerobes	≥ 21	20-17	< 17	≤ 4/2	> 8/2
Ticarcillin+Clavulanate	75+15 µg	TI+CL						
			Anaerobes	≥ 24	23-22	< 22	≤ 8/2	> 16/2
Piperacillin+Tazobactam	100+10µg*PI+TZ							
			Anaerobes	≥ 21	20-19	< 19	≤ 8/4	> 16/4
Cephalothin	30 µg	CEP30						
e) Haemophilus spp.				-	-	< 17	-	> 8 (Amp R)
			<i>Campylobacter</i>	≥ 18	17-12	< 12	≤ 8	> 32
Ceftizoxime	30 µg**	ZOX30						
b) S. pneumoniae				≥ 28	-	< 28	≤ 0.5	> 0.5
			(valid for 3rd gen. cephalosporins)					
Cefotaxime	30 µg	CTX30						
/ Ceftriaxone	30 µg	CTR30						
b) S. pneumoniae				(use Ceftizoxime)			≤ 0.5	> 2
			<i>H. influenzae</i>	≥ 32	-	-	≤ 0.12	-
			<i>N. gonorrhoeae</i>	≥ 40	-	-	≤ 0.12	-
			<i>N. meningitidis</i>	≥ 38	-	-	≤ 0.12	-
			<i>Campylobacter</i> spp.	≥ 26	25-23	< 23	≤ 1	> 2
Cefoxitin	30 µg	CFO30						
			Anaerobes	-	-	< 19	-	> 32
Imipenem	10 µg	IMI10						
			<i>S. pneumoniae/streptococci</i>	≥ 24	-	< 24	≤ 2	> 2
			Anaerobes	≥ 24	23-17	< 17	≤ 2	> 8
Doripenem	10 µg	DOR10						
				≥ 23	-	< 23	≤ 1	> 1
Meropenem	10 µg	MRP10						
			<i>S. pneumoniae/streptococci</i>	≥ 24	-	< 24	≤ 2	> 2
			<i>N. meningitidis</i>	≥ 34	33-32	< 32	≤ 0.25	> 0.25
			Anaerobes	≥ 22	21-15	< 15	≤ 2	> 8
Ertapenem	10 µg	ERTAP						
			<i>S. pneumoniae/streptococci</i>	≥ 28	27-24	< 24	≤ 0.5	-
			Anaerobes	≥ 26	-	< 26	≤ 1	> 1
			<i>Haemophilus/Moraxella</i>	≥ 30	-	-	≤ 0.5	-
Azithromycin	15 µg	AZI15						

EUCAST-and CLSI potency NEO-SENSITABS™

Page 2 / 5 - Document: 6.2.0

Interpretation according to MIC
Breakpoints of SFM (France) Using
CLSI Potency Neo-Sensitabs

**Haemophilus spp., S. pneumoniae,
Streptococcus spp., N. gonorrhoeae, N.
meningitidis, Campylobacter spp. and
Anaerobes**

			Zone diameter in mm			Concentrations critiques	
NEO-SENSITABS	POTENCY	CODE	S	I	R	S	R
c)	Haemophilus		≥ 28	27-14	< 14	≤ 0.12	> 4
	N.gonorrhoeae		≥ 27	26-24	< 24	≤ 0.25	> 0.5
	Erythromycin	15 µg* ERY15					
	S. pneumoniae / Streptococcus spp.		≥ 26	25-24	< 24	≤ 0.25	> 0.5
	Campylobacter spp.		≥ 22	21-17	< 17	≤ 1	> 4
n)	Helicobacter pylori		≥ 22	21-17	< 17	≤ 1	> 4
	Clindamycin	2 µg* CLI.2					
	S. pneumoniae / Streptococcus spp.		≥ 19	-	< 19	≤ 0.5	> 0.5
	Anaerobes		≥ 15	-	< 15	≤ 4	> 4
	Linezolid	30 µg LINEZ					
	S. pneumoniae / Streptococcus spp.		≥ 24	-	< 24	≤ 2	> 4
	Anaerobes		≥ 24	-	< 24	≤ 2	> 4
	Telithromycin	15 µg TEL15					
	S. pneumoniae / Streptococcus spp.		≥ 24	23-21	< 21	≤ 0.25	> 0.5
	Tetracyclines	30 µg* TET30					
	Haemophilus spp.		≥ 23	22-20	< 20	≤ 1	> 2
	S. pneumoniae / Streptococcus spp.		≥ 23	22-21	< 21	≤ 1	> 2
	N. gonorrhoeae (tet M)		-	-	< 19	≤ 0.5	> 1
	Campylobacter spp.		≥ 19	18-17	< 17	≤ 4	> 8
	Minocycline	30 µg MIN30			< 19	≤ 0.5	> 1
	Chloramphenicol	30 µg CLR30					
	Haemophilus spp.		≥ 30	29-26	< 26	≤ 1	> 2
	S. pneumoniae / Streptococcus spp.		≥ 23	-	< 23	≤ 8	> 8
	N. gonorrhoeae		≥ 28	27-24	< 24	≤ 4	> 16
	N. meningitidis		≥ 30	-	-	≤ 2	> 4
h)	Campylobacter spp.		≥ 23	22-19	< 19	≤ 8	> 16
	Anaerobes		≥ 23	-	< 23	≤ 8	> 8
	Rifampicin	5 µg RIF.5					
	Haemophilus spp.		≥ 18	-	< 18	≤ 0.5	> 0.5
	Streptococcus spp;Pneumococci		≥ 26	25-22	< 22	≤ 0.06	> 0.5
c ₂)	N. meningitidis		≥ 24	-	-	≤ 0.25	-
	Anaerobes		≥ 14	13-10	< 10	≤ 4	> 16
	Gentamicin	10 µg GEN10					
	Campylobacter spp.		≥ 18	17-16	< 16	≤ 2	> 4
	Tobramycin	10 µg TOB10					
	Campylobacter spp.		≥ 18	17-16	< 16	≤ 2	> 4
	Norfloxacin	10 µg* NORFX					
	S. pneumoniae		-	-	< 12	Reduced susceptibility to quinolones	
	Ciprofloxacin	5 µg CIPR5					
	Haemophilus spp.		use Nalidixic acid			≤ 0.5	-
	N. gonorrhoeae		use Nalidixic acid			≤ 0.03	> 0.06
	Campylobacter spp.		≥ 25	24-22	< 22	≤ 0.5	> 1
	Helicobacter pylori		≥ 20	-	< 20	≤ 1	> 1
	S. pneumoniae		≥ 30	29-19	< 19	≤ 0.12	> 2
	N. meningitidis		≥ 36	35-32	< 32	≤ 0.03	> 0.06
	Ofloxacin	5 µg OFL.5					
	Haemophilus spp.		use Nalidixic acid			≤ 0.5	-
	N. gonorrhoeae or		use Nalidixic acid			≤ 0.12	> 0.25
	Levofloxacin	5 µg LEVOF					
c ₂)	Haemophilus spp.		use Nalidixic acid			≤ 1	-
	S. pneumoniae		≥ 17	-	< 17	≤ 2	> 2
	Streptococcus spp.		≥ 20	19-17	< 17	≤ 1	> 2
	Gatifloxacin	5 µg GATIF					
	Haemophilus spp.		use Nalidixic acid			≤ 1	-
c ₂)	S. pneumoniae / Streptococcus spp.		≥ 20	19-17	< 17	≤ 1	> 2

EUCAST-and CLSI potency NEO-SENSITABS™

Page 3 / 5 - Document: 6.2.0

Interpretation according to MIC Breakpoints of SFM (France) Using CLSI Potency Neo-Sensitabs

**Haemophilus spp., S. pneumoniae,
Streptococcus spp., N. gonorrhoeae, N.
meningitidis, Campylobacter spp. and
Anaerobes**

				Zone diameter in mm			Concentrations critiques	
NEO-SENSITABS		POTENCY	CODE	S	I	R	S	R
c ₂)	Moxifloxacin	5 µg	MOXIF					
	<i>Haemophilus</i> spp.			use Nalidixic acid			≤ 0.5	-
	<i>S. pneumoniae</i>			≥ 24	-	< 24	≤ 0.5	> 0.5
	<i>Streptococcus</i> spp.			≥ 24	23-21	< 21	≤ 0.5	> 1
	Anaerobes			≥ 21	20-18	< 18	≤ 1	> 2
f) k) i)	Nalidixan	30 µg	NAL30					
	<i>Haemophilus</i> spp.			-	-	< 21	decreased susceptibility to quinolones	
	<i>N. gonorrhoeae</i>			-	-	< 25		
	<i>N. meningitidis</i>			-	-	< 21	-	-
	<i>Campylobacter</i> spp.			≥ 20	19-15	< 15	≤ 8	> 16
	Trimethoprim+Sulfa	1.25+23.75µg	SxT25					
	<i>Haemophilus</i> spp.			≥ 24	-	-	≤ 0.5/9.5	> 1/19
	<i>S. pneumoniae</i> / <i>Streptococcus</i> spp.			≥ 19	18-16	< 16	≥ 1/19	> 2/38
	Gentamicin	10 µg*	GEN10					
	<i>Haemophilus</i> spp.			≥ 16	15-14	< 14	≤ 2	> 4
	<i>Campylobacter</i> spp.			≥ 16	15-14	< 14	≤ 2	> 4
	Tobramycin	10 µg	TOB10					
	<i>Campylobacter</i> spp.			≥ 18	17-16	< 16	≤ 2	> 4
	Spectinomycin	200 µg	SPECT					
	<i>N.gonorrhoea</i>			≥ 20	-	< 20	≤ 64	> 64
	Teicoplanin	30 µg	TPN30					
	<i>S. pneumoniae</i> / <i>Streptococcus</i> spp.			≥ 17	-	-	≤ 4	> 8
	Anaerobes			≥ 17	-	-	≤ 4	> 8
	Tigecycline							
	<i>Streptococcus</i> spp.			≥ 22	-	< 22	≤ 0.25	> 0.5
	Anaerobes			≥ 21	-	-	≤ 4	> 8
	Vancomycin	30 µg	VAN30					
	<i>S. pneumoniae</i> / <i>Streptococcus</i> spp.			≥ 17	-	-	≤ 4	> 8
	Anaerobes			≥ 17	-	-	≤ 4	> 8
m)	Metronidazole	16 µg*	MTR16					
	Anaerobes			≥ 21	-	< 21	≤ 4	> 4
	<i>Clostridium difficile</i>			≥ 26	-	< 26	≤ 4	> 4

#) Differences of potency recommended by SFM and CLSI (if any available).

*) There are no potency recommendations from CLSI so far.

**) Special potency Neo-Sensitabs for detection of resistance mechanism.

Pneumococci

- Oxacillin 1 µg is used for the detection of reduced sensitivity to Penicillin in pneumococci. Penicillin resistant isolates from the meninges must be considered resistant to Ampicillin/Amoxycillin, Amox+Clav and first and second generation cephalosporins.
- Cefotaxime and Ceftriaxone must not be tested against pneumococci by the diffusion method. A surrogate test is used instead: Ceftizoxime, Ceftizoxime detects reduced sensitivity to third generation cephalosporins. Strains sensitive to Ceftizoxime show currently MIC < 0.5 µg/ml towards Cefotaxime/Ceftriaxone (susceptible), while isolates resistant to Ceftizoxime should be tested by an MIC method.
- Erythromycin: Interpretation valid for Azithromycin and Clarithromycin.
- c₂) Screening of pneumococci for reduced sensitivity to fluoroquinolones is done using Norfloxacin 10 µg Neo-Sensitabs. If the inhibition zone is < 12 mm (or the MIC is >16 µg/ml) there is a high risk of development of resistant mutants in vivo.

Interpretation according to MIC Breakpoints of SFM (France) Using CLSI Potency Neo-Sensitabs

**Haemophilus spp., S. pneumoniae,
Streptococcus spp., N. gonorrhoeae, N.
meningitidis, Campylobacter spp. and
Anaerobes**

Streptococci

Penicillin resistant strains of Group A streptococci have not yet been recognized. Viridans streptococci isolated from blood or CSF should be tested for Penicillin or Ampicillin susceptibility using an MIC method. Group B streptococci with reduced susceptibility to Penicillin have been isolated.

- d) Oxacillin 1 µg is useful for screening for Penicillin susceptibility in viridans streptococci.
- c) Erythromycin: Interpretation valid for Azithromycin and Clarithromycin.

H. influenzae

- e) Beta-lactamase negative, Ampicillin resistant strains (BLNAR) are best detected using Ampicillin 2 µg Neo-Sensitabs. Cephalothin Neo-Sensitabs is also useful to detect BLNAR strains (zone < 17 mm). BLNAR isolates must be considered resistant to Amoxycillin, Amox+Clav, as well as first and second generation cephalosporins, no matter the size of the inhibition zone.
- f) Strains resistant to Nalidixic acid should be suspected of having reduced susceptibility to quinolones. Strains with reduced susceptibility to ciprofloxacin (MIC ≥ 0.125 µg/ml) show decreased susceptibility to all quinolones.

Meningococci

- g) Oxacillin 1 µg is used routinely for the detection of reduced sensitivity to penicillins, in meningococci (chromosomal resistance).
- h) Rifampicin: Used for prophylaxis only (not treatment).
- i) Nalidixan is useful to screen for strains with reduced susceptibility to quinolones.

Gonococci

- k) Nalidixan is useful to detect strains with reduced susceptibility to quinolones. Ciprofloxacin resistant gonococci should presumably be resistant to all quinolones. A positive beta-lactamase test predicts resistance to Penicillin, Amoxycillin/Ampicillin, Piperacillin and Ticarcillin.
- l) Oxacillin 1 µg Neo-Sensitabs is useful to detect beta lactamase negative gonococci with decreased susceptibility to Penicillin (chromosomal resistance).

Campylobacter

For *Campylobacter* spp. the absence of zone of inhibition around β-lactams, aminoglycosides, macrolides or quinolones indicates high level resistance.

Anaerobes

Vancomycin 5 µg, Kanamycin 500 µg and Colistin 10 µg Neo-Sensitabs are very useful for the identification of the most important gram negative bacilli: *B. fragilis* group are resistant to Vancomycin 5 µg, Kanamycin 500 µg and Colistin 10 µg. *Prevotella* is resistant to Kanamycin 500 µg and Vancomycin 5 µg (zone < 18 mm), while it is variable to Colistin 10 µg. *Porphyromonas* is sensitive to Vancomycin 5 µg (zone > 18 mm) and resistant to Kanamycin 500 µg and Colistin 10µg. *Fusobacterium* is sensitive to Kanamycin 500 µg and Colistin 10 µg and resistant to Vancomycin 5 µg.

For species showing slow growth it may be difficult to establish a correlation between MIC's and zone sizes. Use an MIC method.

- m) Metronidazole: Certain strains may show false resistance to Metronidazole if anaerobiosis is not correct.

Helicobacter pylori

- n) Interpretation valid for Clarithromycin.

References:

- 1) Barbut F. et al: Antimicrobial susceptibilities and serogroups of clinical strains of *Clostridium difficile* isolated in France in 1991 and 1997. Antimicrob. Ag. Chemother., **43**, 2607-11,1999.
- 2) Communiqué 2004 de la Société Française de Microbiologie (CA-SFM).

EUCAST-and CLSI potency NEO-SENSITABS™

Page 5 / 5 - Document: 6.2.0

Interpretation according to MIC Breakpoints of SFM (France) Using CLSI Potency Neo-Sensitabs

**Haemophilus spp., S. pneumoniae,
Streptococcus spp., N. gonorrhoeae, N.
meningitidis, Campylobacter spp. and
Anaerobes**

- 3) Communiqué Janvier 2005 de la Société Française de Microbiologie (CA-SFM).
- 4) Communiqué Janvier 2006 de la Société Française de Microbiologie (CA-SFM).
- 5) Communiqué Janvier 2007 de la Société Française de Microbiologie (CA-SFM).
- 6) Communiqué Janvier 2008 de la Société Française de Microbiologie (CA-SFM).
- 7) Communiqué Janvier 2009 de la Société Française de Microbiologie (CA-SFM).
- 8) Communiqué Janvier 2010 de la Société Française de Microbiologie (CA-SFM).
- 9) Communiqué Janvier 2011 de la Société Française de Microbiologie (CA-SFM).
- 10) Communiqué Janvier 2012 de la Société Française de Microbiologie (CA-SFM).

ORIGINAL ROSCO DOCUMENT

EUCAST-and CLSI potency NEO-SENSITABS™

Interpretation according to MIC
Breakpoints of BSAC (UK and
Ireland) using BSAC and CLSI
potency Neo-Sensitabs

Rapidly growing bacteria

Rapidly growing bacteria

Interpretation according to the MIC break-points adapted to recommendations in the BSAC
Standardized Disc Testing Method (Version 11), UK and Ireland (2012)

Media: Iso Sensitest, Inoculum: Semiconfluent growth (ICS) according to BSAC

				Zone diameter in mm			Break-points MIC µg/ml	
NEO-SENSITABS	POTENCY	CODE		S	I	R	S	R
Beta-lactams:								
	Penicillin	1 unit	PENG1	≥ 25	-	≤ 24		
a)	Staph.						≤ 0.12	penase
b)	Oxacillin	1 µg	OXA.1					
	Staph. aureus (MH Agar)			≥ 15	-	≤ 14	≤ 2	> 2
	Coag. neg. staph. (MH Agar)			≥ 20	-	≤ 19	≤ 0.25	> 0.25
	Ampicillin	10 µg	AMP10	≥ 15	-	≤ 14	≤ 8	> 8
	Enterococci			≥ 20	-	≤ 19	≤ 4	> 8
	Amoxycillin	25 µg	AMX25	≥ 22	-	≤ 21	≤ 8	> 8
	Amoxycillin+Clav.	20+10 µg	AMC30	≥ 21	-	≤ 20	≤ 8	> 8
	Mecillinam (U)	10 µg	MEC10	≥ 14	-	≤ 13	≤ 8	> 8
	Temocillin	30 µg	TEMOC	≥ 20	-	≤ 19	≤ 8	> 8
	Urine			≥ 15	-	≤ 14	≤ 32	> 32
d)	Piperacillin	100 µg	PIPER	≥ 24	-	≤ 23	≤ 8	> 16
	Piperacillin+Tazo.	75+10µg	PTZ85	≥ 23	22-21	≤ 20	≤ 16	> 16
	Pseudomonas spp.			≥ 25	-	≤ 24	<16	> 16
	Ticarcillin	75 µg	TIC75	≥ 21	-	≤ 20	≤ 16	> 16
	Pseudomonas spp.			≥ 20	-	≤ 19	≤ 16	> 16
k)	Ticarcillin+Clav.	75+10 µg	TIM85	≥ 23	-	≤ 22	≤ 8	> 8
	Pseudomonas spp.			≥ 20	-	≤ 19	≤ 16	> 16
o)	Cloxacillin		CLOXA	Detection of plasmid-mediated AmpC				
	Cephalothin	30 µg	CEP30	≥ 27	-	≤ 26	≤ 8	> 8
	Cephalexin	30 µg	CFLEX					
	Urine (<i>E. coli</i>)			≥ 16	-	≤ 15	≤ 16	> 16
	Cefuroxime (oral) UTI	30 µg	CXM30	≥ 20	-	≤ 19	≤ 8	> 8
	Cefuroxime (inj)	30 µg	CXM30	≥ 20	-	≤ 19	≤ 8	> 8
	Cefoxitin	10 µg	CFO10					
	<i>S. aureus</i> (Iso-Sensitest)			≥ 22	-	< 22	Oxa S	Mec A pos.
	<i>S. aureus</i> (Mueller-Hinton)			≥ 18	-	< 18	Oxa S	Mec A pos.
	coag. neg. staph. (Iso-s)			≥ 27	-	≤ 21	Oxa S	Mec A pos.
	Cefoxitin	30 µg	CFO30	≥ 23	-	-	AmpC Screen	
d) j)	Cefpodoxime	10 µg	CPD10	≥ 20	-	≤ 19	≤ 1	> 1
	(Screen ESBL)			-	-	< 24		
	Cefpodoxime+Clav.	10+1 µg	CPD+C	Screening for ESBL			Screen ESBL	
g) d)	Cefotaxime	30 µg	CTX30	≥ 30	29-24	≤ 23	≤ 1	> 2
	Cefotaxime+Clav.	30+10 µg	CTX+C	Detection of ESBL				
g) d)	Ceftriaxone	30 µg	CTR30	≥ 28	29-24	< 23	≤ 1	> 2

EUCAST-and CLSI potency NEO-SENSITABS™

Interpretation according to MIC
Breakpoints of BSAC (UK and
Ireland) using BSAC and CLSI
potency Neo-Sensitabs

Rapidly growing bacteria

			Zone diameter in mm			Break-points MIC µg/ml	
NEO-SENSITABS	POTENCY	CODE	S	I	R	S	R
k) g) d) Ceftazidime	30 µg	CAZ30	≥ 27	26-23	≤ 22	≤ 1	> 4
<i>Pseudomonas</i> spp.			≥ 24	-	≤ 23	≤ 8	> 8
Ceftazidime	10 µg	CAZ10	≥ 24	23-16	≤ 15	≤ 1	> 8
<i>P. aeruginosa</i>			≥ 20	-	≤ 19	≤ 8	> 8
<i>Acinetobacter</i>			≥ 16	-	≤ 15	≤ 8	> 8
Ceftazidime+Clav.	30+10µg	CZ+CL	Detection of ESBL				
k) d) Cefepime	30 µg	FEP30	≥ 32	31-26	< 26	≤ 1	> 4
<i>Pseudomonas</i> spp.			≥ 23	-	≤ 22	≤ 8	> 8
Cefepime+Clav.	30+10µg	CP+CL	Detection of ESBL				
k) g) d) Aztreonam	30 µg	AZT30	≥ 28	27-22	< 22	≤ 1	> 8
<i>Pseudomonas</i> spp.			≥ 36	35-20	≤ 19	≤ 1	> 16
h) l) Imipenem	10 µg	IMI10	≥ 21	20-17	≤ 16	≤ 2	> 8
<i>Ent. faecalis</i>			≥ 19	18-17	≤ 16	≤ 4	> 8
<i>Pseudomonas</i> spp.			≥ 25	24-17	≤ 16	≤ 4	> 8
<i>Acinetobacter</i> spp.			≥ 25	24-14	≤ 13	≤ 2	> 8
p) Imipenem+EDTA	10+750µg	IM10E	Detection of metallo-β-lactamases				
h) l) Meropenem	10 µg	MRP10	≥ 27	26-19	< 19	≤ 2	> 8
h) <i>Pseudomonas</i> spp.			≥ 20	19-16	≤ 15	≤ 2	> 8
<i>Acinetobacter</i> spp.			≥ 25	24-19	≤ 18	≤ 2	> 8
Ertapenem	10 µg	ETP10	≥ 28	27-16	≤ 15	≤ 0.5	> 1
Doripenem	10 µg	DOR10	≥ 24	23-19	≤ 18	≤ 1	> 4
<i>Pseudomonas</i> spp.			≥ 32	31-25	≤ 24	≤ 1	> 4
<i>Acinetobacter</i> spp.			≥ 24	23-18	≤ 17	≤ 1	> 4
<i>Enterococci</i> spp.			≥ 22	21-19	≤ 18	≤ 1	> 4
Aminoglycosides:							
Amikacin	30 µg	AMI30	≥ 19	18-16	≤ 15	≤ 8	> 16
<i>Pseudomonas/Acinetobacter</i> spp			≥ 22	21-16	≤ 15	≤ 8	> 16
Coag. neg. staph.			≥ 25	24-22	≤ 21	≤ 1	> 1
f) Gentamicin	10 µg	GEN10	≥ 20	19-16	< 16	≤ 2	> 4
<i>Pseudomonas</i>			> 18	-	≤ 17	≤ 4	> 4
<i>S.aureus; Acinetobacter</i> spp.			≥ 20	-	≤ 19	≤ 1	> 1
Coagulase neg. staph.			≥ 28	-	≤ 27	≤ 0.25	> 0.25
Gentamicin (HLR)	250 µg	GN250	-	-	< 17	-	HLR
Streptomycin	10 µg (enterob)	STR10	≥ 13	-	< 12	≤ 8	> 8
f) Tobramycin	10 µg	TOB10	≥ 21	20-18	≤ 17	≤ 2	> 4
<i>Pseudomonas/Acinetobacter</i> spp.			≥ 20	-	≤ 19	≤ 4	> 4
<i>Staph. aureus</i>			≥ 21	-	≤ 20	≤ 1	> 1
Coagulase neg. staph.			≥ 30	-	≤ 29	≤ 0.25	> 0.25
Other:							
Erythromycin	15 µg	ERY15	≥ 20	19-17	≤ 16	≤ 1	> 2
Azithromycin	15 µg	AZI15	≥ 20	-	≤ 19	≤ 1	> 2
<i>Salmonella typhi</i>			≥ 19	-	≤ 18	-	-
Clarithromycin	15 µg	CLA15	≥ 20	19-17	≤ 16	≤ 1	> 2
Clindamycin	2 µg	CLI.2	≥ 26	25-23	≤ 22	≤ 0.25	> 0.5
Quinupristin	15 µg	SYN15	≥ 22	21-19	<19	≤ 2	> 2
/Dalfopristin							
<i>Enterococcus</i> (<i>E. faecium</i>)			≥ 20	19-12	≤11	≤ 1	> 4
<i>Staphylococci</i>			≥ 22	21-19	≤18	≤ 1	> 2

EUCAST-and CLSI potency NEO-SENSITABS™

Interpretation according to MIC
Breakpoints of BSAC (UK and
Ireland) using BSAC and CLSI
potency Neo-Sensitabs

Rapidly growing bacteria

NEO-SENSITABS			Zone diameter in mm			Break-points MIC µg/ml			
	POTENCY	CODE	S	I	R	S	R		
n) q)	Linezolid	10 µg	LIZ10	≥ 20	-	≤ 19	≤ 4	> 4	
	Telithromycin	15 µg	TEL15	≥ 27	-	≤ 26	≤ 0.5	> 0.5	
	Fucidin	10 µg	FUC10	≥ 30	-	≤ 29	≤ 1	> 1	
	Doxycycline	30 µg	DOX30	≥ 29	-	≤ 28	≤ 1	> 1	
	Minocycline	30 µg	MIN30						
	Staphylococcus spp.			≥ 28	-	≤ 27	≤ 0.5	> 0.5	
	Tetracyclines	30 µg	TET30	≥ 29	-	≤ 28	≤ 1	> 1	
	Chloramphenicol	30 µg	CLR30	≥ 21	-	≤ 20	≤ 8	> 8	
	Rifampicin (staph)	5 µg	RIF.5	≥ 32	31-26	≤ 25	≤ 0.06	> 0.5	
	Teicoplanin	30 µg	TPN30						
	2+18 h prediffusion (MH agar,McF.0.5)								
	Staphylococcus spp.			≥ 22	-	< 20	≤ 2	> 2 (VISA/ GISA)	
	Enterococcus spp.			≥ 20	-	≤ 19	≤ 2	> 2 (VRE) - -	
	c) n) q) q) i) e) e) q) q)	Vancomycin	30 µg	VAN30					
		2+18 h prediffusion (MH agar, McF 0.5)							
Staphylococcus spp.				≥ 24	-	≤ 22	≤ 2	>2 (VISA/GISA)	
Enterococcus spp.				≥ 16	15-12	< 12	≤ 4	VRE	
Daptomycin		30 µg	DAPCa						
2+18 hours' prediffusion									
Staphylococcus spp.				≥ 22	21-20	< 20	≤ 1	> 2	
Enterococcus spp.				≥ 12	-	≤ 11	≤ 4	> 8	
Nitrofurantoin (U)		200 µg	NI200	≥ 17	-	≤ 16	≤ 64	> 64	
Nalidixic acid (U)		30 µg	NAL30	≥ 18	-	≤ 17	≤ 16	> 16	
Enterobacteriaceae/Acinetobacter spp.				≥ 20	-	≤ 19	Reduced susceptibility to quinolones		
Salmonella spp.				≥ 20	-	≤ 19	Reduced susceptibility to quinolones		
Norfloxacin (syst.)		2 µg	NOR2	≥ 26	25-19	≤ 18	≤ 0.5	> 1	
Norfloxacin (UTI)		2 µg	NOR2	≥ 16	-	≤ 15	≤ 4	> 4	
Ciprofloxacin		1 µg	CIPR1	≥ 20	19-17	≤ 16	≤ 0.5	> 1	
Staphylococcus spp.			≥ 14	-	≤ 13	≤ 1	> 1		
e) e) e) q) q) e)	Pseudomonas spp.			≥ 21	-	≤ 20	≤ 1	> 1	
	Urine			≥ 20	-	≤ 19	≤ 4	> 4	
	Acinetobacter spp.			≥ 21	-	≤ 20	≤ 1	> 1	
	Ciprofloxacin	5 µg	CIPR5						
	Pseudomonas spp.			≥ 30	29-20	≤ 19	≤ 0.5	> 1	
	Ofloxacin	5 µg	OFL.5	≥ 29	28-25	< 25	≤ 0.5	> 1	
	Levofloxacin	5 µg	LEVOF	≥ 22	21-19	≤ 18	≤ 2	> 2	
	Pseudomonas			≥ 22	21-17	≤ 16	≤ 1	> 2	
	Colistin	10 µg	Co.10						
	2+18 hours' prediffusion			≥ 15	14-11	≤ 10	≤ 2	> 2	
	Gatifloxacin	5 µg	GATIF	≥ 28	27-24	≤ 23	≤ 0.5	> 1	
	Moxifloxacin	5 µg	MOXIF	≥ 28	27-24	≤ 23	≤ 0.5	> 1	

EUCAST-and CLSI potency NEO-SENSITABS™

Interpretation according to MIC Breakpoints of BSAC (UK and Ireland) using BSAC and CLSI potency Neo-Sensitabs

Rapidly growing bacteria

			Zone diameter in mm			Break-points MIC µg/ml	
NEO-SENSITABS	POTENCY	CODE	S	I	R	S	R
			≥ 31	30-25	≤ 24	≤ 0.5	> 1
			≥ 20	19-16	≤ 15	≤ 0.5	> 1
	Sulphonamides (U)	240 µg SULFA	≥ 28	-	≤ 27	≤ 32	> 32
	Trimethoprim	2.5 µg TP2.5	≥ 20	19-15	≤ 14	≤ 0.5	> 4
	Enterococci (U)		≥ 33	32-26	≤ 25	≤ 0.03	> 1
	<i>Staph.</i>		≥ 20	-	≤ 19	≤ 1	> 1
	Urine		≥ 17	16-14	≤ 13	≤ 2	> 4
m)	Trimethoprim	1.25+23.75µg SxT25	≥ 16	-	≤ 15	≤ 2	> 4
	+ Sulfa						
	Staphylococci		≥ 17	16-14	≤ 13	≤ 2	> 4
	<i>S. maltophilia</i>		≥ 20	-	≤ 19	≤ 4	> 4
	Fosfomycin (U)	200+50 µg FO200	≥ 25	-	≤ 24	≤ 32	> 32
	<i>Staph.</i>		≥ 34	-	≤ 33	≤ 32	> 32
	Novobiocin (U)	5 µg NOVO5	≥ 16	-	≤ 15	-	-
	Mupirocin (staph)	10 µg MUPIR	≥ 18	17-10 (LLR)	≤ 9 (HLR)	≤ 4	> 128
r)	Tigecycline	15 µg TIG15					
	<i>Enterobacteriaceae/Acinetobacter spp.</i>		≥ 24	23-20	≤ 19	≤ 1	> 2
	<i>Staphylococcus spp.</i>		≥ 26	-	≤ 25	≤ 0.5	> 0.5
	<i>Enterococcus spp.</i>		≥ 21	-	≤ 20	≤ 0.25	> 0.5
Special tests							
o)	Boronic Acid	BORON	Detection of AmpC and KPC beta lactamases				
p)	Dipicolinic Acid	D.P.A.	Detection of metallo-β-lactamases				

Rapidly growing bacteria includes *Enterobacteriaceae.*, *Pseudomonas spp.*, *Acinetobacter spp.*, *Staphylococcus spp.*, *Enterococcus spp.*

Remarks:

- Staphylococci resistant to Penicillin should be reported as resistant to Amoxycillin, Ampicillin, Piperacillin and Ticarcillin.
- Iso-sensitest will not reliably detect resistance of staphylococci towards Oxacillin (test Cefoxitin). Use Mueller-Hinton Agar (BSAC 2001). Staphylococci resistant to Cefoxitin (Oxacillin) should be reported resistant to **all** other beta-lactams including beta-lactamase inhibitor combinations and carbapenems. Cefoxitin testing on Iso-sensitest agar should be performed at a temperature not exceeding 35°C to avoid false susceptibility results (Skov et al. (7))
- When testing enterococci against Vancomycin 5 µg, plates should be incubated for full 24 hours and examined carefully for the presence of a haze or other growth within the zone (indicates resistance).
- Strains of *Klebsiella*, *E. coli* and *Salmonella* that produce ESBL may be clinically resistant to therapy with Penicillin, cephalosporins or Aztreonam, despite apparent in vitro susceptibility. Read in Neo-Sensitabs User's Guide 20th Ed. 2009 on ESBL screening and confirmatory tests.

Interpretation according to MIC Breakpoints of BSAC (UK and Ireland) using BSAC and CLSI potency Neo-Sensitabs

Rapidly growing bacteria

- e) Enterobacteriaceae resistant to Nalidixic acid (zone < 20 mm) show a **decreased** susceptibility to quinolones MIC CIPRO ≥ 0.125 µg/ml). Therefore Nalidixic acid is a good screening for the detection of decreased quinolone susceptibility in *Salmonella* spp. Strains of *E. coli*, Klebsiella and Salmonella with Ciprofloxacin MIC's of 0.25 and 0.5 µg/ml may be reported as resistant to the quinolones. According to Wareham et al. (4) *E. coli* and *K. pneumonia* isolates with plasmid-mediated quinolone resistance had a zone diameter of ≤ 16 mm around Ciprofloxacin 1 µg disks.
- f) Staphylococci that are **resistant** to Gentamicin should be reported as resistant against both Netilmicin and Tobramycin (enzymes APH (2'') + AAC (6°)). Staphylococci resistant to Kanamycin should be reported as resistant to Amikacin.
- h) According to BSAC zones for Meropenem and *P. aeruginosa* should be S ≥ 20 mm, I: 19-16, R ≤ 15 mm. Several carbapenemase producing isolates will be falsely reported as susceptible to Meropenem.
- i) Klebsiella, Enterobacter and *Proteus* spp. should be reported **R** to Nitrofurantoin.
- j) *E. coli*/Klebsiella/Salmonella strains showing zone < 24 mm with Cefpodoxime Neo-Sensitabs should be suspected of producing ESBL. For confirmatory tests use Ceftazidime+Clavulanate and Cefepime+Clavulanate compared to Ceftazidime and Cefepime Neo-Sensitabs (see ESBL).
- k) Synergism between Ticarcillin+Clavulanate and Aztreonam/Ceftazidime/Cefepime permit the detection of ESBL producing strains (*Ps. aeruginosa*).
- l) Carbapenem testing on Iso-sensitest Agar may give falsely susceptibility for isolates that harbour metallo-beta-lactamases, testing on Mueller-Hinton Agar is to be preferred (BSAC 2001) (2)
- m) For *S. maltophilia* and *B. cepacia* use incubation at 30 °C. (3,6) For *B. cepacia* complex should be developed specific criteria.
- n) For detection of VISA, GISA, hVISA strains, see chapter in Neo-Sensitabs User's Guide 20th Edition 2009 on detection of staphylococci with decreased susceptibility to Vancomycin.
- o) Cloxacillin and Phenylboronic acid are used for the detection of plasmid-mediated AmpC beta-lactamases. See description of the procedure in Neo-Sensitabs User's Guide 20th Ed. 2009. For further information on methodology, see leaflet "Screening and detection of AmpC beta lactamases" (www.rosco.dk).
- p) Imipenem+EDTA and Dipicolinic acid are used for the detection of metallo-β-lactamases. See further information in User's Guide "Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™ 2010". For further information on methodology, see leaflet "Screening and detection of carbapenemases" (www.rosco.dk).
- q) Special technique for susceptibility testing of high molecular weight antimicrobials (Colistin, Daptomycin, Teicoplanin, Vancomycin): 2+18 hours' prediffusion method. Description can be found in document **1.5.0**.
- r) Disk diffusion valid only for *E. coli*. For Acinetobacter use MIC method. According to Hope et al. (5) the best advice for Tigecycline testing is to perform broth dilution tests on isolates with borderline results, especially in severe infections and when Tigecycline use is intended.

References:

- 1) BSAC Standardized Disc Testing Method (last update Version 11) (2012).

EUCAST-and CLSI potency NEO-SENSITABS™

Page 6 / 6 - Document: 7.1.0

Interpretation according to MIC Breakpoints of BSAC (UK and Ireland) using BSAC and CLSI potency Neo-Sensitabs

Rapidly growing bacteria

- 2) Wootton M. et al: Examining the effect of media on carbapenem resistance interpretation in *P. aeruginosa*. Abstract presentation D-730, 49th ICAAC, 2009.
- 3) Pitman K. et al: Comparison of methods for susceptibility testing of clinically important *Burkholderia* spp. Abstract presentation D-1443, 49th ICAAC, 2009.
- 4) Amin A., Wareham D.: Ability of Ciprofloxacin 1 µg Discs and BSAC Breakpoints to Detect Plasmid-Mediated Quinolone Resistance in Urinary Enterobacterial isolates. Abstract P586, 20th ECCMID, Vienna, April 2010.
- 5) Hope R. et al: Tigecycline activity: low resistance rates but problematic disc breakpoints revealed by a multicentre sentinel survey in the UK. *J. Antimicrob. Chemother.* 65, 2602-2609, 2010.
- 6) Pitman K et al: Disc susceptibility testing for *Burkholderia cepacia* complex. Presentation D-752, 50th ICAAC, Sept. 2010.
- 7) Skov R et al: Effects of temperature on the detection of methicillin resistance in *S. aureus* using Cefoxitin disc diffusion testing with Iso-sensitest agar. *J. Antimicrob. Chemother.* **63**, 699-703, 2009.

EUCAST-and CLSI potency NEO-SENSITABS™

Page 1 / 5 - Document: 7.2.0

Interpretation according to MIC
Breakpoints of BSAC (UK and
Ireland) using BSAC and CLSI
potency Neo-Sensitabs

Haemophilus spp. *S. pneumoniae*,
Streptococcus spp., *N. gonorrhoeae*, *N.*
meningitidis, *Moraxella catarrhalis*,
Coryneforms, *Campylobacter* spp.,
Pasteurella multocida and anaerobes

Fastidious organisms

Interpretation according to the MIC break-points adapted to recommendations by the BSAC
(Version 11) in UK and Ireland.

Inoculum, media and incubation conditions according to BSAC (2012).

				Zone diameter in mm			Break-points MIC µg/ml	
NEO-SENSITABS	POTENCY	CODE		S	I	R	S	R
e)	Ampicillin 10 µg	AMP10	<i>Haemophilus</i> spp.	Use Amp 2.5 µg			≤ 1	> 1
			<i>Moraxella catarrhalis</i>	≥ 30	-	≤ 29	≤ 1	beta lact.pos
			<i>Pasteurella multocida</i>	≥ 30	-	≤ 29	≤ 1	> 1
e)	Ampicillin 2 µg	AMP.2	<i>Haemophilus</i> spp.	≥ 18	-	≤ 17	≤ 1	> 1
			<i>Streptococci</i>	≥ 24	23-15	≤ 14	≤ 0.5	> 2
			<i>Moraxella catarrhalis</i>	≥ 30	-	≤ 29	≤ 1	> 1
			<i>Pasteurella multocida</i>	≥ 20	-	≤ 19	≤ 1	> 1
			<i>N. meningitidis</i>	≥ 32	-	≤ 31	-	-
a)	Penicillin 1 unit	PENG1	<i>Streptococci</i> haemol.	≥ 20	-	≤ 19	≤ 0.25	> 0.25
			α-haem strep.	≥ 17	16-11	≤ 10	≤ 0.25	> 2
			<i>N. gonorrhoeae</i>	≥ 26	25-18	≤ 17	≤ 0.06	> 0.12
			<i>N. meningitidis</i>	≥ 29	28-15	≤ 14	≤ 0.06	> 0.25
			<i>Coryneforms</i>	≥ 20	-	≤ 19	≥ 0.12	> 0.12
			<i>Moraxella catarrhalis</i>	BL neg.	-	BL pos.	S	R
			<i>Pasteurella multocida</i>	≥ 22	-	≤ 21	≤ 0.12	beta lact.pos
			Anaerobes (<i>Clostridium</i>)	≥ 23	-	≤ 22	≤ 0.25	> 0.5
a)	Oxacillin 1 µg	OXA.1	<i>S. pneumoniae</i>	≥ 20	19-11	≤ 10	≤ 0.06 (pen)	2(MICtest)
d)			<i>Streptococcus</i> spp.	≥ 17	-	≤ 16	≤ 0.12 (pen)	MIC test
l)			<i>N. gonorrhoeae</i>	≥ 15	-	≤ 14	≤ 0.06 (pen)	MIC test
g)			<i>N. meningitidis</i>	≥ 12	-	≤ 11	≤ 0.06 (pen)	MIC test
	Amoxycillin+Clav. 2+1 µg	AMC.3	<i>Haemophilus</i> spp.	≥ 14	-	≤ 13	≤ 2	> 2
			<i>Moraxella catarrhalis</i>	≥ 19	-	≤ 19	≤ 1	> 1
	Amoxycillin+Clav. 20+10 µg	AMC30	<i>Haemophilus</i> spp.	≥ 29	-	≤ 28	≤ 1	> 1
			<i>Moraxella catarrhalis</i>	≥ 29	-	≤ 28	≤ 1	> 1
			Anaerobes	≥ 29	28-21	≤ 20	≤ 4	> 8
			<i>Clostridium</i>	≥ 32	-	≤ 31	≤ 4	> 8
	Piperacillin+Tazobactam 75+10µg	PTZ85	Anaerobes	≥ 27	-	≤ 26	≤ 16	> 16
	Cefaclor 30 µg	CCL30	<i>Haemophilus</i> spp.	≥ 28	-	≤ 27	≤ 0.5	> 0.5
			<i>Moraxella catarrhalis</i>	≥ 38	-	≤ 37	≤ 0.12	> 0.12
	Cefuroxime 30 µg	CXM30						

EUCAST-and CLSI potency NEO-SENSITABS™

Page 2 / 5 - Document: 7.2.0

Interpretation according to MIC
Breakpoints of BSAC (UK and
Ireland) using BSAC and CLSI
potency Neo-Sensitabs

**Haemophilus spp. S. pneumoniae,
Streptococcus spp., N. gonorrhoeae, N.
meningitidis, Moraxella catarrhalis,
Coryneforms, Campylobacter spp.,
Pasteurella multocida and anaerobes**

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
			≥ 32	-	≤ 31	≤ 1	> 2
			≥ 38	-	≤ 37	≤ 1	> 1
			≥ 24	-	≤ 23	≤ 4	> 8
Cefixime	5 µg	CFM.5	≥ 35	-	≤ 34	≤ 0.12	> 0.12
N.gonorrhoeae							
Cefotaxime	5 µg	CTX.5					
			≥ 25	-	≤ 24	≤ 0.12	> 0.12
Haemophilus spp.			≥ 34	-	≤ 33	≥ 1	> 1
Pasteurella multocida			≥ 40	-	≤ 39	≤ 0.12	> 0.12
N. meningitidis			≥ 30	-	≤ 29	≤ 0.12	> 0.12
N. gonorrhoeae			≥ 23	-	≤ 22	≤ 0.5	> 0.5
Streptococcus spp.							
b) S. pneumoniae					use ceftizoxime	-	-
Ceftriaxone	30 µg	CTR30					
			≥ 38	-	≤ 37	≤ 0.12	> 0.12
b) Haemophilus spp.					use ceftizoxime	-	-
S. pneumoniae							
N. gonorrhoeae			≥ 38	-	≤ 37	≤ 0.12	> 0.12
N. meningitidis			≥ 38	-	≤ 37	≤ 0.12	> 0.12
Streptococcus spp.			≥ 35	-	≤ 34	≤ 0.5	> 0.5
Ceftizoxime	30 µg	ZOX30					
b) S. pneumoniae			≥ 30	29-24	≤ 23	≤ 1	MIC
(valid for 3rd gen. cepha.)							
Imipenem	10 µg	IMI10					
			≥ 25	-	≤ 24	≤ 2	> 2
S. pneumoniae			≥ 23	-	≤ 22	≤ 2	> 2
Haemophilus spp.							
Meropenem	10 µg	MRP10					
			≥ 28	-	≤ 27	≤ 2	> 2
S. pneumoniae			≥ 32	31-28	<28	≤ 0.25	> 1
S. pneumoniae (meningitidis)			≥ 23	-	≤ 22	≤ 2	> 2
Haemophilus spp.			≥ 34	33-30	<30	≤ 0.25	> 1
Haemophilus spp.(meningitidis)			≥ 26	25-19	≤ 18	≤ 2	> 8
Anaerobes *							
Ertapenem	10 µg	ETP10					
			≥ 33	-	≤ 32	≤ 0.5	> 0.5
S. pneumoniae			≥ 35	-	≤ 34	≤ 0.5	> 0.5
Streptococci			≥ 35	-	≤ 34	≤ 0.5	> 0.5
M. catarrhalis			≥ 33	-	≤ 32	≤ 0.5	> 0.5
Haemophilus spp.							
Doripenem	10 µg	DOR10					
			≥ 23	-	≤ 22	≤ 1	> 1
Streptococci							
S. pneumoniae							
Erythromycin	15 µg	ERY15					
c) S. pneumoniae/Streptococcus spp.			≥ 24	23-22	≤ 21	≤ 0.25	> 0.5
Moraxella cat			≥ 28	-	≤ 27	≤ 0.25	> 0.5
Campylobacter			≥ 22	-	≤ 21	≤ 4	> 4
Haemophilus spp.			≥ 32	31-16	≤ 15	≤ 0.25	> 32
Azithromycin	15 µg	AZI15					
			≥ 35	34-20	≤ 19	≤ 0.12	> 4
Haemophilus spp.			≥ 28	-	≤ 27	≤ 0.25	> 0.5
m) N. gonorrhoeae			≥ 22	21-20	≤ 19	≤ 0.25	> 0.5
Streptococci, S. pneumoniae							
Clarithromycin	15 µg	CLA15					
			≥ 28	27-26	≤ 25	≤ 0.25	> 0.5
Moraxella catarrhalis			≥ 24	23-22	≤ 21	≤ 0.25	> 0.5
Streptococci, S. pneumoniae			≥ 20	-	≤ 19	≤ 1	> 1
Campylobacter spp.							

EUCAST-and CLSI potency NEO-SENSITABS™

Page 3 / 5 - Document: 7.2.0

Interpretation according to MIC
Breakpoints of BSAC (UK and
Ireland) using BSAC and CLSI
potency Neo-Sensitabs

Haemophilus spp. S. pneumoniae,
Streptococcus spp., N. gonorrhoeae, N.
meningitidis, Moraxella catarrhalis,
Coryneforms, Campylobacter spp.,
Pasteurella multocida and anaerobes

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
Haemophilus spp.			≥ 24	23-10	no zone	≤ 1	> 32
Clindamycin	2 µg	CLI.2					
S. pneumoniae/streptococci			≥ 20	-	≤ 19	≤ 0.5	> 0.5
Anaerobes			≥ 10	-	no zone	≤ 4	> 4
Quinupristin	15 µg	SYN15					
/Dalfopristin							
S. pneumoniae			≥ 20	-	≤ 19	≤ 2	> 2
Linezolid	10 µg	LIZ10					
Moraxella catarrhalis			≥ 20	-	≤ 19	≤ 2	> 4
S. pneumoniae/Streptococcus spp.			≥ 20	-	≤ 19	≤ 4	> 4
Telithromycin	15 µg	TEL15					
S. pneumoniae			≥ 29	-	≤ 28	≤ 0.25	> 0.5
Streptococcus spp.			≥ 29	-	≤ 28	≤ 0.25	> 0.5
M. catarrhalis			≥ 30	-	≤ 29	≤ 0.25	> 0.5
Haemophilus spp.			≥ 31	30-16	≤ 15	≤ 0.12	> 8
Tetracyclines	30 µg	TET30					
S. pneumoniae/Streptococcus spp.			≥ 26	-	≤ 25	≤ 1	> 2
Campylobacter spp.			≥ 22	21-19	≤ 18	≤ 2	> 2
N. gonorrhoeae			≥ 36	35-30	≤ 29	≤ 0.5	> 1
P. multocida			≥ 30	-	≤ 29	≤ 1	> 1
Haemophilus			≥ 26	25-22	≤ 21	≤ 1	> 2
Tigecycline	15 µg	TIG15					
Streptococcus spp.			≥ 25	24-20	≤ 19	≤ 0.25	> 0.5
Haemophilus spp.			≥ 28	-	≤ 27	≤ 0.25	> 0.25
Chloramphenicol	30 µg	CLR30					
S. pneumoniae			≥ 21	-	≤ 20	≤ 8	> 8
N. meningitidis			≥ 26	-	≤ 25	≤ 2	> 4
Haemophilus spp.			≥ 30	29-27	≤ 26	≤ 2	> 2
Moraxella catarrhalis			≥ 32	-	≤ 31	≤ 1	> 2
Rifampicin	5 µg	RIF.5					
S. pneumoniae			≥ 23	22-21	≤ 20	≤ 0.06	> 0.5
N. meningitidis			≥ 34	-	≤ 33	≤ 0.25	> 0.25
h) Vancomycin	30 µg	VAN30					
S. pneumoniae/streptococci			≥ 18	-	≤ 17	≤ 4	> 4
Coryneforms			≥ 20	-	≤ 19	≤ 4	> 8
Cl. difficile			-	-	-	≤ 2	> 2
f) Nalidixic acid (U)	30 µg	NAL30					
Haemophilus spp.			≥ 25	-	no zone	decreased susceptibility to quinolones	
Pasteurella multocida			≥ 28	-	≤ 27		
k) N. gonorrhoeae			≥ 32	31-10	no zone		
Moraxella catarrhalis (screen)			≥ 18	-	≤ 17	decreased susceptibility to quinolones	
Campylobacter spp.			≥ 20	-	≤ 19		
Norfloxacin	2 µg	NOR2					
S. pneumoniae			-	-	< 10	decreased susceptibility to quinolones	

EUCAST-and CLSI potency NEO-SENSITABS™

Page 4 / 5 - Document: 7.2.0

Interpretation according to MIC
Breakpoints of BSAC (UK and
Ireland) using BSAC and CLSI
potency Neo-Sensitabs

**Haemophilus spp. S. pneumoniae,
Streptococcus spp., N. gonorrhoeae, N.
meningitidis, Moraxella catarrhalis,
Coryneforms, Campylobacter spp.,
Pasteurella multocida and anaerobes**

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
f)	Ciprofloxacin 1 µg	CIPR1					
	<i>Haemophilus</i> spp.		≥ 28	-	≤ 27	≤ 0.5	> 0.5
	<i>Campylobacter</i> spp.		≥ 26	-	≤ 25	≤ 1	> 1
	<i>Pasteurella multocida</i>		≥ 29	-	≤ 28	≤ 1	> 1
	<i>S. pneumoniae</i>		≥ 25	24-10	≤ 9	≤ 0.12	> 2
	Coryneforms		≥ 17	16-12	≤ 11	≤ 0.5	> 1
	<i>N. meningitidis</i>		≥ 32	-	≤ 31	≤ 0.03	> 0.06
	<i>N. gonorrhoeae</i>		≥ 29	-	≤ 28	≤ 0.03	> 0.06
	<i>Moraxella catarrhalis</i>		≥ 18	-	≤ 17	≤ 0.5	> 0.5
	Ofloxacin 5 µg	OFL.5					
	<i>S. pneumoniae</i>		≥ 28	27-16	≤ 15	≤ 0.12	> 4
	<i>M. catarrhalis</i> /Haemophilus spp.		≥ 35	-	≤ 34	≤ 0.5	> 0.5
	Levofloxacin 5 µg	LEVOF					
	<i>S. pneumoniae</i> /Streptococcus spp.		≥ 18	-	≤ 17	≤ 2	> 2
	<i>M. catarrhalis</i> / Haemophilus spp.		≥ 28	-	≤ 27	≤ 1	> 1
	Moxifloxacin 5 µg	MOXIF					
	<i>S. pneumoniae</i> /Streptococcus spp.		≥ 22	-	≤ 21	≤ 0.5	> 0.5
	<i>Moraxella catarrhalis</i>		≥ 30	-	≤ 29	≤ 0.5	> 0.5
	<i>Haemophilus</i> spp.		≥ 30	-	≤ 29	≤ 0.5	> 0.5
	<i>Cl. difficile</i>		-	-	-	-	> 4
	Nitrofurantoin 200 µg	NI200					
	Streptococci		≥ 19	-	≤ 18	≤ 64	> 64
	Trimethoprim 2.5 µg	TP2.5					
	Haemophilus spp.		≥ 21	-	≤ 20	≤ 0.5	> 0.5
	Trimethoprim+ Sulfa 1.25+23.75µg	SxT25					
	<i>Haemophilus</i> spp.		≥ 21	20-18	≤ 17	≤ 0.5	> 1
	Streptococci, <i>S. pneumoniae</i>		≥ 20	19-17	≤ 16	≤ 1	> 2
	<i>Moraxella catarrhalis</i>		≥ 17	-	≤ 16	≤ 0.5	> 1
	Spectinomycin 200 µg	SPECT					
	<i>N. gonorrhoeae</i>		≥ 23	22-20	≤ 19	≤ 64	> 64
	Metronidazole 5 µg	MTR.5					
	Anaerobes		≥ 18	-	≤ 17	≤ 4	> 4
	<i>Cl. difficile</i>		-	-	-	≤ 2	> 2
m)	Metronidazole 16 µg	MTR16					
	<i>Campylobacter</i> spp.		≥ 28	-	≤ 27	≤ 4	> 4
	Anaerobes		≥ 24	-	≤ 23	≤ 4	> 4

**B. fragilis* showing inhibition zones < 30 mm show increased expression of *cfiA* and should be reported as resistant to meropenem (1).

Remarks:

Pneumococci

- a) Oxacillin 1 µg is used for the detection of reduced sensitivity to Penicillin in pneumococci. Penicillin resistant isolates from the meninges must be considered resistant to Ampicillin/Amoxycillin, Amox+Clav and first and second generation cephalosporins.

Interpretation according to MIC Breakpoints of BSAC (UK and Ireland) using BSAC and CLSI potency Neo-Sensitabs

**Haemophilus spp. S. pneumoniae,
Streptococcus spp., N. gonorrhoeae, N.
meningitidis, Moraxella catarrhalis,
Coryneforms, Campylobacter spp.,
Pasteurella multocida and anaerobes**

- b) Cefotaxime and Ceftriaxone must not be tested against pneumococci by the diffusion method. A surrogate test is used instead: Ceftizoxime. Ceftizoxime detects reduced sensitivity to third generation cephalosporins. Strains sensitive to Ceftizoxime show currently MIC < 0.5 µg/ml towards Cefotaxime/Ceftriaxone (susceptible), while isolates resistant to Ceftizoxime should be tested by an MIC method.
- c) Erythromycin: Interpretation valid for Azithromycin and Clarithromycin.

Streptococci

Penicillin resistant strains of Group A streptococci have not yet been recognized.

Viridans streptococci isolated from blood or CSF should be tested for Penicillin or Ampicillin susceptibility using an MIC method. Group B streptococci with reduced susceptibility to Penicillin have been isolated.

- d) Oxacillin 1 µg is useful for screening for Penicillin susceptibility in streptococci.
- c) Erythromycin: Interpretation valid for Azithromycin and Clarithromycin.

H. influenzae

- e) Beta-lactamase negative, ampicillin resistant strains (BLNAR) are best detected using Ampicillin 2.5 µg Neo-Sensitabs. BLNAR isolates must be considered resistant to Amoxycillin, Amox+Clav, as well as first and second generation cephalosporins, no matter the size of the inhibition zone.
- f) Strains resistant to nalidixic acid should be suspected of having reduced susceptibility to quinolones. Strains with reduced susceptibility to Ciprofloxacin (MIC ≥ 0.125 µg/ml) show decreased susceptibility to all quinolones.

Meningococci

- g) Oxacillin 1 µg is used routinely for the screening of reduced sensitivity to penicillins, in meningococci (chromosomal resistance).
- h) Rifampicin: Used for prophylaxis only (not treatment).

Gonococci

- k) Nalidixic acid is useful to detect strains with reduced susceptibility to quinolones. Ciprofloxacin resistant gonococci should presumably be resistant to all quinolones. A positive beta-lactamase test predicts resistance to Penicillin, Amoxycillin/Ampicillin, Piperacillin and Ticarcillin.
- l) Oxacillin (1 µg) Neo-Sensitabs are useful to detect beta-lactamase negative gonococci with decreased susceptibility to Penicillin (chromosomal resistance).
- m) Patients treated for Chlamydia with Azithromycin cannot be assumed to have been adequately treated for gonococcal infection. UK Department of Health recommends that Azithromycin should not be used to treat gonorrhoea.

Campylobacter

For *Campylobacter* spp. the absence of zone of inhibition around β-lactams, aminoglycosides, macrolides or quinolones indicates high level resistance.

Anaerobes

Interpretation is valid for *Bacteroides* spp and *Clostridium difficile*.

For species showing slow growth it may be difficult to establish a correlation between MIC's and zone sizes. Use an MIC method.

- m) Metronidazole: Certain strains may show false resistance to Metronidazole if anaerobiosis is not correct.

References:

- 1) Ang Lei et al: Carbapenem resistance in *Bacteroides fragilis*. J.A.C., **59**, 1043-44, 2007.

EUCAST-and CLSI potency NEO-SENSITABS™

Interpretation according to MIC
breakpoints of BSAC (UK and
Ireland) using BSAC potency Neo-
Sensitabs

BSAC urines

BSAC urines

MIC and zone diameter breakpoints for isolates from URINES, according to BSAC 2012

Neo-Sensitabs potency according to BSAC. Iso-sensitest Agar. Inoculum according to BSAC.

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
<u>Amoxycillin</u>	25 µg	AMX25					
Enterobacteriaceae			≥ 22	-	≤ 21	≤ 8	> 8
Enterococci			≥ 20	-	≤ 19		
S. saprophyticus			≥ 26	-	≤ 25		
Group B streptococci			≥ 26	-	≤ 25		
<u>Amoxycillin+Clav.</u>	20+10 µg	AMC30					
Enterobacteriaceae			≥ 13	-	≤ 12	≤ 32	> 32
Enterococci			≥ 21	-	≤ 20		
S. saprophyticus			≥ 28	-	≤ 27		
Group B streptococci			≥ 28	-	≤ 27		
<u>Ampicillin</u>	10 µg	AMP10					
Enterobacteriaceae			≥ 15	-	≤ 14	≤ 8	> 8
<u>Cefepime</u>	30 µg	FEP30					
Enterobacteriaceae			≥ 32	31-26	< 26	≤ 1	> 8
Pseudomonas			≥ 23	-	≤ 22	≤ 8	> 8
<u>Cefepime+Clav.</u>	30+10 µg	FEP+C			Detection ESBL		
<u>Cefotaxime</u>	30 µg	CTX30	≥ 30	29-24	≤ 23	≤ 1	> 2
<u>Cefotaxime+Clav</u>	30+10 µg	CTX+C			Detection ESBL		
<u>Cefpodoxime</u>	30 µg	CPD10	≥ 20	-	≤ 19	≤ 1	Screen ESBL
<u>Cefpodoxime+Clav.</u>	10+1 µg	CPD+C			Screening ESBL		
<u>Ceftazidime</u>	30 µg	CAZ30					
Enterobacteriaceae			≥ 30	29-17	< 17	≤ 1	> 8
Pseudomonas			≥ 26	-	≤ 25	≤ 8	> 8
<u>Cefazidime+Clav.</u>	30+10 µg	CAZ+C			Detection EBSL		
<u>Cefuroxime</u>	30 µg	CXM30	≥ 20	-	≤ 19	≤ 8	> 8
<u>Cephalexin</u>	30 µg	CFLEX					
E.coli/Klebsiella			≥ 16	-	≤ 15	≤ 32	> 32
P. mirabilis			≥ 12	-	≤ 11		
Group B streptococci			≥ 24	-	≤ 23		
<u>Ciprofloxacin</u>	1 µg	CIPR1					
Enterobacteriaceae			≥ 20	19-17	≤ 16	≤ 0.5	> 1
Enterococci			≥ 12	-	≤ 11		
S. saprophyticus			≥ 18	-	≤ 17		
Group B streptococci			≥ 13	-	≤ 12		
<u>Ertapenem</u>	10 µg	ETP10	≥ 28	-	≤ 27	≤ 0.5	> 1
<u>Fosfomycin</u>	200+50 µg	FO200					
E.coli			≥ 25	-	≤ 24	≤ 32	> 32
P. mirabilis			≥ 37	-	≤ 36		
Enterococci			≥ 20	-	≤ 19		
S. saprophyticus			≥ 20	-	≤ 19		
<u>Gentamicin</u>	10 µg	GEN10	≥ 20	19-17	≤ 16	≤ 2	> 4
<u>Imipenem</u>	10 µg	IMI10	≥ 21	20-17	≤ 16	≤ 2	> 8
<u>Mecillinam</u>	10 µg	MEC10					
E.coli/Klebsiella/P.mirabilis			≥ 14	-	≤ 13	≤ 8	> 8
<u>Mecillinam</u>	33 µg	MECIL					
S. saprophyticus			≥ 10	-	no zone	≤ 64	> 64
<u>Meropenem</u>	10 µg	MRP10	≥ 27	26-20	≤ 19	≤ 2	> 8
<u>Nalidixic acid</u>	30 µg	NAL30					
Enterobacteriaceae			≥ 18	-	≤ 17	≤ 16	> 16
<u>Nitrofurantoin</u>	200 µg	NI200					
E.coli			≥ 17	-	≤ 16	≤ 64	> 64

EUCAST-and CLSI potency NEO-SENSITABS™

Interpretation according to MIC
breakpoints of BSAC (UK and
Ireland) using BSAC potency Neo-
Sensitabs

BSAC urines

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm			Break-points MIC µg/ml	
			S	I	R	S	R
Enterococci			≥ 20	-	≤ 19		
S. saprophyticus			≥ 20	-	≤ 19		
Group B streptococci			≥ 19	-	≤ 18	≤ 64	> 64
Norfloxacin	2 µg	NOR.2	≥ 16	-	≤ 15	≤ 4	> 4
Piperacillin+Tazo	(75+10) µg	TIM85	≥ 22	-	< 21	≤ 16	> 16
Temocillin	30 µg	TEMOC					
Enterobacteriaceae			≥ 12	-	≤ 11	≤ 32	> 32
Tetracyclines	30 µg	TET30	≥ 29	-	≤ 28	≤ 1	> 1
Tigecycline	15 µg	TIG15	≥ 24	23-20	≤ 19	≤ 1	> 2
Trimethoprim	2.5 µg	TP2.5					
Enterobacteriaceae			≥ 17	16-14	≤ 13	≤ 2	≥ 4
Enterococci			≥ 50	49-22	≤ 21	≤ 0.03	> 1
S. saprophyticus			≥ 15	14-13	≤ 12	≤ 2	≥ 4
Group B streptococci			≥ 16	-	≤ 15	≤ 2	≥ 4
Trimethoprim+Sulfa	1.25+23.75 µgSxT25		≥ 16	-	< 15	≤ 2	> 2
Special tests							
a)	Cloxacillin	CLOXA	Detection of plasmid mediated AmpC				
a)	Boronic acid	BORON	Detection of AmpC and KPC beta lactamases				
b)	Dipicolinic acid	D.P.A	Detection of metallo-β-lactamases				
b)	Imipenem+EDTA	10+750 µg IM10E	Detection of metallo-β-lactamases				

- a) Description of the methodology can be found in a leaflet "Screening and detection of AmpC beta lactamases" and user's guide "Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™ 2010" (www.rosco.dk).
- b) Description of the methodology can be found in a leaflet "Screening and detection of carbapenemases" It also includes KPC and metallo-β-lactamases (www.rosco.dk).

**Interpretation of the AntibioGramme with Neo-Sensitabs
MIC break-points according to CLSI (M31-A4/M31-S2 2013)
Inoculum according to Kirby-Bauer / confluent colonies.**

Mueller-Hinton agar and McFarland 0.5 inoculum.

NEO-SENSITABS			Zone diameter in mm			Break-points MIC µg/ml	
POTENCY CODE			S	I	R	S	R
	Amoxicillin+Clav.	20+10 µg AMC30	≥ 30	29-28	≤ 27	≤ 0.25	≥ 1 (dogs, cats)
	<i>Staphylococcus</i> spp.			Use cefoxitin		≤ 4/2	≥ 8/4
	Other organisms (UTI only)		≥ 18	17-14	≤ 13	≤ 8/4	≥ 32/16
	Amikacin 30 µg	AMI30	≥ 17	16-15	≤ 14	≤ 16	≥ 32
	Amoxycillin 30 µg	AMOXY					
	<i>Enterococcus</i> spp.		≥ 20	-	≤ 19	≤ 8	≥ 16
	Other organisms		≥ 20	19-17	≤ 16	≤ 8	≥ 32
	Ampicillin 10 µg	AMP10					
	Enterobacteriaceae		≥ 17	16-14	≤ 13	≤ 8	≥ 32
	<i>Staphylococcus</i> spp.			Use penicillin			
	<i>Enterococcus</i> spp.		≥ 17	-	≤ 16	≤ 8	≥ 16
	<i>Streptococcus</i> spp. (not <i>S. pneumoniae</i>)		≥ 24	23-17	≤ 16	≤ 0.25	≥ 8
	<i>Listeria monocytogenes</i>		≥ 20	-	≤ 19	≤ 2	≥ 4
	<i>Mannheimia haemolytica</i>		≥ 26	25-20	≤ 19	≤ 0.25	≥ 1
c)	Apramycin 40 µg	APRAM	≥ 23	22-20	≤ 19	≤ 4	≥ 16
c)e) h)	Cefadroxil 30 µg	CDX30	≥ 23	22-20	≤ 19	≤ 8	≥ 32
e)h)	Cefazolin 30 µg	CFZ30	≥ 23	22-20	≤ 19	≤ 2	≥ 8
f)	Cefoxitin 30 µg	CFO30					
	<i>S. aureus</i>		≥ 22	-	≤ 21	OxaS	MecA pos
	Coag. neg. staph.		≥ 25	-	≤ 24	OxaS	MecA pos
a)c) h)	Ceftriaxone 30 µg	CTR30	≥ 23	22-20	≤ 19	≤ 1	≥ 4
	(Cefoperazone, Cefovecin)						
h)	Cefpodoxime (dogs)	10 µg CPD10	≥ 21	20-18	≤ 17	≤ 2	≥ 8
c)h)	Cefquinome 30 µg	CFQUI	≥ 23	22-20	≤ 19	≤ 2	≥ 8
h)	Ceftiofur 30 µg	CFTIO	≥ 21	20-18	≤ 17	≤ 2	≥ 8
c)h)	Cefuroxime 30 µg	CXM30	≥ 18	17-15	≤ 14	≤ 8	≥ 32
c)e) h)	Cephalexin 30 µg	CFLEX	≥ 20	19-17	≤ 16	≤ 8	≥ 32
c)e) h)	Cephalothin 30 µg	CEP30	≥ 23	22-20	≤ 19	≤ 2	≥ 8
	Cephapirin 30 µg	CFP30	≥ 18	17-15	≤ 14	≤ 8	≥ 32
	Chloramphenicol 30 µg	CLR30					
	<i>S. pneumoniae</i>		≥ 21	-	≤ 20	≤ 4	≥ 8
	<i>Streptococcus</i> spp.		≥ 21	20-18	≤ 17	≤ 4	≥ 16
	Other organisms		≥ 18	17-13	≤ 12	≤ 8	≥ 32
j)	Ciprofloxacin 1 µg	CIPR1					
	<i>Salmonella</i> spp.		≥ 30	-	≤ 30	≤ 0.06	≥ 0.12
i)	Clindamycin 2 µg	CLIN2					
	<i>Staphylococcus</i> spp.		≥ 21	20-15	≤ 14	≤ 0.5	≥ 4
	Colistin 10 µg	CO.10	≥ 15	-	≤ 14	≤ 2	≥ 2
c)	Doxycycline 30 µg	DOX30	≥ 14	13-11	≤ 10	≤ 4	≥ 16
	Enrofloxacin 5 µg	ENR.5	≥ 23	22-19	≤ 18	≤ 0.25	≥ 1
	Enrofloxacin 10 µg	ENROF	≥ 25	24-21	≤ 20	≤ 0.25	≥ 1
i)	Erythromycin 15 µg	ERY15					
	<i>Streptococcus</i> spp.		≥ 21	30-16	≤ 15	≤ 0.25	≥ 1
	<i>Staph./Enterococci</i>		≥ 23	22-14	≤ 13	≤ 0.5	≥ 8
	Florfenicol 30 µg	FFC30					
	Cattle		≥ 19	18-15	≤ 14	≤ 2	≥ 8
	Swine		≥ 22	21-19	≤ 18	≤ 2	≥ 8
	Fish pathogens		≥ 28	27-24	≤ 23	≤ 2	≥ 8
c)	Flumequine 30 µg	FLUME	≥ 20	19-17	≤ 16	≤ 2	≥ 4
	Fosfomycin (U) 200 µg	FO200	≥ 18	-	≤ 18	≤ 32	> 32

EUCAST-and CLSI potency NEO-SENSITABS™

Page 2 / 8 - Document: 9.1.0

Veterinary practice according to
CLSI breakpoints

Veterinary practice

NEO-SENSITABS				Zone diameter in mm			Break-points MIC µg/ml	
	POTENCY	CODE		S	I	R	S	R
c)	Fucidin	100 µg	FUCID	≥ 28	27-24	≤ 23	≤ 1	≥ 4
	Plain agar			≥ 26	25-23	≤ 22	≤ 1	≥ 4
	Blood agar			≥ 21	-	≤ 20	≤ 1	≥ 1
c)	Furazolidone	50 µg	FURAZ	≥ 23	22-20	≤ 19	≤ 4	≥ 8
	Gentamicin	10 µg	GEN10	≥ 16	15-13	≤ 12	≤ 2	≥ 8
l)	Gentamicin	250 µg	GN250	-	-	< 14	-	HLR
	Enterococci			≥ 18	17-14	≤ 13	≤ 16	≥ 64
c) i)	Kanamycin	30µg	KAN30	≥ 20	19-17	≤ 16	≤ 2/4	≥ 5/16
	Lincomycin/Neomycin	15+60 µg	LIN+N					
	(Albionic Forte)							
i)	Linco-spectin	15+200µg	LI+SP	≥ 20	19-17	≤ 16	≤ 4/32	≥ 16/64
	Imipenem	10 µg	IMI10	≥ 16	15-14	≤ 13	≤ 4	≥ 16
c)	Marbofloxacin	5 µg	MAR.5	≥ 20	19-15	≤ 14	≤ 1	≥ 4
b) c)	Metronidazole	16 µg	MTR16	≥ 28	27-24	≤ 23	≤ 4	≥ 8
	(anaerobes)							
c)	Naf-Pen-Strep	5+2+20 µg	N+P+S	≥ 20	19-17	≤ 16	≤ 1/1/4	≥ 2/2/16
c)	Neomycin	120 µg	NEOMY	≥ 23	22-20	≤ 19	≤ 6	≥ 25
	Neomycin/Oxacillin(Cloxacillin)	60/5 µg		≥ 20	19-17	≤ 16	≤ 4/2	≥ 16/4
c)	Nitrofurantoin	300 µg	NI300	≥ 17	16-15	≤ 14	≤ 32	≥ 128
	Nitrofurantoin	100 µg	NI100	≥ 15	14-13	≤ 12	≤ 32	≥ 128
c)	Novobiocin	5 µg	NOVO5					
	Blood agar			≥ 13	12-11	≤ 10	≤ 2	≥ 4
	Plain agar			≥ 16	15-14	≤ 13	≤ 2	≥ 4
f)	Oxacillin	1 µg	OXA.1					
	Coag. neg. staph.			≥ 18	-	≤ 17	≤ 0.25	≥ 0.5
	S. intermedius (pseud)			-	-	-	-	-
	S. pneumoniae (penicillin)			≥ 20	-	-	≤ 0.06	≥ 2
c)	Oxolinic acid	10 µg	OXOLI	≥ 20	19-17	≤ 16	≤ 4	≥ 8
	Penicillin/Novo	10 U + 30 µg	PEN+N					
	Mastitis			≥ 18	17-15	≤ 14	≤ 1/2	≥ 4/8
	Other organisms			≥ 17	-	≤ 16	≤ 1/2	≥ 4/8
	Penicillin	10 units	PEN10					
k)	Staphylococcus spp.			≥ 29	-	≤ 28	≤ 0.12	≥ 0.25
	Streptococcus viridans			≥ 26	25-13	≤ 12	≤ 0.12	≥ 4
	Beta haemolytic			≥ 24	-	-	≤ 0.12	-
	Enterococci			≥ 15	-	≤ 14	≤ 8	≥ 16
	Other organisms			≥ 18	17-11	≤ 10	≤ 1	≥ 4
	Pirlimycin	10 µg	PIRLI	≥ 18	-	≤ 17	≤ 2	≥ 4
	Rifampicin	5 µg	RIF5					
	S. pneumoniae			≥ 20	19-17	≤ 16	≤ 1	≥ 4
	Enterococci			≥ 20	19-17	≤ 16	≤ 1	≥ 4
	Spectinomycin	200 µg	SPECT	≥ 18	17-15	≤ 14	≤ 32	≥ 128
	(Pasteur., Haemoph.)							
	Spiramycin	200 µg	SPIRA	≥ 26	25-23	≤ 22	≤ 2	≥ 8
c)	Streptomycins	10 µg	STR10	≥ 15	14-12	≤ 11	≤ 6	≥ 25
	Streptomycin (ent. cocci)	500 µg	ST500	-	-	< 14	-	HLR
	Sulphonamides (U)	240 µg	SULFA	≥ 17	16-13	≤ 12	≤ 256	≥ 512
	Tetracyclines	30 µg	TET30					
	Staphylococci			≥ 19	18-15	≤ 14	≤ 4	≥ 16
	Cattle (BRD)			≥ 23	22-19	≤ 18	≤ 2	≥ 8
	Swine (BRD)			≥ 27	26-24	≤ 23	≤ 0.5	≥ 2
	S. pneumoniae			≥ 15	14-12	≤ 11	≤ 4	≥ 16
	Tiamulin	30 µg	TIAMU					
	Spirochaetae			≥ 28	27-24	≤ 23	≤ 1	≥ 4
	Actinobacillus			≥ 11	-	no zone	≤ 16	≥ 32
	Streptococci			≥ 25	24-19	< 18	≤ 1	≥ 4
c)	Tiamulin+Tetra	30+15 µg	-	≥ 16	15-14	≤ 13	≤ 8/4	≥ 16/8
	Ticarcillin	75 µg	TIC75					
	Ps. aeruginosa			≥ 24	23-16	≤ 15	≤ 16	≥ 128
	Gram neg. ent.			≥ 20	19-15	≤ 14	≤ 16	≥ 128

NEO-SENSITABS				Zone diameter in mm			Break-points MIC µg/ml	
POTENCY CODE				S	I	R	S	R
c) g)	Ticarcillin+Clavulanate	75+10 µg	TIM85	≥ 24	23-15	≤ 14	≤ 16/2	≥ 128/2
	<i>Ps. aeruginosa</i>							
	Tilmicosin	80 µg	TILMI	≥ 18	17-15	≤ 14	≤ 8	≥ 32
	Bovine RD			≥ 15	-	≤ 14	≤ 16	≥ 32
	Swine RD			≥ 16	15-11	≤ 10	≤ 4	≥ 16
c) g)	Trimethoprim	5 µg	TRIM5	≥ 19	18-16	≤ 15	≤ 0.5/9.5	≥ 4/76
	Trimethoprim+Sulfa	1.25+23.75µ	SxT25	≥ 16	15-11	≤ 10	≤ 0.5/9.5	≥ 4/76
	<i>S. pneumonia, Haemophilus</i>			≥ 16	15	≤ 14	≤ 2/38	≥ 4/76
	Systemic infection			≥ 26	25-23	≤ 22	≤ 4	≥ 16
	Urine			≥ 17	16-15	≤ 14	≤ 4	≥ 32
c)	Tylosin	150 µg	TYLOS					
	Vancomycin	30 µg	VAN30					
<u>Special tests</u>								
Cloxacillin				Detection of plasmid mediated AmpC				
Phenylboronic acid				Detection of AmpC and KPC beta lactamases				
Imipenem+EDTA				Detection of metallo-β-lactamases				
Dipicolinic acid				Detection of metallo-β-lactamases				
<u>Kits</u>								
ESBL, AmpC, ESBL + AmpC Confirm kit				ESBL + AmpC				
ESBL + AmpC Screen Kit				ESBL + AmpC				
KPC, MBL, OXA-48 Confirm Kit				Carbapenemases				

Remarks:

- Results with Ceftriaxone are valid for Cefoperazone and Cefovecin.
- Metronidazole 16 µg is the representative of the Nitroimidazole-group, including Ronidazole, Ornidazole, Iprnidazole, and Moxnidazole. Results obtained with Metronidazole are applicable to the others.
- MIC break-points have not yet been given by the CLSI.
- From August 2005, the FDA no longer allows the use of Enrofloxacin for treating infections in poultry (to avoid development of resistance in *Campylobacter* spp.).
- Results of Cephalothin susceptibility tests are used to predict susceptibility to the first generation cephalosporins, such as Cephadroxil and Cephalexin.
- Results of Cefoxitin and Oxacillin with staphylococci are used to predict susceptibility to Cloxacillin. Cefoxitin resistant staphylococci should be reported as resistant to all beta-lactams. In case of discordant results between Cefoxitin and Oxacillin, report the strains as resistant (R). Use Cefoxitin for testing *S. aureus*. For non *S. aureus* staphylococci use Oxacillin 1 µg: R ≤ 17 mm (MIC ≥ 0.5 µg/ml)
- The results of Trimethoprim+Sulfa can be used to predict the susceptibility of other potentiated sulphomanides with Trimethoprim.

- h) For detection of ESBL (CTX-M) and AmpC beta-lactamases (CMY) in *Salmonella* spp. see user's guide "**Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™**" (www.rosco.dk) on "Detection of Beta-Lactamases" (7,9).
- i) Routine screening of Clindamycin inducible resistance in staphylococci/streptococci should be performed (double disk/induction test. Results are also valid for Lincomycin.
- j) *Salmonella* spp. Resistant to Ciprofloxacin 1 µg should be reported as resistant to other fluoroquinolones (Enrofloxacin, Marbofloxacin etc.).
- k) With *S. aureus* and Penicillin 10 U: Sharp zone edge (cliff) indicates beta-lactamase positive. Fussy zone edge (beach) indicates beta-lactamase negative.
- l) High level Aminoglycoside resistance if zone is < 14 mm.

Staphylococci in animals

- a) Screening of beta-lactamase production in *S. aureus*: Look at the inhibition zone around Penicillin 10 U Neo-Sensitabs. If the zone edge is sharp (cliff) the isolate produces beta lactamase. If the zone edge is fuzzy (beach) the isolate is beta-lactamase negative. This test is only valid for *Staphylococcus aureus*.
- b) Screening for Methicillin (Oxacillin) resistance: Test *S. aureus* using both Cefoxitin 30 µg and Oxacillin 1 µg Neo-Sensitabs. In case of discordant results between Cefoxitin and Oxacillin, report the strain as resistant. For non- *S. aureus* staphylococci use Oxacillin 1 µg Neo-Sensitabs. R ≤ 17 mm (MIC ≥ 0.5 µg/ml).
- c) Test for Clindamycin induction: Place Clindamycin 2 µg and Erythromycin 15 µg Neo-sensitabs 20 mm (edge to edge) apart from each other (for staphylococci). For streptococci the distance should be 12 mm. Flattening of the Clindamycin zone adjacent to Erythromycin 15 µg indicates inducible Clindamycin resistance. Result is reported as Clindamycin R.
- d) MRSA and similar:

Bemis et al. (11) found that the Cefoxitin disk diffusion test had low sensitivity for detection of Oxacillin resistance in members of the *S. intermedius* (*S. intermedius*, *S. schleiferi subsp schleiferi* and *S. schleiferi subsp coagulans*). Oxacillin disk diffusion had a high sensitivity and specificity for detecting *mecA* mediated Oxacillin resistance.

Perreten et al. (12) describe the spread of Methicillin resistant *S. pseudointermedius* in Europe and North America. Besides MR, the isolates showed resistance to Trimethoprim 90.3%, Gentamicin 88.3%, Streptomycin 90.3%, Macrolides 89.3%, Fluoroquinolones 87.4%.

High occurrence of MRSA and Methicillin resistant *S. pseudointermedius* (15,16) was found in horses and small animals respectively.

Petersen et al (24) mention that the recent discovery of human and bovine MRSA isolates carrying a new *mecA* gene (now designated *mecC*) caused concern, because they are not detected by conventional confirmatory tests for MRSA. They conclude, that that *mecC* carrying MRSA can be exchanged between humans and ruminants.

Walther et al (25) found that the new MRSA variant (*mec C*) is not restricted to ruminants or humans, but it was also found in companion animals (dogs and cats). The isolates were resistant to 6 µg/ml Cefoxitin, but were not detected by the *mecA* test.

Vanderhaegen et al (26) isolated MR non-aureus staphylococci from bovines. They were *S. sciuri*, *S. lentus* and *S. fleuretti*. Both *mecA* and *mecC* genes were detected. Cefoxitin MIC was a poor indicator

for mecA mediated resistance in non-aureus staphylococci from animals. Both Oxacillin and Cefoxitin should be tested. Resistance to one or both indicates Methicillin resistance.

Enterococci

Detection of High Level resistance:

Use Gentamicin 250 µg Gentamicin and Streptomycin 500 µg Neo-Sensitabs.

Enterococci showing zones of inhibition < 14 mm around Gentamicin 250 µg Neo-sensitabs should be reported as High Level Resistant to Gentamicin (no synergism with penicillins).

Enterococci showing zones of inhibition < 14 mm with Streptomycin 500 µg Neo-Sensitabs should be reported as High level resistant to Streptomycin (no synergism with penicillins).

ESBL, AmpC and carbapenemases

Vo et al. (13) found that 17 Ceftiofur resistant isolates from horses (4 *E. coli* and 3 *K. pneumoniae*) were multidrug resistant. 5 produced ESBL and 1 produced AmpC with integrons in 6 isolates.

Hopkins et al. (14) report a marked increase of *Salmonella enterica* serovar resistant to Ampicillin, Streptomycin, Sulfonamides and Tetracycline, has been noted in food-borne infections, and in pigs, pig meat in several European countries. To prevent a global epidemic of these newly emerging clones or strains, as occurred with *S. tiphimurium* DT104, intervention strategies are needed as soon as possible.

CTX-M , ESBL enzymes have been found in egg, bovine mastitis, raw chicken, broiler chickens and turkeys in the UK and Europe (17,18,19,20).

Madec et al (21) found multidrug resistant *Salmonella enterica* serovar tiphimurium DT104 in cattle in France. They were Ceftiofur resistant and Cefoxitin susceptible, showing a typical Cephalosporin/Clavulanate synergy, indicating ESBL production.

Leverstein et al (22) conclude that intestinal carriage with ESBL-producing bacteria in food-producing animals and contamination of retail meat may contribute to increased incidences of infections with ESBL-producing bacteria in humans. Transmission of ESBL genes from poultry to humans takes place most likely through the food chain.

Fischer et al (23) isolated an *E. coli* producing VIM-1 (MBL) carbapenemase on a pig farm and conclude that the prevalence of carbapenemases in bacteria from livestock is probably underestimated.

Dierikx et al (27) found a high prevalence (6.8%) of ESBL/AmpC producing *E. coli* at Dutel broiler farms and a high prevalence of ESBL/AmpC in *E. coli* from farmers. Most common genes were CTX-M1, SHV-12 (ESBLs) and CMY-2 (AmpC).

Sunde et al (28) reports the first isolate (*E. coli* 1248) of animal origin (broiler) detected in Norway, with reduced susceptibility to cephalosporins. The ESBL was detected using Ceftazidime and Cefepime Neo-sensitabs with and without Clavulanate. The detection was unexpected and, probably, it has been part of the bacterial flora of animals taken to Norway for breeding purposes.

Cottell et al (29) detected ESBL, CTX-M32 in 29 of 88 steers over a 26 day period. Besides, CTX-M positive bacteria were found in feces in greater numbers than previously reported in the United States. The authors fear that the CTX-M genes may spread among animals.

Detection of ESBL, AmpC and carbapenemases in animals

Screening:

Enterobacteriaceae showing the following zones of inhibition:

Cefpodoxime 10 µg Neo-Sensitabs : ≤ 17 mm and/or

Cefotaxime 30 µg or Cefquinome 30 µg Neo-Sensitabs : ≤ 27 mm and/or

Ceftazidime 30 µg and/or Ceftiofur 30 µg Neo-Sensitabs : ≤ 22 mm

Should be suspected of possessing an **ESBL** and should be tested to confirm it.

Isolates producing zones of inhibition ≤ 16 around Cefoxitin 30 µg Neo-Sensitabs should be suspected of possessing an **AmpC** beta-lactamase and should be tested to confirm it.

Use the ESBL + AmpC Screen kit or the Total ESBL + AmpC Confirm Kit from Rosco Diagnostica.

Isolates producing zones of inhibition < 23 mm around Imipenem 10 µg Neo-Sensitabs should be suspected of possessing a **Carbapenemase** and should be tested to confirm it. Use the KPC, MBL OXA-48 Confirm Kit from Rosco Diagnostica.

Quality Control of Veterinary Antibiotics

			Zone diameter in mm			
NEO-SENSITABS	POTENCY	CODE	<i>E. coli</i> ATCC 25922	<i>S. aureus</i> ATCC 25923	<i>Ps.</i> <i>aeruginosa</i> ATCC 27853	<i>S.</i> <i>pneumoniae</i> ATCC 29212
Apramycin	40 µg	APRAM	21-28	22-30	18-24	-
Cefquinome	30 µg	CFQUI	28-36	25-33	-	30-38
Ceftiofur	30 µg	CFTIO	26-31	27-31	14-18	30-36
Enrofloxacin	5 µg	ENR.5	32-40	27-31	15-19	-
Florfenicol	30 µg	FFC30	22-28	22-29	-	24-31
Marbofloxacin	5 µg	MAR.5	29-37	24-30	20-25	-
Penicillin/Novo	10/30 µg	PEN+N	-	30-36	-	24-30
Pirlimycin	10 µg	PIRLI	-	25-32	-	-
Tiamulin	30 µg	TIAMU	-	22-27	-	-

Acceptable Q.C ranges for *Campylobacter jejuni* ATCC 33560.

Incubation at 36-37° C for 24 hours:

NEO-SENSITABS	POTENCY	CODE	Zone diameter in mm
Nalidixan	30 µg	NAL30	25-34
Ciprofloxacin	5 µg	CIPR5	32-45
Erythromycin	15 µg	ERY15	26-38

References:

- 1) M37-A2 Development of in vitro susceptibility testing criteria and quality control parameters for veterinary antimicrobial agents; approved guideline 2002.
- 2) M31-A4 Performance standards for antimicrobial disk and dilution susceptibility tests for bacteria isolated from animals; draft, 2013, and M31-S2 (2nd Informational Supplement 2013).
- 3) Jones R. N. et al: Tiamulin activity against fastidious and non-fastidious veterinary and human bacterial isolates: initial development of in vitro susceptibility test methods. J. Clin. Microbiol., **40**, 461-5, 2002.
- 4) Petersen A et al: Harmonization of antimicrobial susceptibility testing among veterinary diagnostic laboratories in the five nordic countries. Microbial Drug Resistance, **9**, 381-388, 2003.
- 5) Gray J.T. et al: Antibiotic susceptibility testing of bacteria isolated from animals. Clin. Microbiol. Newsletter, **27**, 131-5, 2005.
- 6) Rich M. et al: Clindamycin-resistance in MRSA isolated from animals. Veterinary Microbiology, **111**, 237-40, 2005.
- 7) Xian-Zhi Li: Beta lactam resistance and beta-lactamases in bacteria of animal origin. Vet. Microbiol., **121**, 197-214, 2007.
- 8) Lüthje P. et al: Molecular basis of resistance to macrolides and lincosamines among staphylococci and streptococci from various animal sources collected in the resistance monitoring program B of T-Germ. Vet. Int. J. Antimicrob. Agents, **29**, 528-535, 2007.
- 9) Smet A. et al: Diversity of ESBL and Class C beta-lactamases among cloacal E. coli isolates in Belgian broiler farms. Antimicrob. Ag. Chemother **52**, 1238-43, 2008.
- 10) Hunter P.A. et al: Antimicrobial-resistant pathogens in animals and man: prescribing, practices and policies. J. Antimicrob. Chemother. **65**, suppl1, 3-17, 2010.
- 11) Bemis et al.: Comparison of Test to Detect Oxacillin Resistance in *S. intermedius*, *S. schleiferi* and *S. aureus* isolates from canine hosts. J. Clin. Microbiol., **44**, 3374-3376, 2006.
- 12) Perreter V. et al.: Clonal spread of MR *S. pseudointermedius* in Europe and North America: An International Multicenter Study. J. Antimicrob. Chemother., **65**, 1145-1154, 2010.
- 13) Vo ATT et al.: Characteristics of Extended-Spectrum Cephalosporin Resistant *E. coli* and *K. pneumonia* isolates from horses. Vet. Microbiol., **124**, 248-255, 2007.
- 14) Hopkins K. L. et al.: Multiresistant *Salmonella enterica* serovar 4, in Europe: a new pandemic strain?. Eurosurveillance, **15**, 3, June, 2010
- 15) Van den Eede A et al: High occurrence of MRSA ST398 in equine nasal samples. Accepted manuscript. Vet. Microbiology 2008.
- 16) Ruscher C et al: Widespread rapid emergence of distinct Methicillin and multidrug-resistant *S. pseudointermedius* (MRSP) genetic linkage in Europe. Vet. Microbiol. **144**, 340-346, 2010.
- 17) Bortolaia V et al: E. coli producing CTX-M-1,-2 and -9 group beta-lactamases in organic chicken egg production. Antimicrob. Ag. Chemother. **54**, 3527-3528, 2010.
- 18) Locatelli C et al: CTX-M1 ESBL-producing *K. pneumonia* isolated from cases of bovine mastitis. J. Clin. Microbiol. **48**, 3822-3823, 2010.
- 19) Dhanji H et al: Cephalosporin resistance mechanisms in *E. coli* isolated from raw chicken imported into the UK. J. Antimicrob. Chemother. **65**, 2534-2537, 2010.
- 20) Randall LP et al: Prevalence of E.coli carrying ESBLs (CTX-M and TEM-52) from broiler chickens and turkeys in Great Britain between 2006 and 2009. J. Antimicrob. Chemother. **66**, 86-95, 2011.
- 21) Madec JY et al: ESBL's CTX-M1 gene carried on an IncI1 plasmid in multidrug resistant *S. enterica* serovar Typhimurium DT104 in cattle in France. J. Antimicrob. Chemother. **66**, 942-943, 2011.
- 22) Leverstein MA et al: Dutch patients, retail chicken, meat and poultry share the same ESBL genes, plasmids and strains. Clin. Microbiol, Infect. Febr. 2011 (accepted article)
- 23) Fischer J et al: E. coli producing VIM-1 carbapenemase isolated on a pig farm. J. Antimicrob. Chemother. **67**, 1793-1795, 2012.
- 24) Petersen A et al: Epidemiology of MRSA carrying the novel mecC gene in Denmark, corroborates a zoonotic reservoir with transmission to humans. Clin. Microbiol. Infect. **19**, E16-E22, 2013.
- 25) Walther B. et al: MRSA variant in companion animals. Emerg. Infect. Dis. **18**, 2017-2019, 2012.
- 26) Vanderhaegen W. et al: Characterization of MR non-*S. aureus* staphylococci carriage isolates from different bovine populations. J. Antimicrob. Chemother. **68**, 300-307, 2013.
- 27) Dierckx C et al: ESBL and AmpC producing *E. coli* in Dutch broilers and broiler farmers. J. Antimicrob. Chemother. **68**, 60-67, 2013.

- 28) Sunde M et al: E. coli of animal origin in Norway contains a TEM-20 carrying plasmid closely related to TEM-20 and TEM-52 plasmids from other European countries. J. Antimicrob. Chemother. **63**, 215-216, 2009.
- 29) Cottell et al: CTX-M32 an IncN plasmid in E. coli from beef cattle in the United States. Antimicrob. Ag. Chemother. **57**, 1096-1097, 2013.

ORIGINAL ROSCO DOCUMENT